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Basically, DNA extraction in potato usually can obtain  $100\text{ng} / \mu\text{L}$

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## Testing for Foods Produced by Recombinant DNA Techniques

### 1. Sampling Methods

#### 1.1 Sampling of Foods Produced by Recombinant DNA Techniques

##### 1.1.1 Sampling of Corn Kernels and Soybeans

On the assumption that corn kernels or soybeans produced by recombinant DNA techniques exist inhomogeneously, obtain a representative sample from a lot in accordance with lot size, type of packing (bulk or sacked) and packaging style of the lot. Take notice that corn kernels or soybeans from other lots do not commingle when you obtain a sample. To obtain a sample, use utensils, containers and packaging materials that are disposable or cleaned thoroughly.

Mix a sample well until it becomes homogenous, and weigh out certain amount of grains\* that you need for tests. Then, grind them until they become homogeneous with equipments (e.g. grinding machine).

\*Approximately 500g of corn kernels or soybeans are needed for a quantitative test.

##### 1.1.1.1. Sampling Sacked Grains

Obtain samples according to the following table.

Lot Size	Number of Sacks to be Unpacked	Amount of Sample (Kg)	Number of Samples	
			Qualitative	Quantitative
□ 15	2	1	1	3
16 ~ 25	3	1	1	3
26 ~ 90	5	1	1	3
91 ~ 150	8	1	1	3

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151 ~ 280	13	1	1	3
281 ~ 500	20	1	1	3
501 ~ 1,200	32	1	1	3
1,201 ~ 3,200	50	1	1	3
3,201 ~ 10,000	80	1	1	3
10,001 ~ 35,000	125	1	1	3
35,001 ~ 150,000	200	1	1	3
150,001 ~ 500,000	315	1	1	3
□500,001	500	1	1	3

### 1.1.1.2. Sampling Bulk Grains

#### 1.1.1.2.1 Sampling Silos During Loading

When foods are loaded into silos, regard a silo as a lot and obtain a representative sample of a lot with tools (e.g. auto-sampler). Draw corn kernels or soybeans for 15 times at appropriate intervals in order that the total amount weighs more than 10 Kg. Then, divide the grains homogeneously and use a portion (more than 1 Kg) as a sample per silo.

Regarding corn kernels or soybeans that have already been loaded into silos, obtain a sample in the same way when they are reloaded into another silo.

#### 1.1.1.2.2 Sampling Barges During Loading

When foods are loaded into barges (including domestic vessels), regard a barge as a lot and obtain a representative sample of a lot with tools (e.g. auto-sampler). Draw grains for 15 times at appropriate intervals in order that the total amount weighs more than 10 Kg. Then, divide the grains homogeneously and use a portion (more than 1 Kg) as a sample per barge.

#### 1.1.1.2.3 Sampling Barges After Loading

When foods have already been loaded into barges, regard a barge as a lot and obtain a representative sample of a lot by drawing corn kernels or soybeans from 15 points (every 5 points from the high, middle, and lower layer of the barge) in order that the total amount

---

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weighs more than 10 Kg. Then, divide the grains homogeneously and use a portion (more than 1 Kg) as a sample per barge.

### 1.1.2. Sampling of Papaya

Regarding sampling of papayas, obtain a sample according to the following table.

Lot Size	Number of Packages to be Unpacked	Amount of Sample (Number of Pieces)	Number of Sample
□ 50	2	2	1
51 ~ 500	3	3	1
501 ~ 35,000	5	5	1
□ 35,001	8	8	1

### 1.2 Sampling of Processed Foods

Regarding sampling of processed foods, obtain a sample according to the following table.

#### 1.2.1. Ground Corn and Soybean (e.g. corn grits, corn flour, corn meal)

Obtain samples according to the table in "1.1.1.1. Sampling Sacked Grains".

#### 1.2.2. Processed Foods other than 1.2.1.

Obtain samples according to the following table.

Lot Size	Number of Packages to be Unpacked	Amount of Sample (g)	Number of Samples	
			Qualitative	Quantitative
□ 15	2	120	1	3
16 ~ 50	3	120	1	3
51 ~ 150	5	120	1	3
151 ~ 500	8	120	1	3
501 ~ 3,200	13	120	1	3
3,201 ~ 35,000	20	120	1	3
35,001 ~ 500,000	32	120	1	3
□ 500,001	50	120	1	3

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2. Tests for foods produced by recombinant-DNA techniques whose safety assessments are not reviewed

2.1 Test methods

2.1.1 Test of corn (CBH351)

The lateral flow method is used in the test of corn kernels. This method is also used in the test of ground food products such as corn grits, corn flour and corn meal (hereinafter referred to as "corn semi-products"), in which proteins newly expressed by genetic modification are not subjected to physical and chemical changes.

The qualitative PCR method is used for other corn processed food products.


Corn semi-products are tested by the lateral flow method and then subjected to the identification test using the qualitative PCR method.

2.1.1.1 Detection of CBH351 in corn grains

2.1.1.1.1 Lateral flow method

The commercially available kit, Trait  Bt9 Corn Grain Test Kit (Part 7000003) manufactured by Strategic Diagnostics Inc.(SDI), is used. The procedures described below are basically the same as directed in the kit's explanatory leaflet. When an experiment is performed in a laboratory, RO water purified using a reverse osmotic membrane or distilled water should be used unless otherwise specified.

2.1.1.1.1.1 Testing procedures

Grind kernels randomly selected from a sample, and place the ground corn\* in an open-mouthed jar with a lid (a jar holds approximately 500 mL). Add 288 mL of water and shake the jar for 10-20 seconds to thoroughly wet all the mixture. If there is no free liquid at this stage, add a small amount of water, shake the jar, and check for the presence of the free liquid. Repeat this procedure until a few mL of free liquid develop. Transfer 0.5 mL of the free liquid to a 1.5-mL sample tube provided with the test kit, and place a Trait  Bt9 Test Strip vertically into the sample tube.

\* Normally, grind 230 g of corn kernels (if the 230 g of kernels are made up of fewer than 800 kernels, grind 800 kernels).

2.1.1.1.1.2 Judgment of results

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10 minutes after placing the test strip\*, observe the result window of the test strip. If two red lines appear on the result window, the result is positive; if only a control line appears, the result is negative; and if no line appears, the test is invalid

\* If the strip is placed in the tube for more than 10 minutes, the red lines may deepen and hinder an accurate judgment.

### 2.1.1.2 Detection of CBH351 in processed corn products

In accordance with "Extraction and purification of DNA from processed food products" (2.2.3), extract one sample in duplicate. Subject the obtained DNA solution to qualitative PCR under the conditions specified below.

#### 2.1.1.2.1 Qualitative PCR

Qualitative PCR is the method to amplify some parts of the extracted DNA by use of primers and to detect the amplified DNA by electrophoresis.

\* As a very small amount of template DNA is amplified in PCR reaction, it is very important to avoid contamination of DNA other than the target DNA (particularly PCR byproducts). Moreover, DNase, secreted from the human skin surface, must be prevented from entering the PCR reaction mixture as it decomposes DNA. Therefore, disposable tubes and tips must be autoclaved at 121°C for                    or more just prior to use. As for the water used in qualitative PCR, unless otherwise specified, use RO water purified using a reverse osmotic membrane or super-purified water prepared by deionizing distilled water to 17 MW/cm by Milli-Q etc., and then autoclaved at 121°C for 20 minutes or more.

##### 2.1.1.2.1.1 PCR amplification

The PCR reaction mixture is prepared as follows in a PCR reaction tube: add 2.5  $\mu$ L of DNA sample solution adjusted at                    in an ice bath to a mixture including PCR buffer\*1, 0.20 mmol/L dNTP, 3 mmol/L MgCl<sub>2</sub>, 0.2  $\mu$ mol/L 5' and 3' primers\*2 and 0.625-units Taq DNA polymerase\*3, to bring the total volume to 25  $\mu$ L. Place the reaction tube in a PCR amplifier\*4. After starting the reaction by keeping the mixture at 95°C for 10 minutes, maintain the temperature of the mixture at 95°C for 0.5 minutes, at 60°C for 0.5 minutes, and then at 72°C for 0.5 minutes as one cycle, and repeat this cycle for 40 times. Stop the reaction by keeping the mixture at 72°C for 7 minutes and store the mixture at 4°C as the PCR amplified reactant. As the blank reactant, both of the PCR reaction mixtures without primers or the DNA sample solution should be subjected to PCR reaction. To confirm the extraction of DNA, subject a mixture including a

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positive-control (to detect reference gene) primer pair\*5 instead of the primer pairs for detection of CBH351 to PCR reaction for each DNA sample solution.

\*1 PCR buffer

Use a PCR Buffer II (manufactured by PE Biosystems) or a buffer producing an equivalent result.

\*2 Primer pair for detection of CBH351

F-primer (CaM03-5'): 5'-CCT TCG CAA GAC CCT TCC TCT ATA-3'

R-primer (CBH02-3'): 5'-GTA GCT GTC GGT GTA GTC CTC GT-3'

\*3 Taq DNA polymerase

Use an AmpliTaq Gold DNA Polymerase (manufactured by PE Biosystems) or a polymerase producing an equivalent result.

\*4 PCR amplifier

Use a GeneAmp PCR System 9600 (manufactured by Perkin Elmer Co., Ltd.) or a system producing an equivalent result.

\*5 Primer pair for positive control (to detect reference gene)

F-primer (Zein n-5'): 5'-CCT ATA GCT TCC CTT CTT CC-3'

R-primer (Zein n-3'): 5'-TGC TGT AAT AGG GCT GAT GA-3'

### 2.1.1.2.1.2. Agarose gel electrophoresis

The PCR amplified reactant is separated by agarose gel electrophoresis to identify the amplified DNA band.

#### 2.1.1.2.1.2.1 Preparation of agarose gel

Add TAE buffer\*1 to the required quantity of agarose, and heat the mixture to dissolve the agarose. Then, add 5  $\mu$ L of ethidium bromide (10 mg/mL)\*2 per 100 mL of the above solution, cool the mixture to approximately 50°C, pour it into a gel maker, and cool and solidify it at room temperature to make a gel\*3. Although it is desirable to use the gel

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immediately after preparation, it can be preserved in the buffer for several days. The concentration of the gel (1.0% - 4.0%) should be decided based on the band length of the goal product, as it depends on the length of the DNA to be electrophoresed.

### \*1 TAE buffer

Prepare a solution having final concentrations for 40 mmol/L Tris-acetic acid and 1 mmol/L EDTA with distilled water, and use the solution as the TAE buffer.

### \*2 Ethidium bromide

Ethidium bromide is a fluorescent agent that enters the interstices of the double-stranded DNA and has a strong carcinogenic and toxic effect. Wear gloves and a mask during handling.

### \*3 Pre-staining

Pre-staining is used here. Without the additional ethidium bromide here, post-staining of gel following completion of electrophoresis can be used as directed in 2.1.1.2.1.2.3.

#### 2.1.1.2.1.2.2 Electrophoresis

Set the prepared gel in an electrophoresis bath filled with TAE buffer. Add 7.5  $\mu$ L of PCR amplified reactant to a proper quantity of gel-loading buffer, and inject the mixture into the well of gel. If sample injection into the gel takes too long, the DNA will diffuse, resulting in a vague pattern. Next, perform electrophoresis at a constant voltage of 100 V until the BPB in the gel-loading buffer advances to the 1/2-2/3 position of the gel.

#### 2.1.1.2.1.2.3 Staining of gel (post-staining)

When pre-staining is used, the procedures specified in this section are not necessary. Transfer the electrophoresed gel to a container containing a sufficient amount of TAE buffer to soak the gel thoroughly. Add 5  $\mu$ L of ethidium bromide (10 mg/mL) per 100 mL buffer, place the container on a shaker, and perform staining for 30 minutes while shaking the container mildly.

#### 2.1.1.2.1.3 Gel image analysis

Place the stained electrophoretic gel on a piece of \* spread on the stage of a gel image analyzer, and radiate ultraviolet rays (312 nm). Identify the electrophoretic pattern on the display of the gel image analyzer. Judge the presence of the goal band by

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Sample No. 9 will be subjected to a third extraction.

\* "+" means "positive", "-" means "negative", and "/" means "no test."

\* The CBH identification primer pair is as follows:

F-primer (Cry9C-5'): 5'-TAC TAC ATC GAC CGC ATC GA-3'

R-primer (35Ster-3'): 5'-CCT AAT TCC CTT ATC TGG GA-3'

### 2.1.1.3 Detection of CBH351 in corn semi-products (corn grits, corn flour, corn meal etc.)

Test the sample by the lateral flow method as directed in 2.1.1.1.1, except using 230 g of the sample without grinding. If the result is positive, extract DNA from the sample in duplicate as directed in 2.2.1 and subject the DNA extract to the qualitative PCR described in 2.1.1.2.1. When the PCR amplified reactant from one extract displays a 157-bp PCR amplified band in the positive-control (detecting reference gene) primer pair lane, and a 170-bp PCR amplified band is found in the CBH351 detection primer lane, judge the sample to be CBH351-positive.

### 2.1.2 Detection of papaya (55-1)

For raw papayas and their processed food, perform qualitative PCR as described in 2.1.1.2.1, except use the 55-1 detection primer pair (NosC-5', CaMVN-3') that produces (detects) a 207-bp amplified band, the papain primer pair (papain-5', papain-3') that produces (detects) a 211-bp amplified band for positive control (detecting reference gene), and the 55-1 identification primer (for 2nd screening) pair (CaM3-5', GUSn-3') that produces (detects) a 250-bp amplified band.

#### 55-1 Detection primer pair

F-primer (NosC-5'): 5'-TTA CGG CGA GTT CTG TTA GG-3'

R-primer (CaMVN-3'): 5'-CAT GTG CCT GAG AAA TAG GC-3'

#### Papain-gene detection primer pair

F-primer (papain-5'): 5'-GGG CAT TCT CAG CTG TTG TA-3'

R-primer (papain-3'): 5'-CGA CAA TAA CGT TGC ACT CC-3'

#### 55-1 Identification primer pair

F-primer (CaM 3-5'): 5'-CCT TCG CAA GAC CCT TCC TCT ATA-3'

R-primer (GUS n-3'): 5'-TCG TTA AAA CTG CCT GGC AC-3'

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### 2.1.3 Detection of New Leaf Y potato (including processed food)

Perform Qualitative PCR as directed in 2.1.1.2.1, using a detection primer pair that produces (detects) a 225-bp amplified band for detection of New Leaf Y and a Pss primer pair (detects the Patatin gene universally existing in potato) that produces (detects) a 216-bp amplified band for positive control (to detect reference gene). For the confirmation test of New Leaf Y presence, use a confirmation primer pair that produces (detects) a 161-bp amplified band for detection of PVY-cp gene.

New Leaf Y detection primer pair

F-primer (p-FMV05-5') : 5'-AAA AGA GCT GTC CTG ACA GC-3'

R-primer (PVY02-3') : 5'-TCC TCC TGC ATC AAT TGT GT-3'

Pss primer pair (for detection of the patatin gene)

F-primer (Pss01n-5') : 5'-TGA CCT GGA CAC CAC AGT TAT-3'

R-primer (Pss01n-3') : 5'-GTG GAT TTC AGG AGT TCT TCG A-3'

PVY-cp gene detection primer pair

F-primer (PVY01-5') : 5'-GAA TCA AGG CTA TCA CGT CC-3'

R-primer (PVY01-3') : 5'-CAT CCG CAC TGC CTC ATA CC-3'

### 2.2 Extraction and purification of DNA

Unless otherwise specified, the water used in the extraction and purification of DNA is RO water purified using a reverse osmotic membrane, or ultra-purified water produced by deionizing distilled water to 17 MW /cm using a Milli-Q and autoclaving at 121°C for 20 minutes or more.

#### 2.2.1 Extraction and purification of DNA from corn kernels and soybeans

The CTAB method, which uses a mixture of the surfactant cetyl trimethyl ammonium bromide (CTAB) and phenol/chloroform in extraction and purification, is applicable to various areas, leaves virtually no PCR inhibitors, and provides high-purity DNA. However, it has the disadvantages of using harmful reagents such as phenol and chloroform and requiring complicated purification procedures. Currently marketed DNA extraction kits eliminate such disadvantages. They include a silica-gel membrane type, a silica-base resin type, an ion-exchange resin type and a magnet-adsorptive bead type. All of these allow the extraction and purification of DNA usable in PCR from corn kernels and soybeans. This

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section describes purification using the CTAB method, a silica-gel membrane-type kit (QIAGEN DNeasy Planet Mini) or a silica-base resin-type kit (Promega Wizard DNA Clean-Up System).

### 2.2.1.1 CTAB method

Place 2 g of uniformly ground sample into a polypropylene centrifuge tube (50 mL), add 15 mL of CTAB buffer\*1, and homogenize the mixture. Add 30 mL of CTAB buffer while washing the edge of the tube and the tip of the homogenizer, mix them by overturning the tube, and leave the mixture at 55°C for 30 minutes. Stir the mixture, place 600  $\mu$ L of the homogeneous suspension in a micro centrifuge tube (1.5 mL), add 500  $\mu$ L of phenol/chloroform mixtures\*2, mix them by overturning, suspend the mixture mildly using a mixer, centrifuge it at 7,500 x g for 15 minutes, and transfer the water layer (upper) to another micro centrifuge tube. Here, be careful not to disturb the medium layer. Add 500  $\mu$ L of a chloroform/isoamyl alcohol mixtures\*3, mix by overturning, suspend the mixture mildly using a mixer, centrifuge it at 7,500 x g and room temperature for 15 minutes, and transfer the water layer (upper) to another micro centrifuge tube. Add the same volume of isopropyl alcohol (room temperature), mix by overturning, perform centrifugation at 7,500 x g and room temperature for 10 minutes, and throw away the supernatant by decantation. Gently add 500  $\mu$ L of 70% ethanol along the wall of the tube, perform centrifugation at 7,500 x g and room temperature for 1 minute, and suck out as much ethanol as possible without touching the precipitate. Then, vacuum-dry the precipitate for 2-3 minutes, being careful not to dry it completely. Add 50  $\mu$ L of TE buffer\*4, mix fully, and leave the mixture at room temperature for 15 minutes. During this process, overturn the tube occasionally to dissolve the precipitate completely. Add 5  $\mu$ L of RNase A and leave the mixture at 37°C for 30 minutes. After adding 200  $\mu$ L of CTAB buffer, add 250  $\mu$ L of chloroform/isoamyl alcohol mixture, mix them by overturning, suspend the mixture mildly using a mixer, centrifuge it at 7,500 x g and room temperature for 15 minutes, and transfer the water layer (upper) to another micro centrifuge tube without disturbing the medium layer. Add 200  $\mu$ L of isopropyl alcohol, mix by overturning, centrifuge the mixture at 7,500 x g and room temperature for 10 minutes, and throw away the supernatant by decantation. Next, gently add 200  $\mu$ L of 70% ethanol along the wall, perform centrifugation at 7,500 x g and room temperature for 1 minute, and suck out as much ethanol as possible, without touching the precipitate. Vacuum-dry the precipitate for 2-3 minutes, being careful not to dry it completely. Add 50  $\mu$ L of water to the precipitate, mix, and leave the mixture at room temperature for 15 minutes. During this process, overturn the tube for mixing occasionally to dissolve the precipitate completely. Use this solution as the DNA sample stock solution.

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### \*1 CTAB buffer

Place 8 mL of 0.5 mol/L EDTA (pH 8.0), 20 mL of 1 mol/L Tris/HCl (pH 8.0), and 56 mL of 5 mol/L salt solution in a beaker, add enough water to bring the volume to approximately 150 mL, add 4 g of cetyl trimethyl ammonium bromide (CTAB) while stirring, and dissolve it completely. Add enough water to bring the total volume to 200 mL, and autoclave the solution. Use the solution as the CTAB buffer.

### \*2 Phenol/chloroform mixtures

Mix 1 mol/L Tris/HCl (pH 8.0) saturated phenol with chloroform/isoamyl alcohol at 1:1 (v/v), and use this as the phenol/chloroform mixture.

### \*3 Chloroform/isoamyl alcohol mixtures

Mix chloroform and isoamyl alcohol at 24:1 (v/v), and use this as the chloroform/isoamyl alcohol mixture.

### \*4 TE buffer

Prepare a solution having final concentrations of 10 mmol/L Tris/HCl (pH 8.0) and 1 mmol/L EDTA (pH 8.0) with water, and use this solution as the TE buffer.

#### 2.2.1.2 Silica-gel membrane-type kit method

Place 2 g of uniformly ground sample in a polypropylene centrifuge tube (50 mL), add 10 mL of AP1 buffer previously warmed at 65°C and 20  $\mu$ L of RNase A, mix vigorously with a vortex mixer to pulverize lumps of the sample, and leave the mixture at 65°C for 15 minutes. During this process, shake the sample by overturning the centrifugal tube 2-3 times. Add 3250  $\mu$ L of AP2 buffer, leave the mixture on ice for 5 minutes, and centrifuge it at 3,000 x g or more for 5 minutes. Then, apply 500  $\mu$ L of the supernatant to a QI Ashredder spin column, and perform centrifugation at 10,000 x g or more for 4 minutes and place the eluate in a micro centrifuge tube (15 mL). After repeating this procedure, add 1.5 times the eluate volume of AP3 buffer/ethanol mixture\*1. Apply 500  $\mu$ L of this solution to a Mini spin column, and perform centrifugation at 10,000 x g or more for 1 minute. Apply 500  $\mu$ L of the remaining solution to the same Mini spin column, perform centrifugation under the same conditions, and dispose of the eluate. Repeat the same procedures until all of the solution is used. Apply 500  $\mu$ L of AW buffer to the column, perform centrifugation at 10,000 x g or more for 1 minute, and dispose of the eluate. Apply the AW buffer again and repeat the same procedures. After disposing of the eluate, subject the Mini spin column to centrifugation at 10,000 x g or more for 15 minutes to dry it.

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Transfer the Mini spin column to a centrifuge tube of the kit, add 70  $\mu$ L of previously warmed water, leave the mixture for 5 minutes, and centrifuge it at 10,000 x g or more for 1 minute to elute DNA. Add water again and repeat the same procedures. Put the obtained eluates together, and use this as the DNA sample stock solution.

\*1 AP3 buffer/ethanol mixture

Mix AP3 buffer and ethanol (96%-100%) at 1:2 and use this as the AP3 buffer/ethanol mixture.

### 2.2.1.3. Silica-based resin-type kit method

Place 2 g of uniformly ground sample in a polypropylene centrifuge tube (50 mL), add 17.2 mL of extraction buffer\*1, 2 mL of 5 mol/L guanidine-HCl, and 0.8 mL of 20 mg/mL Proteinase K, mix them vigorously using a vortex mixer, and keep the mixture at 55°C - 60°C for 3 hours while shaking the tube. Cool the mixture to room temperature and centrifuge it at 3,000 x g for 10 minutes. If the supernatant is turbid, transfer a portion of it to a micro centrifuge tube (1.5 mL) and further centrifuge it at 14,000 x g for 10 minutes. Place 500  $\mu$ L of the obtained clear supernatant and 1 mL of DNA Clean-Up Resin in a micro centrifuge tube (1.5 mL), and mix them by overturning. Connect an injection syringe to the upper end of a Mini column, and attach them to a manifold (aspirator\*2). After confirming that the cock of the manifold is closed, apply the mixture to the Mini column through the injection syringe. Open the cock, remove the solvent completely by aspiration, and add 2 mL of 80% isopropyl alcohol through the syringe to wash the column. Remove the syringe, attach the Mini column to a micro centrifuge tube (1.5 mL), and centrifuge at 1,000 x g and room temperature for 2 minutes to dry the column. Transfer the Mini column to another micro centrifuge tube (1.5 mL), and drip it with 50  $\mu$ L of water previously warmed to 65°C - 70°C. After 1 minute, centrifuge it at 1,000 x g or more and room temperature for 1 minute to elute the DNA, and use the obtained eluate as the DNA sample stock solution.

\*1 Extraction buffer

10 mmol/L Tris-HCl buffer solution (pH 7.5) containing 150 mM sodium chloride (NaCl); 2 mmol/L EDTA; and 1% SDS.

\*2 Aspirator

If an aspirator is not available, centrifugation can be substituted.

### 2.2.2 Extraction and purification of DNA

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Remove the seeds from the papayas, cut the flesh into 10 mm cubes, and freeze-dry them. Grind the frozen flesh using a mixing mill. Perform extraction and purification of DNA on the sample, in accordance with the CTAB or silica-gel membrane-type kit (QIAGEN DNeasy Plant Mini) methods.

### 2.2.2.1 CTAB method

Place 20 mg of the ground sample in a micro centrifuge tube (1.5 mL), add 150  $\mu$ L of CTAB buffer, and homogenize the mixture using a micro mixer. Add 450  $\mu$ L of CTAB buffer while washing the inside of the tube, mix the solution by overturning, and leave it at 55°C for 30 minutes. Add 500  $\mu$ L of phenol/chloroform mixture, mix the solution by overturning, suspend the mixture mildly using a mixer, centrifuge it at 7,500 x g and room temperature for 15 minutes, and transfer the water layer (upper) to another micro centrifuge tube without disturbing the medium layer. Add 500  $\mu$ L of chloroform/isoamyl alcohol mixture, mix by overturning, suspend the mixture mildly using a mixer, centrifuge it at 7,500 x g and room temperature for 15 minutes, and transfer the water layer (upper) to another micro centrifuge tube. Add the same quantity of isopropyl alcohol (room temperature), mix the solution by overturning, centrifuge it at 7,500 x g and room temperature for 10 minutes, and discard the supernatant by decantation. Gently add 500  $\mu$ L of 70% ethanol along the wall of the tube, centrifuge the precipitate at 7,500 x g and room temperature for 1 minute, and suck out as much of the ethanol as possible without touching the precipitate. Then, vacuum-dry the precipitate for 2-3 minutes, being careful not to dry the precipitate completely. Add 50  $\mu$ L of TE buffer, mix the solution well, and leave it at room temperature for 15 minutes. During this process, overturn the tube occasionally to completely dissolve the precipitate. Add 5  $\mu$ L of RNase A, and leave the mixture at 37°C for 30 minutes. After adding 200  $\mu$ L of CTAB buffer, add 250  $\mu$ L of chloroform/isoamyl alcohol mixture, mix the solution by overturning the tube, suspend the mixture mildly using a mixer, centrifuge it at 7,500 x g and room temperature for 15 minutes, and transfer the water layer (upper) to another micro centrifuge tube. During this process, be careful not to disturb the medium layer. Add 200  $\mu$ L of isopropyl alcohol, mix the solution by overturning the tube, perform centrifugation at 7,500 x g and room temperature for 10 minutes, discard the supernatant by decantation, and remove the isopropyl alcohol using a pipette. Gently add 200  $\mu$ L of 70% ethanol along the tube wall, perform centrifugation at 7,500 x g and room temperature for 1 minute, and suck out as much ethanol as possible without touching the precipitate. Then, vacuum-dry the precipitate for 2-3 minutes, being careful not to dry the precipitate completely. Add 50  $\mu$ L of water, mix the solution, and leave it at room temperature for 15 minutes. During this process, overturn the tube occasionally to dissolve the precipitate completely. Use the solution as the DNA sample stock solution.

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### 2.2.2.2 Silica-gel membrane-type kit method

Place 80 mg of ground sample in a micro centrifuge tube (2 mL), and perform extraction and purification on the sample using a silica-gel membrane kit (QIAGEN DNeasy Plant Mini) as directed below. Add 400  $\mu$ L of AP1 buffer previously warmed to 65°C and 4  $\mu$ L of RNase A of the kit, use a homogenizer to mix them fully and pulverize lumps of the sample, and leave it at 65°C for 15 minutes. During this process, overturn the tube several times to shake it. Then, add 130  $\mu$ L of AP2 buffer, leave the mixture on ice for 5 minutes, and perform centrifugation at 10,000 x g and room temperature for 5 minutes. Apply the supernatant to a QIAshredder spin column, perform centrifugation at 10,000 x g and room temperature for 2 minutes, and transfer the eluate to a micro centrifuge tube (2 mL). Add 1.5 times the eluate volume of AP3 buffer/ethanol mixture to the centrifuge tube and stir the mixture for 10 seconds using a vortex mixer. Apply 500  $\mu$ L of the mixture to a Mini spin column, perform centrifugation at 10,000 x g and room temperature for 5 minutes,\* and discard the eluate. Apply 500  $\mu$ L of the remaining mixture to the same Mini spin column, perform centrifugation under the same conditions, and discard the eluate. Repeat the same procedures until all of the mixture is used up. Add 500  $\mu$ L of AW buffer to the column, perform centrifugation at 10,000 x g and room temperature for 5 minutes, and discard the eluate. Add the AW buffer again and repeat the same procedures. After throwing away the eluate, perform centrifugation at 10,000 x g or more for 15 minutes to dry the Mini spin column. Transfer the Mini spin column to the kit's centrifuge tube, add 50  $\mu$ L of previously warmed water, leave the tube for 5 minutes, and perform centrifugation at 10,000 x g for 1 minute to elute DNA. Add water again and repeat the same procedures. Put the obtained eluates together, and use this mixture as the DNA sample stock solution.

\* If there are deposits in the solution, the column may become clogged. In such a case, prolong the centrifugation time to 10 minutes for complete elution.

### 2.2.3 Extraction and purification of DNA from processed food

The extraction and purification of DNA from processed food products prescribed in Attached Table 2, which are provided in Article 8 of Notification No. 517 of the Ministry of Agriculture, Forestry and Fisheries (2000), is performed as directed in the JAS analytical test handbook, "Genetically Modified Food Test and Analysis Manual for Individual Products," issued by IAI Center for Food Quality, Labeling, and Consumer Services. Extraction and purification of DNA from tacos and tortillas (restricted to products processed using heat), potato (including processed food) and canned papaya are performed as specified below.

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### 2.2.3.1 Extraction and purification of DNA from tacos and tortillas (restricted to products processed using heat)

Freeze-dry the sample and grind it using a mixing mill etc. Place 1 g of the ground sample in a polypropylene centrifuge tube (50 mL), and perform extraction and purification of DNA using an ion-exchange resin type of DNA extraction and purification kit

, as specified below. Add 4 mL of G2 buffer\*1 to the sample, mix the solution vigorously using a vortex mixer, add a further 4 mL of G2 buffer, 100  $\mu$ L of Proteinase K\*2, and 10  $\mu$ L of RNase A, mix the solution by shaking, and leave it at 50°C for 2 hours. During this process, overturn the centrifuge tube 2-3 times to mix. Then, perform centrifugation at 3,000 x g or more and low temperature (4°C) for 15 minutes, transfer the supernatant to a polypropylene centrifuge tube (15 mL), and centrifuge it mildly again. Then, apply the supernatant by dividing it into 2-mL portions in a QIAGEN Genomic-tip 20/G equilibrated with 1 mL of QBT buffer\*1. Wash the tip 3 times with 2 mL of QC buffer\*1, transfer the tip to another centrifuge tube, add 1 mL of QF buffer\*1 previously warmed to 50°C two times, and elute the DNA. Transfer the eluate to a centrifuge tube, add 0.7 times the eluate volume of isopropyl alcohol, mix the solution fully, centrifuge the mixture at 10,000 x g or more and low temperature (4°C) for 15 minutes, and discard the supernatant. Add 1 mL of 70% ethanol and perform centrifugation at 10,000 x g or more and low temperature (4°C) for 5 minutes. After throwing away the supernatant, dry the precipitate using an aspirator, add 100  $\mu$ L of water, leave it at 65°C for 5 minutes, dissolve the DNA by pipetting, and use the solution as the DNA sample stock solution.

\*1 G2, QBT, QC, and QF buffers are included in the kit. If they are not enough, they can be prepared as directed in the kit's explanatory leaflet.

\*2 Use QIAGEN's proteinase or that having an equivalent effect.

Grind the sample using a mixing mill etc. and place 200mg of the ground sample in a polypropylene centrifuge tube (15 mL). If the water content of the sample is high, perform centrifugation at 8,000 x g for 15 minutes and discard the supernatant. Perform extraction and purification of DNA on the sample using a silica-gel membrane-type kit (QIAGEN DNeasy Plant Mini), as directed below.

Add 1.5 mL of AP1 buffer previously warmed to 65°C and 10  $\mu$ L of RNaseA to the sample, mix them vigorously using a vortex mixer to pulverize lumps of sample, and leave the mixture at 65°C for 15 minutes. During this process, shake the sample by overturning the tube for several times. Then, add 400  $\mu$ L of AP2 buffer, leave the tube on ice for 5

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minutes, and perform centrifugation at 10,000 x g and room temperature for 5 minutes. Transfer the supernatant to another centrifuge tube, apply 500  $\mu$ L of the supernatant to a QIAshredder spin column, perform centrifugation at 10,000 x g and room temperature for 2 minutes, and transfer the eluate to the kit's centrifuge tube. Repeat these procedures until all of the supernatant has been applied. Divide the obtained eluate into 2 equal parts, transfer them to separate centrifuge tubes (2 mL), add 1.5 times the eluate quantity of AP3 buffer/ethanol mixture to each tube, and stir them for 10 seconds using a vortex mixer to obtain a solution. Apply 500  $\mu$ L of the obtained solution to a Mini spin column, perform centrifugation at 10,000 x g and room temperature for 1 minute, and discard the eluate. Of the remaining solution, apply a further 500  $\mu$ L to the same Mini spin column, perform centrifugation under the same conditions, and discard the eluate. Repeat the same procedures until all of the solution is used up. Next, add 500  $\mu$ L of AW buffer to the column, perform centrifugation at 10,000 x g and room temperature for 1 minute, and discard the eluate. Add the AW buffer again and repeat the same procedures. After discarding the eluate, dry the Mini spin column by centrifugation at 10,000 x g or more for 15 minutes. Transfer the Mini spin column to the kit's centrifuge tube, add 50  $\mu$ L of previously warmed water, leave the column for 5 minutes, and perform centrifugation at 10,000 x g and room temperature for 1 minute to eluate the DNA. Add water again, perform the same procedures, put the obtained eluates together, and use the eluate as the DNA sample stock solution.

### 2.2.3.3 Extraction and purification of DNA from canned papaya

Wash the sample fully with water, freeze-dry it, and grind it using a mixing mill. Place 2 g of the ground sample in a polypropylene centrifuge tube (50 mL) and perform extraction and purification of DNA on it using an ion-exchange resin-type DNA extraction and purification kit (QIAGEN Genomic-tip), as directed below.

Add 7.5 mL of G2 buffer to the sample, mix them vigorously using a vortex mixer, add further 7.5 mL of G2 buffer, 200  $\mu$ L of QIAGEN Proteinase K, and 20  $\mu$ L of RNase A, mix them by shaking, and leave the mixture at 50°C for 2 hours. During this process, overturn the tube 2-3 times to mix. Perform centrifugation at 3,000 x g or more and low temperature (4°C) for 15 minutes, transfer the supernatant to a polypropylene centrifuge tube (15 mL), and centrifuge it mildly again. Apply the supernatant by dividing it into 2-mL portions in a QIAGEN Genomic-tip 20/G equilibrated with 1 mL of QBT buffer. Wash the tip with 2 mL of QC buffer 3 times, transfer it to another centrifuge tube, and add 1 mL of QF buffer previously warmed to 50°C twice, and eluate the DNA. Transfer the eluate to a centrifuge tube, add a 0.7 times the eluate volume of isopropyl alcohol, mix them fully, perform centrifugation at 10,000 x g or more and low temperature (4°C) for 15 minutes, discard the supernatant, add 2 mL of 70% ethanol, and perform centrifugation

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again at 10,000 x g or more and low temperature (4°C) for 5 minutes. Discard the supernatant, dry the precipitate with an aspirator, add 100  $\mu$ L of water, leave the solution at 65°C for 5 minutes, dissolve the DNA by pipetting, and use the solution as the DNA sample stock solution.

### 2.2.4 Determination of the purity of DNA in DNA sample stock solutions and storage of DNA sample solution

Dilute a proper amount of DNA sample by 10 times with TE buffer, determine the UV absorption spectrum of this dilution at 200-300 nm, and record its absorbencies at 260 nm and 280 nm (O.D. 260 and O.D. 280\*). Using the value for O.D. 260 as 50 ng/ $\mu$ L DNA, calculate the concentration of DNA and O.D. 260/O.D. 280. If the ratio is 1.7-2.0, the DNA is fully purified. Based on the determined DNA concentration, dilute the DNA stock solution with water to a concentration required for subsequent PCR reaction, dispense the solution in 20  $\mu$ L portions into micro tubes, freeze it, and store it at -20°C for use as the DNA sample solution. Use the sample solution immediately after melting. The remaining sample should not be used and must be discarded. Unless the concentration of the DNA sample stock solution reaches that required for PCR, use it as the DNA sample solution without dilution.

\* Regard the value for O.D. 260 as the absorbance derived from DNA, and the value for O.D. 280 as the absorbance derived from impurities such as proteins.

## 3. A Method for the foods produced by recombinant-DNA techniques whose safety assessments were reviewed

### 3.1. Soybeans

#### 3.1.1 ELISA-Method

The ELISA Method is used for detecting CP4EPSPS protein in a sample. Take 0.1 g of a ground sample sieved through a 100-mesh sieve (the mesh length is 150  $\mu$ m), subject them to the GMO Soya RUR Test Kit manufactured by SDI whose method is explained in the leaflet. The method is described as follows:

Measure 0.1 g of both a sample and a standard sample into centrifuging tubes made of polypropylen (50 mL). Add 20 mL of Soya Extraction buffer and, mix for at least 1 minute with a vortex mixer, and then let the mixture settle for five minutes until a clear supernatant appears. Use the supernatant as a sample solution. As the scope of the standard curve is 0-2.5% for SDI Kit, it is necessary to prepare also 1/10 dilution (diluted with Soya Extraction Buffer) in order to interpolate the quantitative value within the scope for a

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sample solution. Put 100  $\mu$ L of each sample solution into Wells and incubate them at room temperature for one hour. Then wash the Wells four times with Wash buffer, add 100  $\mu$ L of the RUR Conjugate 1 and incubate them at room temperature for 30 minutes. Again wash the Wells four times with Wash buffer. Then add 100  $\mu$ L of the RUR Conjugate 2, and incubate it at room temperature for 30 minutes, then wash the Wells four times. Add 100  $\mu$ L of Color Reagent and incubate them at room temperature for 30 minutes. Add 100  $\mu$ L of the Stop Solution to stop color development. Read the absorbance of the developed color at 450 nm using a microtiter plate reader, and determine the GMO content by using the standard curve of the reference standard sample\*. Duplicate the test and average the two data.

\*extra-cost options

### 3.1.2. Quantitative PCR Method

Here the quantitative PCR method is used with Taqman Chemistry. In this method a luminescent origo-nucleotide probe is used. In this method, a luminescent origo-nucleotide probe, which anneals in the DNA sequence between the both primer pairs and has two dyes (reporter and quencher) is used. When this probe is hydrolysed by the 5'-nuclease activity of the DNA polymerase, the reporter dye is separated from the quencher dye and emits luminescence. The strength of luminescence shows exponential increase for the PCR cycle. Therefore, by comparing the PCR cycle number which has reached a certain luminescence, the amount of original DNA can be obtained. The quantitative analysis of GM products is to determine the amount of recombinant gene compared to the internal standard which is an universally present gene in the non-GMOs. In the case of soybeans, use the lectin-gene that is ubiquitous in soybeans as the internal standard. The number of copies of the lectin-gene in the DNA-sample is obtained by quantitative PCR which uses a primer pair (Le1-n02) detecting the lectin-gene and a probe (Le1-Taq). At the same time, the number of copies of the GM DNA recombinant gene in a sample is obtained by another quantitative PCR using a primer pair (RRS-01) and probe (RRS-Taq). The percentage of GMO in a sample is obtained from the following way.

- Divide the number of recombinant gene copies by the number of the lectin-gene copies.
- Divide this value again by the Internal standard ratio\*(which is determined previously), and multiply the result by 100.

In the following paragraphs we will describe the quantitative PCR. As for the water used in Quantitative PCR for 5.1.2., 5.2.2., unless otherwise specified, use water, which is purified

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RO water purified using a reverse osmotic membrane or deionizing distilled water to 17 MΩ/cm by Milli-Q and is autoclaved at 121°C for 20 minutes or more.

\*Internal standard ratio: the ratio of the number of target recombinant gene to the number of internal standard gene (e.g. in the case of Soybeans, lectin gene) in seeds detected by a certain quantitative PCR. The internal standard ratio is provided separately for each primer pair and probe.

### 3.1.2.1 Making the standard solution for the standard curve

Prepare a GM plasmid DNA solution of which the number of copies has already been determined, using salmon sperm DNA solution (5ng/μL), and dilute GM plasmid DNA solution so that plasmid DNA in 2.5 μL makes standard curve (e.g. 10, 125, 1,500, 20,000, 250,000 copies)\*. Use the dilution as the standard solution for the standard curve.

\* Dilution of the DNA solution

### 3.1.2.2 Preparation of the solution for the PCR reaction

The composition of 25 μL solution for the PCR reaction is the following:

- 12.5 μL of Universal PCR Master Mix \*1
- 0.5 μmol/L (12.5 pmol) of Primer
- 0.2 μmol/L (5pmol) of Probe
- 2.5 μL (50ng) of 20ng/μL DNA sample solution or 2.5 μL of the standard solution for the standard curve, or 2.5 μL of 5ng/μL the salmon sperm DNA solution (blank experimental solution)

In order to minimize the differences between the PCR reactions, prepare PCR reaction solution for 3wells per each DNA sample solution\*2. Firstly prepare the mixture (Master Mix) of the Universal PCR Master Mix with the primer pair and the probe, which are a little more than quantity required, then mix up Master Mix with the each DNA sample and/or the blank sample in different centrifuge tubes. Divide this mix solution of 25 μL into each well on a plate and put the lid of the plate on. At this time, to be impartial, close alternately from both sides. Then, using a special roller, seal up the wells. Lastly observe the bottom of the wells, and if there are bubbles, knock lightly a rim of the plate to dissolve the bubbles.

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### \*1 Universal Master Mix

This solution has high viscosity. Therefore, when adding to and mixing with other solutions, it is necessary to make sure that they are thoroughly mixed. If not, the PCR reaction possibly might not work well. Therefore after mixing it with a vortex for about 3 seconds just before usage, centrifuge lightly, and use after the solution has amassed at the bottom of the tube. When dividing the mixture into the wells, remember stirring and centrifuging will be difficult after this, put it into the bottom of the well accurately.

### \*2 Preparation of the solution for the quantitative PCR reaction

Take the reagents from the refrigerator, and if necessary, melt them at room temperature, then keep them on ice. If use the same tip successively for dividing, the usual handling will not work from the second time onward because the air is cooled in the pipette. It is necessary to read up on how to use the pipette at low temperatures in the manual of the pipette before starting the procedure.

#### 3.1.2.3 Quantitative PCR

ABI PRISM7700 or 5700 is used for the quantitative PCR. Set the plate onto the device. Check that the Cover temperature is around 105°C, begin the reaction and the reading of the data. The conditions for the reaction are as follows: After maintaining it at 50. C for 2 minutes, keep it at 95. C for 10 minutes. After starting the reaction with a hot start method, conduct 40 cycles of amplification reaction (30 seconds at 95. C and 1 minute at 59. C for 1 cycle). After the end of the reaction, analyse the results.

#### 3.1.2.4. Drawing the standard curve

The standard curve for the internal standard gene and the recombinant gene can be obtained from the following way:

On the amplification curve plotting the increasing amount of luminescent signals (DRn), select the DRn part where luminescent signals from both of the standard solution for the standard curve and DNA sample solution amplify exponentially, and draw a threshold line (Th). At this time, a threshold line (Th) must not cross the non-specific amplified curve from the blank sample solution. The exact drawing of a threshold line (Th) is performed as directed in the JAS analytical test handbook, "Genetically Modified Food Test and Analysis Manual for Individual Products," issued by IAI Center for Food Quality, Labeling, and Consumer Services of MAFF (the Ministry of Agriculture, Forestry and Fisheries). Use the point, where the Th and the luminescent signal of the standard solution for the standard curve cross, as the "threshold cycle" (Ct). Next, plot the Ct value (y-axis) for the

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logarithm of copies of the (each) standard solutions for standard curve (x-axis). Use the approximate straight line obtained from each Ct as the standard curve.

### 3.1.2.5. Calculation of the content of the foods produced by recombinant-DNA techniques in samples

Obtain the Ct for unknown DNA sample, by using the Th of the standard curve. For each of the internal standard gene and recombinant gene, calculate the numbers of original genes' copies for all 3 Wells from the standard curve. Use the average of these values is to be the number of copies of the original internal standard gene and original recombinant gene. Calculate the amount of GMO contained in the products from following formula:

The GMO content in the products (%) = [(number of copies of the original recombinant gene) / (number of copies of the original internal standard genes x internal standard ratio)] x 100

In the case of soybeans, relevant GM variety is only "Round Up Ready Soybean", the value calculated from the number of copies of Le1 and RRS genes shows the percentage GMO in the product.

### 3.1.3 Judgment of result

Extract 3 times for each IP handling samples. If the average of the amount of GMO by ELISA or Quantitative PCR is more than 5%, the IP handling has possibly failed.

## 3.2 Corn

In case of corn, use only the quantitative PCR, not ELISA to quantify. Because not only there are various GM varieties which have different objective proteins, but also the expression amounts of the proteins in each GM variety are different even if a single kind of protein manifests for each GM variety.

### 3.2.1 Quantitative PCR method

As there are some GM varieties for corn, firstly subject to the screening, and quantify each GM variety separately based on the screening results, and judge from the result of total of the quantitative value for each GM variety. In the case of corn, use the starch synthase IIb (SSIb) gene as the universally present 'internal standard' gene. Calculate the amount of GMO content using following numbers formula as shown in 3.1.2.5.

- The number of copies of starch synthase IIb (SSIb) gene (obtained from the

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primer pair SSIIb and probe SSIIb-Taq, which detects starch synthase IIb (SSIIb) gene).

The number of copies of target recombinant gene (obtained from the primer pair and probe, which detects target recombinant gene,)

### 3.2.1.1 Screening

#### 3.2.1.1.1. The Quantification of GM varieties inserted 35S promoter derived from the Cauliflower mosaic virus

As the GM varieties, Event 176, Bt11, T25 and Mon810, are inserted the 35S promoter (CaM) derived from the Cauliflower mosaic virus, concerning the mixtures in these varieties, it is possible to infer the approximate amount of GM contained from an index of these DNA contents. The method for the quantification is the same method as the quantitative PCR method described for the soybeans apart from the primer and the probe. Use starch synthase IIb (SSIIb) gene as internal standard gene, and the primer pair SSIIb and the probe SSIIb-Taq, which detect starch synthase IIb (SSIIb) gene. The primer pair and the probe for target gene is P35S-1 and P35S-Taq\*, respectively. The amount of GM inserted the 35S promoter (CaM) is calculated using a internal standard ratio that has been defined in the Annex.

\*The internal standard ratio in case of using P35S-1 and P35S-Taq, refer to calculated value for Mon810. Mon810 is the GM variety which is produced widest in the United States, there is only one copy of the 35S promoter in a recombinant DNA of Mon810, then there is little possibility in underestimating the amount of GMO in a sample. The concentration of P35S-Taq is 0.1  $\mu$  mol/L (2.5 pmol).

#### 3.2.1.1.2 Determination of GA21

GA21 GM variety is not inserted CaM gene. Therefore, it is necessary to quantify the content of this variety using the same DNA sample, which is used for the determination of CaM gene. The amount of GA21 is obtained from calculating the number of copies of GA21 with same method as described in 3.2.1.1.1., using GA21-3 and GA21-Taq (GA21 specific primer and probe, respectively), for the same sample as analyzing CaM.

#### 3.2.1.1.3 Judgment of the Results

Each sample is extracted once and subject to quantitative PCR. If the sums of the amount of GMO inserted CaM gene sequence and the amount of GA21 is higher than 4.5%, extract more twice. In this case, perform Specific determination of GM corn varieties for each DNA sample solution (from 3 times extraction).

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### 3.2.1.2. Specific determination of GM corn varieties

#### 3.2.1.2.1. Determination of Event 176, Bt 11, T25 and Mon 810

As quantitative primer and probe, use E-176-2 and E 176-Taq, Bt 11-3 and Bt 11-Taq, and T25-1 and T25-Taq, and M810-2 and M810-Taq, for GM varieties of Event 176, Bt 11, T25 and Mon 810, respectively. The number of copies of the varieties Event 176, Bt 11, T25 and Mon 810 are obtained by the same method as described in 3.2.1.1. Then determine the content of Event 176, Bt 11, T25 and Mon 810.

#### 3.2.2 Judgement of the Results

Calculate the total amount for each DNA sample solution of GA21, Event 176, Bt 11, T25 and Mon 810 obtained under 3.2.1.1.2. If the average of each total for a sample is more than 5%, the IP handling of the sample has possibly failed.

(Annex) Internal Ratio.

Food	Varieties	Internal Ratio	Standard	Remarks
Soybeans	Roundup Ready Soybean	0.95		Le1-n02 and Le1-Taq, RRS-01 and RRS-Taq
Corn	Not specified (Screening)	0.39		SSIIb and SSIIb-Taq, P35S-1 and P35S-Taq
Corn	GA21	1.40		SSIIb and SSIIb-Taq, GA21-3 and GA21-Taq
Corn	Event176	2.05		SSIIb and SSIIb-Taq, E176-2 and E176-Taq
Corn	Bt11	0.50		SSIIb and SSIIb-Taq, Bt11-3 and Bt11-Taq
Corn	T25	0.34		SSIIb and SSIIb-Taq, T25-1 and T25-Taq
Corn	Mon810	0.38		SSIIb and SSIIb-Taq, Mon810-2 and Mon810-Taq

(For your information)

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1 You can purchase AP1/AP2 buffer and RNase A separately used in "2.2.1.2 Silica-gel membrane-type kit method" from QIAGEN (Forefront Tower II, 3-13-1 Kachidoki Chuo-Ward Tokyo 104-0054 Japan, Tel: +81-3-5547-0811, Fax: +81-3-5547-0818).

2 You can purchase GM standard plasmid DNA used in "3.1.2.1 Making the standard solution for the standard curve" from NIPPON Gene (1-29 Tonya-cho Toyama-City Toyama 930-0982 Japan, Tel: +81-76-451-6548, Fax: +81-76-451-6547).

3 You can purchase the primer pairs and the probes used in "3.1.2.2 Preparation of the solution for the PCR reaction" from NIPPON GENE and FASMAC (5-1-3 Midorigaoka, Atsugi-City, Kanagawa 243-0041 Japan, Tel: +81- 46-295-8787, Fax: +81- 46-294-3738).

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


Dilute each extracted DNA by H<sub>2</sub>O to obtain ~100  $\mu$  L 10 ng /  $\mu$  L DNA solutions.

Store the original DNA solutions at -20°C.

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<sup>5</sup> Primers for detection of endogenous gene, PSS (potato sucrose synthetase)、Primers for detection of the region between exogenous regulatory gene and exogenous structure gene

<sup>6</sup> Primers for detection of exogenous structure gene

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目的：Validation of inspection for GM potato by qualitative PCR methods<sup>1</sup>

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1. Non-GM：DANSHAKU 男爵<sup>2</sup>
2. GM：0.1% NewLeaf plus、1% NewLeaf plus、0.1% NewLeaf Y, and 1% NewLeaf plus<sup>3,4</sup>

DNA extraction (2 copies)：

DNA conc. of each extraction：

	230 nm	260 nm	280nm	320 nm	260 nm / 280 nm	Conc. (ng / $\mu$ L)
DANSYAKU 男爵 No.1	0.097	0.144	0.096	0.021	1.532	144
DANSYAKU 男爵 No.2	0.091	0.141	0.094	0.020	1.5	141
0.1% NewLeaf plus No.1	0.081	0.119	0.081	0.019	1.469	119
0.1% NewLeaf plus No.2	0.074	0.117	0.079	0.016	1.481	117
1% NewLeaf plus No.1	0.064	0.096	0.063	0.011	1.524	96
1% NewLeaf plus No.2	0.074	0.108	0.074	0.018	1.459	108
0.1% NewLeaf Y No.1	0.068	0.105	0.069	0.012	1.522	105
0.1% NewLeaf Y No.2	0.061	0.089	0.060	0.011	1.483	89
1% NewLeaf Y No.1	0.062	0.095	0.063	0.012	1.508	95
1% NewLeaf Y No.2	0.054	0.086	0.055	0.007	1.564	86

Dilute each extracted DNA by H<sub>2</sub>O to obtain ~100  $\mu$ L 10 ng /  $\mu$ L DNA solutions.  
Store the original DNA solutions at -20°C.

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<sup>2</sup> The major strain of potato eaten in Japan.

<sup>3</sup> Standard are obtained from MONSATO Co.

<sup>4</sup> Add 100% GM to Non-GM DANSHAKU potato to make 0.1% and 1% GM samples.

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### 定性試驗 (Qualitative PCR)

#### The compositions of reaction solution<sup>5</sup>

	Final conc.
10 × PCR buffer	1 ×
MgCl <sub>2</sub>	3.0 mM
dNTP	2.0 mM
F-primer	0.2 μM
R-primer	0.2 μM
Taq polymerase <sup>6</sup>	0.625-units
Template DNA	25 ng
H <sub>2</sub> O	Fill to bring a 25 μL final reactant

#### Preparation

	Volume
10 × PCR buffer	2.5 μL
MgCl <sub>2</sub>	3 μL
dNTP	2.5 μL
F-primer	1 μL
R-primer	1 μL
Taq polymerase	0.125 μL
Template DNA (10 ng/μL)	2.5 μL
H <sub>2</sub> O	12.375 μL
Total volume	25 μL

#### The PCR program<sup>7</sup>

95°C 10 min  
 95°C 30 sec   
 60°C 30 sec 40 cycles  
 72°C 30 sec   
 72°C 7 min

<sup>5</sup> See "Testing for Foods Produced by Recombinant DNA Techniques in Japan", 2.1.1.2.1.

<sup>6</sup> AmpliTaq Gold, 5 units/μL

<sup>7</sup> See "Testing for Foods Produced by Recombinant DNA Techniques in Japan", 2.1.1.2.1

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### Detection of New Leaf Y potato (including processed food)<sup>8</sup>

Using a detection primer pair that produces (detects) a 100 bp amplified band for detection of New Leaf Y and a Pss<sup>9</sup> primer pair (detects the Patatin gene universally existing in potato) that produces (detects) a 100 bp amplified band for positive control (to detect reference gene). For the confirmation test of New Leaf Y presence, use a confirmation primer pair that produces (detects) a 100 bp amplified band for detection of PVY-cp<sup>10</sup> gene.

#### New Leaf Y detection primer pair

F-primer (p-FMV05-5'): 5'-AAA AGA GCT GTC CTG ACA GC-3'

R-primer (PVY02-3'): 5'-TCC TCC TGC ATC AAT TGT GT-3'

#### Pss primer pair (for detection of the patatin gene)

F-primer (Pss01n-5'): 5'-TGA CCT GGA CAC CAC AGT TAT-3'

R-primer (Pss01n-3'): 5'-GTG GAT TTC AGG AGT TCT TCG A-3'

#### PVY-cp gene detection primer pair

F-primer (PVY01-5') : 5'-GAA TCA AGG CTA TCA CGT CC-3'

R-primer (PVY01-3') : 5'-CAT CCG CAC TGC CTC ATA CC-3'

### Detection of New Leaf plus potato

Using a detection primer pair that produces (detects) a 100 bp amplified band for detection of New Leaf plus and a Pss primer pair for positive control. For the confirmation test of New Leaf plus presence, use a confirmation primer pair that produces (detects) a 100 bp amplified band for detection of PLRV-erp<sup>11</sup> gene.

#### New Leaf plus detection primer pair

F -primer (p-FMV02-5'): 5'-AAA TAA CGT GGA AAA GAG CTG TCC TGA-3'

R-primer (PLRV01-3'): 5'-AAA AGA GCG GCA TAT GCG GTA AAT CTG-3'

Pss primer pair (for detection of the patatin gene)--- The same to the New Leaf Y

---

<sup>8</sup> See "Testing for Foods Produced by Recombinant DNA Techniques in Japan", 2.1.3.

<sup>9</sup> Potato sucrose synthetase

<sup>10</sup> Potato virus V core protein

<sup>11</sup> Potato leaf roll virus replicase

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PLRV-rep gene detection primer pair

F-primer (PLRV-rep1-5') : 5'-CTT CTT TCA CGG AGT TCC AG-3'

R-primer (PLRV-rep1-3') : 5'-TCG TCA TTA AAC TTG ACG AC-3'

Perform PCR reaction in ABI PCR system 9700

Electrophoresis

Observe the results under UV light


Results :

	PSS 216 bp band	Detection of New Leaf Y 225 bp	Detection of New Leaf plus 234 bp	Confirmation of New Leaf Y 161 bp	Confirmation of New Leaf plus 172 bp
DANSYAKU 男爵 No 1	+	-	-	-	-
DANSYAKU 男爵 No.2	+	-	-	-	-
0.1% NewLeaf plus No.1	+	-	+	-	+
0.1% NewLeaf plus No 2	+	-	+	-	+
1% NewLeaf plus No.1	+	-	+	-	+
1% NewLeaf plus No 2	+	-	+	-	+
0.1% New Leaf Y No 1	+	+	-	+	-
0.1% New Leaf Y No 2	+	+	-	+	-
1% New Leaf Y No 1	+	+	-	+	-
1% New Leaf Y No.2	+	+	-	+	-
H <sub>2</sub> O negative control	-	-	-	-	-

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
### Detection PCR

#### 1. PSS detection

From lane 1 to lane 13 : Marker 、 Non-GM1 、 Non-GM2 、 0.1% NLP1 、 0.1% NLP2 、 1% NLP1 、 1% NLP2 、 0.1% NLY1 、 0.1% NLY2 、 1% NLY1 、 1% NLY2 、 NTC and marker. All PCR reaction produces a 216 bp amplicon, except for NTC.


#### 2. Detection of New leaf plus

From lane 1 to lane 13 : Marker 、 Non-GM1 、 Non-GM2 、 0.1% NLP1 、 0.1% NLP2 、 1% NLP1 、 1% NLP2 、 0.1% NLY1 、 0.1% NLY2 、 1% NLY1 、 1% NLY2 、 NTC and marker. Lane 4 、 5 、 6 and 7 have 225 bp products. The other PCR reactions have no product.

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### 3. Detection of New leaf Y

From lane 1 to lane 13 : Marker 、 Non-GM1 、 Non-GM2 、 0.1% NLP1 、 0.1% NLP2 、 1% NLP1 、 1% NLP2 、 0.1% NLY1 、 0.1% NLY2 、 1% NLY1 、 1% NLY2 、 NTC and marker. Lane 8 、 9 、 10 and 11 have 234 bp products. The other PCR reactions have no product.


### 4. Confirmation of New leaf plus

From lane 1 to lane 13 : Marker 、 Non-GM1 、 Non-GM2 、 0.1% NLP1 、 0.1% NLP2 、 1% NLP1 、 1% NLP2 、 0.1% NLY1 、 0.1% NLY2 、 1% NLY1 、 1% NLY2 、 NTC and marker. Lane 4 、 5 、 6 and 7 have 172 bp products. The other PCR reactions have no product.

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### 5. Confirmation of New leaf Y


From lane 1 to lane 13 : Marker 、 Non-GM1 、 Non-GM2 、 0.1% NLP1 、 0.1% NLP2 、 1% NLP1 、 1% NLP2 、 0.1% NLY1 、 0.1% NLY2 、 1% NLY1 、 1% NLY2 、 NTC and marker. Lane 8 、 9 、 10 and 11 have 161 bp products. The other PCR reactions have no product.

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
**Extraction Methods Using QIAGEN Plant Maxi Kit**

**For 2 gram samples ( for soybean and rice ) :**

Weight 2 g of sample powder, put into a 50 mL tube

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

**For 1 gram samples ( for maize, wheat, barley and buckwheat ) :**

Weight 1 g of sample powder, put into a 50 mL tube

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Translated from 「Protocol of DNA Extraction – Using for Blinder test」, Version 1.0 (Jan. 26. 2001), Japan.

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### Use QIAGEN DNeasy Plant Maxi Kit to Extract DNA of Cereals

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#### Grinding and sampling :

1. Use Retch MM 200 to mill an appropriate amount of sample granules.
2. Weight 1 g or 2 g sample powder :

	Sampling weight
Barley	1 g
Rice	2 g
Buckwheat	1 g
Corn Bt11	1 g
Corn T25	1 g
Corn MON 810	1 g
Corn Event 176	1 g
Soybean Round up Ready	2 g

#### Extraction of DNA :

##### Method :

1. Use QIAGEN DNeasy Plant Maxi Kit
2. Use「Protocol of DNA Extraction – Using for Blinder test」, Version 1.0( Jan. 26. 2001 ) , Japan.

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The conc. and quality of extracted DNAs :

	230 nm	260 nm	280nm	320 nm	260 nm / 280 nm	Conc. ( ng / $\mu$ L
Barley	0.138	0.179	0.113	0.020	1.584	179 <sup>1</sup>
Rice	0.052	0.044	0.027	0.005	1.630	44 <sup>2</sup>
Buckwheat	0.136	0.068	0.046	0.000	1.478	68 <sup>3</sup>
Corn Bt11	0.234	0.444	0.261	-0.000	1.701	444
Corn T25	0.436	0.899	0.530	0.003	1.692	899
Corn MON 810	0.431	0.884	0.525	0.013	1.684	884
Corn Event 176	0.233	0.432	0.252	-0.006	1.714	432
Soybean Round up Ready	0.656	0.124	0.133	0.048	0.932	124 <sup>4</sup>

<sup>1</sup> Extraction of DNA from 1 g of barley obtained DNA in low conc.

<sup>2</sup> Previous extraction of DNA from 1 g of rice obtained very rare DNA. In this time, extraction of DNA from 2 g of rice still obtained only a few of DNA.

<sup>3</sup> Extraction of DNA from 1 g of buckwheat also obtained a few of DNA.

<sup>4</sup> The DNA quality in this extraction is not good, and the DNA conc. is low

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### Quantitative PCR Performed by ABI 7000

**Purpose :** Blinder tests for validation of coefficient values in quantitative PCR performed by ABI 7000

**Introduction :** Although the coefficient values of 5 GM maize and 1 GM soybean<sup>1</sup> in ABI 5700 and ABI 7700 machines were established<sup>2</sup>, Japan government still plan to establish the coefficient values of these GM cereals in ABI 7000 and in ABI 7900 machines. The preliminary works have obtained the coefficient values of these GM cereals in ABI 7000 and in ABI 7900 machines<sup>3</sup>. Now the blinder test are undergoing to validation these new coefficient values, and to confirm methods of GMO inspections that used ABI 7000 and in ABI 7900 machines. In this experiment, 3 detection of MON 810 GM corns is performed. Sample A is judged as Non-GM. Sample B contains 5% MON 810. Sample C is 100% MON 810 corn.

**Reference materials :**

Plasmid DNA in 20、125、1,500、20,000 and 250,000 copies/ $\mu$ L

**Unknown samples :**

Unknown sample A、Unknown sample B and Unknown sample C

**Detection<sup>4</sup> :**

1. SSIIb gene ( endogenous gene )
2. MON 810 specific region
3. 35S promoter region

Each performs 3 repeats

---

<sup>1</sup> GM maize : Bt11、Event 176、T25、MON 810 and GA21. GM soybean : Round-up Ready.

<sup>2</sup> Testing for Foods Produced by Recombinant DNA Techniques, version 2, 2002/4/30

<sup>3</sup> Data are still closed

<sup>4</sup> The information of primers and of probes are closed

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PCR program :

50°C 2min  
 95°C 10 min  
 95°C 30 sec }  
 40 cycles  
 59°C 1 min }

Results :

1. Determine the threshold line of SSIb runs<sup>5</sup> :

m	2 <sup>m</sup>	Th. Line	Corr.	Slope	y-intercept	增 (A)	ΔA  (%)	ΔA  条	備考	採用値
0	1	0.001	0.009	-0.323	18.229	1247.95				
1	2	0.002	0.003	-0.182	18.275	10758.63	24801.61			
2	4	0.004	0.013	-0.406	20.101	289.39	99.91			
3	8	0.008	0.416	-3.216	33.305	2.05	99.29			
4	16	0.016	0.810	-4.047	39.383	77	13.68			
5	32	0.032	0.987	-3.816	40.282	1.83	3.51			
6	64	0.064	0.994	-3.698	40.703	1.86	1.93			
7	128	0.128	0.996	-3.647	41.364	1.88	0.87			
8	256	0.256	0.996	-3.610	42.204	1.89	0.66			*
9	512	0.512	0.995	-3.612	43.414	1.89	0.03			
10	1024	1.024					0.71		Plot out	
11	2048	2.048								

<sup>5</sup> The rule to determine the threshold line is formulated by specialists of statistics in NFRI. The principles will not be introduced here.

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The threshold line of SSIIb runs is determined as 0.128.

2. Determine the threshold line of 35S runs :

m	2 <sup>m</sup>	Th. Line	Corr.	Slope	y-intercept	增 (A)	ΔA  (%)	ΔA  条	備考	採用値
0	1	0.001	0.103	-0.909	19.209	12.58	2.48			
1	2	0.002	0.105	-0.918	19.264	12.27	4.76			
2	4	0.004	0.199	-0.937	19.373	11.68	57.30			
3	8	0.008	0.165	-1.433	21.944	4.99	12.77			
4	16	0.016	0.141	-1.333	22.446	5.63	36.60			
5	32	0.032	0.193	-1.811	27.969	3.57	45.15			
6	64	0.064	0.979	-3.431	37.622	1.96	0.12	2		
7	128	0.128	0.995	-3.425	38.616	1.96	1.06	2		
8	256	0.256	0.994	-3.480	39.758	1.94	1.02	2		*
9	512	0.512	0.994	-3.535	40.873	1.92	0.05	2		
10	1024	1.024	0.994	-3.552	42.128	1.92				
11	2048	2.048							Plot out	

The threshold line of 35S runs is determined as 0.256.

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3. Determine the threshold line of MON 810 specific runs :

m	2 <sup>m</sup>	Th. Line	Corr.	Slope	y-intercept	增 (A)	ΔA  (%)	ΔA  条	備考	採用值
0	1	0.001	0.491	-0.737	18.556	22.78				
							35.58			
1	2	0.002	0.205	-0.857	19.541	14.67				
							49.14			
2	4	0.004	0.193	-1.146	20.880	7.46				
							63.80			
3	8	0.008	0.400	-2.317	27.369	2.70				
							10.61			
4	16	0.016	0.508	-2.612	31.016	2.41				
							19.38			
5	32	0.032	0.986	-3.457	36.570	1.95				
							1.43			
6	64	0.064	0.997	-3.384	37.409	1.97				
							0.60	1		
7	128	0.128	0.997	-3.355	38.188	1.99				
							0.64	1		
8	256	0.256	0.998	-3.387	39.262	1.97				*
							0.19	1		
9	512	0.512	0.998	-3.377	40.387	1.98				
							0.24	1		
10	1024	1.024	0.998	-3.365	41.869	1.98				
11	2048	2.048							Plot out	

The threshold line of MON 810 specific runs is determined as 0.256.

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Key in each threshold line determined to obtain the raw data of each run.

	SSIib	35S	
Unknown A	11833.24	0	
Unknown B	12346.86	264.36	
Unknown C	19164.82	9856.39	

The values are average values of 3 repeats of each run.

Use raw data to calculate the MON 810 percentages of unknown samples.<sup>6, 7, 8</sup>

	35S	MON 810 specific
Unknown A	0 %	0 %
Unknown B	5.490029 %	5.341439 %
Unknown C	131.8708 %	109.14695 %

**Conclusion :**

1. Sample unknown A is Non-GM sample.
2. Sample unknown B contains 5% MON 810.
3. Sample unknown C is 100% MON 810.

<sup>6</sup> Number of copies of the original recombinant gene / [ (number of copies of the original internal standard genes) × internal standard ratio] × 100

<sup>7</sup> See Testing for Foods Produced by Recombinant DNA Techniques, version 2, 2002/4/30








<sup>8</sup> The internal standard ratios of MON 810 specific and of 35S, which are used in these calculations are in compliance to 「Testing for Foods Produced by Recombinant DNA Techniques, version 2, 2002/4/30」. Actually, the internal standard ratios are specified for ABI 5700 and for ABI 7700, not for ABI 7000. The internal standard ratios of MON 810 and of 35S for ABI 7000 are already decided. However, they still are closed. So, the calculations used the internal standard ratios for ABI 5700/ABI 7700. Since the internal standard ratios for ABI 7000 are more or less the same to those for ABI 5700/ABI 7700, the judgments are correct.

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1. Eggs	Albumin
2. Milk	Casein
3. Wheat	Gliadin
4. Peanut	Arah 1
5. Buckwheat <sup>3</sup>	Small molecular proteins <sup>4</sup>

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- A. ELISA methods are available for detection of the major allergen in these five kinds of food.
- B. PCR methods have been developed for detection of the major allergen in wheat、peanut and buckwheat<sup>6</sup>.
- C. Western blotting methods are also available for detection of the major allergen in these five kinds of food.

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
---




<sup>5</sup> Methods for detection of major allergenic materials have been established. The ELISA methods should be carried on in first while performing inspections. The ELISA methods are not only methods of detection but also quantitative methods. The PCR methods and the Western blotting methods are used for confirmation.

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


## Extraction Methods Using QIAGEN Genomic-tip Kit


**For Tacos, Tortillas, Corn Chips and Cornflake, and processed potato food (restricted to products processed using heat) :**

Weight 1 g of frozen-dried sample powder, put into a 50 mL tube (polypropylene centrifuge tube)

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

### For canned papaya:

Weight 2 g of frozen-dried sample powder, put into a 50 mL tube ( polypropylene centrifuge tube )

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Translated from 「Testing for Foods Produce by Recombinant DNA Techniques」,  
Version 2.0 ( Apr. 30. 2002 ) , Japan.

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
目的：To detect if this processed potato food contains New Leaf Y<sup>3</sup>

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1. Grinding : Use a miller usually used to mill coffee beans.

<sup>1</sup> A stick shape snack made of potato. There are more than 4 kinds of flavors, サラダ (salad), うすしお味 (salty), チーズ (cheese)、ジャーマンポテト (Germany style) and so on.

<sup>2</sup> ジャがりも is potato in Japanese.

<sup>3</sup> Since New Leaf and New Leaf plus GM potato were pass the safety assessment in Japan, these two strains do not need to perform qualitative test now. The safety assessment of New Leaf Y is undergoing. Therefore, this GM potato is still illegal. However, New Leaf Y was detectable in a processed potato food previously. So, the Japan government wants to check if there are other processed potato products contain New Leaf Y.

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
2. Sampling : 2 g
3. DNA extraction ( 2 copies ) : Use QUAGEN Genomic tips kits<sup>4</sup>

Note-- add 100  $\mu$ L of  $\alpha$ -Amylase( 1 mg/mL ) after adding 7.5 mL of G2 buffer, and incubate at 37°C for 30 minutes before adding 200  $\mu$ L of Proteinase K and 20  $\mu$ L of RNase.

DNA conc. of each extraction :

	230 nm	260 nm	280nm	320 nm	260 nm / 280 nm	Conc. ( ng / $\mu$ L
No.1	0.044	0.099	0,070	0.005	1.414	99
No.2	0.048	0.107	0.077	0.008	1.39	107

定性試験 ( Qualitative PCR )<sup>5</sup> :

Detection of GM potato ( New Leaf Y and New Leaf plus ) in processed potato food, ジャがりこ---サラダ. Lane1 and lane 17 are markers. Lane 2, lane7 and lane 12 are sample 1. Lane 3, lane 8 and lane 13 are sample 2. Lane 4, lane 9 and lane 14 are Non GM DANSYAKU. Lane 5, lane 10 and lane 15 are 1% New Leaf plus. Lane 6, lane 11 and lane 16 are 1% New Leaf Y. Lane 2 to lane 6 detected 216 bp amplicons by Pss primer pair. Lane 7 to lane 11 shows the PCR result of New Leaf Y detection primer pair. Only 1% New Leaf Y ( lane 11 ) obtains a 225 bp band. Lane 12 to lane 16 shows the PCR result of New Leaf plus detection primer pair. Only 1% New Leaf plus ( lane 15 ) obtains a 234 bp band.

<sup>4</sup> Testing for Foods Produced by Recombinant DNA Techniques, version 2, 2002/4/30

<sup>5</sup> Testing for Foods Produced by Recombinant DNA Techniques, version 2, 2002/4/30

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

(A)

(B)

## Designated Processed Foods in GMO Labeling

1. Tofu
2. Frozen-tofu, yuba and other tofu
3. Natto
4. Soy milk
5. Miso paste
6. Boiled soybeans
7. Canned or bottled soybeans
8. Soy power
9. Roasted soybeans
10. Processed food made from 1~9
11. Processed food made from soybeans
12. Processed food made from soy powder
13. Processed food made from soy protein
14. Processed food made from green soybean
15. Processed food made from soybean sprout
16. Corn snack
17. Corn starch
18. Pop-corn
19. Frozen-corn
20. Canned or bottled corn
21. Processed food made from corn flour
22. Processed food made from corn grits
23. Processed food made from corn
24. Processed food made from 16~20
25. Frozen potato
26. Dried potato flakes
27. Potato starch
28. Potato snack
29. Processed food made from 25~28
30. Processed food made from potato

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## Protocol for Isolation of DNA from Plant Tissue with the DNeasy Plant Mini Kit

### Important notes before starting

- If using DNeasy Plant Mini Kits for the first time please read "Technical Information" (page 11).
- Buffers AP1 and AP3/E concentrate may form precipitates upon storage. If necessary, warm to 65°C to redissolve (before adding ethanol to Buffer AP3/E) Do not heat Buffer AP3/E after ethanol has been added
- Buffer AP1 may develop a yellow color upon storage This does not affect the procedure
- Buffers AW and AP3/E are supplied as concentrates. Before using for the first time, add the appropriate amount of ethanol (96–100%) as indicated on the bottle to obtain a working solution.
- Preheat a water bath or heating block to 65°C
- Preheat Buffer AE to 65°C
- All centrifugation steps are carried out at room temperature in a microcentrifuge

1. Grind plant or fungal tissue under liquid nitrogen to a fine powder using a mortar and pestle. Transfer the tissue powder and liquid nitrogen to an appropriately sized tube and allow the liquid nitrogen to evaporate. Do not allow the sample to thaw. Continue immediately with step 2.

Note: See "Disruption of plant material" (page 11)

2. Add 400 µl of Buffer AP1 and 4 µl of RNase A stock solution (100 mg/ml) to a maximum of 100 mg of ground (wet weight) or 20 mg (dried) plant or fungal tissue and vortex vigorously.

No tissue clumps should be visible Vortex or pipet further to remove any clumps Clumped tissue will not lyse properly and will therefore result in a lower yield of DNA In the rare case where clumps cannot be removed by pipetting and vortexing, a disposable micropestle may be used.

Note: Do not mix Buffer AP1 and RNase A prior to use

3. Incubate the mixture for 10 min at 65°C. Mix 2–3 times during incubation by inverting tube.

This step lyses the cells

4. **Add 130  $\mu$ l of Buffer AP2 to the lysate, mix, and incubate for 5 min on ice.**

This step precipitates detergent, proteins, and polysaccharides.

**(Optional) Centrifuge the lysate for 5 min at full speed.**

Some plant materials can generate very viscous lysates and large amounts of precipitates during this step resulting in shearing of the DNA in the next step (see "Lysate filtration with QIAshredder", page 12). In this case optimal results are obtained if the majority of these precipitates are removed by centrifugation for 5 min at maximum speed. After centrifugation, apply supernatant to QIAshredder spin column and continue with step 5.

5. **Apply the lysate to the QIAshredder spin column (lilac) sitting in a 2 ml collection tube and centrifuge for 2 min at maximum speed.**

It may be necessary to cut the end off the pipette tip to apply the lysate to the QIAshredder column. QIAshredder removes most precipitates and cell debris, but a small amount will pass through and form a pellet in the collection tube. Be careful not to disturb this pellet in step 6.

6. **Transfer flow-through fraction from step 5 to a new tube (not supplied) without disturbing the cell-debris pellet.**

Typically 450  $\mu$ l of lysate is recovered. For some plant species less lysate is recovered. In this case determine volume for the next step.

7. **Add 1.5 volumes of Buffer AP3/E to the cleared lysate and mix by pipetting.**

Example: To 450  $\mu$ l lysate add 675  $\mu$ l Buffer AP3/E. Reduce the amount of Buffer AP3/E accordingly if less lysate is recovered. A precipitate may form after the addition of ethanol but this will not affect the DNeasy procedure.

**Note:** Ensure ethanol has been added to Buffer AP3/E (see "Important notes before starting").

**Note:** It is important to pipet Buffer AP3/E directly onto the cleared lysate and to mix immediately.

8. **Apply 650  $\mu$ l of the mixture from step 7, including any precipitate which may have formed, to the DNeasy mini spin column sitting in a 2 ml collection tube (supplied). Centrifuge for 1 min at  $\geq 6000 \times g$  (corresponds to  $\geq 8000$  rpm for most microcentrifuges) and discard flow-through.\***

Reuse the collection tube in step 9.

9. **Repeat step 8 with remaining sample. Discard flow-through\* and collection tube.**

10. **Place DNeasy column in a new 2 ml collection tube (supplied), add 500  $\mu$ l Buffer AW to the DNeasy column and centrifuge for 1 min at  $\geq 6000 \times g$  ( $\geq 8000$  rpm). Discard flow-through and reuse the collection tube in step 11.**

**Note:** Ensure ethanol is added to Buffer AW (see page 15).

\* Flow-through fractions contain Buffer AP3/E, and are therefore not compatible with bleach

11. Add 500  $\mu$ l Buffer AW to the DNeasy column and centrifuge for 2 min at maximum speed to dry the membrane.

It is important to dry the membrane of the DNeasy column since residual ethanol may interfere with subsequent reactions. This spin ensures that no residual ethanol will be carried over during elution. Discard flow-through and collection tube.

After washing with Buffer AW, the DNeasy mini spin column membrane is usually only slightly colored. In the rare case that the membrane remains significantly colored after washing with Buffer AW, refer to "Darkly colored membrane" in the Troubleshooting Guide on page 21.

**Note:** Following the spin, remove the DNeasy column from the collection tube carefully so the column does not contact the flow-through as this will result in carryover of ethanol.

12. Transfer the DNeasy column to a 1.5 ml or 2 ml microcentrifuge tube (not supplied) and pipet 100  $\mu$ l of preheated (65°C) Buffer AE directly onto the DNeasy membrane. Incubate for 5 min at room temperature and then centrifuge for 1 min at  $\geq 6000 \times g$  ( $\geq 8000$  rpm) to elute.

Elution with 50  $\mu$ l (instead of 100  $\mu$ l) increases the final DNA concentration in the eluate significantly, but also reduces overall DNA yield. If larger amounts of DNA ( $>20$   $\mu$ g) are loaded, eluting with 200  $\mu$ l (instead of 100  $\mu$ l) increases yield. See "Elution" on page 12.

13. Repeat elution (step 12) once as described.

A new microcentrifuge tube can be used for the second elution step to prevent dilution of the first eluate. Alternatively, the microcentrifuge tube can be reused for the second elution step to combine the eluates.

**Note:** More than 200  $\mu$ l should not be eluted into a 1.5 ml microcentrifuge tube because the DNeasy column will contact the eluate.

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
# National Food Research Institute

# MISSION

National Food Research Institute (NFRI) is one of independent administrative institutions (IAI) under the supervision of the Ministry of Agriculture, Forestry and Fisheries (MAFF). The mission of NFRI is to provide the society with healthy and enriched life and secure supply of safe food through conducting research. Its food research projects cover a wide range of scientific and technical fields and are aimed at developing technological systems that strengthen the agriculture, forestry, fisheries, and food processing industries.

These research projects include: the elucidation of the relation of various foods to health, namely, functional evaluation of food components; improvement of technology for analysis; development of technology to ensure food safety; improvement of food distribution and processing methods; and exploration and utilization of functions of living organisms.

## Medium Term Plan for Research (2001-2005)

Elucidation of functional properties of food and development of technologies to utilize these properties

Development of methods of analysis in support of the food labeling

Development of technologies to ensure the safety and high quality of food

Development of technologies for utilizing various food materials, and advancement of food manufacturing techniques which meet the environmental needs of the society

Development of technologies for effectively using microorganisms and enzymes

Elucidation of biological functions of molecules through the latest biochemical methods and development of technologies for utilizing these molecules

## Cooperation, Research Planning and Organization Division for Collaboration with the Industry


- ▷ Joint research with private enterprises, various public corporations, and local governmental agencies;
- ▷ Research entrusted to the Institute by private enterprises, various public corporations, and local governmental agencies, and research entrusted to other organizations or bodies by the Institute on specific topics of the Institute's projects;
- ▷ Acceptance of trainees from private enterprises, various public organizations, and local governmental agencies, and of researchers with special status and post-doctorate fellows;
- ▷ Cooperation with international organizations, such as the United Nations University, and exchange of researchers through fellowship programs and bilateral cooperation arrangements;
- ▷ International/national seminars, workshops and symposia on a variety of topics including food-related research; and
- ▷ Tour in the Institute and explanation of its research projects for visitors upon request

# Research Activities

## Research Teams

Five research teams in the division of planning and coordination for needed research relating to foods

- ➊ Development of measures for food protection
- ➋ Food function and quality
- ➌ Promotion of genome information analysis of useful microorganisms for fermentation industry
- ➍ Creation of functionally-strengthening macromolecules based upon protein engineering, 3D analyses and bioinformatics
- ➎ Development of micro-channel arrays and their applications including blood rheology measurement


*Staphylococcus aureus*

*Enterohaemorrhagic Escherichia coli O157:H7*

## Food Function Division

The nutritional, gustatory, physical and physiological functionalities of foods and food components

- ➊ The physiological activity of food components to reduce serum lipid and prevent obesity
- ➋ Mechanism of taste response and its physiological regulation
- ➌ Physicochemical properties of food and the physiological effects on humans
- ➍ Research into functional food factors in vegetables and fruits and their effects
- ➎ Mechanism of digestion and absorption, and improving effects on brain function


Lipoprotein production and hepatic fatty acid metabolism

## Food Safety and Quality Division

Food safety and preservation of quality during distribution of food resources and foods

- ➊ Evaluation and method development of safety of food resources and foods
- ➋ Control of microorganisms and microbial toxins contaminating food resources and foods
- ➌ Quality control of vegetables and fruits based on physiological properties
- ➍ Food packaging and quality maintenance of packaged foods
- ➎ Ecology, physiology and control of food insect pests

*Aspergillus flavus*: aflatoxin-producing fungus


*Trichogramma chilonis* Ishii, the egg parasitoid of *Coreva cephalonica* Stanton

## Food Material Division

Biochemical and structural properties of food ingredients and development of new technologies for processing of agricultural products

- ➊ Development of new food carbohydrates and related enzymes
- ➋ Evaluation of quality and function of food proteins
- ➌ Elucidation of structure and function of hydrophilic macromolecules
- ➍ Clarification of factors affecting rice quality and development of technology for quality evaluation and utilization
- ➎ Development of new food manufacturing technology of cereals


The differences between *Japonica* and *Indica* rice grains

## Food Engineering Division

Food processing, distribution, nano-measurement, information technology and other technologies

- Development of advanced food processing and waste treatment system
- Development of novel dispersion, separation and reactor system
- Development of novel analyzing methodology for nano-structure and function of foods and biomaterials
- Development of technology such as soft electron, quality measurement and network application system
- Development of environmentally-friendly system for high quality food distribution


Microchannel(MC) Emulsification Mechanism

## Applied Microbiology Division

Molecular analysis of food microorganisms and application of microbial functions in food processing

- Molecular biological analysis of the genes expressed in Koji mold used in Miso, Soy sauce and Sake production
- Molecular and cellular biology of baker's yeast and lactic acid bacteria
- Molecular analysis of the bacterial  $\gamma$ -polyglutamic acid (the sticky polymer of Natto) production
- Development of enzymatic system for the production of useful amino-sugars


Metabolism of  $\gamma$ -polyglutamic acid, the sticky polymer of Natto

## Biological Function Division

Engineering and biotechnology of biomolecules and cells

- Macromolecular structure and catalytic mechanism
- Improvement of enzyme function by gene manipulation
- Secondary metabolism, and ribosome function of actinomycetes and *Bacillus subtilis*
- Cell metabolism and regulation


Tertiary structure of a protease exhibiting improved heat stability, constructed by gene replacement

## Analytical Science Division

Studies on food analysis and quality evaluation

- Development of sensitive and precise analytical methods and procedures to determine food components
- Analysis of structure and dynamics of food components at molecular level in relation to their function
- Development of nondestructive methods for bio-measurement and food analysis
- Development of spectral analysis and image analysis methods for quality evaluation of foods
- Collection and preparation of data for "New Standard Tables of Food Composition in Japan"

High performance instrument for the analysis of molecular structure and dynamics

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