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
**piezosystem jena**

**product  
news**

**NANOSTEP WILL BE THE FUTURE - WALK WITH US**

## Piezoelectric Micro Objective Positioning System **MIPOS 5**

- piezo focus adjustment
- very compact design
- 500  $\mu\text{m}$  adjustment range
- parallelogram construction
- high resonant frequency
- easy to mount on microscopes

The MIPOS 5 was developed for fine adjustment of micro objectives. The maximum motion is 500  $\mu\text{m}$ . The resolution of MIPOS 5 is very high and, in practice, only limited by the voltage noise of the power supply.

A special integrated pre-load design results in the following advantages:

- very small angular deviation
- highly parallel large range of motion
- high resonant frequency

On the top and bottom side of MIPOS 5 all standard threads are available from W 0.8 x 1/36" up to M27 x 0.75. Mounting this system on the microscope is very easy. You screw a thread-ring into the microscope and mount the MIPOS 5 on this ring with a clamping-screw.

recommended amplifier control system: NV 40 / 1 CL (see page 11)

### technical data


motion ( $\pm 10\%$ )	500 / 400 $\mu\text{m}$ (open / closed loop)
operating voltage	-10... +150 V
capacitance ( $\pm 20\%$ )	21 $\mu\text{F}$
resonant frequency ( $\pm 20\%$ )	200 Hz (with micro objective $m = 80 \text{ g}^*$ ) 180 Hz (with micro objective $m = 100 \text{ g}^*$ )
connector	LEMO
length of cable	1 m
option	strain gauge measurement system

\* small signal value

Please note: additional information for dynamical use after request

**Piezo Actuator series PA / T  
with M14 outside threading  
PA35/T14; PA50/T14; PA80/T14**

**Piezo Actuator series P / S  
hermetic insulated  
P 50/S13**

- outside thread for mounting
- versatile use
- excellent dynamical behavior by integrated pre-load
- very small and compact design
- high force generation
- high stiffness

These actuators are mechanically pre-loaded to make them ideal for dynamical applications. They can handle loads up to 1000 N. The housing of the actuator in combination with the springs allows a motion without mechanical play. The M14 outside threading is used to mount the actuator. It also helps to adjust the position precisely. Depending on the chosen model maximum motions from 35  $\mu\text{m}$  up to 80  $\mu\text{m}$  are available. The operating voltage is -10 V...+150 V.

- universal use
- independent from environmental conditions
- very small and compact design
- stable against high pressure forces

The piezo actuators of the series P /S13 have a hermetically sealed housing. The special housing allows its use in aggressive and humid environments. The actuators are not pre-loaded. The construction of the housing tightly encloses the stack, thereby avoiding mechanical play. They can handle a load up to 1000 N. The M3-thread in the bottom plate allows an easy adaptation in existing systems and guaranties an exact positioning. Depending on the type, actuators with a stroke of 8  $\mu\text{m}$  up to 50  $\mu\text{m}$  are available. The operating voltage is -10V...+150 V.


**technical data**

motion ( $\pm 10\%$ )	35 / 50 / 80 $\mu\text{m}$
voltage	-10...+150 V
max. force	1000 N
stiffness ( $\pm 20\%$ )	30 / 20 / 12.5 N/ $\mu\text{m}$
capacitance ( $\pm 20\%$ )	3.6 / 5.4 / 7.2 $\mu\text{F}$
length	53 / 71 / 89 mm
diameter	14 mm
thread	M 14
spanner flats	13 mm
pre-load ( $\pm 10\%$ )	300 N

motion ( $\pm 10\%$ )	50 $\mu\text{m}$
voltage	-10...+150 V
max. force	1000 N
stiffness ( $\pm 20\%$ )	20 N/ $\mu\text{m}$
capacitance ( $\pm 20\%$ )	4.5 $\mu\text{F}$
length	71 mm
diameter	max. 13.3 mm
thread (bottom plate)	M 3
spanner flats	12 mm


## Piezo-Positionierer PXZ 400

- motion 400  $\mu\text{m}$  in x-z axis
- compact and stiff design
- high resonant frequency

## Piezo-Positionierer PZ 400

- motion: 400  $\mu\text{m}$  in z-axis
- compact design
- high resonant frequency

These new positioning systems allow a parallel motion of 400  $\mu\text{m}$  in the X- and Z- axes. With proper control electronics, a resolution in the nm range can be reached. An innovative new design combines a larger motion with a high resonant frequency. The compact size allows an easy installation in existing applications. When strain gauge (SG) position sensors are built in a drift- and hysteresis-free positioning is guaranteed. Vacuum and low temperature options are available.

recommended amplifier control system: NV 40 / 1 CL (see page 11); ENV 19" modular amplifier system

### technical data

part-no.	T-209-00 / T-209-01*	T-118-00 / T-118-01*
motion per axis ( $\pm 10\%$ )	400 / 320 $\mu\text{m}$	400 / 320 $\mu\text{m}$
capacitance per axis ( $\pm 20\%$ )	14 $\mu\text{F}$	4.6 $\mu\text{F}$
voltage	-10 V...+150 V	-10 V...+150 V
stiffness x/z axis ( $\pm 20\%$ )	0.35 / 0.3 N/ $\mu\text{m}$	0.25 N/ $\mu\text{m}$
resonant frequency x/z axis ( $\pm 20\%$ )**	300 / 180 Hz	280 Hz
connector	LEMO	LEMO
cable length	1 m	1m
size (l x w x h)	80 x 80 x 35,5 mm	48 x 52 x 28,5 mm
weight ( $\pm 10\%$ )	685 g	190 / 195 g


\* strain gauge version

\*\* small signal value

Please note: additional information for dynamical use after request.

## Nano Positioning System **PENTOR - 5 axis**

- combination of a three axes translation system and a two axis tilting system
- free center hole with 17 mm diameter
- integrated preload
- translation each axis: 100  $\mu\text{m}$
- tilt each direction: +/- 2.5 mrad
- solid state flexure hinges in parallelogram design without mechanical play
- integrated measurement system

For fiber alignment optical stages often require an open central space (e.g. for passing light). For such applications, the new **PENTOR** (T-405-00) translation and tilting system was developed with a 17 mm free center hole. It offers a motion in XYZ of 100  $\mu\text{m}$  and a tilt of +/-2.5 mrad on two orthogonal axes.

Flexure hinges for the three translation axes guarantee the highest degree of parallelism over the travel range of the stage. Each axis has a built in integrated mechanical preload, making the new **PENTOR** very well suited for dynamic applications.

The tilting part is designed for "plus-minus tilting". The construction is temperature compensated: changes in the surrounding temperature do not affect the tilting angle. As an option, the **PENTOR** can be offered with an integrated strain gauge measurement system part.no. T-450-01.

### technical data

#### **PENTOR**

##### **Whole Actuator:**

operating voltage  
measurement system  
weight ( $\pm 15$  g)  
dimensions (length x width x height)  
free center space  
thread pattern  
cable length  
connector

##### **Translation:**

motion x,y,z ( $\pm 10\%$ )  
capacitance each axis ( $\pm 20\%$ )  
resonant frequency x/y/z-axis ( $\pm 20\%$ )\*\*  
max. load  
stiffness each axis ( $\pm 20\%$ )  
temperature range

##### **Tilt:**

number of axes  
tilt ( $\pm 10\%$ )  
pivot points  
capacitance each axis ( $\pm 20\%$ )

##### **Vibration options:**

#### **T-450-00**

-10...+150 V  
-  
450 g  
68 x 69 x 40 mm  
17 mm  
50 mm M4  
1m (39 inches)  
LEMO

100  $\mu\text{m}$   
1.7  $\mu\text{F}$   
680 / 750 / 580 Hz  
5 N  
0.5 N/ $\mu\text{m}$   
-20°C to 80°C

2  
 $\pm 2.5$  mrad  
center of the stage, 4.5 mm below top plate  
0.85  $\mu\text{F}$   
vacuum, low temperature option

#### **T-450-01\***


-10...+150 V  
strain gauge  
450 g  
68 x 69 x 40 mm  
17 mm  
50 mm M4  
1m (39 inches)  
LEMO

100  $\mu\text{m}$  (80  $\mu\text{m}$ \*)  
1.7  $\mu\text{F}$   
680 / 750 / 580 Hz  
5 N  
0.5 N/ $\mu\text{m}$   
-20°C to 80°C

\* closed loop version; \*\* small signal value  
Please note: additional information for dynamical use after request.

## Nano Positioning System TRITOR 400 / 400 CAP

- motion up to 400  $\mu\text{m}$  per axis
- high resonant frequency
- max. load 50 N
- free hole 19 mm

New in the family of well-known compact, multi-dimensional translation stages TRITOR is the serie TRITOR 400 and TRITOR 400 CAP. The TRITOR 400 was designed for large, highly parallel motion (up to 400  $\mu\text{m}$ ) in all three axes. Potential applications for this system include: optics and laser tuning, fiber positioning, micro-positioning and scanning. The TRITOR 400 incorporates a highly compact design, parallel motion due to its parallelogram design, and high reliability because of the solid state hinges. Due to its capacitive measurement system (CAP), the PX 400 CAP translation stage offers resolutions in the low nm range. As an option, all TRITOR 400 may be specially prepared for vacuum and/or low temperature applications.

All TRITOR stages consist of actuators integrated into a housing with an internal lever transmission. This housing is constructed of stainless steel, in combination with anodized aluminum, which results in a rugged yet highly compact design.

Selected components can be mounted on the top plate by four tapped holes and can be accurately located by using the precision pin holes.

### technical data

part. no.:	TRITOR 400 T-406-00	TRITOR 400 CAP T-406-06
motion x,y,z ( $\pm 10\%$ )	400 $\mu\text{m}$	400 / 320 $\mu\text{m}$ (open/closed loop)
integrated measurement system	-	capacitive
operating voltage	-10...+150 V	-10...+150 V
resolution	0.53 nm*	1 nm
capacitance each axis ( $\pm 20\%$ )	14 $\mu\text{F}$	14 $\mu\text{F}$
max. load (statical use)	10 N	10 N
temperature range	-20°C bis 80°C	-20°C bis 80°C
connector	LEMO	LEMO
cable length	1 m	1 m
dimensions (length x width x height)	116 x 106 x 40 mm	116 x 106 x 40 mm
free center space	19 mm	19 mm
thread pattern	32 mm x 32 mm M3x6	32 mm x 32 mm M3x6

#### options

vacuum option, low temperature option


\*measured with E-103-18 amplifier

## Fiberswitch Multimode 1x1 to 1x2; 1x3 to 1x9

The simplest way for the creation of optical fiber switches is to achieve a direct movement of the optical fiber ends with respect to each other. Piezoelectric actuators are ideal candidates for this purpose due to their typical movement ranges and actuation speed. In order to exclude the influences of hysteresis and creep on the accuracy of such a switch, the actuator end positions can be fixed by special stopping configurations. In the case of multimode fibers, the required accuracy (10 microns or more depending on the fiber core diameter) can be achieved by different approaches of precision engineering.

Multimode fiber switches are used mainly in optical metrology applications and optical quality control systems. They are ideally suited to combine up to 9 measurement points with only one spectrometer. Therefore, the end-user derives a cost benefit, and is able to directly compare different optical channels using only one measurement equipment.

- fast switching time
- low insertion loss
- high optical isolation
- compact design
- no additional wavelength dependence

### technical data

#### OPTICAL CHARACTERISTICS

insertion loss	typ.* 0.7 dB (max. 1.2 dB) for 1 x 1 to 1 x 3 switch
	typ.* 1.4 dB (max. 1.8 dB) for 1 x 4 to 1 x 9 switch
cross talk	typ.* -60 dB (max. -55 dB)
repeatability (1000 cycles over 2h)	typ.* 0.02 dB (max. 0.03 dB)

#### TECHNICAL CHARACTERISTICS

switching time	typ. 2 ms (max. 3 ms)
lifetime	10.000.000 cycles
operating temperature	0°C to 60°C
storage temperature range	-40 to +80°C
humidity	55 % RH, >55 % RH-special version

#### ELECTRICAL CHARACTERISTICS

operating voltage	5 V DC / 7-12 V DC**
power consumption	100 mA (f = 0 Hz 50 mA ; f = 50 Hz 300 mA)
control signal	binary code (BCD) / 5 V - TTL

\* depends on fiber core diameter  
 \*\* 1x1, 1x2, 1x3 to 1x9

### types of switches

no. of channels	part no.	housing L/B/H	fiber core diameter
1 x 1	F-101-04	85/105/44 mm	< 100/140 µm
	F-101-04	85/105/44 mm	200/220 µm
	F-101-05	175/105/44 mm	400/440 µm
	F-101-05	175/105/44 mm	600/660 µm
1 x 2	F-102-04	85/105/44 mm	< 100/140 µm
	F-102-05	175/105/44 mm	200/220 µm
	F-102-05	175/105/44 mm	400/440 µm
	F-102-05	175/105/44 mm	600/660 µm
1 x 3	F-103-05	175/105/44 mm	< 100/140 µm
	F-103-03	225/105/44 mm	200/220 µm
1 x 4	F-104-05	175/105/44 mm	< 100/140 µm
	F-104-03	225/105/44 mm	200/220 µm
1 x 6	F-106-05	175/105/44 mm	< 100/140 µm
	F-106-03	225/105/44 mm	200/220 µm
1 x 9	F-109-05	175/105/44 mm	< 100/140 µm
	F-109-03	225/105/44 mm	200/220 µm

### ordering instructions

Please note: when ordering you will need the part number of the switch and also the part no. of the fiber (which includes the optical connector). For each switch you need one input fiber and various numbers of output fibers depending on the switching channels. All fibers must be assembled with a connector, otherwise piezosystem jena can not determine the optical parameters of insertion loss and back reflection. For custom configuration, piezosystem jena retains the right to cancel an order if unforeseen additional efforts are needed which may cause the cost of the systems to exceed the quoted cost. This is related especially to special fiber types piezosystem jena has never worked with or which are supplied by our customers. Standard fiber length will be 1m (±15%). Cost for additional fiber length will be added.

piezosystem jena is using the following types of fibers and connectors for all standard applications:

fiber [µm]		index-profile	numerical aperture	wavelength range	part no. for fiber (without optical connector)
core	cladding				
50	125	graded-index	0.22	850/1300 nm	C-319-*
62,5	125	graded-index	0.28	850/1300 nm	C-329-*
100	140	graded-index	0.29	850/1300 nm	C-339-*
100	110	step-index	0.22	180-1100 nm	C-230-*
100	140	step-index	0.22	350-2600 nm	C-130-*
200	220	step-index	0.22	350-2600 nm	C-140-*
200	220	step-index	0.22	180-1100 nm	C-240-*
400	440	step-index	0.22	350-2600 nm	C-160-*
400	440	step-index	0.22	180-1100 nm	C-260-*
600	660	step-index	0.22	350-2600 nm	C-180-*
600	660	step-index	0.22	180-1100 nm	C-280-*


\* when ordering please use the end number of the optical connectors as follows:

ST = end number -10; SMA = end number -20; FC/PC = end number -30; FC/APC = end number -50

Additional fiber extension or a reduction of fiber length available on request.

## MICI 180 + Micrometer Stages (application sample)

- motion: 180  $\mu\text{m}$
- parallelogram construction allows to increase the resolution of micrometer screw driven stages
- integrated strain gauge measurement system

picture: MICI 180 mounted on a standard Newport stage. Additional modifications of the stage are not necessary.

This element was developed for quality control in the optical industry.  
The piezostage is joined with micrometer drives to be mounted into a micrometer positioning system.  
At first the position of the actuator can be adjusted by the micrometer drive and then this position can be fine-tuned with the aid of the piezostage.  
The resolution is in the nm and sub-nm range.

### technical data

part no.	S-704-00	S-704-01
motion ( $\pm 10\%$ )	180 $\mu\text{m}$	150 $\mu\text{m}$
operating voltage	-10 ... +150 V	-10 ... +150 V
capacitance ( $\pm 20\%$ )	1.7 $\mu\text{F}$	1.7 $\mu\text{F}$
resolution	0.1 nm (open loop)*	2 nm (closed loop)**
temperature range:	-20 to 80 °C	-20 to 80 °C
connector:	LEMO	LEMO
length of cable:	1 m	1 m



\* measured with E-100-18 amplifier  
\*\* measured with E-240-100 amplifier

**NV 40 / 1 CL**  
*voltage amplifier with integrated  
 strain gauge measurement  
 amplifier and PID controller*

**12 V 40**  
*OEM voltage amplifier  
 input voltage 12V / 24V*

- stand alone casing
- compact dimensions
- excellent cost effectiveness
- low noise amplifier (0.25 mV RMS)
- RS 232 interface

- compact size
- 12 V or 24 V power supply
- excellent cost effectiveness
- small dimensions

The module **NV40 / 1 CL** was especially developed for low voltage piezoelements equipped with strain gauge sensors. By using the analog input the amplifier can control by an analog voltage of 0...+10V. There is also a BNC socket to control the output voltage. Special circuits suppress voltage peaks to avoid damage of the actuator. The system has very low noise characteristics of less than 0.25 mV RMS and thus provides resolutions in the sub-nm range. The closed loop system will be realized by an integrated PID controller. Please note the compact dimensions of the amplifier!

The voltage amplifier **12V40** is well suited for low voltage piezo elements. Two different types of casing are available. The amplifier requires only 12 V DC or 24 V DC for supply. The amplifier can be controlled via MOD input (front- or backside). The output voltage can be supervised on the monitor output. Special circuits are integrated to protect the piezo element from voltage peaks and excessive voltages. Due to very low voltage noise of only 0.3 mV RMS this amplifier is well suited for sub-nm positioning tasks. All sockets are also available on the backside so this amplifier can be easily integrated as an OEM-product.

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
## Long-Travel and High Load Motorized Actuators

**NEW**  
FROM NEWPORT

### Key Features

- Up to 50 mm travel in a space saving design
- Non-rotating tip improves motion smoothness and has no wear
- Exceptional position sensitivity provided by low-friction design and ultra-high resolution encoder feedback
- Adjustable limit switch prevents damage from over-travel
- Convenient manual positioning knob

The LTA actuators are the successor of our popular 850G and VM series and provide up to 50 mm motorized travel in a compact package. LTAs are designed to fit into your existing manual stages and other opto-mechanical components as a direct replacement for manual micrometers.

The LTA series features a space-saving design with the motor and lead screw side-by-side. This reduces the actuator length by 50% and minimizes the negative effects of long cantilever loads on micro-positioning equipment. The non-rotating tip prevents wear and avoids periodic motion variations which can be caused by variable contact made between a rotating tip and the mating surface.

Precision motion is accomplished by a miniature DC servo motor with an optimized output torque. This allows for faster motion with higher load capacity. The LTA-HL is recommended for heavy load applications.

It features a strong, 8 mm diameter rod and a M12-0.5 mounting bezel that is compatible with our (M-)UMR8 and (M-)MVN80 linear stages and SK and SL series optical mounts. The LTA-HS is optimized for high-speed applications and provides the longer travel range. The mounting interface of the LTA-HS is compatible with a large number of Newport and others manual components.

A movable limit switch prevents equipment from over-travel. The switches position can be changed in minutes to adjust the maximum travel position. A manual adjustment knob permits quick positioning of the actuator while the motor is off. Convenient laser marked scales indicate coarse actuator position in both millimeters and inches.

  
**Newport.**  
Experience | Solutions

## Design Details

Base Material	Stainless steel body with aluminum cover
Drive mechanism	Non-rotating lead screw (rotating nut)
Drive screw pitch (mm)	1.0
Reduction gear	LTA-HS: 1:14 LTA-HL: 1:66
Feedback	Motor mounted rotary encoder, 2,048 cts/rev.
Limit switches	Mechanical switches, both ends, max. travel limit is adjustable
Origin	Uses min. travel limit for homing, typically <20 µm repeatability
Motor	DC-servo motor UE1724SR
Cable length (m)	3
Vacuum compatibility	Vacuum compatible versions are available up to 10 <sup>-6</sup> Torr
MTBF	10,000 h at 25 N load and a 10% duty cycle
Weight	LTA-HS: 0.38 kg LTA-HL: 0.41 kg

## Motion Controller Options

For optimum performance and seamless compatibility, we recommend using one of the following Motion Controllers/Drivers:

XPS

ESP6000 with UNIDRIV6000

ESP300

ESP100

## Specifications

Caution: Values in the table below are indicated for actuators operating with the rod in -Cx direction.

	LTA-HS	LTA-HL
Travel (mm)	50	25
Resolution (µm)	0.035	0.0074
Minimum Incremental Motion (µm)	0.1	0.05
Uni-directional Repeatability (µm) guaranteed	0.5	0.5
Bi-directional Repeatability (µm) <sup>(1)</sup> guaranteed	2	2
On-Axis Accuracy (µm) guaranteed	15	10
Maximum Speed (mm/s)	5	1

<sup>(1)</sup> After backlash compensation: value provided with each actuator.

## Load Characteristics

+Cx Direct Axial Load Capacity (N)  
-Cx Inverse Axial Load Capacity (N)

Caution: To reach the specifications stated for LTA actuators, a minimum axial load must be applied at the end of the rod.



	LTA-HS	LTA-HL
Direct Axial Load Capacity +Cx (N) <sup>(2)</sup>	40	100
Inverse Axial Load Capacity -Cx (N) <sup>(2)</sup>	50	120
Minimum Axial Load Capacity (N)	2.5	5
Side Load Capacity (N) <sup>(3)</sup>	5	20

<sup>(2)</sup> Do not apply alternatively +Cx and -Cx loads during an operating cycle.



<sup>(3)</sup> No side loads during motion.

## Dimensions


### LTA-HL & LTA-HS Actuators

### LTA-M12 Adapter

SL mount with LTA actuators.

The LTA Series actuator is compatible with most standard Newport manual stages and opto-mechanical components like the M-436 Series translation stage.

## Ordering Information

Description	Model
Motorized actuator, high speed version, 0-50 mm travel	LTA-HS
Motorized actuator, high load version, 0-25 mm travel	LTA-HL
LTA-HS Actuator, vacuum compatible to 10 <sup>-6</sup> Torr *	LTA-HSV6
LTA-HL actuator, vacuum compatible to 10 <sup>-6</sup> Torr *	LTA-HLV6
M12-0.5 mounting adapter for LTA-HS to M-UMR80, M-MVN80, SL and SK products	LTA-M12

\* Available February 2004.

ISO 9001  
No. 2727

Newport Corporation, Irvine, California, has been certified compliant with ISO 9001 by the British Standards Institution.


Belgium	France	Germany	Italy
Newport BV filiaal België	MICRO-CONTROLE	Newport GmbH	Newport/Micro-Controle Italia
Tel +32-(0)11 6402927	Tel +33-(0)1-60-91-68-68	Tel +49-(0)6151 / 3621-0	Tel +39-(0)2-929-0921
Fax +32-(0)11 6402227	Fax +33-(0)1-60-91-68-69	Fax +49-(0)6151 / 3621-52	Fax +39-(0)2-923-2448
<b>Netherlands</b>	<b>Sweden/Denmark</b>	<b>United Kingdom/Ireland</b>	
Newport B.V.	NPT Instruments AB	Newport Ltd.	
Tel +31-(0)30 6592111	Tel +46-(0) 8591 18844	Tel +44-(0)1635 521757	
Fax +31-(0)30 6592130	Fax +46-(0) 859440423	Fax +44-(0)1635 521348	

Visit Newport Online at [www.newport.com](http://www.newport.com)

LTA-DS (10/03)

## Integrated Motion Controller/Driver System

Preliminary  
Data Sheet

Newport's XPS is a high-performance, integrated motion controller/driver system offering outstanding trajectory accuracy and exceptional programming functionality. It combines operational simplicity with advanced features to precisely control the most diverse motion sequences and synchronize them via measurement, command, or external data acquisition tools. Furthermore, sporting a real-time, multi-tasking functionality, XPS can also operate as a process controller to run your application.

XPS supports any Newport positioner and many other motion devices, whether they use stepper, DC brush, DC brushless rotary or linear, voice coil or piezoelectric actuation.

A unique motion profiler auto-adapts the jerk time depending on move distance and greatly reduces the excitation of the system's mechanical resonances and stress, resulting in faster settling, more accurate trajectories and a longer life time—all that without sacrificing the move time.


XPS is based on a 10/100 Base-T Ethernet communication link with TCP/IP protocol and uses a web site approach for all software tools and FTP for data transfer. This makes it almost independent from the operating system of the user. When networked, Unix-, Linux-, or Windows-users can access the same controller from any place in the world for remote control, code development, data transfer or diagnostics.

XPS also excels as a process controller for running your applications. Based on a real-time, multi-tasking functionality, XPS is able to execute user-defined applications in real time on the motion controller using TCL scripts.


For applications that require analog data analysis in exact relation to the position, XPS's gathering mode captures all important axis information and I/O's with a time jitter of less than 50 ns, storing data in a custom-configured table.

### Key Features


- High performance, 1-8 axes motion controller/driver
- High-speed 10/100 Base-T Ethernet communication interface, typ. 0.3 ms command execution
- TCP/IP communication protocol allows for multi-tasking and multi-user applications from most operating systems: e.g. Unix, Linux, and Windows
- Powerful, intuitive and object-oriented command language with multi-parameter API's
- Highly developed motion trajectories including continuous path contouring, splines and PVT mode
- Advanced servo loop with variable PID's, low-pass and notch filters
- Broad range of compensation features including 3D mapping
- Hardware position latch and position compare output for exact synchronization with less than 50 ns latency
- Extensive I/O functionalities including four uncommitted analog inputs, four uncommitted analog outputs and 30 TTL in- and outputs
- High-speed data acquisition at up to 10 kHz rate
- Complete library of LabView VI's
- P4 2 GHz motion processor manages various important tasks

XPS controller configuration is supported by a convenient software tool. When used with Newport stages, all settings can be made in just a few clicks.

The command screen lists all APIs by groups, including all necessary or available parameters. It is a convenient tool to learn the XPS syntax, allowing for simple programming, code testing and debugging.

ISO 9001  
PM 07207

Newport Corporation, Irvine, California, has been certified compliant with ISO 9001 by the British Standards Institution.

**Belgium**  
Newport BV filiaal België  
Tel: +32-(0)11 6402927  
Fax: +32-(0)11 6402227

**Netherlands**  
Newport B.V.  
Tel: +31-(0)30 6592111  
Fax: +31-(0)30 6592120

**France**  
MICRO-CONTROLE  
Tel.: +33-(0)1-60-91-68-68  
Fax: +33-(0)1-60-91-68-69

**Sweden/Denmark**  
NPT Instruments AB  
Tel.: +46-(0) 859118844  
Fax: +46-(0) 859440423

**Germany**  
Newport GmbH  
Tel.: +49-(0)6151 / 3621-0  
Fax: +49-(0)6151 / 3621-52

**United Kingdom/Ireland**  
Newport Ltd.  
Tel.: +44-(0)1635 521757  
Fax: +44-(0)1635 521348

**Italy**  
Newport/Micro-Controle Italia  
Tel.: +39-(0)2-929-0921  
Fax: +39-(0)2-923-2448

Visit Newport Online at [www.newport.com](http://www.newport.com)  
XPS-DS (10/03)

# NAIS

Quality control by...

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Imagechecker

**Machine Vision Systems in the  
packing and food industries**

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## **About this booklet ...**

This booklet from Matsushita shows you just a few possibilities to check the quality of products by machine vision systems.

Matsushita offers a large range of machine vision systems adapted to different applications : the systems A100 and A200 with their compact size and multifunctional performance, the M100 and M200 series for special applications such as character verification and high speed measurement and the P400 system for very complex applications with up to nine cameras.

The philosophy of parametrical systems is valid for all Matsushita machine vision systems.

### **What are Imagecheckers ?**

It is a camera based on image checking system. This means each product manufacture, assembled, packed or printed on a production line will be photographed and the image will be analysed. Depending on parameters defined with a „good“ product, the calculator can determine if a product is conform or not. In the case of no conformity the Imagechecker will send a signal that permits to eject automatically all no conform products. The time between a picture taken and the indication of a no conform product can go down to 20 ms. In other words, even 3000 products per minute or 180,000 products per hour can be checked. Of course the more criteria you will observe on your product, the more time the calculator needs.

For your future „IMAGECHECKER“ applications we added a questionnaire on the last page which will help us to carry out feasibility tests with your samples.

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Matsushita offers a large range of machine vision systems which can respond to your different applications. Our smallest units are the A100 and A200 with a size of 100 mm x 50 mm. They only use a space of 14 mm on the fixing rail in the electrical cabinet and can store up to 96 different types of programs. All parameters can be set with a small keypad.

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This series offer different program packages for special applications. The package OCV (character verification) for example has been developed to inspect the printing of dates, logos or serial numbers. It can detect missing, incorrect and unreadable printings. The mould watcher package, another special development, prevents a double injection of plastic inside a mould.

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This is our PC based machine vision system, the most powerful in our range. Up to 9 cameras can be handled on the same system and the number of programs is not limited. You can generate automatically files, indicating the quality of your products and use also high resolution cameras. The NAI's vision software runs under Windows NT on an industrial PC.

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
Due to the high importance of the light source in most of the machine vision system applications, Matsushita also offers a wide range of different lighting systems. Depending on your application and the test done in our laboratory, we offer back lights or ring lights with high frequency fluorescence and LED lighting in different sizes, colors and forms. Colour and polarised filters are used for getting the best contrast to carry out the optical control.

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The world's leading CMOS and CCD image sensors, cameras, and components.

# Vision for Machines

**DALSA**  
technology with vision

# We're passionate about innovation.

Being passionate about innovation has been our formula for success for over 20 years. It's transformed us into one of the largest and most successful companies in the imaging industry. Yet, innovation applies to more than just our products and technology. It's also an important part of how we think and manage our company, and has helped us grow to over 650 employees across seven international locations.

## We push the limits of performance


DALSA has a reputation for leading the industry. For the performance factors that matter most in image capture — speed, resolution, dynamic range, and responsivity — few companies can match our long list of industry firsts:

- 30MHz operation (1985)
- Commercial Time Delay and Integration (TDI) products (1988)
- 25 million pixel devices (1991)
- 1GHz total throughput (1992)
- Scientific quality, mega-pixel camera with real time video rates (1995)
- TDI with antiblooming (1998)
- 8 million pixel 60fps digital cinema device (1999)
- 1200fps VGA CMOS device with non-rolling shutter (2000)
- 19.2Gigabit/sec Camera (2001)

With design groups in Waterloo, Colorado Springs, Tucson and Eindhoven, and a specialized wafer foundry in Bromont, we continue to build on our leadership in innovation. We own dozens of patents for the design of CCD and CMOS image sensor chips and digital camera as well as dozens more for fabricating the semiconductors that drive them. The net result is digital imaging solutions that meet the needs of the imaging industry's most challenging applications.

## We're here for the long term

At DALSA, we pride ourselves on having a strong and resilient foundation and an approach that is focused and deliberate. As a team, we will continue to build on these attributes as we grow in scope, revenues and profitability over the coming years. If you're looking for a strong partner for the long run, insist on a DALSA solution.

# Line Scan Imaging

**D**ALSA line scan cameras have a strong reputation because they're fast and they're highly responsive. Whether you're inspecting printed circuit boards, scanning documents or reading bar codes, the ability to run at the highest possible line rate without running out of light and still maintaining dynamic range defines the success of your application.

## Faster is definitely better

We recognize your need for speed by offering the industry's fastest line scan cameras. From single output models running at 40MHz to multiple output models running at up to 200MHz, you can count on a DALSA line scan camera to provide the line speed your application needs. Take our **Piranha2** cameras; at a combined data rate of 160MHz, Piranha2 is the fastest standard line scan camera around. It also requires less light than our original Piranha camera, even at these high speeds.

## A brighter idea

Our line scan cameras are the most responsive in the industry. Our **Advanced** line scan and **High Sensitivity** products, for example, are designed to provide up to 100 times more responsivity than comparable standard line scan products — all at the highest possible line rates.

If you're looking for the highest performance with the greatest selection of models, look to DALSA for your next line scan application.

*Piranha2 is the fastest standard line scan camera around with a total data rate of 80/160MHz (2 or 4 outputs). It's also extremely responsive, actually requiring less light than our original Piranha camera, even at a 60% greater data rate.*

## Standard Line Scan

- Data rates to 160MHz
- Line rates up to 79kHz
- Resolutions up to 8192 pixels
- 1 to 4 outputs
- 8 cameras (one with color) and 4 sensor products to choose from

## Advanced Line Scan

- Ultra-high responsivity (up to 100 times that of standard line scan)
- Pixel-by-pixel flat field correction on some models
- Data rates up to 40MHz
- Line rates up to 64.1kHz
- Resolutions up to 2048 pixels
- Single output

## High Sensitivity Line Scan

- Ultra-high responsivity
- Data rates up to 200MHz
- Line rates up to 83kHz
- Resolutions up to 4096 pixels
- 3, 4 and 8 outputs
- 4 cameras (one in color) and 1 sensor product to choose from

# Area Scan Imaging

If you're looking for high performance area scan, we've got it covered — from feature-rich, fast-shuttering megapixel models to scientific quality cameras with true 12-bit dynamic range.

## “Stop Action” for smear-free capture of fast moving objects

For the growing number of applications that require mega+ resolution combined with electronic shuttering capabilities, we're proud to offer the new DALSTAR “Stop Action” family of area cameras. These new cameras are designed to “stop the action” and provide smear-free capture of fast moving objects.

Drawing on the best available CMOS and CCD technologies, DALSTAR “Stop Action” means best in class performance for 100% PCB inspection, industrial metrology, traffic management and robotics. Features include 1024 x 1024, 1600 x 1200 or 1920 x 1080 pixel resolutions with 8 or 10 bit dynamic range and CameraLink™ interface. In select cameras we've included exclusive features such as LINLOG™ response which provides tremendously wide dynamic range (up to 120dB) — perfect for imaging dark scenes with very bright elements, such as in night-time license plate and other traffic management imaging.

## Scientific-grade for superior resolution and ultra-high image fidelity

While the majority of machine vision applications use RS-170/CCIR resolution CCD imagers, a growing number require the superior spatial resolution and higher image fidelity of scientific-grade cameras. Higher resolution, higher speed and higher dynamic range are the keys to success in numerous mission-critical electronics, semiconductor, flat panel, and medical imaging applications. DALSTAR cameras meet this need with mega and multi-mega-pixel resolutions, ultra high speed and true 12-bit dynamic range. DALSTAR cameras truly set the new standard for high performance.

*The DALSTAR 2M30-SA is our newest “Stop Action” area scan camera. At 30 frames per second, it features a two-megapixel interline transfer CCD with either 1600 x 1200 or 1920 x 1080 pixels — twice the resolution of “standard” 1k x 1k cameras.*

## “Stop Action” Electronic Shuttering


- Fast electronic shutter for smear-free image capture
- 1024 x 1024, 1920 x 1080 or 1600 x 1200 pixels
- Up to 75 fps
- 8 or 10 bit
- LINLOG™ and windowing on CMOS models
- 3 cameras (one with color option) and 1 sensor product to choose from

## 8-bit Frame Transfer

- 128 x 128 to 1024 x 1024 pixels
- Up to 955 frames per second
- High Responsivity
- 3 cameras and 1 sensor product to choose from

## 12-bit Frame Transfer/ Full Frame

- 1024 x 1024 to 3072 x 2048 pixels
- True 12 bit dynamic range with CDS
- Up to 30 fps
- 1 to 4 outputs
- 6 cameras (2 with color) to choose from

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Fax: +81 3 5960 6354


E-mail: asia@dalsa.com

e2v technologies continues to lead the world market, supplying high performance Charge Coupled Devices (CCDs) for a broad range of demanding applications. Operating across the electro-magnetic spectrum from X-rays, through ultraviolet and visible light to infrared, our solid-state sensor range includes very large area (VLA) and very high speed (VHS) arrays, alongside a comprehensive range of very low noise (VLN) spectroscopic devices.

Delivered to highly customized or standard designs, our sensors serve space & astronomy, dental, scientific & medical, and industrial markets. Our CCD technology is also extensively qualified in the defense sector.

Over two decades of continued investment, research and development has enabled e2v to break several imaging boundaries to date. Leading-edge technologies, such as L3Vision™, are ensuring that the company's high-end imaging products remain world leaders.

CCD47-20 with custom filters attached

#### **MEDICAL**

As an established supplier of detectors and sub-assemblies, e2v technologies continues to make significant contributions to medical imaging, principally to x-ray and biopsy analysis. Our unique capabilities include:

- Large area sensors
- Proprietary scintillator and phosphor deposition for excellent image quality
- Fiber-optic hard coupling with proven reliability
- High full well, low dark current for high dynamic range
- Readily available applications support

Further developments of our medical imaging capabilities include the introduction of our extreme low light level CCD technology, L3Vision™ (see page 4). L3Vision™ is advancing pre-emptive diagnostic and minimally invasive surgical techniques worldwide.

#### DENTAL

e2v's intra-oral X-ray imaging sensors are revolutionizing the way X-ray images are captured, stored and retrieved in dental offices. Our encapsulated package sensors offer high performance imaging from low radiation doses. Special-to-purpose packaging and X-ray conversion coatings have resulted in a discrete, easy to use sensor with exceptional performance.


For OEMs, we can also offer a driver board for both high and low resolution sensors. e2v technologies' DXDB01 dental driver board has passed certification as a Class 1 device under medical electrical equipment IEC60601 and provides a convenient interface between PC and e2v intra-oral dental sensors. In addition, we can offer a suite of software and detailed documentation to simplify the integration into the OEM application.

#### SCIENTIFIC

Ensuring our customers have the correct CCD for the specific task is a central philosophy at e2v technologies. Our range of standard and custom devices encompasses many state-of-the-art processes and technologies, which results in the scientific community being well served across the application spectrum, including spectroscopy, microscopy and general scientific imaging.

Our CCDs are well suited to biosciences, bioluminescence, chemical fluorescence and drug discovery applications, such as high throughput screening. Several imaging arrays in our portfolio offer low noise readout, high and medium resolution, and formats of 256x256, 512x512, 1kx1k and 2kx2k pixels (please see e2v's CCD Selection Guide for more details). Bioscience applications, such as single molecule fluorescence, have benefited greatly from the improved sensitivity at high speeds that our L3Vision™ technology offers.


High quantum efficiency on back illuminated devices (typically 95% peak), very low noise (2e-rms noise equivalent signal) or high speed (above 35 MHz) are key features of our CCDs. When combined with the options of Inverted Mode Operation/Multi-phase Pinned (IMO/MPP), back thinning and a wide selection of conversion and anti-reflective coatings, we have the capability to match demanding customer requirements.

L3Vision™ is a unique, groundbreaking CCD technology from e2v technologies that offers the best low light CCD performance available today. With performance in excess of the best ICCDs, L3Vision™ finds its main areas of application in photon-starved situations, such as night-time surveillance or the detection of fluorescent and luminescent markers in lifesciences. Additionally, the sensor can detect single photons, opening the possibility of photon counting with a solid state imaging sensor. However, it is the ability to operate across the full range of light levels, with the consequent increase in system flexibility, which really makes this technology stand out.


In lifesciences, fast, high sensitivity operation can be combined with quantitative, slow scan and time integrated imaging measurements. For surveillance, true 24-hour capability is now possible in one sensor. L3Vision™ CCDs achieve their performance at very low light levels by amplifying the image signal on the chip while it is still in the charge domain. This allows the sensor to operate in real time with sub-electron equivalent read-out noise, enabling very dim sources and very dark scenes to be imaged.

L3Vision™'s low light bio- and chemi-luminescence applications

e2v's L3Vision™ sensors are capable of TV imaging at faceplate illumination levels up to 1 lux. In bright sunlight, this is achieved with an auto iris lens. At the other extreme, the sensor will image at faceplate illumination levels down to 0.1 mlux, equivalent to overcast starlight with an F1.4 lens. In overcast starlight conditions, the sensor is producing an image from only a few photons per pixel. The sensor achieves this low light performance with a simple optical path and minimizes scintillations, halo, shading and lag, creating an image that is easy on the eye and provides superb resolution and detail. The sensor exhibits excellent intrascene contrast (improved mid-band MTF) and benefits from a wide spectral response (400-1060nm). Undamaged by sudden or prolonged light overloads, L3Vision's™ ruggedness has a profound impact on system-lifecycle costs.

L3Vision™ image of scene with complex lighting taken at night

#### LOW LIGHT CAMERA SOLUTIONS

Building on the success of our low light video applications, we offer our own specialist CCD camera solutions within the L3Vision™ range. The L3C65 series camera, housing e2v's new integral Peltier sensor, offers RS170- or CCIR-compatible video imaging from bright daylight down to overcast starlight conditions. Our L3C60 camera, available with the back-thinned CCD60 at its core, utilizes a similar design approach to enable frame rate imaging up to a maximum of 500fps.

Both cameras are available as fully working modules or as unboxed sub-assemblies, which can be easily integrated into alternate camera systems. Subject to order quantity, both cameras can be specially configured to customer specification.


### SPACE & ASTRONOMY

Building partnerships coupled with exemplary technology has made us one of the leading imaging suppliers to the space industry. Headline space organizations such as NASA, ESA and ISAS have selected e2v technologies' CCDs for major research missions, including ENVISAT, the Hubble Space Telescope, XMM-Newton and Solar-B. We have also successfully participated in many missions with Universities and Research Institutes from around the world.

Mostly custom devices, these highly complex CCDs have been developed for a wide range of applications including spectroscopy, atmospheric monitoring and planetary surveys using x-ray, ultraviolet, visible light and near-infrared.

Our understanding of the effect radiation has on CCDs has led to the development of specialist techniques to reduce the impact of radiation without degrading performance. Sophisticated packaging techniques are used to optimize the thermal, mechanical, optical and practical aspects of our devices ensuring effective, reliable operation.

ENVISAT environmental monitoring satellite with an image taken of Sicily, courtesy of ESA.

### MILITARY & SURVEILLANCE


Long established in the defense and surveillance sectors, our pre-eminent CCDs continue to offer outstanding performance in exacting applications. The parameters where we excel in the space market are those that also make serving the military market natural for us: quality, performance and security of supply. e2v technologies supplies a range of area CCDs for military video rate applications, as well as area and linear devices for high resolution surveillance purposes. Our CCDs are often used in conjunction with image intensifiers for night-time operation.

Increasingly, an L3Vision™ technology CCD is selected over an intensifier tube/CCD combination. L3Vision's™ unique capability is revolutionizing the 24-hour surveillance market. The ability to image a wide range of light levels, from overcast starlight to bright sunlight, offers a single surveillance solution that can be adapted for use on land, at sea or in the air.

S e l e c t i o n   G u i d e

e2v technologies

This guide describes ezv technologies' standard CCD (Charge Coupled Device) product line. ezv is known as a leader in sensitivity through the use of high quantum efficiency, back-illumination and very low read noise output amplifiers. Other options, including custom CCD designs, are also available to extend the wavelength range and sensitivity of devices across the electro-magnetic spectrum, from X-rays, through ultraviolet and visible light, to near infrared. Technical support is readily available and it is highly recommended that all new applications be discussed with ezv's sales and applications staff.

#### **PIXEL SIZE**

The selection of pixel size is a balance between basic performance parameters. Smaller pixels will give increased resolution, assuming that resolution is not limited elsewhere in the system. As the pixel size decreases, well capacity (the amount of charge that can be stored in a pixel) and the overall dynamic range decrease. As well as having a higher dynamic range, larger pixels gather more photons and hence may improve overall system sensitivity in light-starved applications. However, larger pixels also accumulate more dark current per pixel and a lower operating temperature may be necessary to offset this effect.

#### **NUMBER OF PIXELS**

This is determined by the active area and the pixel size. The number of pixels required for a given application may be limited by the readout time needed. Increasing the readout frequency will increase the read noise and may limit low signal performance. By reading out through multiple outputs, the frequency per output and hence the read noise may be reduced while maintaining the readout time.

#### **OUTPUT AMPLIFIER TYPE**

Several different types of output amplifier designs are used on ezv CCDs. The Very Low Noise (VLN) amplifier is optimized for the lowest possible read noise performance at slow scan speeds. The VLN amplifier can have a noise floor as low as 2 electrons. The Large Signal (LS) amplifier is optimized for high output node capacity to maximize dynamic range at slow scan speeds. It typically has a somewhat higher noise floor than the VLN amplifier. The High Speed (HS) amplifier is optimized for low noise at speeds of 10 MHz or higher. The dental CCDs use a moderate speed output amplifier with an integral sample and hold (S+H) to minimize the circuitry needed outside the CCD and to drive the long connection cable.

The L3Vision™ devices utilize a HS amplifier and a unique gain register that multiplies signal electrons so that they may be read at high speeds over the on-chip amplifier read noise, giving an effective noise of less than 1 electron.

#### **BACK-THINNING**

This is the process by which a front-illuminated CCD is mounted upside down on a support substrate and the bulk silicon is removed to allow the incident photons to reach the sensitive pixel wells unimpeded by intervening gate structures and insulators. Back-thinned devices have much higher quantum efficiency than standard (front-illuminated) devices. At UV and soft X-ray wavelengths, back-thinning will provide the highest quantum efficiency of any technology available today. Back-thinned devices are generally used in light-starved applications where a high signal to noise ratio is needed at very low light levels.

Back-thinned devices are available with either the basic thinning process or for better UV quantum efficiency the enhanced (or "astro") process. For most applications, the basic process is the best type to use, and combines the highest quantum efficiency with the lowest cost.

The enhanced process in combination with AIMO (Advanced Inverted Mode Operation) is only available on smaller sensors so that high row shift rates may be maintained. Larger sensors are available with the enhanced process, but only of the NIMO (Non-IMO) type.

Back-thinned devices are supplied with an anti-reflection coating, optimized for a particular wavelength, as shown in the Quantum Efficiency (QE) curves. Not all coatings are available for all types of chip, except to special order. Other coatings are available upon request.

Dental	Number of active pixels		Pixel pitch (µm)		Image size (mm)		Output amplifier type	Maximum readout rate (MHz)
	V	H	V	H	V	H		
CCD38-20	684	456	44	44	30.1	20.1	S+H	1
CCD76-40	1640	1250	22	22	36.1	27.5	S+H	1
CCD78-20	1368	912	22	22	30.1	20.1	S+H	1

Full-frame spectroscopic and scientific sensors	Number of active pixels		Pixel pitch (µm)		Image size (mm)		Output amplifier type	Maximum readout rate (MHz)
	V	H	V	H	V	H		
CCD30-11	256	1024	26	26	6.7	26.7	LS	5
CCD42-10	512	2048	13.5	13.5	6.9	27.6	2 x VLN	3
CCD42-20	1024	2048	13.5	13.5	13.8	27.6	2 x VLN	3
CCD42-40	2048	2048	13.5	13.5	27.6	27.6	2 x VLN	3
CCD42-80	4096	2048	13.5	13.5	55.2	27.6	2 x VLN	3
CCD42-90	4608	2048	13.5	13.5	62.2	27.6	2 x VLN	3
CCD44-82	4096	2048	15	15	61.4	30.7	2 x VLN	3
CCD47-10	1024	1024	13	13	13.3	13.3	2 x VLN	5
CCD55-20	1152	770	22.5	22.5	25.9	17.3	VLN + LS	7
CCD55-30	1152	1252	22.5	22.5	25.9	27.9	VLN + LS	7
CCD62-06	578	385	22	22	12.7	8.5	HS	12
CCD77-00	512	512	24	24	12.3	12.3	2 X LS	7

Frame transfer scientific sensors	Image section Number of active pixels		Store section Total pixels		Pixel pitch (µm)		Image size (mm)		Output amplifier type	Maximum readout rate (MHz)
	V	H	V	H	V	H	V	H		
CCD39-01	80	80	80	80	24	24	1.9	1.9	4 x VLN	3
CCD39-02	80	80	80	80	24	24	1.9	1.9	1 x VLN	3
CCD44-82	2048	2048	2052	2098	15	15	30.7	30.7	2 x VLN	3
CCD47-20	1024	1024	1033	1056	13	13	13.3	13.3	2 x VLN	5
CCD48-20	1024	1024	1033	1056	13	13	13.3	13.3	2 x HS	10
CCD55-20	576	770	576	804	22.5	22.5	13.0	17.3	VLN + LS	7
CCD55-30	576	1252	576	1276	22.5	22.5	13.0	27.9	VLN + LS	7
CCD57-10	512	512	528	536	13	13	6.7	6.7	2 x VLN	3
CCD62-06	288	385	290	395	22	22	6.3	8.5	LS	12
CCD67	256	256	264	268	26	26	6.7	6.7	2 x VLN	5

L3Vision™	Image section Number of active pixels		Store section Total pixels		Pixel pitch (µm)		Image size (mm)		Output amplifier type	Maximum readout rate (MHz)
	V	H	V	H	V	H	V	H		
CCD60	128	128	130	130	24	24	3.0	3.0	Low noise high capacity	18
CCD65 (525 line)	244	576	250	591	35.5	20	8.6	11.5	Low noise high capacity	16
CCD65 (625 line)	288	576	296	591	30	20	8.6	11.5	Low noise high capacity	16
CCD87	512	512	528	536	16	16	8.2	8.2	Low noise high capacity & VLN	13

**NOTE:**

Whilst ezv technologies has taken care to ensure the accuracy of the information contained herein, it accepts no responsibility for the consequences of any use thereof and also reserves the right to change the specification of goods without notice. ezv accepts no liability beyond that set out in its standard conditions of sale in respect of

infringement of third party patents arising from the use of devices in accordance with information contained herein. Not all combinations of variants are available for a single device. Users are advised to contact ezv to confirm if their particular requirements are available with a standard device before designing their system.

Readout noise (e-)	Package type	IMO	Coatings
			X-ray phosphor
150 (RNA)	Sealed plastic with cable & plug	■ (A)	■
150 (RNA)	Sealed plastic with cable & plug	■ (A)	■
150 (RNA)	Sealed plastic with cable & plug	■ (A)	■





Readout noise (e-)	Package type	Back-thinned	IMO	TE cooled peltier package	Anti-blooming	Coatings	Fiber optic	Open electrode	Hi-res deep depletion silicon
						Lumogen X-ray phosphor			
4 (RNB)	20-pin DIL ceramic	■	■ (A)			■ □	■	■	□
2 (RNB)	20-pin DIL ceramic	■	■ (A)		□	■ □	■		□
2 (RNB)	Ceramic PGA	□	□ (A)		□	□ □	□		□
2 (RNB)	24-pin DIL ceramic	■	■ (A)		□	■ □	■		□
2 (RNB)	3-side metal buttable	■	□ (A)		□				□
2 (RNB)	3-side metal buttable	■	□ (A)		□				□
2 (RNB)	3-side metal buttable	■	□						□
2 (RNB)	24-pin DIL ceramic	■	■ (A)	□	□	■ □	■		□
3 (RNB)	44-pin flat PGA	■	■ (A)		□	■ □	■		□
3 (RNB)	44-pin flat PGA	■	■ (A)		□	■ □	■		□
4 (RNB)	20-pin DIL ceramic	□	■ (A)		□	■ □	■		□
3 (RNB)	24-pin DIL ceramic	■	■	□		■ □	■		□


Readout noise (e-)	Package type	Back-thinned	IMO	TE cooled peltier package	Anti-blooming	Coatings	Fiber optic	Open electrode	Hi-res deep depletion silicon
						Lumogen X-ray phosphor			
3 (RNB)	24-pin DIL ceramic	■	■ (A)	■					
3 (RNB)	24-pin DIL ceramic	■	□ (A)	■					
2 (RNB)	3-side metal buttable	□	□						□
2 (RNB)	32-pin DIL Ceramic	■	■ (A)	■	■ (N)	■ □	■		□
35 (RNA)	32-pin DIL Ceramic	■	□	□	■ (N)	■ □	□		□
3 (RNB)	44-pin flat PGA	■	■ (A)		□	■ □	■		□
3 (RNB)	44-pin flat PGA	■	■ (A)		□	■ □	■		□
2 (RNB)	24-pin DIL ceramic	■	■ (A)	■	■ (N)	■ □	■	□	□
4 (RNB)	20-pin DIL ceramic	□	■ (A)		□	■ □	■		□
4 (RNB)	24-pin DIL ceramic	■	□ (A)	■	□	■ □	□		□

Readout noise (e-)	Package type	Back-thinned	IMO	TE cooled peltier package	Anti-blooming	Coatings	Fiber optic
						Lumogen X-ray phosphor	
< 1 (RNA)	24-pin DIL ceramic	■	■	□		■ □	□
< 1 (RNA)	36-pin PGA		■	■	■	■ □	□
< 1 (RNA)	36-pin PGA		■	■	■	■ □	□
< 1 (RNA)	30-pin DIL ceramic	■	■	□		■ □	□


- Key**
- Available
  - Subject to tooling and/or order quantity
  - (A) Advanced IMO
  - (N) Non-IMO only
  - IMO Inverted Mode Operation
  - VLN Very Low Noise scientific amplifier
  - LS Large Signal scientific amplifier
  - S+H Sample and Hold integral to output
  - HS High Speed scientific amplifier
  - DIL Dual in line
  - PGA Pin grid array
  - V Vertical
  - H Horizontal
  - TE Thermo-electric
  - RNA At maximum readout rate, FI device
  - RNB At 20kHz readout rate, FI device

# Quantum Efficiency(QE) curves - BACK-THINNED


Typical spectral response - Deep depletion, basic -20°C



Typical spectral response - Deep depletion, basic -100°C

Quantum Efficiency(QE) curve - FRONT-ILLUMINATED

Typical spectral response +20°C


**Back Thinned CCD with Improved Etaloning Interference**

Back illuminated CCD's are thin devices which become semitransparent in the near infrared. Reflections between the nearly parallel front and back surfaces of these devices cause them to act as etalons. This etalon-like behavior leads to unwanted fringes of constructive and destructive interference which artificially modulate a spectrum. The extent of modulation can be very significant (over 20%) and the spectral spacing of fringes, typically 5 nm, is close enough to make them troublesome for almost all NIR spectroscopy.

In order to reduce the etalon effect, Hamamatsu presents the new Back-Thinned Area Image Sensors S9060 and S9061 series. These products are high sensitivity, back thinned area image sensors, ideal for high precision measurement applications in the near IR region.

With a pixel size of 24 µm x 24 µm, these CCD are available in two pixels formats 512 x 256 (suffix -0908) and 1024 x 256 (suffix -1008). The S9061 has a one-stage thermoelectric cooler. The cooler is hermetically sealed with the chip in the same package. Therefore no dry air is required, allowing easy handling. Hamamatsu have a fully comprehensive range of back-thinned CCD chips and digital CCD cameras, for all demanding applications.

Author: Alberto Verga  
Hamamatsu Photonics Italy



S9060, S9061

**Features:**

- Improved etaloning (interference) type
- Pixel size: 24 µm x 24 µm
- High IR sensitivity
- Line/pixel binning operation

**Applications:**

- Fluorescence spectroscopy
- Raman spectroscopy
- Semiconductor inspection
- DNA sequencer
- Low-light level detection

Etaloning graph for back thinned CCD

**Specifications:**

Type no	Cooling	Number of active pixel	Active area H (mm) x V (mm)	Spectral response range (nm)
S9060-0908	Non-cooled	512 x 250	12.288 x 6.000	200-1100
S9060-1008		1024 x 250	24.576 x 6.000	
S9061-0908	One-stage TE cooled	512 x 250	12.288 x 6.000	
S9061-1008		1024 x 250	24.576 x 6.000	

**CCD Multichannel Detector Head**

The C7027 is a multichannel detector head designed for front-illuminated CCD image sensors, that offers high UV sensitivity and quantum efficiency.

This new multichannel detector head incorporates a low noise driver/amplifier circuit, that provides reliable operation from a simple external signal. A highly stable temperature controller allows cooling down to a temperature level of -50°C.

The housing configuration is designed for good heat dissipation. Four mounting holes on the front panel allow connection to other devices, such as monochromators for low light level detection applications.

Hamamatsu offers two series of CCD Image Sensors suitable for use with the C7027, the S7017 and S7847 series.

Author: Federica Gornati  
Hamamatsu Photonics Italy

**Specifications:**

Type no	Cooling	Number active pixel
S7017-1007	Four stage TE cooled	1024 x 1
S7017-1008		1024 x 2
S7847-1007		1024 x 1
S7847-1008		1024 x 2

**Time Line**

● Represents Management / Capital

Production / Technology


**BT-CCD Multichannel Detector Head****Features:**

- Designed for front-illuminated CCD image sensor
- Driver / amplifier circuit for low noise CCD operation
- High stability temperature controller
- Line binning and area scanning operation

**Applications:**

- Fluorescence spectroscopy
- Raman spectroscopy
- Low-light level detection

C7027

Active area (mm) x V (mm)	Spectral response range (nm)
24.576 x 2.976	400 - 1100
24.576 x 6.048	400 - 1100
24.576 x 2.976	400 - 1100
24.576 x 6.048	High IR sensitivity type

The CCD area scan camera C8677 series employs a high-speed, high resolution CCD (S8677-1010) and can be used in a wide range of applications including high-speed imaging, bio-photon imaging and semiconductor inspection. The CCD used in the C8677 series is sensitive in the ultraviolet region down to 130 nm (when CCD-faceplate is removed). The C8677 series is ideal for applications involving UV Light sources and UV Lasers. As soon as power is turned on, the C8677 series begins to operate in free running mode with a pre-set integration time and outputs continuous image data to an external unit. Image data is output in 12-bit parallel along with various timing signals for image acquisition, allowing easy connection to commercially available frame grabbers (C8677: RS-422, C8677-01: LVDS). All camera settings can be made through serial communication ports (RS-232C).

In addition to various external trigger functions, the C8677 series also provides one digital output. The C8677 series has sub-array binning to enhance the frame rate, computer-controlled gain and offset adjustments.

The CCD mounted in the C8677 series operates at high speeds and produces considerable heat so dark current increases in the CCD. To minimise dark current caused by this heat, a temperature control module is also built into C8677 series to maintain the CCD temperature near room temperatures ( $T_c = +20^\circ\text{C}$ ).

The camera body is compact yet designed to radiate heat efficiently away from the internal circuit and sensor. The CCD focal plane is flush with the front panel of the camera body making optical design easier.


Author: Arie van Gool  
Hamamatsu Photonics Deutschland

**Features:**

- Digital output
- High UV sensitivity
- High resolution
- Fast frame rate: up to 10 frames per second
- Sub-array binning: 2 x 2, 4 x 4, 8 x 8, 16 x 16
- Sensor temperature control:  $+20^\circ\text{C} \pm 0.1^\circ\text{C}$
- Easy handling (ie controllable from PC through RS-232C)

**Applications:**

- Spectrophotometry
- High-speed UV imaging
- Bio-photon observation
- Semiconductor inspection

C8677

**Specifications:**

The table below shows CCD area image sensor applicable for the C8677 series. The C8677 series does not include a CCD sensor.

Type no.	Number of pixels	Number of active pixels	Active area H (mm) x V (mm)
S8677-1010	1044 x 1032	1024 x 1024	18.43 x 18.43

**1953 - 1960**

- Hamamatsu TV Co Ltd (former name) established with a capital of 500,000 Yen
- Production of phototubes began
- Trial production of photomultiplier tubes began

Image pickup tubes were put on the market  
Successful production of night vision tubes

CdS cells were put on the market  
Side-on photomultiplier tubes were put on the market

- Capital increased to 8 million Yen

**1961 - 1969**

- Tokyo office opened
- Production of PbS detectors started
- Infrared Vidicons were put on the market

**BT-CCD Multichannel Detector Head**

The CCD multichannel detector head C9265 series employs a 2/3-inch frame transfer CCD (S7987-01) for NTSG B/W and can be used in a wide range of applications including high-speed imaging, bio-photon imaging and semiconductor inspection. The CCD used in C9265 series is sensitive in the ultraviolet region (down to 130 nm after removing CCD faceplate) and is ideal for applications involving UV light sources and UV lasers.

As soon as power is turned on, the C9265 series begins to operate in free running mode with pre-set integration time and outputs continuous image data to an external unit. Image data is output in 12-bit parallel format (C9265: RS-422 or C9265-01: LVDS) along with various timing signals for image acquisition, allowing easy connection to commercially available frame grabbers.

All camera settings can be made through serial communication ports (RS-232C). In addition to various external trigger functions, C9265 series also provides one digital output. The C9265 series has a sub-array binning to enhance the frame rate, computer-controlled gain and offset adjustments.

The CCD mounted in C9265 series operates at high speeds and produces considerable heat, so dark current increases in the CCD. To minimise dark current caused by this heat, a temperature control module is also built into C9265 series to maintain the CCD temperature near room temperatures ( $T_c = +10^\circ\text{C}$ ).

**Specifications:**

The table below shows CCD area image sensor applicable for the C9265 series. The C9265 series does not include a CCD sensor.

Type no.	Number of pixels	Number of active pixels	Active area H (mm) x V (mm)
S7987-01	680 x 500	658 x 490	9.212 x 6.860

The camera body is compact yet designed to radiate heat efficiently away from the internal circuit and sensor. The CCD focal plane is flush with the front panel of the camera body making optical design easier.


Author: Arie van Gool  
Hamamatsu Photonics Deutschland

**Features:**

- High UV sensitivity: sensitive down to 130 nm
- High resolution
- Fast frame rate: up to 25 fps (frames per second)
- Sub-array binning: 1 x 1, 2 x 2
- Various trigger modes
- Sensor temperature control:  $+10^\circ\text{C} \pm 0.1^\circ\text{C}$  ( $T_a = +25^\circ\text{C}$ )
- Easy handling
- Controllable from PC (RS-232C)

**Applications:**

- High-speed UV imaging
- Semiconductor inspection
- Microscope

C9265

**CCD Multichannel Detector Head**

C9047 series is an industrial camera, designed for our high speed and high resolution back-thinned CCD Image Sensors S9037-0902 and S9037-1002 with 512 or 1024 pixel. It is suitable for a wide range of applications including production line inspection and surveying. The CCD Sensors are sensitive in the UV region down to 130 nm (windowless version) which makes it ideal for UV laser applications.

As soon as the power is turned on, C9047 series enters standby mode ready to output image data (12 bit digital output) in synchronization with an external trigger signal. All camera settings can be made through serial communication ports (RS-232C).

The camera body is compact yet designed to radiate heat efficiently away from the internal circuit and sensor. The CCD focal plane is flush with the front panel of the camera body making optical design easier.

Author: Andreas Hafner  
Hamamatsu Photonics Deutschland

**Features:**

- 12 bit digital output
- High UV sensitivity down to 130 nm (windowless version)
- Fast line rate  
S9037-0902: 16 kHz  
S9037-1002: 9 kHz
- Easy handling

**Applications:**

- Spectrophotometry
- High-speed UV imaging
- Bio-photon observation
- Semiconductor inspection

Silver-bismuth photomultiplier tubes put on the market  
Successful production of 1/2 inch Vidicon  
Sale of infrared video cameras began  
● Ichino factory started operation

UV-TRON and solar blind photomultiplier tubes put on the market  
Head-on type photomultiplier tubes and storage tubes put on the market

● New York office opened  
Experimental production began of hollow cathode lamps and deuterium lamps

1/2 inch side-on type photomultiplier tubes, electron multiplier tubes, pick-up storage tubes and glow modulator tubes were put on the market.

## New Photodiode Array with Amplifier

Hamamatsu have increased their range of S8865 series photodiode arrays. The S8865-256 is a photodiode array with a signal processing circuit chip. This signal processing circuit chip is formed by a C-MOS process and incorporates timing generator, shift register, charge amplifier array and clamp circuit, making the external circuit configuration simple. A long, narrow image sensor can also be configured by arranging multiple arrays in a row since the S8865 device can be stacked "end to end". For X-ray detection applications, types with fluorescent phosphor affixed to the active area are also available.

Author: Oliver Roesler  
Hamamatsu Photonics Deutschland

### Features:

- S8865-256: 0.2 mm pitch; 256 ch; 51.2 mm
- 5 V single power supply operation
- Sequential readout with a shift register (Data rate: 1 MHz max.)
- Integrated clamp circuit allows low noise and wide dynamic range
- Integrated timing generator allows operation at two different pulse timings

### Applications:

- Long line sensors
- Line sensors for X-ray detection

## New Compact C-MOS Driver Circuit

The C9118 series C-MOS driver circuit is designed for the S8865 series photodiode array with amplifier.

The C9118 series operates a linear image sensor by just inputting two signals (M-CLK and M-RESET) and a single +5 V supply voltage. It is assembled on a compact board, allowing downsizing of the measuring equipment. This driver circuit is intended for single use or parallel connections, while the C9118-01 is suitable for cascade connections.


Author: Oliver Roesler  
Hamamatsu Photonics Deutschland

### Features:

- Single power supply
- Operation with two input signals
- Compact size, 48 mm x 56 mm

### Applications:

- Image acquisition
- Optical detection equipment

C9118

## Si Photodiodes with Precision Human Eye Spectral Response

For precision measuring instruments demanding a very accurate human eye spectral response the S9219 series is the optimal choice. These photodiodes have an integrated visible compensation filter which closely matches human eye response.

The S9219 series is ideal for demanding photometry and luxmeter applications. The S9219 series filter has been improved. There are two different active areas and types of package available as standard.

Author: Dieter Habben  
Hamamatsu Photonics Deutschland

### Specifications:

Type no.	S9219	S9219-01
Spectral response range (nm)	380 ... 780	380 ... 780
Peak sensitivity wavelength (nm)	555	555
Sensitivity at peak sensitivity wavelength (A/W)	0.24	0.22
Active area (mm)	dia. 11.3	3.6 x 3.6
Package	BNC	TO-5

For less demanding applications, Hamamatsu also offers a broad range of human eye response photodiodes.

S5493-01, S5627-01	filterless, low cost
CdS cells	low cost, resistor
S7565	linear, Photo-IC, low cost
S7686	fast speed, sensitivity close to the human eye
S7123-01, S7123-02	wider operating temperature range, filterless
S9066, S9067	Photo-IC, without filter
S9219	high precision

The S8865 series photodiode array with amplifier was put on the market in 1970. Hamamatsu Photonics Deutschland

Image dissectors, Irisorders and multi-alkali photomultiplier tubes put on the market

### 1970 - 1979

Side-on type multi-alkali photomultiplier tubes, hollow cathode lamps and deuterium lamps were put on the market. 3" photomultiplier tubes launched for the scintillation counting market

1/2 inch head-on photomultiplier tubes, photon-counting photomultiplier tubes, GaAs photomultiplier tubes, Si photodiodes and X-ray Vidicons were put on the market

**HAMAMATSU**

### Low Cost C-MOS Linear Image Sensors

The S9226/S9227 are a new kind of small, easy to use low cost linear image sensors from Hamamatsu. These, low cost devices, are ideal for applications such as position detection and low cost spectroscopy.

In our Euronews 2-2002 we presented our new integrated C-MOS linear image sensors. They allow easy operation, by just inputting 2 signals and a single +5 V supply voltage.

The S9226/S9227 are new (small) members of this family offering the same simple operation.

When high linearity is needed the smaller pixel type S9226 is the choice. It has a built-in signal processing circuit. It can be used for the low current 3.3 V technology. The S9227 allows high speed up to 5 MHz video data rate. Simultaneous integration and built in shutter allows the correct data transfer even when the light input changes.

Author: Dieter Habben  
Hamamatsu Photonics Deutschland

#### Applications:

- Position detection
- Image reading
- Analytical equipment

#### Applications:

- Position detection
- Image reading

#### Specifications:

Type no.	S9226	S9227
Pixel pitch (µm)	7.8	12.5
Pixel height	125	250
Number of pixels	1024	512
Video data rate (kHz)	max. 200	max. 5000
Supply voltage (V)	3.3, 5	5
Spectral response range (nm)	400 - 1000	
Housing	8-pin DIP, 16-pin SMD available	
Timing	Built in timing generator - only start and clock pulse necessary	
Specialities	On-chip charge amplifier	Shutter function

### InGaAs Linear Image Sensors

Hamamatsu has released an improved range of InGaAs linear image sensors.

Compared to the previous types, our G80xx or G81xx series, the Si C-MOS multiplexer chip has been changed. The key point of this new development is to improve the characteristics whilst making it compatible to the previous types.

The input and output (analog and digital) signals, the basic circuit connection, the dimensional outline, the pin connection and the cooling efficiency are the same, compared to former types.

The dynamic range has however been increased from 15000 to 18700 (G9204-512S) and the dark current decreased by a factor of 6 (G9204-512S). The linearity error at a level of 50% of Vsat decreased by 3 times. These improvements enable the user to operate our 2.6 µm cutoff wavelength array up to 100 ms integration time.



Slight changes in existing circuits are only in timing chart (10 extra clocks) and the reference voltage, which has to be adjusted to 1.26 V.

Customers can choose between 2 different packages, ceramic DIL and metal package. The metal package is hermetically sealed and

can be used with internal Peltier cooling. The ceramic package is a lower cost, but it can only be operated without internal cooling.

Hamamatsu offers readout circuits for all different types of detectors.

Author: Andreas Hafner  
Hamamatsu Photonics Deutschland

#### Specifications:

λc	Pixel pitch, number of pixel	Pixel height (µm)	Cooling	Dead pixel (max.)	Package	Type number	
1.70 µm	50 / 256	250	No	0%	Metal	G9201-256R	
1.67 µm			1TEC			G9201-256S	
1.70 µm	25 / 512	250	No		Ceramic DIL	G9202-512R	
1.67 µm			1TEC			G9202-512S	
1.70 µm	50 / 256	500	No		Metal	G9203-256D	
1.70 µm			1TEC			G9203-256R	
1.70 µm	25 / 512	500	No		Ceramic DIL	G9204-512D	
1.67 µm			1TEC			G9204-512R	
1.67 µm	50 / 256	250	1TEC		Metal	G9211-256S	
1.67 µm						500	G9212-512S
1.85 µm	50 / 256	250	2 TEC		5%	Metal	G9213-256S
2.05 µm							G9214-512S
2.25 µm				G9205-256W			
2.55 µm				G9206-256W			

- A joint company, Hamamatsu Television Europe GmbH, was established in Germany
- Capital increased to 90 million yen

Photomultiplier tubes for X-ray CT scanners and Imagers were put on the market

Computer compatible video cameras and streak camera systems put on the market

Micro channel plates, GaAsP and Si PIN photodiodes were put on the market

- Osaka sales office opened

PSDs and high temperature photomultiplier tubes were put on the market

**Profile Sensor**

The S9132 profile sensor has many advantages when compared to conventional C-MOS area image sensors and 2D position sensitive detectors (PSDs), including for example, improved output linearity and the ability to detect multiple light spots. The sensor is designed for use in inspection equipment and medical measurement devices. The S9132 acquires a projection profile in the X and Y axis to produce high-speed position and moving object detection. It has a 1.9968 by 1.9968 mm<sup>2</sup> active area with a 256- by 256-pixel array. Spectral response range is 350 to 1150 nm. At 780 nm, peak sensitivity of 40 V/nJ (high gain) or 8 V/nJ (low gain) is achieved. This high performance C-MOS sensor, integrated on the same chip with a timing and bias voltage generator and a 10-bit AD converter circuit, allows operation with very simple external signal processing circuits.


Author: Alberto Verga  
Hamamatsu Photonics Italy

**Features:**

- 10 bit / 8 bit switchable AD converter
- High speed frame rate
- Multiple light spots detection

**Applications:**

- Light spot position detection
- 3D measurement
- Moving object detection

S9132

**Specifications:**


Type no.	Number of pixels	Pixel pitch (μm)	Active area (mm <sup>2</sup> )	Spectral response (nm)	Resolution
S9132	256 x 256	7,8	1.9968 x 1.9968	350 - 1150	10 or 8 bit

**830 nm 1W CW Laser Diode**

Hamamatsu offers a variety of laser diodes emitting at 808 nm. However, some applications need very specific wavelength emission. To meet these specific needs, Hamamatsu has the skill to adapt the LD emission to your required wavelength range. This can be either a customised product or one from our existing large range of standard laser diode products. We have newly released the 830 nm 1 W CW laser diode L9277-42. Other LD with wavelength at 850 nm, 920 nm and 980 nm will be released soon.

The output power of L9277 is far over 1 W having an emission area of 1 μm x 50 μm for excellent focussing characteristics. It comes in a sealed 9.0 mm CD package ideal for single laser applications. OHS package with side-out (-62) or head-out (-72) laser diode emission directions are also available. The OHS package is taking account of the heat dissipation by a special heat sink used as needed e.g. for compact or multiple laser integration. For most applications such as printing or measuring instruments, the peak emission variation within ±10 nm is not relevant and decreases cost. A smaller value of variation can be selected if required. The L9277 series is an excellent small-band NIR emitter used in many applications such as printing, optical measuring instruments, medical instruments and pumping source for solid state lasers. They are known for their excellent differential energy efficiency. The high stability and long life of this compact laser diode will provide reliable operation for your optical instruments and measurements.

Author: Marius Metzger  
Hamamatsu Photonics France

L9277-42

**1980 - 1984**

Hamamatsu Corporation began production of photomultiplier tubes

Image Intensifiers and Night Viewers were put on the market

- Tenno Glass Works completed

Digital radiography video cameras, 20 inch photomultiplier tubes and 2D photon counting imager tubes were put on the market

- Company name changed to Hamamatsu Photonics K.K
- Capital increased to 495 million Yen


# Light Platform

HASO's optomechanical add-on

## PRINCIPLE

- Increase in HASO wavefront sensor's functionalities by adding a light source and an optical system.
- Now able to measure wavefront generated by the system to characterise and then its shape imperfections. HASO still keeps its classical functions: MTF and PSF processing.

A deported laser diode (630 nm), a beam splitter and lenses, with the HASO, allow to measure and process data on convergent systems which numerical aperture is between 4 and 30, that allows the direct measurement of concave mirrors (N= 2 to N= 15).

**Legend:**  
 HASO's elements  
 LIP's elements  
 Optical path


## ADVANTAGES

### Light Platform advantages:

- Cost reduction for setting up
- Time gain (setting up and measurement)
- Adaptable on HASO II 32 and 64
- Direct characterisation of opened mirrors, till N=2.

### HASO advantages:

- Insensibility to vibrations and atmospheric perturbations.
- Real time display and processing (average possible for stabilised measures.)
- Use HASO with or without Light Platform :
  - Simple and repeatable fixation on HASO II sensor.
  - Compatible with classical HASO's functionalities.
  - Only one software to run both HASO and Light Platform.

## APPLICATIONS

Concave

Convergent

Aspherical

### Main application:

- Direct characterisation of convergent optical systems.
- Real time controls and adjustments.
- WFE (Wave Front Error) and concave mirrors defaults characterisation

### And also:

- Specific engineering  
Custom developments for dedicated applications (mirror N=1 ...).
- MTF Measurement  
MTF processing and displaying.


# Light Platform

## TECHNICAL DATA:

Modules Available in standard version  
(numerical apertures):  
6.2 / 10 / 16 / 24

Modules available on request (numerical apertures):  
4 / 30


With the 5 axis stage (optional), the optical axis' height is adjustable between 168 and 242 mm.

## SPECIFICATIONS (HASO II PLUS LIGHT PLATFORM)

	HASO 32	HASO 64
Analysing aperture dimension	5 x 5 mm <sup>2</sup>	12 x 12 mm <sup>2</sup>
Number of sub-apertures	32 x 32	64 x 64
Wavefront measurement accuracy (rms)	~ λ/100 <sup>(1)</sup>	
Pupil Conjugation <sup>(2)</sup>	1 m	
Source wavelength	630 nm	
Working Temperature	15 - 30 °C	
Power supply	5,-12,12 V / 6 W	+5,-5,-15,28 V / 20 W

- (1) Depending on reference optic quality  
(2) System design for a 1 meter pupil conjugation.

## SET UP :

The minimal configuration for convergent systems characterisation require the following elements:

- 1 Light Platform
- 1 HASO II
- 1 five axis stage (Optional)
- 1 refrence mirror (Optional)

Development of specific configurations on request.

## REFERENCES:

HASO and Light platform are already in use at:

ALCATEL SPACE INDUSTRIE, CYBERNETIX, NASA, NATIONAL OBSERVATORY OF CHINA ...

Imagine Optic - 18, rue Charles de Gaulle - 91400 Orsay - France  
tel: + 33 1 64 86 15 60 - fax: + 33 1 64 86 15 61 - e-mail: contact@imagine-optic.com

## HASO II

*new generation wavefront sensor*

### PRINCIPLE

*More Performant*

### APPLICATIONS

*Beam Analysis*

*Phase and Intensity*

*MTF*

*Characterization*

*Adaptive Optics*

### ADVANTAGES

*Real Time*

*High Dynamic*

*Interactivity*

➤ **Optical wavefront sensor**

These wavefront sensors (Shack-Hartmann type) allow to obtain a map of local slopes of the analysed wavefront with the help of a microlenses array. This principle used for a long time in astronomy was amply improved as well on hardware as on software. The CCD detectors reached the maturity and the necessary quality/price ratio to certify an irreproachable metrology with moderate prices; the fabrication of microlenses arrays has become reliable. The software innovations of Imagine Optic allow to use trough optimal way this simple and direct measurement.

➤ **Beams analysis in phase and intensity**

Laser, laser diode and all the coherent or incoherent light sources: the HASO allows to measure divergence, position and size of the waist,  $M^2$  parameter, Strehl ratio, astigmatism values, ...

➤ **Components and systems characterization**

Especially adjusted to middle size components (aperture of 5 to 12 mm), the HASO dynamic allows the analysis of small and large components too with divergent or convergent beam.

➤ **Focal length, MTF and point spread function measurement**

➤ **Adjustment, mounting, alignment**

The operator visualizes in real time the wavefront and some measured parameters, which show adjustments quality and allow to improve them (tilts and curvature values, space source position, third order aberrations which can move during an adjustment).

➤ **Active and adaptive optics**

Wavefront correction in real time with  $\lambda/10$ - $\lambda/20$  pv.

➤ **Achromatic and non sensitive to the coherence**

Particularly, monochromatic wavefront analysis.

➤ **Display and calculation in real time**

Adjustment and on-line control.

➤ **High spatial resolution**

Up to 4096 points measured in a standard configuration on the pupil.

➤ **Strong divergent and convergent beams analysis**

➤ **Wavefronts measurement in absolute or relative mode**

The HASO is factory calibrated. The measure is absolute except if a reference wavefront is specified.

# HASO II

## SPECIFICATIONS

	HASO 16	HASO 32	HASO 64
Aperture dimension	2,5 x 2,5 mm <sup>2</sup>	5 x 5 mm <sup>2</sup>	12 x 12 mm <sup>2</sup>
Number of sub-apertures dedicated for analysis (1)	16 x 16	32 x 32	64 x 64
Tilt dynamic range	> ± 3 ° (260 λ)	> ± 3 ° (520 λ)	> ± 3 ° (1200 λ)
Focus dynamic range	± 0,025 m to ± ∞ (50 λ)	± 0,025 m to ± ∞ (200 λ)	± 0,04 m to ± ∞ (800 λ)
Repeatability (rms)	< λ/200		
Wavefront measurement accuracy in absolute mode rms (2)	λ/100		
Wavefront measurement accuracy in relative mode rms (3)	λ/150		
Tilt measurement sensitivity	6 μrad	3 μrad	1 μrad
Focus measurement sensitivity	5.10 <sup>-3</sup> m <sup>-1</sup>	10 <sup>-3</sup> m <sup>-1</sup>	5.10 <sup>-4</sup> m <sup>-1</sup>
Spatial resolution	~160 μm	~160 μm	~190 μm
Max acquisition frequency	955 Hz	77 Hz	20 Hz (opt. 40 Hz)
Processing frequency	8-40 Hz	4-20 Hz	2-10 Hz
Working wavelength	350-1100 nm		
Calibrated wavelength band	400 - 600 nm / 500 - 700 nm / 650 - 900 nm / 800 - 1100 nm		
Extended calibrated wavelength band	400 - 700 nm / 500 - 900 nm / 650 - 1100 nm		
Working temperature	15 - 30 °C		
Dimensions / weight	90x100x130 mm / 1125 g	80x90x100 mm / 980 g	90x100x130 mm / 1200 g
Power supply	5,-12,12 V / 6 W	5,-12,12 V / 6 W	+5,-5,-15,28 V / 20 W

these data are given for the central wavelength of the band defined

(1) other configurations are available: HASO 50 for instance.

(2) wavefront seen by the analyser. performance kept on the whole band

(3) difference between the real wavefront and a reference wavefront obtained in similar conditions (10 λ of shift max).

## SOFTWARE

### Wavefront acquisition

Continuous or ext. trigger mode – absolute (factory calibration) or referenced (user calibration) mode – one-shot or flying acquisition – camera signal – background signal : save and subtract in real time – for the sensor adjustment : tilt visualization ; manual adjustment of integration time.

### Wavefront reconstruction

Wavefront display (unit in microns or lambda) : 3D, 2D projection, fines mode, xy profile, save function (text or excel format, comments) – display in real time without tilt and/or focus and/or 3<sup>rd</sup> order aberrations – intensity display : 3D, 2D projection, xy profile, save function (text or excel format, comments) – reconstruction : zonal or modal (Zernike or Legendre) – pupil shape : with or without central occultation – location (x y z) of the focus point – polynomial coefficient display (unit in microns or lambda), historic function : display of the time changes of 4 coef. – local slopes display : with or without tilt and focus, zoom, rms and pv of the slopes - pupil calculation : automatic or manual (save and load functions) – spot diagram in real time for various observation positions.

### Extended wavefront reconstruction

Zernike polynom coefficient processing in real time for any pupil shape – display of the calculated wavefront surface and display of the residual – manual or automatic choice of the pupil - save functions: text or excel format, comments

### Point Spread Function measurement

PSF display : 3D, 2D projection, xy profile – parameters : best focus plane, user specified plane, focal plane of ideal lens – zoom (1 to 8) – Strehl ratio – numerical value in table - save and load functions (text or excel format, comments).

### MTF measurement

MTF display : 2D, x profile (user choice) – parameters : best focus plane, user specified plane, focal plane of ideal lens – zoom (1 to 8) – historic function : display of the time change of a specified frequency - save and load functions (text or excel format, comments).

### Laser characterisation

Beam location, dimension and position of the waist – display of the beam shape – M<sup>2</sup> calculation – calculation of intensity shape in any plane - save and load functions (text or excel format, comments).

### Dynamic library

dll in C or VI for Labview : all the software functions to build specific soft applications.

Imagine Optic - 18, rue Charles de Gaulle - 91400 Orsay - France  
Tel : +33 1 64 86 15 60 - fax: +33 1 64 86 15 61 - e-mail: contact@imagine-optic.com

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## HASO NIR

*new wavefront sensor for telecom applications*

### PRINCIPLE

*More Performant*

*Low Cost Solution!*

### APPLICATIONS

*Beam Analysis*

*Phase and Intensity*

*MTF*

*Characterization*

*Fiber Coupling*

### ADVANTAGES

*Real Time*

*High Dynamic*

*Interactivity*

#### > **Optical wavefront sensor**

These wavefront sensors (Shack-Hartmann type) allow to obtain a map of local slopes of the analysed wavefront with the help of a microlenses array. This principle used for a long time in astronomy was amply improved as well on hardware as on software. But the CCD technology does not allow measurement in the telecom range, and especially around 1,5  $\mu\text{m}$ . The InGaAs photodiodes technology is extremely expensive and does not answer to quality criteria to guaranty a good Shack-Hartmann metrology (reliability, dead pixels, etc...). Some detectors (based on a phosphor coating technology on a CCD) reached the maturity and the necessary quality/price ratio to certify an irreproachable metrology with moderate prices; the fabrication of microlenses arrays has become reliable. The software innovations of Imagine Optic allow to use trough optimal way this simple and direct measurement.

#### > **Beams analysis in phase and intensity**

1500 nm laser source: the HASO allows to measure optical quality, divergence, position and size of the waist,  $M^2$  parameter, Strehl ratio, astigmatism values, ...

#### > **Components and systems characterization**

Especially adjusted to middle size components (aperture of 5 mm), the HASO dynamic allows also the characterization of small and large components with divergent or convergent beam.

#### > **Adjustment, mounting, alignment**

The operator visualizes in real time the wavefront and some measured parameters, which show adjustments quality and allow to improve them (tilts and curvature values, space source position, third order aberrations which can move during an adjustment). When the user ends operation, he can get the optical quality (and especially the MTF).

#### > **Point spread function**

The phase and intensity measurement lead to the PSF (shape and location of the focus spot). Then it is possible to know accurately the fiber coupling efficiency and to understand some bad performances.

#### > **Achromatic and non sensitive to the coherence**

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#### > **Display and calculation in real time**

Adjustment and on-line control..

#### > **High spatial resolution**

Up to 1200 points measured in a standard configuration on the pupil.

#### > **Strong divergent and convergent beams analysis**

Compatible with fiber and laser source.

#### > **Wavefronts measurement in absolute or relative mode**

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\*: all specifications given with an average operating mode of 10 acquisitions.

Note 1: For the central wavelength: 120 λ for the HASO NIR 32, 180 λ for the HASO NIR 40.

Note 2: wavefront directly measurement by the wavefront sensor (no added lens).

Note 3: difference between a referenced wavefront and the measured wavefront, in a range < 10 λ.

## SOFTWARE

### Wavefront acquisition (standard software)

Absolute (factory calibration) or referenced (user calibration) mode – one-shot or flying acquisition – camera signal – background signal : save and subtract in real time – for the sensor adjustment : tilt visualization ; manual adjustment of integration time.

### Wavefront reconstruction (standard software)

Wavefront display (unit in microns or lambda) : 3D, 2D projection, fines mode, xy profile, save function (text or excel format, comments) – display in real time with or without tilt and/or focus and/or 3<sup>rd</sup> order aberrations – intensity display : 3D, 2D projection, xy profile, save function (text or excel format, comments) – reconstruction : zonal or modal (Zernike or Legendre) – pupil shape : with or without central occultation.

### Wavefront : other tools (standard software)

Location (x y z) of the focus point – polynom coefficient display (unit in microns or lambda), historic function : display of the time changes of 4 coef. – local slopes display : with or without tilt and focus, zoom, rms and pv of the slopes - pupil calculation : automatic or manual (save and load functions) – spot diagram in real time for various observation positions.

### Point Spread Function measurement (PSF software)

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### Laser characterisation (M square software)

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
### Coupling efficiency (Preliminary)

Beam location, dimension and position of the waist – display of the beam shape – calculation of intensity shape in any plane, and comparison with dimension and aperture of the fiber - save and load functions (text or excel format, comments).

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## HASO UV

*new wavefront sensor for UV applications*

### PRINCIPLE

*More Performant*

➤ **Optical wavefront sensor**

These wavefront sensors (Shack-Hartmann type) allow to obtain a map of local slopes of the analysed wavefront with the help of a microlenses array. This principle used for a long time in astronomy was amply improved as well on hardware as on software.

The CCD detectors reached the maturity and the necessary quality/price ratio to certify an irreproachable metrology with moderate prices; the fabrication of microlenses arrays has become reliable. The software innovations of Imagine Optic allow to use through optimal way this simple and direct measurement.

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From 244 to 350 nm laser source: the HASO allows to measure optical quality, divergence, position and size of the waist,  $M^2$  parameter, Strehl ratio, astigmatism values, ...

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Adjustment and on-line control.

➤ **High spatial resolution**

Up to 2500 points measured in a standard configuration on the pupil.

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Note 1: For the central wavelength: 150 λ for the HASO UV 25, 470 λ for the HASO UV 50.

Note 2: wavefront directly measurement by the wavefront sensor (no added lens).

Note 3: difference between a referenced wavefront and the measured wavefront, in a range < 10 λ.

## SOFTWARE

### Wavefront acquisition

Continuous or ext. triggered mode – absolute (factory calibration) or referenced (user calibration) mode – one-shot or flying acquisition – camera signal – background signal : save and subtract in real time – for the sensor adjustment : tilt visualization ; manual adjustment of integration time.

### Wavefront reconstruction

Wavefront display (unit in microns or lambda) : 3D, 2D projection, fringes mode, xy profile, save function (text or excel format, comments) – display in real time without tilt and/or focus and/or 3<sup>rd</sup> order aberrations – intensity display : 3D, 2D projection, xy profile, save function (text or excel format, comments) – reconstruction : zonal or modal (Zemike or Legendre) – pupil shape : with or without central occultation.

### Wavefront : other tools

Location (x y z) of the focus point – polynom coefficient display (unit in microns or lambda), historic function : display of the time changes of 4 coef. – local slopes display : with or without tilt and focus, zoom, rms and pv of the slopes - pupil calculation : automatic or manual (save and load functions) – spot diagram in real time for various observation positions.

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
## 10" SCHMIDT-CASSEGRAIN TELESCOPE CHARACTERISATION THANKS TO A SHACK HARTMANN WAVEFRONT SENSOR

**Charles Brossollet, Jérôme Ballesta, Laurent Couvet , Guillaume Dovillaire**  
*Imagine Optic, 18 rue Charles de gaulle, 91400 Orsay, France*

### *Abstract*


*HASO (Shack-Hartmann wavefront sensor manufactured by Imagine Optic) is in part dedicated to optical system metrology. Its wide dynamic and its easy to use and to integrate approach allow this sensor to characterise a wide range of dioptricals and catadioptricals systems. This document will present how to make a wide aperture optical system characterization with an HASO+LIP. Results will be also presented and effects of misalignment on the WFE will be appreciated.*

*For more technical or commercials information please contact us*  
[contact@imagine-optic.com](mailto:contact@imagine-optic.com)

## 1. Set up description

This application shows how to measure Schmidt-Cassegrain WFE with a HASO plus LIP upgrading system. The illuminating beam is adapted to the numerical aperture of the telescope thanks to a lens ( $f/10$ ). As the telescope has got 10" aperture, a 8"  $\lambda/150$  rms flat mirror has been used for the autocollimation.

## 2. Measurement description

First of all the measurement the set up must be calibrated and. Here a better than  $\lambda/150$  rms flat mirror<sup>1</sup> has been used for the autocollimation and the whole process has been done vertically. A perfect sphere has been used for the HASO+LIP calibration. In this configuration, WFE measurement accuracy was better than  $\lambda/50$  rms. Please find below the procedure for the Telescope characterisation:

1. **Calibration of LIP thanks to a perfect sphere ( $\lambda/150$  rms). Saving of the LIP's WFE file**
2. **Alignment of the LIP in front of the telescope (rotation of the flat mirror).**
3. **Measurement of the telescope's WFE in referenced mode (using the data file described in 1).**


---

<sup>1</sup> This mirror will limit the accuracy of the measurement. Several tests allow to check its flatness.

### 3. Results

#### 3.1 Before Alignment


HASO.exe displays shown us that a field aberration was. In fact a misalignment of the secondary mirror introduced coma when the Lip was illuminated the system aligned on the mechanical axis.


#### 3.2 After Alignment

Real time measurement allow us to process the alignment of the whole system. The secondary mirror was aligns on the primary in order in reducing to 0 field aberrations.

WFE both rms and PV value allow us to appreciate the effect of the alignment. With a such accurate wavefront measurement a lot of processing are reachable, and then it's possible to calculate MTF of the system for example.

1999

# HIGH-TECH IN GLASS.


## KNOW-HOW AT THE CUTTING EDGE.


We're an innovative, modern company with outstanding know-how like our advanced ion beam polishing process. Our technological resources are complemented by CNC machines specially developed for precision glass processing, modern thin film sputtering technology, and precision measuring techniques.

Quality, innovation, reliability, and flexibility: they're not mere mottos but the yardsticks by which we measure our success.

In close co-operation with our customers, we develop individual solutions for special technological challenges – and accompany every project from start to finish with our expertise and advice.

In every phase of cooperation we maintain an intensive dialogue and teamwork with our customers. This ensures maximum flexibility throughout development, production, assembly, adjustment, and certification.


## SAFEGUARDING PRECISION.


Seeing, protecting, measuring, guiding, converging, dispersing, reflecting, recording – accurately, precisely, reliably – with glass and comparable materials. The Berliner Glas workforce in Germany, Switzerland, and the United States are the guarantee for success.

From the idea through prototype development to series production, 630 highly skilled men and women give their best.

Committed teams of experienced experts from a wide range of specialized fields work on customer-orientated optical solutions – flat, cylindrical, spherical and aspherical – in the spectral ranges EUV, DUV, UV, VIS, NIR.

Glass components are produced and developed in optical and mechanical design, manufacturing, site of optical coating, and adjustment of components, assemblies and the finishing of glass for technical applications. We also produce individual measuring systems.




The precision shaping and surface finishing of glass and similar materials for demanding technological uses is our core business.

Production and measuring processes are constantly being developed and refined to ensure a wide range of cutting-edge applications for our products.

Safeguarding precision is our profession.

To keep it that way, our staff receives ongoing training in comprehensive programs. We have a network with research institutions, suppliers and customers.

## FORM FOLLOWS FUNCTION.

The application and purpose of a glass component alone determine its form – the finishing of shaped glass optimizes its functionality.

We offer technical glass for any conceivable application – mirrors or two-sided antireflection-coated glass for pictures, scanners, auxiliary TV screens or large-screen projectors.

We shape glass and glass ceramics with computer numerical control (CNC) processing centers. Parts can be created for entirely new functions with a high degree of precision and reproducibility.

Our screen printing process permits highest quality printed images on glass and ceramics, metal and composite materials for measuring instruments or design elements.



We can produce glass surfaces in various degrees of roughness and/or gloss structures both chemically and mechanically.

For protective glasses or pressure-loaded display devices we offer flat or curved, single thickness or laminated safety glass. The glass is bonded as required under pressure and heat with transparent high-performance films – which can perform additional functions.

We update our production processes continuously and we're proud of the results.


For standard and custom optical lenses, line-to-standard, printing, precision optics and peak-to-valley heights, we also produce customized etching structures to meet the transmission and reflection requirements of a given glass.

Optimized, product-specific measuring processes subject all work to constant scrutiny.





PRODUCTION:

- Edge processing
- Chemical strengthening
- Screen printing
- Laminating glass and various materials
- Sand blasting
- Etching
- Bending
- Thermal moulding
- CNC shaping, grinding and drilling

INSPECTION:

- Thickness and pinhole measurement
- Interferometer
- Coordinate measuring machines
- 3D surface measurement
- Roughness measurement
- Glass measurement systems

## THE INTELLIGENCE OF LIGHT.

Optical systems determine the path of light. They bundle and disperse, direct and reflect, reduce and enlarge, record and reverse.

Optic solutions are our specialty – from consultation to design, development and production. And for a vast range of applications whether flat or round, cylindrical, spherical or aspherical, complex modules, multifunctional or unmounted systems. Many different materials are used: glass, ceramics,  $\text{CaF}_2$ ,  $\text{MgF}_2$ , germanium, silicon and diverse special materials.


Precision is our forte: we achieve flatnesses of  $\lambda/30$  pv, owing especially to optimized ion beam processing. And we achieve surface roughnesses of  $\text{RMS} \approx 0.4$  nm and less.

Optical substrate coating makes light into a tool. Reflection is reduced, rays directed, mirrored or filtered.

Thin film technology in state-of-the-art processes ensures impeccable results in silvering or antireflection coating for various materials and layer thicknesses – according to customer requirements. This applies also to the manufacture of graduated or non-graduated filters.

We make the greatest demands on mounting and assembly adjustment, which is predominantly restricted to the clean room. This is where, for example,  $\text{CaF}_2$  assemblies for excimer laser are produced.

High quality requires constant control. The individual development of appropriate measurement setups and software adaptation are therefore a matter of course. With test certificates, documented measurement results and thus guaranteed measured precision.

DEVELOPMENT AND DESIGN:



- Optical design
- Mechanical design
- Product development
- Design of measurement set-ups

PRODUCTION:

- Grinding and polishing of
  - optical glass
  - glass ceramics
  - synthetic quartz
  - natural quartz
  - CaF<sub>2</sub>, MgF<sub>2</sub>
  - silicon
  - germanium
  - special materials
- Fine correction techniques
  - ion beam
  - mechanical
- CNC shaping
- CNC centering
- Thin film coating
- Micro-structuring
- Clean-room mounting and adjustment of
  - electro-optical and optomechanical
  - instruments

MEASUREMENT:

- 4" to 24"-interferometer
- 3D measurement (CNC-controlled)
- Sub-nanometer surface roughness test
- AFM power microscope
- Goniometer
- Shack-Hartmann wave front sensor
- Spectral photometer
- Spherometer
- Computer-supported MTF measurement
- Focus measuring device "Optomatic"
- Video precision measuring technology (Smart Scope)

# Enter into the fascinating World of High Performance Optics



Solutions in optics

**and discover the special capability of SwissOptic ...**



 swissoptic

R & D and Production  
Technology

# Peak Performances require advanced experience and know-how.

Your highly demanding tasks will be translated into a full process chain where consultancy, design, engineering, production, testing and documentation can be offered to you. A know-how gathered during several decades and most modern technology allow us to be your all-round partner for optical solutions.

**R & D and Production Technology are areas which offer many novel ideas to you. Please discuss your project with us!**

#### **R & D**

Development of optical components, subassemblies and systems in the spectral ranges UV, VIS, NIR and FIR

#### **Optical components**


- Transmissive and reflective components in the diameter range of 2 to 500 mm
- Multifunctional and monolithic components; several mechanical and optical functionalities are integrated into one glass component.
- R & D and production of thin-film coatings

#### **Subassemblies and Systems**

- Multifunctional optical sensor systems
- Optical Systems of highest precision for a large variety of applications.



#### **Capabilities**

- Optical design and tolerancing
- Opto-mechanical design
- R & D of thin-film coatings
- Finite element analyses
- Environmental simulations
- Prototyping
- Measuring Technologies
- Understanding of International Standards such as MIL, DIN, ISO
- System optimisations based on efficient production technology
- Quality testing devices




# Innovative solutions are made possible by complex shaping.

Plano Optics

Modern CNC-Technology adapted to the needs of modern optical production is our answer to new and unique demands in the field of high-tech optical engineering. One single glass component can integrate several optical and mechanical functions. These fascinating opportunities are part of our daily work.

**Multifunctional and monolithical components make the most novel designs possible, especially in plano optics.**

**Some technical data**

- Processing of different optical materials, such as glass, fused silica, ceramics, germanium, CaF<sub>2</sub>,...
- Angular accuracy up to 1"
- Parallelism up to 1"
- Flatness up to  $\lambda/30$  at  $\lambda = 633$  nm
- Typical dimensions 2 to 500 mm
- High-precision cementing
- Micro-roughness up to 2 Å RMS


**Novel design possibilities due to CNC-manufacture of components**

- Manufacture of complex shapes in glass, fused silica, ceramics
- Several mechanical functionalities integrated into one optical component (multifunctional optics)
- Several optical functionalities integrated into one optical component (monolithical optics)
- Design of components adapted to assembly requirements
- Reduction of number of components in a system

- Integration of fixture devices
- Uniformity of expansion coefficients
- Reduction in volume and weight

**Some typical plano-optical components**

- Filters
- Plano-parallel plates
- Wedges
- Dove prisms
- Pentaprisms
- Polygons
- Roof prisms
- Schmidt/Pechan prisms
- Windows
- Corner cubes
- Beamsplitting plates & cubes
- Inverting systems


# Enthusiasm for dimensions beyond the common.

Spherical Optics

A knowledge beyond pure optics and an understanding for the final application are qualifications required for the development and manufacture of high-tech spherical optics. Our specialised staff and latest CNC manufacturing technology assure extraordinary performances.


**Top quality spherical optics for outstanding insights and outlooks.**

**Some technical data**

- Lens diameters from 2 to 500 mm
- Radii from 1.2 mm to  $\infty$
- Irregularity up to  $\lambda/30$  PV at  $\lambda = 633$  nm
- Centering accuracy  $10''$
- Center thickness tolerance  $\pm 5 \mu\text{m}$
- Diameter tolerance  $\pm 5 \mu\text{m}$




**Additional information**

- Various materials being processed, such as Glass, Fused Silica, Germanium, Silicon, other special materials
- Single and compound lenses
- Special know-how in optical cementing
- Lacquered diaphragms, lacquered circumferences
- Special shaping of contours
- High-precision aspherical components, on and off-axis



# Finest precision created by most modern technologies.


Microstructured patterns and thinfilm coatings

Complete solutions integrating microstructured patterns and coatings are very sought-after. Our offer to you includes also consultancy, R & D and documented measuring results. Very sophisticated manufacturing procedures deliver also a high-volume output with constant quality and extremely high precision.

**Optical components with microstructured patterns and coatings at the limits of feasibility.**

**Reticles, masks, diaphragms**

- Diameters from 2 mm to 140 mm
- Line width  $> = 1 \mu\text{m}$
- Line precision up to  $0.5 \mu\text{m}$
- Black, etched or luminescent

**Circles**

- Encoded circles
- Graduated circles

**Gradient filters**

- Continuous or customer-specific behaviour curve

**Coated patterns**

- Mosaic filters
- Customer specific patterns

**Spectral ranges**

- DUV, UV, VIS, NIR, IR

**Antireflection Coatings**


- Broadband AR
- V-coatings; Reflection  $< 0.1 \%$
- Dual-Band-Coatings

**Filters and Mirrors**

- Long pass filters
- Short pass filters
- Narrowbandpass filters
- Heat protection filters
- Cold mirrors
- Metallic or dielectric mirrors
- Laser mirrors



**Beamsplitting Coatings**

- B/S cubes
- B/S plates
- Polarisers
- Unpolarising B/S coatings





# Seeking for intelligent solutions with new combinations ...

Subassemblies  
and Systems

The combination of components, assemblies and modules, leading to high-tech optical systems, is one of our strengths. Carefully finished and tested components, an inventive spirit, technical know-how and an excellent infrastructure are our advantages leading to convincing optical solutions.

**Complex optical systems for the most demanding applications result from high-tech assembly of individual components.**

#### **Complete Solutions**

- Assembly of optical components, mechanical mounts and electrical components
- Assembly of systems
- System testing and documentation


#### **Assembly-related know-how**

- Gluing
- Clamping
- Lacquering
- Adjustment/Alignment
- Ball chuck centering; centering precision up to 2  $\mu\text{m}$
- Assembly of lens systems; diameter range 2 to 300 mm
- Cleanroom assembly (Flow-boxes)
- Ultrasonic welding of plastics to glass
- Ultrasonic cleaning of optics and mechanics

#### **Measuring Technology**

- Computer based MTF-measurements
- Interferometrical measurements
- Microscopical image resolution tests
- Wavefront measurements by means of Shack-Hartmann sensor
- Measurements by means of goniometers
- Testing of gastightness of assemblies

**Please come to test our Efficiency.**

**R & D**

**Plano Optics**


**Spherical Optics**

**Microstructures**

**Coatings**

**Assemblies**

**Systems**

**Products that integrate optics from SwissOptic are our best references. Your demand will be welcome for detailed discussions beyond this brochure ...**

**[www.swissoptic.com](http://www.swissoptic.com)**  
**Solutions in optics**

SwissOptic AG . Heinrich-Wild-Strasse . CH-9435 Heerbrugg (Switzerland)  
Tel +41 71 727 30 74 . Fax +41 71 727 46 86 . e-mail: [swissoptic@swissoptic.com](mailto:swissoptic@swissoptic.com)

LEISTER®


LEISTER Microsystems

micro-optics


imagine the future of light

**LEISTER**  
Microsystems


## From components

*Refractive lens arrays are provided in a range of diameters and focal lengths.*

*Hexagonal lens arrays provide very high fill factors with minimal space between the lenses.*

*Diffractive fan-out structures are used for multiple beam generation.*

LEISTER Microsystems

LEISTER Microsystems provides a range of standard and customer defined optical microstructures. These include:

- refractive microlens arrays and
- diffractive structures

which are used either as stand-alone items or as enabling components for advanced optical systems. Micromechanical structures such as V-grooves complement these structures as building blocks for micro-optical systems.

### Optical functions

- Collimation
- Fiber coupling
- Beam steering
- Beam shaping
- Diffusing
- Beam homogenizing
- Fan-out / beam splitting

### Materials

To provide an optimal solution for your needs our micro-optical components are available in a range of materials, including:

- Fused silica
- Silicon
- Plastics
- Polymer on glass / silicon / GaAs

## for applications in:

- Telecom
- Displays
- Data storage
- Imaging

## → to solutions

We start with your early ideas to design optimized micro-optical components and systems to suit your needs.

Whether you need rapid prototyping or large volumes of micro-optics our varied in-house fabrication technologies, combined with external partnerships, allow us to offer you a reliable production solution for your micro-optical devices.

### Design

- Optical design and modelling


### Fabrication

- Prototyping and small volume manufacturing
- Large volume fabrication
- Replication on custom specific substrates and parts


### Systems

- Stacked arrays
- MOEMS capabilities
- Custom specific developments
- Assembled subsystems


- Laser material processing
- Metrology
- Optical sensors
- Biomedical

*The ability to measure micro-optics is crucial to achieve a quality product. The interferogram shows a microlens with diffraction limited quality.*

*Good optical design reduces overall development time and ensures optimized final systems.*

*An important application of microlens arrays is fill factor enhancement in CCD cameras.*

LEISTER Microsystems  
ISO 9001:2000

# Microlens Arrays

## Applications Include:

- Laser beam shaping
- Imaging systems
- Displays
- Sensors
- Fiber optics
- Lighting / illumination

## Microlens Array (MLA) Parameters

Array Size	up to 130 mm
Lens Aperture Size	10 $\mu\text{m}$ to >1.5 mm
Aperture Shape	Circular, Hexagonal, Elliptical, Square, Other
F-number	F/0.66 to F/ $\infty$ (NA $\leq$ 0.60)
Materials	Fused Silica, Silicon, Plastics

LEISTER Process Technologies offers high quality refractive microlens arrays in a wide range of focal lengths, aperture sizes, aperture shapes, lens spacing, and arrangements.

Our experienced optics design experts cooperate closely with customers to find optimal standard or OEM solutions for their requests.

Well-controlled ISO 9001:2000 certified production processes, with dedicated equipment and specialized characterization tools, combined with 50 years of manufacturing experience, make us a reliable supplier.

## Contact Information

### Headquarters

LEISTER Process Technologies  
Microsystems Division

Riedstrasse

CH-6060 Sarnen/Switzerland

Phone +41 41 662 74 74

Fax +41 41 660 20 61

e-mail [microsystems@leister.com](mailto:microsystems@leister.com)  
[www.leister.com](http://www.leister.com)

### USA

LEISTER Technologies, LLC

1253 Hamilton Parkway

Itasca, IL 60143/USA

Phone (630) 760 1000

Fax (630) 760 1001

e-mail [sales@leisterusa.com](mailto:sales@leisterusa.com)  
[www.leisterusa.com](http://www.leisterusa.com)

LEISTER Microsystems Exclusive Distributor



LEISTER Microsystems  
ISO 9001:2000

# Fiber Collimating/Coupling Lens Arrays

*Microlens arrays for collimating and  
coupling of single mode optical fibers*

## **Benefits**

- Standard 250  $\mu\text{m}$  pitch fiber array compatible
- Focal length uniformity  $< \pm 0.45\%$  lens-lens
- Uniform pointing performance
- High surface quality, very low roughness
- Fused Silica material, Telcordia compliant
- AR coating
- 100 % inspection

## **Applications**

Collimated optical fiber arrays for various multi-channel communications devices, including DWDM MUX, DeMUX and OADM, optical amplifiers, VOAs, optical isolators, as well as waveguide to fiber coupling.

LEISTER Microsystems  
ISO 9001:2000

### **FCA 250 Fiber Collimating Lens Array Preliminary Specifications**

Array pitch	250 $\mu$ m	Tolerance: $\pm$ 0.25 $\mu$ m
Number of lenses per array	FCA 250-4/720: 4 $\times$ 1 FCA 250-8/720: 8 $\times$ 1 FCA 250-16/720: 16 $\times$ 1	Larger arrays available upon request
Lens array dimensions	FCA 250-4/720: 1.5 mm $\times$ 1.0 mm $\times$ 1.0 mm FCA 250-8/720: 2.5 mm $\times$ 1.0 mm $\times$ 1.0 mm (length $\times$ width $\times$ thickness)	Tolerances: L, W $\pm$ 0.05 mm, T $\pm$ 0.05 mm
Lens diameter	240 $\mu$ m	
Lens surface profile	Optimized aspheric design for peak performance	
Surface figure	< 35 nm RMS	
Numerical aperture	N.A. = 0.16	Custom designs available, e.g. for fiber / waveguide coupling
Effective focal length	720 $\mu$ m (FCA 250-x/720)	Tolerances: < $\pm$ 0.45 % within array
Material	Fused Silica (Schott Lithosil® Q1)	
AR coating	Double sided, for $\lambda$ = 1280-1620 nm	R < 0.5 % per surface
Surface roughness	Ra < 2 nm	
Insertion loss	< 1 dB collimator to collimator	

Customized microlens arrays are available upon request, optimized for your specific element dimensions, alignment marks, as well as pitch, focal length and throw distance.

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
# SAC

## Slow Axis Collimator Lens Array

*Array of cylindrical microlenses for collimating the «slow axis» of diode laser beams. The optimal solution for high brightness diode laser sources.*


### Benefits

- Compensation of beam astigmatism
- High power resistant  
Fused Silica material
- High surface quality and  
very low surface roughness  
minimize loss
- Compatible with typical  
high power diode laser bars
- Compatible with standard  
fast axis collimators
- Stackable

### Applications

Enhances fiber coupled laser diodes, diode sources for solid-state laser pumping, diode laser direct welding systems, and other high power diode lasers.


LEISTER Microsystems  
ISO 9001:2000

### SAC 500 Specifications

SAC type	SAC 500/1.5/x	SAC 500/1.1/x	
Array pitch (emitter spacing)	500 $\mu\text{m}$	500 $\mu\text{m}$	Other pitches available upon request
Number of lenses per array	23 $\times$ 1	21 $\times$ 1	
Lens array substrate size	12.0 mm $\times$ 1.5 mm $\times$ 0.6 mm (length $\times$ width $\times$ thickness)	11.0 mm $\times$ 1.1 mm $\times$ 0.6 mm (length $\times$ width $\times$ thickness)	Tolerances: L,W $\pm$ 0.05 mm, T $\pm$ 0.05 mm
Lens surface profile	Cylindrical	Cylindrical	
Effective focal length	1.8 mm (SAC 500/1.5/1.8) 2.0 mm (SAC 500/1.5/2.0) Other values available upon request	1.8 mm (SAC 500/1.1/1.8) 2.0 mm (SAC 500/1.1/2.0) Other values available upon request	Tolerances: < $\pm$ 0.45 % within array < $\pm$ 5 % array to array
Material	Fused Silica (Schott Lithosil® Q1)	Fused Silica (Schott Lithosil®Q1)	
AR coating	Double sided, for $\lambda =$ 770-1080 nm	Double sided, for $\lambda =$ 770-1080 nm	R < 0.5 % per surface
Surface roughness	Ra < 2 nm	Ra < 2 nm	

### Typical surface profile

### Contact Information

#### Headquarters

LEISTER Process Technologies  
Microsystems Division  
Riedstrasse  
CH-6060 Sarnen/Switzerland  
Phone +41 41 662 74 74  
Fax +41 41 660 20 61  
e-mail [microsystems@leister.com](mailto:microsystems@leister.com)  
[www.leister.com](http://www.leister.com)

#### USA

LEISTER Technologies, LLC  
1253 Hamilton Parkway  
Itasca, IL 60143/USA  
Phone (630) 760 1000  
Fax (630) 760 1001  
e-mail [sales@leisterusa.com](mailto:sales@leisterusa.com)  
[www.leisterusa.com](http://www.leisterusa.com)

LEISTER Microsystems Exclusive Distributor