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**Japan Agricultural Standards
Testing and Analysis Handbook Series**

**Instruction Manual for
Testing and Analyzing Genetically Modified Foods
— Part 2 —
Basic Procedures**

June 20, 2002

The Center for Food Quality, Labeling and Consumer Services

I. Introduction

Based on Article 7, Section 1 of the Standards for the Labeling of Processed Foods Derived from Genetically Modified Crops promulgated under the Law for the Standardization and Labeling of Agricultural Products (1950 Law No. 175), so-called the “JAS Law,” and the Standards (MAFF Notice No. 517 of March 31, 2000) published by the Minister of Agriculture, Forestry and Fisheries under Article 7, Section 1 of the Standards for the Labeling of Fresh Food, food labeling for the contents of recombinant materials was mandated in Japan on April 1, 2001. Labeling food products as “non-GMO” is made possible through use of the identity preserved (IP) handling system, in which GM and non-GM crops are strictly segregated and documented throughout the life of the crops at issue from production through transportation, importation and processing to retail distribution. The Center for Food Quality, Labeling and Consumer Services (hereinafter referred to as the “Center”) is responsible for, as part of the JAS regulations compliance monitoring program, the testing of food products for the contents of recombinant materials. The Instruction Manual for Testing and Analyzing Genetically Modified Foods (hereinafter referred to as the “Manual”) were developed to establish the standardized methods for this purpose.

This part of the Manual describes the basic procedures involved in the *qualitative* PCR screening to determine the presence or absence of recombinant materials in food crops and processed foods derived from such crops, according to the stages of the screening process.

1. Purchase and documentation of food samples
2. Sample preparation
3. DNA extraction
4. PCR amplification and gel electrophoresis
5. Determination as either GM or non-GM food.

II. General precautions

PCR is so powerful that even a trace amount of a template DNA can be amplified to a detectable level. It is thus critical to prevent the contamination of samples by undesired DNA, especially carryover PCR products, and by the DNases excreted from human skin that enzymatically degrade PCR samples. While *Part 5—Prevention of Contamination* of the Manual describes in detail the proper procedures to prevent contamination, the following general precautions may be useful:

- (1) Autoclave all solutions other than those that are sensitive to heat. For the purpose of this Manual, “pure water” means deionized water with a conductivity of 0.0056 mS/m or below at 25 °C, and “sterilized water” means the pure water autoclaved at 121 °C for no less than 15 minutes.
- (2) Use disposable micropipette tips and sample tubes to avoid cross-contamination.
- (3) Use micropipettes, tips and 1.5- and 0.2-ml sample tubes only after they have been autoclaved and completely dried in a dryer or, if possible, dry heat-sterilized.
- (4) Prior to handling DNA samples, disinfect the surface of the laboratory bench with ethanol. Always wear rubber gloves also disinfected with ethanol. Gloves should

be powder-free. If your gloves are powdered, wash your hands thoroughly after wearing the gloves. Gloves should be ethanol-disinfected frequently during the procedure. Use a laminar flow clean bench if necessary.

- (5) It is recommended that each analyst re-prepare his/her own sterilized water and TE buffer at least once every month. Note that the quality of sterilized water and TE buffer is critical to successful PCR. DNase contamination, among others, can have extensive and serious effects on PCR including failed reactions and false results.

III. Handling and management of reagents

Some of the chemicals and compounds used in laboratory activities involving genes are powerful mutagens. Proper management of reagents and effluents is thus imperative for laboratories engaged in GM food testing. Handle and manage them properly in accordance with the instructions of reagent suppliers' manuals and *Part 4—Preparation of Reagents* of the Manual.

IV. Basic procedures

1 Purchase and documentation of food samples

See *Part 1—Handling of Test Samples* of the Manual for the documentation of food samples. For a food product to be tested, purchase three packages each on the retail market.

2 Sample preparation

See *Part 3—Specific Procedures by Food Product* of the Manual for the preparation of test samples.

3 DNA extraction

One DNA sample is extracted from each package of the food product to be tested. In other words, three DNA samples will normally be taken for a food product from three different packages.

DNA samples may be extracted from prepared food samples using spin columns with silica matrices, ion exchange columns or cetyltrimethylammonium bromide (CTAB). Of these options, the spin and ion exchange column methods are recommended in view of obtained DNA sizes that fit for PCR as well as health and environmental considerations. Described below are the extraction procedures using Qiagen's DNeasy Plant Maxi kit as an example of the spin column method and those using Qiagen's Genomic tip 20/G as an example of the ion exchange column method. While these kits and procedures are given only as examples, any alternative method to be selected should be able to provide the same quality^(Note 1) of DNA extracts.

^(Note 1) For the purpose of Stage 3 "DNA extraction" of this part of the Manual, the "same quality" means that not only the obtained DNA solution has the same purity as determined by the $OD_{260\text{ nm}}/OD_{230\text{ nm}}$ ratio, but also the selected method successfully extracts DNA samples suitable for PCR from each of the long array of food products listed in *Part 3—Specific Procedures by Food Product* of the Manual. The Center verifies the suitability of a method by comparing the results from the quantification of the DNA sequence specific to an endogenous gene of a crop plant,

Precautions

- Use sterilized laboratory utensils and equipment.
- Take all necessary measures to prevent contamination, such as wearing disinfected gloves.

3.1 DNA extraction using DNeasy Plant Maxi kit

3.1.1 Brief description of the procedure

This procedure involves the lysis of the ground food sample in a lysis buffer, digestion by RNase, and removal of proteins and other impurities, followed by adsorption of DNA to a silica-gel membrane in the presence of high concentrations of chaotropic salt, and finally elution of purified DNA in a low-salt buffer or sterilized water.

3.1.2 Authority

Kuribara H., Shindo Y., Matsuoka T., Takubo K., Futo S., Aoki N., Hirao T., Akiyama H., Goda Y., Toyoda M., Hino A.: Quantification Methods using Novel Reference Molecules for Detection of Genetically Modified Maize and Soybean, *Journal of AOAC International* (in press).

3.1.3 Applicability of the procedure to different food products

See Part 3—Specific Procedures by Food Product of the Manual

3.1.4 Equipment

- Centrifuge with a swing-bucket rotor and capable of spinning 50 ml polypropylene tube at 3,000 x g;
- Refrigerated centrifuge capable of spinning 2 ml microtubes;
- Micropipettes for volume ranges of 0.5-10, 10-100, 100-1,000 and 1,000-5,000 $\mu\text{l}^{(\text{Note } 2)}$;
- Water bath; and
- Vortex test tube mixer.

3.1.5 Reagents

Qiagen's DNeasy Plant Maxi kit.

using the method described in Part 6 - Quantitative PCR of the Manual, that are present in the DNA extracts from the quoted and alternative extraction methods.

^(Note 2) These volume ranges are only examples (??)

3.1.6 Extraction procedure

3.1.6.1 Procedure A: DNA extraction using DNeasy Plant Maxi kit

- (1) Measure an appropriate amount of the ground sample into a 50 ml tube and add 20 μ l RNase (supplied in the kit) using a 10-100 μ l micropipette and 10 ml Buffer AP1 preheated to 65°C using a 1,000-5,000 μ l micropipette. After mixing the contents of the tube by repeatedly turning the tube upside down until no visible sample remains on the bottom of the tube, mix the solution on a vortex tube mixer.
- (2) Incubate the solution for one hour in a 65°C water bath. Every 15 minutes during incubation or 15, 30 and 45 minutes after incubation started, mix the solution by vigorously inverting the tube several times and then mix it at maximum speed on a vortex tube mixer.
- (3) Spin the solution at 3,000 x g for 10 minutes at room temperature in a swing-out rotor.
- (4) Transfer 7 ml of the supernatant into a new 15 ml tube using a 1,000-5,000 μ l micropipette, carefully avoiding the precipitate on the bottom of the tube or film on the liquid surface.
- (5) Add 2.5 ml AP2 buffer to the 15 ml tube using a 1,000-5,000 μ l micropipette, mix the solution at maximum speed on a vortex tube mixer for 10 seconds, and incubate in an ice-water bath for 15 minutes.
- (6) Spin the solution at 3,000 x g for 35 minutes at room temperature in a swing-out rotor.
- (7) Aspirate 8 ml of the supernatant using a 1,000-5,000 μ l micropipette, carefully avoiding the precipitate on the bottom of the tube or film on the liquid surface, to load it onto the QIA shredder spin column (lilac).
- (8) Spin the solution at 3,000 x g for 5 minutes at room temperature in a swing-out rotor.
- (9) Transfer 7.5 ml of the supernatant into a new 50 ml tube using a 1,000-5,000 μ l micropipette, carefully avoiding the precipitate on the bottom of the column.
- (10) After mixing on a vortex tube mixer at maximum speed for 10 seconds, transfer 6.8 ml of the solution into a new 50 ml tube using a 1,000-5,000 μ l micropipette.
- (11) Add 10.2 ml AP3/Et-OH buffer using a 1,000-5,000 μ l micropipette, mix the solution on a vortex tube mixer at maximum speed for 10 seconds, and load the entire solution onto the DNeasy spin column (colorless) by decanting.
- (12) Spin the column at 3,000 x g for 15 minutes at room temperature in a swing-out rotor and discard the eluate.
- (13) Add 12 ml Buffer AW to the spin column using a 1,000-5,000 μ l micropipette and spin the column at 3,000 x g for 15 minutes at room temperature in a swing-out rotor.
- (14) Place the spin column in a new 50 ml tube and add 1 ml sterilized water preheated to 65°C in a separate spin column.

- (15) After leaving at room temperature for 5 minutes, spin at 3,000 x g for 10 minutes at room temperature in a swing-out rotor.
- (16) Using a 100-1,000 µl micropipette, measure the volume of the eluate and transfer it to a 2 ml sample tube. Add an equal volume of isopropanol to the tube using a 100-1,000 µl micropipette.
- (17) After mixing by slowly turning the tube upside down 10 times, leave the solution at room temperature for five minutes.
- (18) After spinning the tube at 12,000 x g at 4°C for 15 minutes in a swing-out rotor, discard the supernatant using a 100-1,000 µl micropipette.
- (19) Add 500 µl 70% ethanol using a 100-1,000 µl micropipette and tap the end of the tube with a finger until the precipitate on the bottom of the tube is completely resuspended in the solution.
- (20) After spinning at 12,000 x g for 3 minutes at 4°C in a swing-out rotor, discard the entire supernatant using a 10-100 µl micropipette and dry the precipitate.
- (21) Add 100 µl TE buffer, pH 8.0, using a 10-100 µl micropipette to dissolve the precipitate. See 3.1.10 *Additional notes* for further information.
- (22) Repeat tapping the tube with a finger and spinning it until the liquid drops on the tube wall completely disappear. Incubate the solution overnight (12-24 hours) in a refrigerator.
- (23) A DNA lysate is finally produced when the solution is free from any visible insoluble matter. If any visible insoluble matter is still present even after 24 hours of incubation, use the supernatant obtained by spinning the solution at 12,000 x g for 3 minutes at 4°C and transferred to a new tube as a DNA lysate. Store the DNA lysate as well as the precipitate at -20°C or below.

3.1.6.2 Procedure B: DNA extraction using DNeasy Plant Maxi kit

- (1) Measure an appropriate amount of the ground sample into a 50 ml tube and add 10 µl RNase (supplied in the kit) using a 0.5-10 µl micropipette and 5 ml Buffer AP1 preheated to 65°C using a 1,000-5,000 µl micropipette. After mixing the contents of the tube by repeatedly turning the tube upside down until no visible sample remains on the bottom of the tube, mix the solution on a vortex tube mixer.
- (2) Incubate the solution for one hour in a 65°C water bath. Every 15 minutes during incubation or 15, 30 and 45 minutes after incubation started, mix the solution by vigorously inverting the tube several times and then mix it at maximum speed on a vortex tube mixer.
- (3) After adding 1.8 ml Buffer AP2 to the tube using a 1,000-5,000 µl micropipette and mixing the solution at maximum speed for 10 seconds on a vortex tube mixer, incubate the tube in an ice-water bath for 15 minutes.
- (4) Spin the solution at 3,000 x g for 15 minutes at room temperature in a swing-out rotor.
- (5) Aspire 4.2 ml of the supernatant using a 1,000-5,000 µl micropipette, carefully avoiding the precipitate on the bottom of the tube or film on the liquid surface, to

load it onto the QIA shredder spin column (lilac).

- (6) Spin the solution at 3,000 x g for 5 minutes at room temperature in a swing-out rotor.
- (7) Transfer 4 ml of the supernatant into a new 50 ml tube using a 1,000-5,000 μ l micropipette, carefully avoiding the precipitate on the bottom of the column.
- (8) After mixing on a vortex tube mixer at maximum speed for 10 seconds, transfer 3.4 ml of the solution into a new 50 ml tube using a 1,000-5,000 μ l micropipette.
- (9) Add 5.1 ml AP3/Et-OH buffer using a 1,000-5,000 μ l micropipette, mix the solution on a vortex tube mixer at maximum speed for 10 seconds, and load the entire solution onto the DNeasy spin column (colorless) by decanting.
- (10) Spin the column at 3,000 x g for 5 minutes at room temperature in a swing-out rotor and discard the eluate.

For the subsequent steps (11) through (21) under 3.1.6.2 *Procedure B: DNA extraction using DNeasy Plant Maxi kit*, refer steps (13) through (23), respectively, under 3.1.6.1 *Procedure A: DNA extraction using DNeasy Plant Maxi kit*.

3.1.7 Determination of the purity and yield of extracted DNA

Aspirate 5 μ l of the DNA lysate solution using a 0.5-10 μ l micropipette and add TE buffer to give a total volume of 50 μ l. Scan the UV absorbance in the wavelength range from 200 to 300 nm to measure the absorbance at 230, 260 and 280 nm^(Note 3) ^(Note 4). Calculate the concentration of DNA by converting an absorbance (O.D.) of 1.0 at 260 nm to 50 ng of DNA per μ l^(Note 5). Prepare a PCR solution with a DNA concentration of 10ng/ μ l using sterilized water. If the concentration of the extracted DNA was too low for use as a PCR solution, perform the DNA extraction procedure again. If the concentration of the re-extracted DNA was still too low for use as a PCR solution, use the undiluted DNA solution as a PCR solution.

3.1.8 Purity of extracted DNA

The purities of the soybean and maize DNA extracted by these procedures as determined by the OD_{260 nm}/OD_{280 nm} ratios will be approximately 1.7-2.0.

3.1.9 Record keeping

Record the dilution factor, absorbance readings, and OD_{260 nm}/OD_{280 nm} ratio. Also record in detail any deviation from the standard procedures you used to obtain the DNA solution for PCR.

3.1.10 Additional notes

In Step 21 under 3.1.6.1 *Procedure A: DNA extraction using DNeasy Plant Maxi kit*,

^(Note 3) The absorbance spectrum for DNA has a dip at 230 nm and a peak at 260 nm whereas impurities such as proteins have highest absorbances at around 280 nm.

^(Note 4) Dilution should be adjusted depending on the DNA yield of an extraction procedure.

^(Note 5) Sambrook, J., Russel, D. W. *Molecular cloning: a laboratory manual*, 3rd Ed. (volume 3), Cold Spring Harbor, N Y. Cold Spring Harbor Laboratory Press, 2001, A8.20. (ISBN 0-87969-577-3 (pbk), ISBN 0-87969-576-5 (cloth)).

adjust the dilution factor depending on the yield of extracted DNA. We recommend using 100 µl TE buffer for soybeans and 50 µl TE buffer for maize and maize products.

3.1.10 Additional notes

In Step 21 under 3.1.6.1 *Procedure A: DNA extraction using DNeasy Plant Maxi kit*, adjust the dilution factor depending on the yield of extracted DNA. We recommend using 100 µl TE buffer for soybeans and 50 µl TE buffer for maize and maize products.

If the concentration of the DNA solution was too low for use in PCR, try the following operations:

1. Concentrate the DNA solution by ethanol precipitation.
2. Perform the DNA extraction procedure from the beginning in accordance with this protocol except that you add 20 µl TE buffer rather than 100 µl in Step 21 under 3.1.6.1 *Procedure A: DNA extraction using DNeasy Plant Maxi kit*.

If the above operations still do not provide the concentration of DNA solution the PCR procedure normally requires, use the DNA solution finally obtained as a starting DNA solution for PCR. Record the volume of DNA in the solution used in PCR.

3.2 DNA extraction using QIAGEN Genomic-tip 20/G

3.2.1 Brief description of the procedure

This procedure involves the lysis of the ground food sample in a lysis buffer, digestion by RNase, and removal of proteins and other impurities, followed by adsorption of DNA to an ion-exchange column in the presence of low concentrations of chaotropic salt, and finally elution of purified DNA in a high-salt buffer.

3.2.2 Authority

The QIAGEN Genomic-tip protocol for DNA isolation has been modified to meet the requirements of the Ministry of Health, Labor and Welfare (MHLW) guidelines, *Testing for Foods Produced by Recombinant DNA Techniques* (MHLW Food Announcement No. 110 of 27 March 2001), and the amendment to the guidelines (MHLW Food Announcement No. 158 of 25 May 2001).

3.2.3 Applicability of the procedure to different food products

See *Part 3—Specific Procedures by Food Product* of the Manual

3.2.4 Equipment

- Centrifuge with a swing-out rotor and capable of spinning 50 ml polypropylene tube at 3,000 x g;
- Refrigerated centrifuge capable of spinning 2 ml microtubes;
- Micropipettes for volume ranges of 0.5-10, 10-100, 100-1,000 and 1,000-5,000

$\mu\text{l}^{(\text{Note } 6)}$;

- Water bath; and
- Vortex test tube mixer.

3.2.5 Reagents

QIAGEN Genomic-tip 20/G (Cat. No. 10223) ion-exchange column;

QIAGEN RNase A (Cat. No. 19101);

QIAGEN Proteinase K (Cat. No. 19131 or 19133);

Buffers G2, QBT, QC and QF^(Note 7).

3.2.6 Extraction procedure: DNA extraction using QIAGEN Genomic-tip 20/G

- (1) Measure an appropriate amount of the ground sample into a 50 ml tube, add 7.5 ml Buffer G2 using a 1,000-5,000 μl micropipette, and vigorously mix the solution on a vortex tube mixer.
- (2) Add another 7.5 ml of Buffer G2 using a 1,000-5,000 μl micropipette, 200 μl QIAGEN Proteinase K using a 100-1,000 μl micropipette, and 20 μl RNase A using a 10-100 μl micropipette to the tube. After mixing the contents of the tube by repeatedly turning the tube upside down until no visible sample remains on the bottom of the tube, mix the solution on a vortex tube mixer.
- (3) Incubate the solution for one hour in a 50°C water bath. Every 15 minutes during incubation or 15, 30 and 45 minutes after incubation started, mix the solution by vigorously inverting the tube several times and then mix it at maximum speed on a vortex tube mixer.
- (4) Spin the solution at 3,000 x g for 15 minutes at 4°C in a swing-out rotor.
- (5) Transfer the entire supernatant into either a new 15 ml or 50 ml tube using a 1,000-5,000 μl micropipette, carefully avoiding the precipitate on the bottom of the tube or film on the liquid surface.
- (6) Flash centrifuge the tube.
- (7) Equilibrate QIAGEN Genomic-tip 20/G with 1 ml Buffer QBT.
- (8) Load the entire supernatant onto the QIAGEN Genomic-tip 20/G in 2 ml aliquots using either a 100-1,000 or 1,000-5,000 μl micropipette, carefully avoiding the precipitate on the bottom of the tube. Allow the supernatant to move through the QIAGEN Genomic-tip 20/G by gravity flow.
- (9) Wash the column by loading 2 ml Buffer QC onto the tip using either a 100-1,000 or 1,000-5,000 μl micropipette and allowing the buffer to move through the QIAGEN Genomic-tip 20/G by gravity flow.

^(Note 6) These volume ranges are only examples. (??)

^(Note 7) Prepare these buffers according to the instructions of the supplier's protocol. You may also purchase either Qiagen's Blood & Cell Culture DNA Mini Kit (Cat. No. 13323) or Qiagen's Genomic DNA buffer set (Cat. No. 19060) in which these buffers are supplied.

- (10) Repeat Step (9) washing operation additional two times.
- (11) Place the tip over a 1.5 ml tube and add 750 μ l Buffer QF prewarmed to 50°C using a 100-1,000 μ l micropipette to elute DNA (Eluate 1).
- (12) Place the tip over a new 1.5 ml tube and add again 750 μ l Buffer QF prewarmed to 50°C using a 100-1,000 μ l micropipette to elute DNA (Eluate 2).
- (13) Measure the volumes of Eluates 1 and 2 and add to each eluate an equal volume of isopropanol using a 100-1,000 μ l micropipette. After mixing the eluates by gently inverting the tubes 10 times, leave them at room temperature for 5 minutes.
- (14) After spinning at 12,000 x g at 4°C for 15 minutes in a fixed angle rotor, discard the supernatant using a 200-1,000 μ l micropipette.
- (15) Add 1,000 μ l 70% ethanol using a 100-1,000 μ l micropipette and mix by slowly turning the tube upside down 10 times.
- (16) After spinning at 12,000 x g at 4°C for 3 minutes in a fixed angle rotor, discard the entire supernatant using either a 10-100 or 100-1,000 μ l micropipette and dry the precipitate.
- (17) Add 50 μ l TE buffer (pH 8.0) to the tube containing Eluate 2 using a 10-100 μ l micropipette and resuspend the precipitate by shaking at 65°C for 15 minutes on a shaker.
- (18) Transfer the entire contents of the tube containing Eluate 2 to the one containing Eluate 1 using a 10-100 μ l micropipette and resuspend the DNA by shaking at 65°C for 15 minutes on a shaker.
- (19) Tap the tube with a finger and place it in a refrigerator for incubation for 12-24 hours.
- (20) A DNA lysate is finally produced when the solution is free from any visible insoluble matter. If any visible insoluble matter is still present even after 24 hours of incubation, use the supernatant obtained by spinning the solution at 12,000 x g for 3 minutes at 4°C and transferred to a new tube as a DNA lysate. Store the DNA lysate as well as the precipitate at -20°C or below.

3.2.7 Determination of the purity and yield of extracted DNA

Same as 3.1.7.

3.2.8 Purity of extracted DNA

The purities of the soybean and maize DNA extracted by this procedure as determined by the $OD_{260\text{ nm}}/OD_{280\text{ nm}}$ ratios will be approximately 1.7-2.0.

3.2.9 Record keeping

Same as 3.1.9.

3.2.10 Additional notes

Same as 3.1.10, except that Step 21 under 3.1.6.1 *Procedure A: DNA extraction using DNeasy Plant Maxi kit* should be read as Step 17 under 3.2.6 *Extraction procedure: DNA extraction using QIAGEN Genomic-tip 20/G*.

3.3 DNA extraction using CTAB

3.2.1 Brief description of the procedure

This procedure involves the lysis of the ground food sample in a lysis buffer, removal of proteins and other impurities, and digestion by RNase to extract DNA. Use of CTAB as the eluent effectively prevents polysaccharides from being extracted.

3.3.2 Authority

Murray, M.G., Thompson, W., *Rapid isolation of high molecular weight plant DNA*, *Nucleic Acids Res.*, 8, 4321-4325(1980), with partial modification.

3.3.3 Applicability of the procedure to different food products

See Part 2—Basic Procedures, Chapter 1 of the Manual

3.3.4 Equipment

- Refrigerated centrifuge capable of spinning 2 ml microtubes;
- Micropipettes for volume ranges of 0.5-10, 10-100, 100-1,000 and 1,000-5,000 μl ^(Note 8);
- Water bath; and
- Vortex test tube mixer.

3.3.5 Reagents

See _____ of the Manual.

3.3.6 Extraction procedure: DNA extraction using CTAB

(1) Cell lysis of the sample

Weigh out an appropriate amount of the sample into a mortar^{(Note 9)(Note 10)}. Add a small amount of silica sand and 2 ml CTAB extraction buffer to the mortar and

^(Note 8) These volume ranges are only examples. (??)

^(Note 9) An excess amount of the sample will lead to an increased volume of the intermediate layer during phenol-mediated removal of proteins and make the subsequent operations more difficult.

^(Note 10) In weighing the sample, you can use the aluminum foil that has been used to wrap sterilized mortar as a replace for a glassine weighing paper. Do not pick up the samples directly with your fingers. Always use a sterilized spatula to scoop and transfer the sample.

grind them with a pestle. Transfer the ground mixture into 1.5 ml tubes^(Note 11).

After incubating at 60°C for 30 minutes, centrifuge at 14,000 rpm for three minutes^(Note 12). (离心 1~2 次)

Transfer 700 µl of the supernatant to a new tube.

(2) PCI-mediated removal of proteins

Add an equal volume of PCI to the tube containing the sample solution, vigorously shake the tube for two minutes, and then centrifuge at 14,000 rpm for 15 minutes at room temperature^(Note 13).

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To remove phenol from the aqueous layer, add an equal volume of CIA to the tube containing the sample solution, vigorously shake the tube for two minutes, and then centrifuge at 14,000 rpm for three minutes at room temperature.

Remove the upper layer into a new tube. Pipette carefully to avoid drawing the intermediate layer. 将上层吸出, 转入 / new tube. 避免中间层

(4) Alcohol precipitation

To precipitate DNA, add an equal volume of isopropanol to the tube containing the sample solution, invert the tube for 30 seconds, and then centrifuge for three minutes at 12,000 rpm or about 13,000 x g^(Note 14).

Discard the supernatant.

Wash the precipitate by adding 800 µl 70% ethanol, mix by inverting, incubate for three minutes, and centrifuge for three minutes at 12,000 rpm.

Discard the supernatant^(Note 15) and vacuum dry for five minutes using a centrifugal vacuum concentrator or a desktop desiccator. Visually confirm the dehydration of the sample.

(5) Dissolving DNA in TE buffer and removing RNA with CTAB eluent

Add 100 µl TE and 2 µl RNase A (10 mg/mL) to dissolve DNA.

After incubating the solution for 30 minutes at room temperature or 37°C, add 400 µl CTAB eluent. (CTAB extraction bu (f.c.e.))

^(Note 11) Proteinase digestion: If the sample is expected to have a high protein content and produce a thick intermediate layer during the PCI treatment, you can add 20 µl of a 20 mg/ml proteinase K solution to each tube to decrease the volume of the protein layer.

^(Note 12) Here, the centrifuge is assumed to be Eppendorf Centrifuge 5417R at 16,000 x g but use of the maximum speed with your centrifuge will generally produce satisfactory results.

^(Note 13) Record in detail the changes in the contents of the tube. If the solution does not successfully separate into layers, repeat either the spinning operation or the PCI-mediated removal of proteins. Perform all spinning operations at room temperature; lower temperatures will cause CTAB to precipitate and the solution does not separate properly.

^(Note 14) These conditions may need to be modified depending on salt concentrations and polysaccharide contents.

^(Note 15) After the supernatant is discarded, flash centrifuge the tube again at 5,000-12,000 rpm to thoroughly remove the supernatant. If the resulting precipitate is gel-like, wash the precipitate with alcohol several times to dehydrate it to some extent.

(6) Re-extraction with CIA

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Add 500 μ l CIA and mix briefly. Centrifuge for 15 minutes at 12,000 rpm and transfer the top layer to a new tube.

(7) Alcohol precipitation

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To remove phenol from the aqueous layer, add an equal volume of isopropanol to the tube containing the sample solution, invert the tube gently for 30 seconds, and then centrifuge for three minutes at 12,000 rpm or about 13,000 \times g^(Note 14).

Discard the supernatant^(Note 15) and vacuum dry for five minutes using a centrifugal vacuum concentrator or a desktop desiccator. Visually confirm the dehydration of the sample.

(8) Dissolving DNA in sterilized water

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Add 100 μ l sterilized water^(Note 16) to make a DNA solution. Divide the solution into several single-use aliquots for frozen storage at -20°C or below because repeated freezing and thawing is not advisable.

3.3.7 Determination of the purity and yield of extracted DNA

Same as 3.1.7.

3.3.8 Purity of extracted DNA

The purities of the soybean and maize DNA extracted by these procedures as determined by the OD_{260 nm}/OD_{280 nm} ratios will be approximately 1.8-2.0.

3.3.9 Record keeping

- (1) See 3.1.9. Also record the following:
- (2) The changes in the contents of the tube in Step (2) *PCI-mediated removal of proteins*; and
- (3) The manner in which precipitation proceeds in Step (4) *Alcohol precipitation* when an equal volume of isopropanol is added to the sample solution, the mixed solution is centrifuged and the supernatant is discarded.

3.3.10 Additional notes

Not applicable.

^(Note 16) In making the DNA solution, use of distilled water rather than TE buffer is recommended because EDTA in TE buffer can attract magnesium ion which can negatively affect PCR reactions.

4. PCR amplification and gel electrophoresis

4.1 Brief description of the procedure

PCR is performed using the extracted DNA as a template, a DNA polymerase and primers that specifically bind to the target sequences of genetically modified crops. The PCR products are then loaded on electrophoresis and visualized under UV light to determine whether the sample contains genetically modified crop varieties. At the same time, to determine that the target PCR products may be obtained, PCR amplification tests are performed using the primer pairs designed to amplify crop-specific endogenous genes.

Safety considerations for ethidium bromide

Personal protection:

The fluorescent dye ethidium bromide is a potent mutagen that can intercalate into a DNA strand. Always wear gloves when working with this dye and a mask when weighing the powder form.

Waste disposal:

For ethidium bromide solutions of lower concentrations, treat them using commercially available decontamination devices. For higher-concentration solutions, dispose of them through specialized hazardous waste treatment service providers.

4.2 Authority

- (1) For PCR, and the primers for soybean and maize samples, Kuribara H., Shindo Y., Matsuoka T., Takubo K., Futo S., Aoki N., Hirao T., Akiyama H., Goda Y., Toyoda M., Hino A.: *Quantification Methods using Novel Reference Molecules for Detection of Genetically Modified Maize and Soybean*, Journal of AOAC International (in press).
- (2) For the primers for potato samples, the Ministry of Health, Labor and Welfare (MHLW) guidelines, *Testing for Foods Produced by Recombinant DNA Techniques* (MHLW Food Announcement No. 110 of 27 March 2001), and the amendment to the guidelines (MHLW Food Announcement No. 158 of 25 May 2001).
- (3) For the electrophoresis procedures, the QIAGEN Genomic-tip protocol for DNA isolation.

4.3 Applicability of the procedure to different food products

This procedure may be used to:

- (1) specifically detect the recombinant Roundup Ready soybean event 40-3-2;
- (2) specifically detect and screen for the five maize recombinants Bt11, Event 176, T25, MON810 and GA21; and
- (3) specifically detect the potato recombinants NewLeaf events Bt6 and SPBT02-05 and NewLeaf Plus events RBMT21-129, RBMT21-350 and RBMT22-82;

Note that this procedure screens maize samples for the 35S promoter from cauliflower mosaic virus and thus is not applicable to those processed food products that contain both maize and other plants that may be transformed using the promoter.

4.4 Equipment

4.4.1 PCR

- Thermal cycler:
Optional thermal cyclers include MJ Research's PTC-200 DNA Engine, Takara Bio's TaKaRa PCR Thermal Cycler MP, Applied Biosystems's GeneAmp system 9700 with the MAX heating mode, and others that provide the same PCR results as the former three.
- Micropipettes in the volume ranges of, for example, 0.5-10, 10-100, 100-1,000 and 1,000-5,000 μ ls.

4.4.2 Gel electrophoresis

- Gel electrophoresis apparatus:
Advance's Mupid[®]-2 gel electrophoretic system, or equivalent;
- Gel maker specified by your gel electrophoresis apparatus;
- Transilluminator;
- Photographic system;
- Micropipettes in the volume ranges of, for example, 0.5-10, 10-100, 100-1,000 and 1,000-5,000 μ ls.

4.5 Reagents

4.5.1 PCR

- Applied Biosystems' AmpliTaq Gold[®] DNA Polymerase, or equivalent;
- 2 mmol/l deoxynucleoside triphosphate (dNTP) solution
(The AmpliTaq Gold[®] DNA Polymerase kit includes the solution.);
- 25 mmol/L magnesium chloride (MgCl₂) solution
(The AmpliTaq Gold[®] DNA Polymerase kit includes the solution.);
- 10 X PCR Buffer II
- Primer pairs
According to the target sequences you are detecting (see Tables 4.1(1)-(3)), you can either purchase your primer pairs or prepare them yourself. The sizes of target sequences to be amplified are shown in Tables 4.1(1)-(3). Primer pairs are available from Nippon Gene Co., Ltd. and FASMAC Co., Ltd. and they will meet most of your amplification needs.

Primer pairs to detect GM soybeans:

- Endogenous *le1* gene (Nippon Gene Cat # 313-05501, FASMAC Cat # S1-1M, or their bulk products); and
- Roundup Resistant soybeans-specific sequence (Nippon Gene Cat # 310-05511, FASMAC Cat # S2-1M, or their bulk products).

Primer pairs to screen for GM maize:

- Endogenous *SSIb* gene (Nippon Gene Cat # 315-05441, FASMAC Cat # M1-1M, or their bulk products);
- CaMV 35S promoter (Nippon Gene Cat # 317-05521, FASMAC Cat # C1-1M, or their bulk products); and
- GA21-specific sequence (Nippon Gene Cat # 312-05451, FASMAC Cat # M2-1M, or their bulk products).

Primer pairs to detect specific GM maize lines:

- Endogenous *SSIb* gene (Nippon Gene Cat # 315-05441, FASMAC Cat # M1-1M, or their bulk products);
- GA21-specific sequence (Nippon Gene Cat # 312-05451, FASMAC Cat # M2-1M, or their bulk products);
- Bt11-specific sequence (Nippon Gene Cat # 319-05461, FASMAC Cat # M3-1M, or their bulk products);
- Event 176-specific sequence (Nippon Gene Cat # 316-05471, FASMAC Cat # M4-1M, or their bulk products); and
- MON810-specific sequence (Nippon Gene Cat # 310-05491, FASMAC Cat # M6-1M, or their bulk products).

Primer pairs to detect GM potatoes:

- Endogenous *Pss* gene (Nippon Gene Cat # 311-05501, FASMAC Cat # G5-1, their bulk products, or self-prepared);
- NewLeaf-specific sequence (Note that no commercial product of primer pairs seems to be available at present that can be used to detect this GM line. For the time being, the GM line may be detected by combination of the primer pairs for the CaMV 35S promoter and the positive control plasmid for either Roundup Resistant soybeans or maize. However, this method should be used with utmost care because the target sequence of PCR is also present in plants infected with cytomegalovirus (CMV) and many GM crops.);
- NewLeaf Plus-specific sequence (Nippon Gene Cat # 311-05511, FASMAC Cat # G1-1, their bulk products, or self-prepared);

Reference plasmid solution^(Note 17) is available from Nippon Gene and FASMAC.

Reference plasmid for soybeans:

^(Note 17) The reference plasmid solution is used as a positive control for both PCR reactions and GM crop varieties. You can also use the DNAs of recombinants if they are available.

- Positive control plasmid for the detection of Roundup-resistant GM soybeans (Nippon Gene Cat # 311-04941, FASMAC Cat # PS-1, or their bulk products)

Reference plasmid for maize:

- Positive control plasmid for the detection of GM maize (Nippon Gene Cat # 314-04811, FASMAC Cat # PM-1, or their bulk products)

Reference plasmid for potatoes (see the Note in the second item under “*Primer pairs to detect GM potatoes*”):

- Positive control plasmid for the detection of NewLeaf Plus GM potatoes (Nippon Gene Cat # 311-05301, FASMAC Cat # PP-1, or their bulk products)

4.5.2 Gel electrophoresis

- Agarose gel
- TBE or TAE buffer
- Ethidium bromide
- Gel loading buffer (blue juice)
- DNA molecular weight marker:
Choose a marker with good resolution for your PCR products. A 100-bp ladder marker should be the best choice.

4.6 Procedure

4.6.1 PCR procedure

Precautions

- Use sterilized laboratory utensils and equipment.
- Take all necessary measures to prevent contamination, such as wearing disinfected gloves.

(1) Number of tubes needed for PCR

Perform PCR for each solution obtained by DNA extraction of a sample. Include the primer pairs that can detect the endogenous gene specific to the crop to be tested to confirm that the PCR successfully yields the product of expected length. Use of such primer pairs should be accompanied by both negative and positive controls. A negative control detects contamination of reagents by sequences that can serve as PCR templates. A positive control with a reference plasmid for the crop, signals that the primer pairs are of good quality, that the PCR solution is thoroughly mixed, and that the PCR reaction proceeds successfully. Set up a “no template” negative control when the primer pairs for the detection of GM varieties are included. Also include in every PCR at least one negative control that has a reference plasmid template but does not have primer pairs that can bind to any GM variety. This control detects contamination of reagents by such primer pairs. If you have the non-GM variety of the crop at issue, you can set up a

negative control with a template prepared from the DNA extracted from the non-GM variety. This control detects contaminants such the primer pairs that bind to the endogenous gene present in those that are designed to bind to GM varieties.

(2) Preparation of master mix

The composition of the master mix is shown in Table 4.2 *Composition of PCR Solution*. Note that primer pairs are added after all other components have been mixed. The number of sterilized 0.2 ml tubes is determined by that of the samples for which PCR is performed. The total volume of PCR solution is then determined on the basis of the number of tubes. Add an extra volume to allow for losses due to, for example, droplets adhering to the inside wall of the tube. The volumes of individual components other than template DNA solution are calculated based on their relative volumes shown in Table 4.2. Mixed solution is then transferred in 22.5 µl aliquots to each reaction tube. Choose your micropipette considering its volume range, especially the minimum volume.

For the “no primer” negative control, prepare a separate master mix in which no primer pairs are added.

(3) Addition of samples

After 2.5 µl aliquots of template DNA solution have been added to the tubes containing the master mix, add the samples to the tube in the order of the extracted DNA, negative control and positive control solutions.

(4) PCR amplification

After all solutions are added, perform PCR on the thermal conditions specified in Table 4.3 *Thermal Cycles*.

(5) Post-PCR procedure

When PCR is complete, store the PCR product in a refrigerator or a freezer, or otherwise immediately run electrophoresis.

4.6.2 Gel electrophoresis

Run gel electrophoresis according to 4.6.2.1 *Pre-staining gel electrophoresis* or 4.6.2.2 *Post-staining gel electrophoresis*. These and other gel electrophoresis procedures described in this manual are based on the Mupid[®]-2 electrophoresis apparatus (Advance Co., Ltd., Japan).

Precautions

- There is no need at this stage for the use of sterilized laboratory utensils and equipment.
- Wear gloves to protect you from contact with hazardous materials.

4.6.2.1 Pre-staining gel electrophoresis

(1) Preparing the gel ready

Assemble the gel maker module and prepare a 2-3% agarose gel. First, weigh out

an appropriate amount of agarose and add it to the TBE buffer^{(Note 18)(Note 19)}. Heat the gel solution to melt agarose^(Note 20). When agarose is completely melted, add the ethidium bromide solution to the gel solution so that each 100 ml gel solution contains 50 µg ethidium bromide. Pour the thoroughly mixed solution into the gel maker module and place the comb carefully so that no air bubble is formed around it. Leave the gel to cool at ambient temperature for about 30 minutes. When the gel has cooled and set, remove the comb carefully from the gel so that you don't make a hole through the well from which the sample solution can leak.

(2) Preparing the gel tank ready

Place the gel in the gel tank. Check to ensure that the wells are on the side of the negative electrode. While it is advisable that you immediately run the gel, you can store it in the buffer for a few days.

Flood the gel with the TBE buffer^(Note 18) to barely cover the gel surface. Note that excess buffer results in longer migration time due to reduced current density.

(3) Running electrophoresis

Add 1 µl gel loading buffer to 5 µl of the post-PCR sample DNA solution^(Note 21) and gently pour the mixture into the wells. Also load the molecular weight marker in one or two lanes on the same gel. Remember that taking too much time in pouring a sample solution into wells can diffuse the DNAs into the gel and produce smeared bands.

When the sample solution has properly been poured into the wells, turn on the power and run the gel at 100 volts. Check if bubbles are generated at the electrodes by the electrolysis of water to ensure that electrophoresis is occurring. Stop running the gel when the BPB marker dye added to the gel loading buffer has migrated half of the way down the gel. Although "half of the way down the gel" is only an indication, it is not advisable to run the gel over a longer distance because the size of the PCR product is relatively small.

Photograph the gel immediately after electrophoresis to prevent diffusion of DNAs resulting in lower resolution of bands.

(4) Photographing the gel

Cover the surface of a UV transilluminator with plastic wrap^(Note 22). Place the gel

^(Note 18) You may use TAE for TBE if it works. In that case, also use TAE in the electrophoresis buffer.

^(Note 19) A 100 ml gel solution is sufficient for two large and one small gels.

^(Note 20) The gel solution can conveniently be melted in the microwave. Handle the flask(?) carefully to avoid burn by sudden boiling.

^(Note 21) You need not to strictly adhere to the 1:5 ratio of the gel loading buffer to the sample DNA solution. Thus, you can make the mixture simply by: (1) aspirating five volumes of the sample solution into the micropipette tip; (2) placing one volume of the buffer on a laboratory film; and (3) mixing the two solutions by drawing and expelling them several times with the micropipette. If you wish to increase the volume of the mixture because of, for example, the size of the wells, you can add one volume of the gel loading buffer per six volumes of the mixture.

^(Note 22) Some food wrap films absorb the UV light depending on the wavelength and the band image may not appear. In that case, use a polyvinylidene chloride film.

on the wrap and turn on the transilluminator^(Note 23). Photograph the gel with a CCD camera and examine the pattern of bands against that of the DNA molecular weight marker to determine if the bands of interest are present.

4.6.2.2 Post-staining gel electrophoresis

(1) Preparing the gel ready

Assemble the gel maker module and prepare a 2-3% agarose gel. First, weigh out an appropriate amount of agarose and add it to the TBE buffer^{(Note 18)(Note 19)}. Heat the gel solution to melt agarose^(Note 20). Pour the thoroughly mixed solution into the gel maker module and place the comb carefully so that no air bubble is formed around it. Leave the gel to cool at ambient temperature for about 30 minutes. When the gel has cooled and set, remove the comb carefully from the gel so that you don't make a hole through the well from which the sample solution can leak.

(2) Preparing the gel tank ready

Place the gel in the gel tank. Check to ensure that the wells are on the side of the negative electrode. While it is advisable that you immediately run the gel, you can store it in the buffer for a few days.

Flood the gel with the TBE buffer^(Note 18) to barely cover the gel surface. Note that excess buffer results in longer migration time due to reduced current density.

(3) Running electrophoresis

Add 1 µl gel loading buffer to 5 µl of the post-PCR sample DNA solution^(Note 21) and gently pour the mixture into the wells. Also load the molecular weight marker in one or two lanes on the same gel. Remember that taking too much time in pouring a sample solution into wells can diffuse the DNAs into the gel and produce smeared bands.

When the sample solution has properly been poured into the wells, turn on the power and run the gel at 100 volts. Check if bubbles are generated at the electrodes by the electrolysis of water to ensure that electrophoresis is occurring. Stop running the gel when the BPB marker dye added to the gel loading buffer has migrated half of the way down the gel. Although "half of the way down the gel" is only an indication, it is not advisable to run the gel over a longer distance because the size of the PCR product is relatively small.

^(Note 23) UV transilluminators are available in a variety of wavelengths. Considering health risk and sensitivity, a UV light of about 312 nm is best recommended.

^(Note 18) You may use TAE for TBE if it works. In that case, also use TAE in the electrophoresis buffer.

^(Note 19) A 100 ml gel solution is sufficient for two large and one small gels.

^(Note 20) The gel solution can conveniently be melted in the microwave. Handle the flask(?) carefully to avoid burn by sudden boiling.

^(Note 21) You need not to strictly adhere to the 1:5 ratio of the gel loading buffer to the sample DNA solution. Thus, you can make the mixture simply by: (1) aspirating five volumes of the sample solution into the micropipette tip; (2) placing one volume of the buffer on a laboratory film; and (3) mixing the two solutions by drawing and expelling them several times with the micropipette. If you wish to increase the volume of the mixture because of, for example, the size of the wells, you can add one volume of the gel loading buffer per six volumes of the mixture.

Stain the gel immediately after electrophoresis to prevent diffusion of DNAs resulting in lower resolution of bands.

(4) Staining the gel

Transfer the electrophoresed gel to a plastic container filled with an amount of new electrophoresis buffer sufficient to flood the gel. Add 50 µg of ethidium bromide per 100 ml of buffer. Place the container on a shaker and gently shake the solution for about 30 minutes to stain the gel.

(5) Photographing the gel

Cover the surface of a UV transilluminator with plastic wrap^(Note 22). Place the gel on the wrap and turn on the transilluminator^(Note 23). Photograph the gel with a CCD camera and examine the pattern of bands against that of the DNA molecular weight marker to determine if the bands of interest are present.

4.7 Verification of successful PCR

Successful PCR can be verified on the photograph of the electrophoresed gel by the absence of any band in the lane for the negative control and the presence of the bands of predicted lengths in the lane for the positive control. Once successful PCR is verified, the presence or absence of GMOs in the sample can be determined according to the procedures described in 5. *Determination of the presence or absence of GMOs*. If contamination is found, immediately take necessary actions specified in Part 5 *Prevention of Contamination* and perform PCR again using, if necessary, newly extracted DNA.

4.8 Specificity (and sensitivity) of the procedure

Combined with the DNA templates extracted from soybean and maize seeds and potatoes using the DNeasy Plant Maxi kit, the PCR procedure described in this manual produces the bands only when the test sample contains the target GMOs. This PCR procedure is also plant-specific in that it produces no band when other main crops including rice, wheat and barley are used as test samples.

4.9 Record keeping

Keep the photographs of electrophoresed gels in the image data file.

4.10 Additional notes

Not applicable.

^(Note 22) Some food wrap films absorb the UV light depending on the wavelength and the band image may not appear. In that case, use a polyvinylidene chloride film.

^(Note 23) UV transilluminators are available in a variety of wavelengths. Considering health risk and sensitivity, a UV light of about 312 nm is best recommended.

5. Determination of the presence or absence of GMOs

Presence or absence of GMOs in a test sample is determined by the presence or absence of the target PCR products through electrophoresis. Presence of a GMO is determined when the PCR products contained the lengths of both the plant's endogenous gene and the GMO-specific gene.

If the PCR products did not contain the ones corresponding to the endogenous gene, first perform PCR again using already extracted DNA. If the resulting PCR products did not show the presence of the endogenous gene, then extract DNA again from the sample crop and carry out PCR. If the PCR still did not amplify the sequence of the endogenous gene, identify the sample crop as "undetectable."

Table 4.1(1)
Primer pairs for soybeans*

Target gene	ID	Size	Amplified region	Remarks
Endogenous <i>Le1</i> gene	<i>Le1-n02</i>	118 bp	<i>Le1</i>	
Introduced <i>P-35S</i> gene	<i>P35S-1</i>	101 bp	<i>P-35S</i>	The results should be evaluated carefully because this primer pair also binds to the genes of non-GM plants infected with CMV.
Introduced <i>NOS-ter</i> gene	<i>NOS-ter-2</i>	151 bp	<i>NOS-ter</i>	The results should be evaluated carefully because this primer pair also primes when the soil bacteria <i>Agrobacterium tumefaciens</i> is present.
Introduced RRS-specific gene	<i>RRS-01</i>	121 bp	<i>cp4 epsps</i> from <i>P. hybrida</i>	

* The Japanese patents jointly applied for by the National Food Research Institute, Asahi Breweries, Ltd. and Nippon Flour Mills Co., Ltd. on these primer pairs are pending.

Table 4.1(2)
Primer pairs for maize**

Target gene	ID	Size	Amplified region	Remarks
Endogenous <i>SSIIb</i> gene	SSIIb	151 bp	<i>zSSIIb</i>	
Introduced <i>P-35S</i> gene	P35S-1	101 bp	<i>P-35S</i>	The results should be evaluated carefully because this primer pair also binds to the genes of non-GM plants infected with CMV.
Introduced <i>NOS-ter</i> gene	NOS-ter-2	151 bp	<i>NOS-ter</i>	The results should be evaluated carefully because this primer pair also primes when the soil bacteria <i>Agrobacterium tumefaciens</i> is present.
Introduced Event 176-specific gene	E176-2	100 bp	<i>cryIA(b)</i> - PEPC intron #9	
Introduced Bt11 gene	Bt11-3	127 bp	<i>adh1-1S</i> - <i>cryIA(b)</i>	
Introduced GA21 gene	GA21-3	133 bp	OTP - <i>m-epsps</i>	
Introduced T25 gene	T25-1	149 bp	<i>Pat</i> - 35S- <i>ter</i>	
Introduced MON810 gene	M810-2	113 bp	<i>hsp70</i> - <i>cryIA(b)</i>	

** The Japanese patents jointly applied for by the National Food Research Institute, Asahi Breweries, Ltd. and Nippon Flour Mills Co., Ltd. on these primer pairs are pending.

Table 4.1(3)
Primer pairs for maize

Target gene	ID	Size	Amplified region	Remarks
Endogenous Pss gene	Pss 01n-5'	216 bp	TGACCT GGA CAC CAC AGT TAT	<i>S. tuberosum</i> sucrose synthase/ sense
	Pss 01n-3'		GTG GAT TTC AGG AGT TCT TCGA	<i>S. tuberosum</i> sucrose synthase/ anti-sense
Introduced NewLeaf-specific gene				
Introduced NewLeaf Plus-specific gene	p-FMV02-5'	234 bp	AAATAACGTGGAAAAGAGCTGTCTCTGA	p-FMV/ sense
	PLRV01-3'		AAAAGAGCGGCATATGCCGTTAAATCTG	PLRV/ anti-sense

Table 4.2
Composition of PCR solution.

Components	Volume (µl) / tube	Final concentration
Sterilized water	15.375	
AmpliTaq™ Gold	0.125	0.625 U
10x PCR buffer II	2.5	1 x
dNTP (2 mmol/l each)	2.5	200 µmol/l each
MgCl ₂ (25 mmol/l)	1.5	1.5 mmol/l
Primer pair (25 µmol/l each)	0.5	0.5 µmol/l each
Template DNA (10 ng/µl)	2.5	25 ng
Total volume	25.0	

**Table 4.3
Thermal Cycles and Conditions**

Stage	Temperature (°C)	Duration	Number of cycle repetitions
Denaturation in the first cycle	95	10 min	1 cycle
Denaturation in the second cycle and after	95	30 sec	39 cycles (2nd-40th)
Annealing	60	30 sec	40 cycles
Extension in the first 39 cycles	72	30 sec	39 cycles (1st-39th)
Extension in the final cycle	72	7 min	1 cycle
Storage	4	NA*	NA*

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Designated Processed Foods in GMO Labeling

1. Tofu
2. Frozen-tofu, yuba and other tofu
3. Natto
4. Soy milk
5. Miso paste
6. Boiled soybeans
7. Canned or bottled soybeans
8. Soy power
9. Roasted soybeans
10. Processed food made from 1~9
11. Processed food made from soybeans
12. Processed food made from soy powder
13. Processed food made from soy protein
14. Processed food made from green soybean
15. Processed food made from soybean sprout
16. Corn snack
17. Corn starch
18. Pop-corn
19. Frozen-corn
20. Canned or bottled corn
21. Processed food made from corn flour
22. Processed food made from corn grits
23. Processed food made from corn
24. Processed food made from 16~20
25. Frozen potato
26. Dried potato flakes
27. Potato starch
28. Potato snack
29. Processed food made from 25~28
30. Processed food made from potato

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sample ID	SSIb	5分之1	NK603	5分之1	NK603/SSIb (CV)
1 (20uM)	22249.04	4449.81	11741.10	2348.22	0.53
2 (20uM)	20927.51	4185.50	11212.59	2242.52	0.54
3 (20uM)	25059.59	5011.92	13920.34	2784.07	0.56
4 (20uM)	21148.70	4229.74	11257.40	2251.48	0.53

sample ID	SSIb	5倍	NK603	5倍	NK603/SSIb (CV)
1 (4uM)	3848.63	19243.15	1908.74	9543.70	0.50
2 (4uM)	3606.57	18032.85	1830.42	9152.10	0.51
3 (4uM)	4824.21	24121.05	2388.30	11941.50	0.50
4 (4uM)	3230.48	16152.40	1805.98	9329.90	0.58

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**Study on the Qualitative PCR
to detect the presence of GMO
in Animal-Derived Products**

2003. 8.

(2003. 7. 10 ~ 8. 4)

National Veterinary Research and Quarantine Service
Livestock Products Safety & Inspection Department
Livestock Products Standard Division
Chung, Myeong-Eun

I . Materials and Methods

1. Sample used in this study

No.	Name	Contents	Information
1	Chicken Nugget(C)	with Corn flour	Purchased from Market
2	Baby Ham(H)	with Soy protein	Purchased from Market
3	Sausage(S)	with Corn starch & soy protein	Purchased from Market
4	White Sausage(Fish Sausage W)	with Corn starch	Purchased from Market
5	Sausage(A)	Contents of Corn starch : 7.5%	Prima Ham Company
6	Sliced Ham(CuH)	Contents of Corn starch : 0.9%	Prima Ham Company
7	FrankFruite Sausage(P.S)	Contents of Corn starch : 3.1%	Prima Ham Company
8	Ham(R.H)	Contents of Corn starch : 6.7%	Prima Ham Company
9	Ham(CT)	Contents of Corn starch : 0%	Prima Ham Company (Heating treated at 80°C for 30min)
10	Ham(10%)	Contents of Corn starch : 8.9%	Prima Ham Company (Heating treated at 80°C for 30min)
11	Ham(15%)	Contents of Corn starch : 12.8%	Prima Ham Company (Heating treated at 80°C for 30min)
12	Ham(20%)	Contents of Corn starch : 16.4%	Prima Ham Company (Heating treated at 80°C for 30min)
13	Non-GM Corn Powder	100% Non-GM Corn Powder	(-) control
14	NK Powder	99.9% Non-GM Corn Powder + 0.1% NK powder	(+) control

2. DNA Preparation and Purification

2.1 Methods of DNA Extraction

2.1.1 DNA Extraction using QIAGEN Genomic-tip 20/G

2.1.2 DNA Extraction using QIAGEN DNeasy Maxi Kit plus CTAB (API-CTAB Buffer 0.2%)

2.1.3 DNA Extraction using QIAGEN DNeasy Maxi Kit

2.1.4 DNA Extraction using CTAB

3. Reaction of PCR

3.1 Initial denaturation : 95°C for 10min

3.2 40 cycles of Amplification

- 95°C for 30sec, 60°C for 30sec, 72°C for 30sec

3.3 Completed at 72°C for 7min

4. PCR Analysis

II. Results

1. Concentration and Purify of Extracted DNA

Name of Sample	QIAGEN DNeasy Mini Kit			API-CTAB Buffer 0.2%			QIAGEN Genomic-tip 20/G			CTAB		
	(DNA) ng/ul	OD ₂₆₀ /OD ₂₈₀	OD ₂₆₀ /OD ₂₈₀	(DNA) ng/ul	OD ₂₆₀ /OD ₂₈₀	OD ₂₆₀ /OD ₂₈₀	(DNA) ng/ul	OD ₂₆₀ /OD ₂₈₀	OD ₂₆₀ /OD ₂₈₀	(DNA) ng/ul	OD ₂₆₀ /OD ₂₈₀	OD ₂₆₀ /OD ₂₈₀
Chicken Nugget(C)	135.49	1.71	0	833.74	1.75	0.96	1680.53	1.82	2.33	49.22	2.04	0.26
Baby Ham(H)	658.42	1.71	1.12	834.23	1.78	1.10	2521.15	1.83	2.39	56.64	2.12	0.28
Sausage(S)	394.60	1.74	1.78	1023.61	1.75	1.14	2856.03	1.82	2.39	53.67	2.13	0.27
White Sausage(W)	0	2.69	0.08	174.14	1.77	0.75	680.41	1.83	2.24	0	0.49	0
Sausage(A)	58.84	1.90	3.60	125.51	1.84	1.27	2538.94	1.84	2.34	59.83	2.36	0.34
Sliced Ham(CuH)	10.7.05	1.73	2.11	1014.07	1.79	1.10	1091.47	1.84	2.34	26.05	3.50	0.15
Frankfruite Sausage(PS)	384.65	1.81	2.23	988.10	1.82	1.53	1319.05	1.83	2.82	69.02	1.96	0.35
Ham(RH)	137.52	1.83	2.42	753.17	1.78	0.97	1257.08	1.82	2.23	22.36	2.21	0.12
Ham(CT)	292.92	1.81	2.49	36.49	1.88	0.45	1143.13	1.82	1.99	0	1.68	1.71
Ham(10%)	194.23	1.80	2.00	189.20	1.85	1.64	1275.29	1.82	2.57	0	1.89	0.38
Ham(15%)	104.03	1.81	2.70	98.60	1.93	1.32	1117.62	1.81	1.89	0	1.71	3.93
Ham(20%)	79.14	1.78	2.32	243.22	1.84	1.69	1122.72	1.82	2.38	0	1.80	1.14
Non-GM Corn	538.65	1.78	7.19	1201.39	1.78	1.87	1174.57	1.79	2.51	24	12.16	0
NK	579.42	1.78	4.95	1383.95	1.80	2.07	1367.44	1.80	2.59	35.76	2.16	1.41

2. Qualitative PCR

2.1 Results of Qualitative PCR with extracted DNA using QIAGEN Genomic-tip 20/G

2.1.1 PCR products with SSIIB primer


Fig. 1. Agarose gel electrophoresis of PCR products amplified from animal-derived products using SSIIB primer
 Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (W), lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(CuH), lane 13 & 14:Frankfruite Sausage(PS), lane 15 & 16:Ham(RH), lane 17 & 18:Ham(CT), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.1.2 PCR products with P35S primer


Fig. 2. Agarose gel electrophoresis of PCR products amplified from animal-derived products using P35S primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (H) lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfruite Sausage(P,S), lane 15 & 16:Ham(R,I,D), lane 17 & 18:Ham(C,T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.1.3 PCR products with GA2I primer


Fig. 3. Agarose gel electrophoresis of PCR products amplified from animal-derived products using GA2I primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (H) lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfruite Sausage(P,S), lane 15 & 16:Ham(R,I,D), lane 17 & 18:Ham(C,T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.2 Results of Qualitative PCR with extracted DNA using QACENDNAasy Mini Kit plus CTAB (API+CTAB Buffer 0.2%)

2.2.1 PCR products with SSIIb primer


Fig. 4. Agarose gel electrophoresis of PCR products amplified from animal-derived products using SSIIb primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (H) lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfruite Sausage(P,S), lane 15 & 16:Ham(R,I,D), lane 17 & 18:Ham(C,T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.2.2 PCR products with P35S-P primer

Fig. 5. Agarose gel electrophoresis of PCR products amplified from animal-derived products using P35S primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (H) lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfruite Sausage(P,S), lane 15 & 16:Ham(R,I,D), lane 17 & 18:Ham(C,T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.3.3 PCR products with GA21 primer


Fig. 6. Agarose gel electrophoresis of PCR products amplified from animal-derived products using GA21 primer
 Lane M : 100bp ladder, lane 1 & 2Chicken Nugget(C), lane 3 & 4Baby Ham(H), lane 5 & 6Sausage(S), lane 7 & 8White Sausage (W), lane 9 & 10Sausage(A), lane 11 & 12Sliced Ham(Cul D), lane 13 & 14Frankfrite Sausage(PS), lane 15 & 16Ham(R.H), lane 17 & 18Ham(CT), lane 19 & 20Ham(10%), lane 21 & 22Ham(15%), lane 23 & 24Ham(20%), lane 25 & 26Non-GM Corn, lane 27 & 28NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.3 Results of Qualitative PCR with extracted DNA using QIAGEN DNeasy Maxi Kit

2.3.1 PCR products with SSIb primer


Fig. 7. Agarose gel electrophoresis of PCR products amplified from animal-derived products using SSIb primer
 Lane M : 100bp ladder, lane 1 & 2Chicken Nugget(C), lane 3 & 4Baby Ham(H), lane 5 & 6Sausage(S), lane 7 & 8White Sausage (W), lane 9 & 10Sausage(A), lane 11 & 12Sliced Ham(Cul D), lane 13 & 14Frankfrite Sausage(PS), lane 15 & 16Ham(R.H), lane 17 & 18Ham(CT), lane 19 & 20Ham(10%), lane 21 & 22Ham(15%), lane 23 & 24Ham(20%), lane 25 & 26Non-GM Corn, lane 27 & 28NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.3.2 PCR products with P35S primer


Fig. 8. Agarose gel electrophoresis of PCR products amplified from animal-derived products using P35S primer
 Lane M : 100bp ladder, lane 1 & 2Chicken Nugget(C), lane 3 & 4Baby Ham(H), lane 5 & 6Sausage(S), lane 7 & 8White Sausage (W), lane 9 & 10Sausage(A), lane 11 & 12Sliced Ham(Cul D), lane 13 & 14Frankfrite Sausage(PS), lane 15 & 16Ham(R.H), lane 17 & 18Ham(CT), lane 19 & 20Ham(10%), lane 21 & 22Ham(15%), lane 23 & 24Ham(20%), lane 25 & 26Non-GM Corn, lane 27 & 28NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.3.3 PCR products with GA21 primer

Fig. 9. Agarose gel electrophoresis of PCR products amplified from animal-derived products using GA21 primer
 Lane M : 100bp ladder, lane 1 & 2Chicken Nugget(C), lane 3 & 4Baby Ham(H), lane 5 & 6Sausage(S), lane 7 & 8White Sausage (W), lane 9 & 10Sausage(A), lane 11 & 12Sliced Ham(Cul D), lane 13 & 14Frankfrite Sausage(PS), lane 15 & 16Ham(R.H), lane 17 & 18Ham(CT), lane 19 & 20Ham(10%), lane 21 & 22Ham(15%), lane 23 & 24Ham(20%), lane 25 & 26Non-GM Corn, lane 27 & 28NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.4 Results of Qualitative PCR with extracted DNA using CTAB

2.4.1 PCR products with SSIIb primer


Fig. 10. Agarose gel electrophoresis of PCR products amplified from animal-derived products using SSIIb primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (S), lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfurter Sausage(PS), lane 15 & 16:Ham(R.H), lane 17 & 18:Ham(C.T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.4.2 PCR products with P35S primer


Fig. 11. Agarose gel electrophoresis of PCR products amplified from animal-derived products using P35S primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (S), lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfurter Sausage(PS), lane 15 & 16:Ham(R.H), lane 17 & 18:Ham(C.T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.4.3 PCR products with GA2I primer


Fig. 12. Agarose gel electrophoresis of PCR products amplified from animal-derived products using GA2I primer

Lane M : 100bp ladder, lane 1 & 2:Chicken Nugget(C), lane 3 & 4:Baby Ham(H), lane 5 & 6:Sausage(S), lane 7 & 8:White Sausage (S), lane 9 & 10:Sausage(A), lane 11 & 12:Sliced Ham(Cul D), lane 13 & 14:Frankfurter Sausage(PS), lane 15 & 16:Ham(R.H), lane 17 & 18:Ham(C.T), lane 19 & 20:Ham(10%), lane 21 & 22:Ham(15%), lane 23 & 24:Ham(20%), lane 25 & 26:Non-GM Corn, lane 27 & 28:NK, lane 29:No DNA, lane 30:No Primer, lane 31: positive Control

2.5 Comparisons for Results of Qualitative PCR with extracted DNA from different DNA extraction methods

Sample	QIAGEN Genomic-tip 20/G			Maxi Kit plus CTAB (AFI-CTAB Buffer 0.2%)			QIAGEN DNeasy Maxi Kit			CTAB		
	SSIIIb	P35S	GA21	SSIIIb	P35S	GA21	SSIIIb	P35S	GA21	SSIIIb	P35S	GA21
Chicken Nugget(C)-1	+	+	-	+	+	-	+	+	-	+	+	-
Chicken Nugget(C)-2	+	+	-	+	Δ	-	+	+	-	+	+	-
Baby Ham(H)-1	-	-	-	Δ	Δ	-	-	-	-	-	+	-
Baby Ham(H)-2	+	-	-	Δ	-	-	-	-	-	-	-	-
Sausage(S)-1	+	-	-	-	-	-	+	-	-	+	+	-
Sausage(S)-2	+	-	-	-	-	-	+	-	-	+	-	-
White Sausage (Fish Sausage(W)-1)	-	-	-	-	-	-	-	-	-	-	-	-
White Sausage (Fish Sausage(W)-2)	-	-	-	-	-	-	-	-	-	-	-	-
Sausage(A)-1	-	-	-	-	-	-	+	+	-	Δ	-	-
Sausage(A)-2	-	-	-	-	-	-	+	-	-	+	-	-
Sliced Ham (CuH)-1	+	+	-	-	-	-	-	-	-	+	+	-
Sliced Ham (CuH)-2	+	+	-	-	-	-	+	+	-	-	+	-
FrankFruite Sausage(PS)-1	Δ	+	-	-	Δ	-	+	+	-	+	+	-
FrankFruite Sausage(PS)-2	Δ	+	-	-	Δ	-	+	+	-	+	+	-
Ham(R.H)-1	+	+	-	-	Δ	-	-	+	-	+	+	-
Ham(R.H)-2	+	+	-	-	Δ	-	+	+	-	+	+	-
Ham(CT)-1	-	-	-	-	-	-	-	-	-	-	-	-
Ham(CT)-2	-	-	-	-	-	-	-	-	-	+	-	-
Ham(10%)-1	+	-	-	-	-	-	+	-	-	+	+	-
Ham(10%)-2	-	-	-	-	-	-	+	+	-	+	+	-
Ham(15%)-1	+	-	-	-	-	-	+	-	-	-	-	-
Ham(15%)-2	-	-	-	-	-	-	+	-	-	Δ	Δ	-
Ham(20%)-1	Δ	+	-	-	-	-	+	-	-	+	-	-
Ham(20%)-2	Δ	-	-	-	-	-	+	+	-	Δ	-	-
Non-GM Corn-1	+	-	-	+	-	-	+	-	-	Δ	-	-
Non-GM Corn-2	+	-	-	+	-	-	+	-	-	+	-	-
NK-1	+	+	-	+	+	-	+	+	-	+	+	-
NK-2	+	+	-	+	+	-	+	+	-	+	+	-

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

No. 5

Sample DNA extraction method :

	CTAB	MAXI	G-Tip
No. 1	OK	OK	OK
No. 2	OK	OK	OK
No. 3	OK	OK	OK
No. 4	OK	OK	OK
No. 5	OK	Loss (I am sorry!)	OK

Sample DNA solution conctratation :



	CTAB	MAXI	G-Tip
No. 1	4.79158	-0.52152	1.00901
No. 2	8.67779	-0.20617	190.12640
No. 3	9.27336	0.38749	3.06020
No. 4	2.05979	0.59363	8.78242
No. 5	31.5500	Loss (I am sorry!)	733.46710

Result :

Extraction method	Sample No.	Sample name	LE 1	RRS
CTAB	1	Sample N0.1-1	—	—
CTAB	2	Sample N0.1-2	—	—
CTAB	3	Sample N0.2-1	—	—
CTAB	4	Sample N0.2-2	+	—
CTAB	5	Sample N0.3-1	—	—
CTAB	6	Sample N0.3-2	—	—
CTAB	7	Sample N0.4-1	+	—
CTAB	8	Sample N0.4-2	+	—
CTAB	9	Sample N0.5-1	+	—
CTAB	10	Sample N0.5-2	+	—
MAXI	11	Sample N0.1-1	—	—
MAXI	12	Sample N0.1-2	—	—
MAXI	13	Sample N0.2-1	—	—
MAXI	14	Sample N0.2-2	—	—
MAXI	15	Sample N0.3-1	—	—
MAXI	16	Sample N0.3-2	—	—
MAXI	17	Sample N0.4-1	—	—
MAXI	18	Sample N0.4-2	—	—
G-tip	19	Sample N0.1-1	+	—
G-tip	20	Sample N0.1-2	—	—
G-tip	21	Sample N0.2-1	++	—
G-tip	22	Sample N0.2-2	+	—
G-tip	23	Sample N0.3-1	+	—
G-tip	24	Sample N0.3-2	—	—
G-tip	25	Sample N0.4-1	++	—
G-tip	26	Sample N0.4-2	—	—
G-tip	27	Sample N0.5-1	+	—
G-tip	28	Sample N0.5-2	+	—
Negative (no-DNA)	29		—	—
Negative (no-primer)	30		—	—
Positive (plasmid)	31		++	+

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CTAB method for extracting DNA from freeze-dried potato powder

Ground sample weight 200 mg sample to put in a 15 mL polypropylene centrifuge tube

- | ← add 5.0 mL CTAB extraction solution*¹ and 30 μ L Proteinase K
- | invert the tube several times to homogenize completely
- | incubate at room temp (RT, 25°C) for 60 min with shaking.
- | after 8,000g centrifugation at RT for 10 min, transfer 4 mL supernatant to a new 15 mL polypropylene centrifuge tube

Phenol and CIA treatment

- | ← add 2.0 mL TE saturated phenol*² into the new tube, shake vigorously for 2 min
- | ← add 2.0 mL CIA, shake vigorously for 2 min
- | after 8,000g centrifugation at RT for 10 min, transfer up-layer solution (~ 4 mL) to a new 15 mL polypropylene centrifuge tube

CIA treatment

- | ← add 4.0 mL CIA*³ into the new tube, shake vigorously for 2 min
- | after 8,000g centrifugation at RT for 10 min, transfer up-layer solution (~ 4 mL) to a new 15 mL polypropylene centrifuge tube

Isopropanol precipitation

- | ← add the same volume isopropanol (~ 4 mL) into the new tube, then invert this tube 10 times to mix well
- | after 8,000g centrifugation at 4°C for 10 min, discard supernatant
- | ← add 4.0 mL 70% alcohol, invert to mix
- | after 8,000g centrifugation at 4°C for 10 min, discard supernatant
- | vacuum dry for 5 min

dissolve DNA and remove RNA

- | add 200 μ L TE to dissolve DNA, then transfer DNA solution to a 2 mL tube
- | ← add 5 μ L RNase (10mg/mL), then incubate at 37°C for 1 hour
- | ← add 800 μ L CTAB extraction solution, mix well

CIA treatment

- | ← add 1 mL CIA, mix gently
- | after 12,000 rpm centrifugation at RT for 10 min, transfer up layer solution to a new tube

Isopropanol precipitation

- | ← add 1 mL isopropanol, invert 10 times to mix well
- | 15,000 rpm centrifugation at RT for 10 min, then discard the supernatant
- | ← add 800 μ L 70% alcohol, invert to mix
- | after 15,000 rpm centrifugation at 4°C for 5 min, discard supernatant

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All samples should be a particle size that passes through a 100 U.S. Mesh Screen

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Extract protein

Weight out 0.5g ground sample to a 15 mL tube

Add 4.5 mL of Soya Extraction Buffer, then vortex 10 sec.

Centrifuge at 2,500 ×g for 15 min

Transfer supernatant to a new 15 mL tube

Dilute the protein extract

Take 20 μ L protein extract, add into 280 μ L Soya assay Buffer, vortex (15× Dilution)

Take 20 μ L protein extract, add into 380 μ L Soya assay Buffer, vortex (300× Dilution)

Take 20 μ L mixed sample 300× Diluent, add into 180 μ L Soya assay Buffer vortex, to get a final 3,000× Diluent

ELISA Test

Add 100 μ L diluted extract to Assay well plate, 2 repeat

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










YOKOHAMA QUARANTINE STATION

History of Examination Center

1991--- Establishment of Examination center

1995--- Construction of Examination center

2001--- Increase of Genetic Modified Organism
Laboratory (March)

2003--- Construction of new Genetic Modified
Organism Laboratory (March)

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
Organization of Yokohama Examination

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

Fig. 1. Schematic diagram of PCR primers designed for detection of NewLeaf Plus[®] potato (NL-P). The two pairs of PCR primers were designed for detection of NewLeaf Plus[®] potato (NL-P). The expected lengths of PCR products are indicated in parentheses.

Table 1. List of PCR Primers

No.	Name	Sequence (5'-3')	Specificity	Amplicon
1	Pss 01n-5' Pss 01n-3'	TGA CCT GGA CAC CAC AGT TAT GTG GAT TTC AGG AGT TCT TCG A	potato sucrose synthetase/sense potato sucrose synthetase/anti-sense	216 bp
2	PLRV-rep1-5' PLRV-rep1-3'	CTT CTT TCA CGG AGT TCC AG TCG TCA TTA AAC TTG ACG AC	PVY/sense PVY/anti-sense	172 bp
3	p-FMV02-5' PLRV01-3'	AAA TAA CGT GGA AAA GAG CTG TCC TGA AAA AGA GCG GCA TAT GCG GTA AAT CTG	p-FMV/sense PVY/anti-sense	234 bp

*3 JRL/ILSI Joint Workshop on Method Development in Relation to Regulatory Requirements for the Detection of GMOs in the Food Chain, Book of Abstracts, p. 49 (11, December, 2000).


Fig. 2. Specificity of potato specific primer (primer pair No. 1)

Arrowheads indicate the expected PCR amplification products.

Lane 1: negative control (no DNA); lane 2: potato DNA; lane 3: wheat DNA; lane 4: barley DNA; lane 5: rice DNA; lane 6: soy DNA; lane 7: maize DNA; M: 100 bp ladder size standard


Fig. 3. PCR products amplified from non-GM and GM potato DNA with positive control primers (primer pair No. 1)

Arrowheads indicate the expected PCR amplification products.

Lane 1: negative control (no DNA); lane 2: NewLeaf[®] potato SPBT02-05 DNA; lane 3: NL-P DNA; lane 4: NewLeaf[®] potato RBBT-06 DNA; lane 5: NewLeaf Y[®] potato DNA; lane 6: non-GM potato (Russet Burbank) DNA; lane 7: non-GM potato (Superior) DNA; M: 100 bp ladder size standard

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
Fig. 6. Agarose gel electrophoresis of PCR products amplified from DNA containing various amounts of NL-P
 Arrowheads indicate the expected PCR amplification products.
 Primer-pair No. 1 (lanes 3, 4), No. 3 (lanes 1, 2, 5~9)
 Lane 1: no potato DNA; lane 2, 3: non-GM DNA; lane 4: NewLeaf YTM potato DNA; lane 5: potato DNA containing 0.01% NL-P; lane 6: potato DNA containing 0.05% NL-P; lane 7: potato DNA containing 0.1% NL-P; lane 8: potato DNA containing 0.5% NL-P; lane 9: potato DNA containing 1.0% NL-P; M: 100 bp ladder size standard

Table 2. Results of NL-P in Processed Foods

Food item	Producing company	Pss 01n-5' p-FMV03-5'	
		Pss 01n-3' PLRV02-3' No. 1	No.
Snack	A	+	-
Snack	B	+	-
Snack	C	+	-
Snack	D	+	+
Snack	D	+	+
Snack	D	+	+
Snack	E	+	-
Snack	F	+	-
Snack	G	+	-
Mashed potato	H	+	-
Mashed potato	I	+	-
Stew or soup	J	-	-
Stew or soup	K	+	-
Stew or soup	D	+	-
Side dish	L	+	-
Frozen food	M	+	-
Crouquette	N	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-
Frozen fried potato	Imported food	+	-

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
Fig. 7. Agarose gel electrophoresis of potato PCR products amplified from potatoes and snack with NL-P specific primer pair No. 3 for detection
 Lane 1: negative control (no primer); lane 2: negative control (no DNA); lane 3: positive control (NL-P); lane 4: snack D1; lane 5: snack D2; lane 6: snack D3.

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A Detection Method for Recombinant DNA from Genetically
Modified Potato (NewLeaf Y[®] Potato)

Hiroshi AKIYAMA^{*1,†}, Takahiro WATANABE^{*1}, Chiseko WAKUI^{*1}, Yoshiko CHIBA^{*2},
Masaaki SHIBUYA^{*3}, Yukihiro GODA^{*1} and Masatake TOYODA^{*1}

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194-8543, Japan; *³Graduate School of Pharmaceutical Sciences, The University of Tokyo:

7-3-1, Hongo, Bunkyo-ku, Tokyo 113-0033, Japan; †Corresponding author)

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A Detection Method for Recombinant DNA from Genetically Modified Potato (NewLeaf Y[®] Potato)

Hiroshi AKIYAMA^{*1,†}, Takahiro WATANABE^{*1}, Chiseko WAKUI^{*1}, Yoshiko CHIBA^{*2},
 Masaaki SHIBUYA^{*3}, Yukihiro GODA^{*1} and Masatake TOYODA^{*1}

(^{*1}National Institute of Health Sciences: 1-18-1, Kamiyoga, Setagaya-ku, Tokyo 158-8501, Japan;
^{*2}Showa Pharmaceutical University: 3-3165, Higashitamagawagakuen, Machida-shi, Tokyo
 194-8543, Japan; ^{*3}Graduate School of Pharmaceutical Sciences, The University of Tokyo:
 7-3-1, Hongo, Bunkyo-ku, Tokyo 113-0033, Japan; [†]Corresponding author)

A detection method using the polymerase chain reaction (PCR) was developed to detect genetically modified (GM) potato (NewLeaf Y[®] potato; NL-Y), of which the mandatory assessment has not yet been completed in Japan. The potato sucrose synthase gene was used as an internal control. We designed a primer pair to specifically detect NL-Y without false-positive results in processed potato foods infected with the potato virus Y (PVY). The DNA introduced into NL-Y using the primer pair could be detected from potato powder samples containing 0.05% NL-Y. In addition, we designed primer pairs for recognizing the *CryIIIa* gene to detect the NewLeaf potato (NL), NewLeaf Plus potato (NL-P) and NL-Y and for recognizing *p-FMV* in order to detect NL-P and NL-Y. The proposed method was applied to the detection of NL-Y in 26 processed potato foods and NL-Y was not detected in any samples.

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
Fig. 1. Schematic diagram of PCR primers designed for detection of NewLeaf Y[®] potato (NL-Y). The expected lengths of PCR products are indicated in parentheses.

*5 JRL/ILSI Joint Workshop on Method Development in Relation to Regulatory Requirements for the Detection of GMOs in the Food Chain, Book of Abstracts, p. 49 (11, December, 2000)

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
Fig. 2. Specificity of GM potato-specific (primer pair No. 2)

Arrow indicates the expected PCR amplification products.

Lane 1: negative control (no DNA); Lane 2: NL SPBT02-05 DNA; Lane 3: NL-P DNA; Lane 4: NL RBBT06 DNA; Lane 5: NL-Y DNA; Lane 6: non-GM potato (Russet Burbank) DNA; Lane 7: non-GM potato (Superior) DNA; M: 100 bp ladder size standard


Fig. 4. Specificity of NL-Y-specific primers (primer pair No. 4)

Arrow indicates the expected PCR amplification products. See footnote for Fig. 2.


Fig. 3. Specificity of NL-P and NL-Y-specific primers (primer pair No. 3)

Arrow indicates the expected PCR amplification products. See footnote for Fig. 2.

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Table 2. Results of NewLeaf Y Monitoring in Processed Foods

Food item	Producing company	PSS 01n-5' p-FMV05-5'	
		PSS 01n-3' No. 1	PVY02-3' No. 5
Snack	A	+	-
Snack	B	+	-
Snack	C	+	-
Snack	D	+	-
Snack	D	-	-
Snack	D	+	-
Snack	E	+	-
Snack	F	+	-
Snack	G	-	-
Mashed potato	H	-	-
Mashed potato	I	+	-
Stew or soup	J	-	-
Stew or soup	K	-	-
Stew or soup	L	-	-
Side dish	M	-	-
Frozen food	N	-	-
Croquette	O	-	-
Starch	P	-	-
Frozen fried potato	Imported	-	-
Frozen fried potato	Imported	-	-
Frozen fried potato	Imported	-	-
Frozen fried potato	Imported	-	-
Frozen fried potato	Imported	-	-
Frozen fried potato	Imported	-	-
Frozen fried potato	Imported	+	-
Frozen fried potato	Imported	+	-
Frozen fried potato	Imported	+	-

+: positive, -: negative


Fig. 5. Specificity of NL-Y-specific primers (primer pair No. 5)
Arrow indicates the expected PCR amplification products. See footnote for Fig. 2.


Fig. 6. Agarose gel electrophoresis of PCR products amplified from potato DNA containing various amounts of NL-Y
Arrowheads indicated the expected PCR amplification products for PSS (216 bp) and NL-Y (225 bp). Primer pair No. 1 (Lanes 3 and 4). Primer pair No. 5 (Lanes 1, 2 and 5~9).
Lane 1: no potato DNA; Lanes 2 and 3: non-GM DNA; Lane 4: NL-Y DNA; Lane 5: potato DNA containing 0.01% NL-Y; Lane 6: potato DNA containing 0.05% NL-Y; Lane 7: potato DNA containing 0.1% NL-Y; Lane 8: potato DNA containing 0.5% NL-Y; Lane 9: potato DNA containing 1.0% NL-Y; M: 100 bp ladder size standard

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Primers for Potato Qualitative PCR

Pss primer pair (for detection of the potato endogenous gene) 216-bp

F-primer : Pss-5' (Pss01n-5'): 5'-TGA CCT GGA CAC CAC AGT TAT-3'

R-primer : Pss-3' (Pss01n-3'): 5'-GTG GAT TTC AGG AGT TCT TCG A-3'

New Leaf Y detection primer pair (cross over FMV and PVC-cp gene) 225-bp

F-primer : NLYd-5' (p-FMV05-5'): 5'-AAA AGA GCT GTC CTG ACA GC-3'

R-primer : NLYd3' (PVY02-3'): 5'-TCC TCC TGC ATC AAT TGT GT-3'

PVY-cp gene detection primer pair (for NLY confirmation) 161-bp

F-primer : NLYc-5' (PVY01-5') : 5'-GAA TCA AGG CTA TCA CGT CC-3'

R-primer : NLYc-3' (PVY01-3') : 5'-CAT CCG CAC TGC CTC ATA CC-3'

New Leaf plus detection primer pair (cross over FMV and PLRV-rep gene) 234-bp

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






2nd antibody binding (60 min)

↓ wash 3 time (15 min)

Alkali phosphatase conjugated Avidin-Biotin incubation(20 min)

↓ wash 3 time (15 min)

Tris/HCl equilibrium (15 min)

↓

Substrate adding, ddH₂O washing and incubate at dark (15 min)

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
Protocol for Mixing GM Sample Powder

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











Protocol for Mixing GM Containing Sample Powder

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<一> DNeasy Plant Mini Kit (QIAGEN) 法

Weight 2 g of ground sample, put into a 50 mL tube (polypropylene centrifuge tube)

↓
Add 10 mL of AP1 buffer (pre-warmed at 65°C)
Add 10 μ L of RNase A

↓
Vortex to mix very well, then incubate at 65°C for 15 min

↓
Invert the tube several times to mix

↓
Add 3.25 mL of AP2 buffer, then incubate at room temperature (RT) for 5 min

↓
3,000 \times g centrifugation at RT for 5 minutes

↓
Transfer supernatant to a new tube quickly (avoid to transfer pellet and upper-layer lipids)

↓
Load 3 mL supernatant (500 μ L per time) to a QIAshredder spin column

↓
10,000 \times g centrifugation at RT for 2 minutes

↓
Transfer flow-through to a new 15 mL tube (finally ~3 mL flow-through is obtained)

↓
Add 1.5 \times volume of AP3/Ethanol solution
AP3 : Ethanol (96~100%) = 1 : 2 (V/V)

↓
Vortex for about 10 seconds with maximum


↓
Apply 500 μ L per time to the Mini spin column

↓
10,000 \times g centrifugation at RT for 1 min

↓
Discard the flow-through


↓
Repeat this step until all the sample solution pass through Mini spin column

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
This method is used to extract DNA from lowly processed 、 low sugar containing and low fat containing food samples mainly.


<二> DNeasy Genomic Tip 20/G Kit (QIAGEN) 法

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This method is mainly used to extract DNA from highly processed including high sugar containing、high fat containing、heat treatment and fermentation food samples.

<三> CTAB 法

Weight 2g of sample, put in a 50 mL polypropylene centrifuge tube

- | ← Add 15 mL CTAB extraction solution *¹
- | To homogenize completely (using a homogenizer is recommend)
- | Incubate at 55°C for 30 min (with shaking)
- | Transfer 600 μL to a new 1.5 mL eppendorf tube

Phenol and Chloroform treatment

- | ← add 500 μL phenol/chloroform *², vortex vigorously to mix well
- | 7,500× g centrifugation at RT for 15 min, transfer up-layer solution (~ 600 μL) to a new 1.5 μL eppendorf tube

Chloroform treatment

- | ← add 500 μL chloroform/Isoamyl alcohol *³, vortex vigorously to mix well
- | 7,500× g centrifugation at RT for 15 min, transfer up-layer solution (~ 600 μL) to a new 1.5 μL eppendorf tube

Isopropanol precipitation

- | ← add the same volume isopropanol (~ 600 μL) into the new tube, then invert this tube several times to mix well
- | 7,500× g centrifugation at RT for 15 min, discard supernatant
- | ← add 500 μL 70% alcohol, wash carefully

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- | 7,500× g centrifugation at RT for 1 min, discard supernatant
- | vacuum dry for 2~3 min

dissolve DNA and remove RNA

- | add 50 μ L TE to dissolve DNA well (place at RT for 15 min and mix several times)
- | ← add 5 μ L RNase (10mg/mL), then incubate at 37°C for 30 min
- | ← add 200 μ L CTAB extraction solution, mix well

Chloroform treatment

- | ← add 250 μ L chloroform/Isoamyl alcohol, mix gently
- | 7,500× g centrifugation at RT for 15 min, transfer up layer solution to a new eppendorf tube

Isopropanol precipitation

- | ← add 200 μ L isopropanol, invert several times to mix well
- | 75,00× g centrifugation at RT for 10 min, then discard the supernatant
- | ← add 200 μ L 70% alcohol , invert to mix
- | after 75,00× g centrifugation at RT for 1 min, discard supernatant
- | vacuum dry for 2~3 min

Dissolve DNA

- | ← add 50 μ L ddH₂O to dissolved DNA

DNA sample solution

*1 CTAB extraction solution

8 mL 0.5mM EDTA (pH 8.0) + 20 mL 1M Tris/HCl (pH 8.0) + 56 mL 5M NaCl + H₂O 150 mL + 4g CTAB stir to dissolve. After dissolve well, add H₂O to a final 200 mL volumn. Autoclave.

*2 Phenol / Chloroform

1M Tris/HCl (pH 8.0) saturated phenol: Chloroform/ Isoamyl alcohol =1:1 (v/v)

*3 Chloroform : Isoamyl Alcohol (CIA)

Chloroform : Isoamyl Alcohol=24:1

This method can be used to extract DNA from both lowly processed and highly processed food samples. However, Genomic Tip 20/G method is recommend when extract DNA from highly processed food samples.

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Reprinted from the Journal of Food Hygienics Society of Japan

Vol. 44, No. 2, February 2003

J. Food Hyg. Soc. Japan
(SHOKUHIN EISEIGAKU ZASSHI)

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Original

Determination of Enzymatic Activity of 5-Enolpyruvylshikimate-3-phosphate Synthase by LC/MS

(Received July 31, 2002)

Haruyo OKUNUKI*¹, Hiroshi AKIYAMA*^{1,†}, Reiko TESHIMA*¹,
Akihiro HINO*², Yukihiko GODA*¹, Jun-ichi SAWADA*¹,
Masatake TOYODA*^{1,*3} and Tamio MAITANI*¹

(*¹National Institute of Health Sciences: 1-18-1, Kamiyoga, Setagaya-ku, Tokyo 158-8501, Japan;

*²National Food Research Institute: 2-1-12, Kannondai, Tsukuba, Ibaraki 305-8642, Japan;

*³Present Address, Jissen Woman's University: 4-1-1, Osakae, Hino 191-8510, Japan;

[†]Corresponding author)

A liquid chromatography-mass spectrometry (LC/MS) method for determining the enzymatic activity of 5-enolpyruvylshikimate-3-phosphate synthase (EPSP synthase), an enzyme of the shikimate pathway, was developed. EPSP synthase catalyzes the formation of 5-enolpyruvylshikimate-3-phosphate (EPSP) from shikimate-3-phosphate (S-3-P) and phosphoenolpyruvate (PEP) in microorganisms and plants. The enzymatic activity of EPSP synthase was assessed by the determination of EPSP after a 30-min incubation with S-3-P and PEP using the LC/MS system. EPSP synthase activity is given in terms of the produced EPSP (pmol/min/mg protein). Glyphosate (*N*-phosphonomethyl glycine)-tolerant EPSP synthase from the *Agrobacterium* sp. strain CP4 (CP4-EPSP synthase) in genetically modified soybeans (GM-soybeans) was found to have an enzymatic activity of 736 EPSP pmol/min/mg protein in the presence of 3 nmol of S-3-P. In contrast, the enzyme activity of non-GM-soybeans was 21 EPSP pmol/min/mg protein. The EPSP synthase activity was markedly decreased in the non-GM-soybeans by the addition of glyphosate, but the enzyme activity of the GM-soybeans was only slightly decreased with this treatment. This LC/MS system could also be applicable to the measurement of EPSP synthase activity in different plant species and the detection of herbicide-tolerant EPSP synthase in GM foods.

Key words: 5-enolpyruvylshikimate-3-phosphate synthase; LC/MS, genetically modified soybeans; shikimate-3-phosphate; glyphosate

Introduction

5-Enolpyruvylshikimate-3-phosphate synthase (EPSP synthase, EC 2.5.1.19), an enzyme of the shikimate pathway, catalyzes the formation of 5-enolpyruvylshikimate-3-phosphate (EPSP) from shikimate-3-phosphate (S-3-P) and phosphoenolpyruvate (PEP)¹⁻³. The shikimate pathway of aromatic amino acid biosynthesis is present in plants and microorganisms, but not in mammals, making it a good target for herbicidal attack⁴. The herbicide glyphosate (*N*-phosphonomethyl glycine) is a specific inhibitor of EPSP synthase and does not affect other PEP-related enzymatic reactions^{5, 6}.

Recently, a number of genetically modified (GM) crops have been introduced into the mainstream cultivation worldwide. Some genetically modified soybeans (GM-soybeans) contain bacterial genes such as epsps from the *Agrobacterium* sp. strain CP4 (CP4-epsps) that have been introduced to confer tolerance to glyphosate¹.

We have already reported a detection method for the recombinant CP4-EPSP synthase DNA in GM-soybeans using a PCR technique⁷. Recently, a variety of geneti-

cally modified soybeans which does not express CP4-EPSP synthase has been found, even though the CP4-EPSP synthase DNA was detected in the variety using the PCR method. Therefore, it is also important to investigate the activity of CP4-EPSP synthase. A radioactive labeling method and an HPLC method with UV detection have been reported for the measurement of EPSP synthase activity^{1, 8}. However, the methods involve time-consuming steps because of many interfering peaks on the HPLC chromatogram or the potential hazard of radioactivity.

In the present report, we describe a novel and simple method based on LC/MS for the determination of the enzymatic activity of EPSP synthase without using complicated extraction steps or radioactive labeling. This LC/MS method was applied to measure the enzymatic activities of CP4-EPSP synthase expressed in the *E. coli* strain BL21 and GM-soybeans, and of EPSP synthase in non-GM-soybeans. This method allows the sensitive measurement of EPSP synthase activity in soybean crude extract, and should be applicable to the determination of EPSP synthase in different species of plants⁹. We also compared the EPSP synthase ac-

tivities in the crude extracts of GM and non-GM soybeans in the presence of glyphosate and confirmed the presence of herbicide-tolerant CP4-EPSP synthase.

Materials and Methods

Preparation of S-3-P and [³²P] S-3-P

S-3-P was prepared the following method of Akiyama *et al.*¹⁰⁾ The reaction buffer contained 10 μ L (protein content 260 μ g) of shikimate kinase II (in destruction buffer of *E. coli* including the *aroL* expression vector), 10 mmol/L shikimic acid (Wako Pure Chemicals Industries, Ltd., Osaka, Japan), 16 mmol/L ATP (Wako Pure Chemicals Industries, Ltd.), 16 mmol/L MgCl₂, 1 mmol/L (NH₄)₆Mo₇O₂₄, 10 mmol/L benzamidine HCl, 50 mmol/L KF and 50 mmol/L HEPES buffer (pH 7.0) in a 100 μ L total volume. The mixture was incubated at 37°C overnight. After deproteinization using a Tosoh Air Press-30 (Tosoh Co., Tokyo, Japan), the reaction mixture was subjected to the HPLC method described below. The S-3-P peak fraction was collected repeatedly and evaporated to dryness at 40°C. The residue was weighed to determine the amount of purified S-3-P.

When [³²P] ATP (2.0 \times 10⁸ dpm/64 nmol) was used for preparing [³²P] S-3-P, the hot reaction mixture was assayed under the radio-HPLC in the conditions described below. For the control experiments, reactions without shikimate kinase II were run in parallel.

HPLC conditions for preparing S-3-P

A Shimadzu LC-10AD pump and an SPD-10AV UV detector (Shimadzu Co., Kyoto, Japan) were used as the HPLC system. The column was a TSK gel CARBON-500 (Tosoh Co.). The column was maintained at 40°C with a flow rate of 1 mL/min. The mobile phase was 0.1% trifluoroacetic acid (TFA), and S-3-P was monitored at 215 nm. The injection volume was 100 μ L. The determination of S-3-P obtained by enzyme reaction (shikimate kinase II) was performed by LC/MS¹⁰⁾.

HPLC conditions for analyzing EPSP

EPSP was analyzed according to the method reported by Washizu with some modifications⁹⁾. A Shimadzu LC-10AD pump and an SPD-10AV UV detector (Shimadzu Co.) were used as the HPLC system. The column was a Senshu Pak NH₂-1251 (4.6 \times 250 mm i.d.; Senshu Scientific Co., Tokyo, Japan). The operating conditions were: column temperature, 40°C; flow rate, 1 mL/min; injection volume, 20 μ L; detector wavelength, 215 nm. The mobile phase was 200 mmol/L potassium dihydrogen phosphate (pH 7.0).

Radio-HPLC conditions for preparing [³²P] S-3-P and analyzing [³²P] EPSP

[³²P] S-3-P were prepared according to the method reported by Washizu with some modifications⁹⁾. A Hitachi L-7100 pump, an L-7400 UV detector and L-7300 column oven (Hitachi Co., Tokyo, Japan) were used as the HPLC system. The column was a Senshu Pak NH₂-1251 (4.6 \times 250 mm i.d.; Senshu Scientific Co., Tokyo,

Japan). The operating conditions were: column temperature, 40°C; flow rate, 1 mL/min; injection volume, 20 μ L; detector wavelength, 215 nm. The mobile phase was 200 mmol/L potassium dihydrogen phosphate (pH 7.0). The eluted fractions (each 1 mL) were collected every minute. A 50 μ L volume of each fraction was added to 9 mL of liquid scintillation fluid (Perkin-Elmer, Inc., Wellesley, Massachusetts, USA) and the [³²P] radioactivity was detected using a liquid scintillation counter, SC-5100 (Aloka Co., Tokyo, Japan).

Preparation of CP4-EPSP synthase in *E. coli*

The entire coding sequence of CP4-EPSP synthase was PCR-amplified using a set of overlapping oligonucleotide primers, PR01-5' (5'-TCCTTCGCAAGACCCTT-CCTCTAT-3') and R1-EPSP synthase primer including *EcoRI* and *SacI* sites (5'-GCTGCCTGATGAGCTCGAA-TTCGA-3'), based on the reported GM-soybean genomic DNA sequence¹¹⁾, using the GeneAmp PCR system 9600 (Applied Biosystems, Foster City, CA, USA). The reaction mixture (25 μ L) contained 0.625 U Taq polymerase, 0.3 μ g GM-soybean DNA, 3.5 μ L of PCR premix (Takara Shuzo Co., Ltd., Tokyo, Japan), and 500 nmol/L primers. The PCR reactions were carried out for 40 cycles (94°C: 30 sec, 62°C: 1 min, 74°C: 2.5 min for each cycle) after incubation at 94°C for 4 min, and final extension was done at 74°C for 10 min. The amplified fragment was purified by gel electrophoresis in 1.5% agarose (GFX PCR DNA and Gel Band Purification Kit, Amersham Pharmacia Biotech Co., Tokyo, Japan) and used as a template for the 2nd PCR using the F1-EPSP synthase primer with an *NdeI* site (5'-TTTCATCATGCTTCAC-GGTGCAAGCAGCCGG-3') and the R1-EPSP synthase primer. The PCR product was extracted with phenol, chloroform, precipitated with ethanol, and digested with restriction enzymes *NdeI* and *SacI* (Toyobo Co., Ltd., Osaka, Japan). The digested fragments were inserted between the *NdeI* and *SacI* sites of the pET23b(+) expression vector (Novagen, Inc., Madison, Wisconsin, USA) and cloned. Briefly, 55.5 ng of the insert DNA and 51 ng of pET23b(+) were ligated at 16°C for 24 hours with 7.5 μ L of Takara Ligation Kit Ver. 2 solution 1 (Takara Shuzo Co., Ltd., Tokyo, Japan). After subcloning of the inserted DNA (1,388 bp), the nucleotide sequence was determined by the dideoxynucleotide chain termination method¹²⁾. *E. coli* strain BL21 was transformed with the CP4-EPSP synthase cDNA expression vector pET23b(+) by the heat shock method. The cloned BL21 cells including CP4-EPSP synthase cDNA-pET23b(+) in 200 mL of NZCYM culture medium (NZ amine 10 g, Bacto yeast extract 5 g, NaCl 5 g, casamino acids 1 g, MgCl₂·7H₂O 2 g/L) (Sigma Chemical Co., St. Louis, MO, USA) were stimulated with 0.5 mmol/L isopropyl β -D-1-thiogalactopyranoside (IPTG; Sigma Chemical Co., St. Louis, MO, USA) to induce the CP4-EPSP synthase protein expression and then suspended in disruption buffer for BL21 (5 mmol/L MgCl₂, 50 mmol/L Tris-HCl (pH 7.0), 0.37% KCl), and sonicated. The supernatant was collected and used as the CP4-


Fig. 1. Shikimic acid pathway and related enzymes

Shikimate kinase II catalyzes the phosphorylation of shikimic acid to shikimate-3-phosphate (S-3-P). 5-Enolpyruvylshikimate-3-phosphate (EPSP) synthase catalyzes the reaction of S-3-P and phosphoenolpyruvate (PEP) to form EPSP and phosphate.

EPSP synthase from *E. coli*.

Extraction of EPSP synthase from GM- and non-GM-soybeans

Powders of 50 g of GM- and non-GM-soybeans, which were products from the state of Ohio, USA, and purchased from a Japanese trading company, respectively, were used. The soybean powder was prepared by the method of Teshima *et al.*¹³ A 200-mL volume of acetone was added to the soybean powder and the mixture was stirred for 30 min at room temperature. The acetone solution was removed by filtration, and 100 mL of acetone was added to the precipitate again. The suspension was stirred for 20 min, then filtered and the precipitates were dried. Next, 250 mL of 50 mmol/L PBS (pH 7.4) with 0.05% sodium azide was added to the dried powder and the mixture was stirred overnight at 4°C. On the next day, it was centrifuged at 40,000 rpm for 1 hour at 4°C, and the supernatant was dialyzed against PBS for 2 days. The dialyzed supernatant was used as the extract for determining the EPSP synthase activity. The protein concentration of the extract was determined by the Bio-Rad protein assay (Bio-Rad Laboratories, Hercules, CA, USA) based on the method of Bradford¹⁴.

Production and determination of EPSP

According to the reports of Padgett *et al.*^{9, 11}, we examined the conditions for the determination of EPSP synthase activity. Among various conditions examined, the following conditions were concluded to be optimal for the determination of EPSP synthase activity using LC/MS. The reaction mixture (100 μ L) containing 500 μ M PEP, 3.2 mmol/L MgCl₂, 10 μ L of 3 nmol S-3-P, 100 mmol/L KCl, 50 mmol/L KF, 5 mmol/L HEPES buffer (pH 7.0) and 1 μ L of crude extract of the *E. coli* expressing CP4-EPSP synthase (protein content 1.5 μ g) or soybeans was incubated at 30°C for 30 min. The reaction tube was put on ice to stop the reaction after incubation since EPSP synthase activity is sensitive to temperature. The reaction mixture was applied to the LC/MS system. In advance, we had prepared an EPSP standard solution (10 ng/mL) to determine the extent of conversion of PEP to EPSP during the CP4-EPSP synthase reaction using the HPLC method for analyzing EPSP. The relationship between PEP reduc-

tion and EPSP production was linear up to about 25% conversion using the HPLC method with UV detection. EPSP was determined with the LC/MS system described below using the standard solution. Control reactions lacking CP4-EPSP synthase were also run in parallel. EPSP synthase activity is given in terms of produced EPSP (pmol/min/mg protein) at 30°C.

LC/MS conditions

The LC/MS system consisted of an Alliance 2790 HPLC system, 996 PDA detector and ZQ MS detector (Waters Co., MA, USA). A Develosil C30 column (4.6 mm i.d. \times 250 mm, Nomura Chemical, Tokyo, Japan) was used as the separation column and was maintained at 50°C. The mobile phase was 3% methanol solution containing 0.1% TFA at the flow rate of 0.7 mL/min. The cone voltage was 30 V. The ionization mode was negative ion mode of electrospray ionization (ESI) was used. The injection volume was 5 μ L.

Results and Discussion

Assay for enzymatic activity of CP4-EPSP synthase expressed in *E. coli*

To develop an assay method for the enzymatic activity of EPSP synthase by LC/MS analysis, CP4-EPSP synthase was over-expressed in *E. coli* for use as the positive control^{11, 15}. Since S-3-P, which is a substrate for measuring the enzymatic activity of EPSP synthase, was not commercially available, shikimate kinase II (S-3-P synthase), which biosynthetically produces S-3-P (Fig. 1), was expressed in *E. coli* as previously described¹⁰. The molecular weights of the expressed shikimate kinase II and expressed CP4-EPSP synthase were determined by SDS-PAGE to be 19 kDa and 44 kDa, respectively (data not shown), similar to the molecular mass of the enzyme predicted from the DNA sequence¹⁶⁻¹⁸. The complete reaction mixture containing CP4-EPSP synthase (protein content 1.5 μ g) and S-3-P (3 nmol) was subjected to HPLC method reported by Washizu with some modifications; an amino-silica column (Senshu Pak NH₂-1251) was eluted with 200 mmol/L KH₂PO₄ (pH 7.0) as shown in Fig. 2⁸. The chromatogram showed an additional peak indicated by the arrow at 13 min (Fig. 2A), which was not obtained from the sample solution without CP4-EPSP synthase (Fig. 2B). Thus, the 13-min peak appears to correspond to EPSP.

Detection of [32 P] S-3-P and [32 P] EPSP on radio-HPLC

To further confirm that the peak detected on HPLC described above was EPSP, we used a radio-HPLC procedure as described in Materials and Methods section⁹⁾. We added [γ - 32 P] ATP (2.0×10^8 dpm/64 nmol) to an enzyme reaction mixture for synthesis of [32 P] S-3-P. A single radioactive peak appeared at the retention time of 6 min (Fig. 3a), whereas the control solution reacted without shikimate kinase II showed no radioactive peak at this position (Fig. 3b). Typically, a 95% conversion

of [32 P] ATP to [32 P] S-3-P was achieved after 24 hours of incubation at 37°C. The peak fraction of [32 P] S-3-P was collected, added to a reaction mixture containing CP4-EPSP synthase from *E. coli* and incubated at 30°C for 30 min. The radio-chromatogram of the incubation mixture is shown in Fig. 3c. After the peak of [32 P] S-3-P was eluted, a new radioactive peak corresponding to [32 P] EPSP appeared at 10 min (Fig. 3c), which was not found in the reaction mixture without CP4-EPSP synthase (Fig. 3d). Though the retention time (10 min) of [32 P] EPSP was different from that (13 min) of EPSP detected above (Fig. 2a), the difference of the retention time was thought to be due to the difference of HPLC systems used.


Fig. 2. Detection of EPSP using the HPLC method. HPLC chromatograms of the reaction mixtures with CP4-EPSP synthase from *E. coli* (A) and without CP4-EPSP synthase (B). The column was a Senshu Pak NH₂-1251 (4.6 × 250 mm i.d.). The operating conditions were as follows: column temperature, 40°C; flow rate, 1 mL/min; injection volume, 20 μ L; detector wavelength, 215 nm. The mobile phase was 200 mmol/L potassium dihydrogen phosphate (pH 7.0).


Fig. 3. Detection of [32 P] S-3-P and [32 P] EPSP using the radio-HPLC method.

HPLC chromatograms of the reaction mixtures with shikimate kinase II (a) and without shikimate kinase II (b). Chromatograms of the reaction mixtures containing [32 P] S-3-P obtained (a) with CP4-EPSP synthase from *E. coli* (c) and without CP4-EPSP synthase (d).


Fig. 4. Total ion chromatogram (a) and selected ion (m/z 323) chromatogram (b) for the LC/MS analysis of EPSP produced by the reaction mixture including GM-soybean crude extract.

The LC/MS conditions were as follows. A devesil C30 column was used as the separation column and was maintained at 50°C. The mobile phase was 3% methanol solution containing 0.1% TFA at the flow rate of 0.7 mL/min. The cone voltage was 30 V. The ionization mode was negative ion electrospray ionization (ESI). The injection volume was 5 μ L.


Fig. 5. ESI mass spectrum for the EPSP peak detected by the LC/MS analysis of EPSP produced by the reaction mixture including GM-soybean crude extract

Determination of EPSP using LC/MS system

To easily determine EPSP without using radio-HPLC, we developed an LC/MS method. Among the several reverse-phase columns examined, Develosil C30 showed the best retention of EPSP in combination with the volatile mobile phase described in the Materials and Methods section. The established HPLC separation system was adapted to the LC/MS system. Negative ion ESI produced deprotonated molecules $[M-H]^-$ (m/z 323) of EPSP without any fragmentation. Furthermore, the selected ion monitoring (SIM) mode unambiguously identified EPSP. Figures 4 and 5 show representative MS chromatograms and an ESI mass spectrum of EPSP, respectively. The relative standard deviation was less than 5% (3 ng/mL of EPSP) after six repeated injections. EPSP could be determined quantitatively in the range of 0.3–5 ng/mL. The detection limit was 0.1 ng/mL. This LC/MS method was used to analyze the EPSP synthase activity.

Determination of EPSP synthase activity in GM-soybean and non-GM-soybean crude extracts by LC/MS

We extracted proteins from GM-soybeans and attempted to detect the EPSP synthase enzymatic activity using the LC/MS method as described above. The enzymatic activity was calculated from the EPSP determination. The crude extract of the GM-soybeans (diluted 10, 50, 100, 500 and 1,000 times) was incubated with reaction mixture containing S-3-P at 30°C for 30 min. We confirmed the correlation between the dilution of the soybean crude extract and the production of EPSP.

Next, we compared the enzyme activity of CP4-EPSP synthase produced by *E. coli* and the soybean crude extracts. The enzymatic activity of EPSP synthase in the *E. coli* and GM- and non-GM-soybean crude extracts was $63,267.7 \pm 4,933.3$ EPSP pmol/min/mg protein, 735.6 ± 67.8 EPSP pmol/min/mg protein and 21.1 ± 10.3 EPSP pmol/min/mg protein, respectively. Each value represents the mean \pm SD from triplicate analyses.

CP4-EPSP synthase in the GM-soybeans is glyphosate-resistant. To measure only the glyphosate-resistant EPSP synthase in the GM-soybeans, we incubated S-3-P,


Fig. 6. EPSP synthase activity of GM- and non-GM-soybean extracts in the presence of glyphosate. EPSP synthase activity of GM- or non-GM-soybean crude extracts was assayed in the presence of 0–500 μ mol/L of glyphosate. EPSP produced by the reaction was measured by LC/MS analysis. The protein concentration of GM- and non-GM-soybean extract was 19 mg/mL. Each value represents the mean \pm SD from triplicate analyses.

PEP and the crude extracts of the GM-soybeans in the presence of glyphosate (0–500 μ mol/L) and analyzed these using the LC/MS method as described above. The crude extract from non-GM-soybeans was examined in parallel. Figure 6 represents the EPSP synthase activity in the GM- and non-GM-soybean crude extracts for several concentrations of glyphosate. The EPSP synthase activity of the GM-soybeans was unchanged or only slightly reduced in the presence of 50–500 μ mol/L of glyphosate, whereas the EPSP synthase activity in the non-GM-soybeans was dose-dependently inhibited by the glyphosate. The concentration of glyphosate which caused 50% inhibition of the EPSP synthase activity of the non-GM-soybeans (IC_{50}) was 258 μ mol/L, and 90% of the EPSP synthase activity was inhibited at 500 μ mol/L of glyphosate. The IC_{50} value obtained with our system was similar to that reported previously^{6, 8, 19}. Therefore, the GM-soybeans show herbicide-resistant EPSP synthase activity.

Padgett *et al.*¹¹ reported a radio-HPLC assay for the determination of EPSP synthase activity using ¹⁴C label. However, this poses a potential hazard. Moreover, Washizu also reported an HPLC method using an

amino-silica column for the determination of EPSP synthase⁹⁾. The method involved time-consuming steps because of many interfering peaks on the HPLC chromatogram. It also seems to be difficult to identify the EPSP peak reliably because of the rapid decline in the quality of the amino-silica column. Therefore, we used a reverse-phase C30 column to retain EPSP, because this should be more durable than the amino-silica column. The separation system of the HPLC was adapted to the LC/MS system for the qualitative and quantitative analysis of EPSP. The established LC/MS method for EPSP analysis proved to be more sensitive and precise than the HPLC method with UV detection. It provides a simple procedure for the determination of CP4-EPSP synthase activity or EPSP synthase activity by LC/MS without using any radioactive label.

Furthermore, we confirmed that the glyphosate-resistant enzymatic activity of CP4-EPSP synthase from enzyme reaction mixtures containing \dot{S} -3-P and PEP could be detected using the LC/MS method in the presence of 500 μ mol/L glyphosate. This method should be applicable to the detection of other glyphosate-resistant EPSP synthase species, which have been introduced into various grains or rapeseeds.

Acknowledgments

We express our thanks to Dr. S. Arami and Dr. H. Miura of the National Food Research Institute, and Ms. S. Tsuzuki of NIHS for their technical assistance and Dr. H. Sasaki of Nihon Waters Co. for advice on LC/MS technique. A part of this study was supported by a grant for Scientific Research Expense for Health and Welfare Programs from the Japanese Government.

References

- 1) Padgette, S. R., Huynh, Q. K., Aykent, S., Sammons, R. D., Sikorski, J. A., Kishore, G. M., Identification of the reactive cysteines of *Escherichia coli* 5-enolpyruvylshikimate-3-phosphate synthase and their nonessentiality for enzymatic catalysis. *J. Biol. Chem.*, **263**, 1,798-1,802 (1988).
- 2) Steinrucken, H. C., Amrhein, N., 5-Enolpyruvylshikimate-3-phosphate synthase of *Klebsiella pneumoniae*. I. Purification and properties. *Eur. J. Biochem.*, **143**, 351-357 (1984).
- 3) Boocock, M. R., Coggins, J. R., Kinetics of 5-enolpyruvylshikimate-3-phosphate synthase inhibition by glyphosate. *FEBS Lett.*, **154**, 127-133 (1983).
- 4) Haslam, E. ed., "The Shikimate Pathway", New York, John Wiley & Sons, 1974, p. 17-19.
- 5) Smart, C. C., Johanning, D., Muller, G., Steinrucken, H. C., Selective overproduction of 5-enol-pyruvylshikimic acid 3-phosphate synthase in a plant cell culture which tolerates high doses of the herbicide glyphosate. *J. Biol. Chem.*, **260**, 16,338-16,346 (1985).
- 6) Steinrucken, H. C., Amrhein, N., The herbicide glyphosate is a potent inhibitor of 5-enolpyruvylshikimic acid-3-phosphate synthase. *Biochem. Biophys. Res. Commun.*, **94**, 1,207-1,212 (1980).
- 7) Matsuoka, T., Kawashima, Y., Akiyama, H., Miura, H., Goda, Y., Sebata, T., Isshiki, K., Toyoda, M., Hino, A., A detection method for recombinant DNA from genetically modified soybeans and processed foods containing them. *J. Food Hyg. Soc. Japan*, **40**, 149-157 (1998).
- 8) Washizu, M., An inhibition assay for EPSP synthase by HPLC. *Mitsubishi Kasei R & D Review*, **5**, 70-75 (1991).
- 9) Padgette, S. R., Huynh, Q. K., Borgmeyer, J., Shah, D. M., Brand, L. A., Re, D. B., Bishop, B. F., Rogers, S. G., Fraley, R. T., Kishore, G. M., Bacterial expression and isolation of *Petunia hybrida* 5-enol-pyruvylshikimate-3-phosphate synthase. *Arch. Biochem. Biophys.*, **258**, 564-573 (1987).
- 10) Akiyama, H., Okunuki, H., Tsuzuki, S., Arami, S., Miura, H., Sakushima, J., Teshima, R., Goda, Y., Hino, A., Sawada, J., Toyoda, M., Construction of the expression system of shikimate kinase II and the preparation of shikimic acid 3-phosphate. *J. Food Hyg. Soc. Japan*, **40**, 438-443 (1999).
- 11) Padgette, S. R., Kolacz, K. H., Delannay, X., Re, D. B., La-Vallee, B. J., Tinius, C. N., Rhodes, W. K., Otero, Y. I., Barry, G. F., Eichholtz, D. A., Peschke, V. M., Nida, D. L., Taylor, N. B., Kishore, G. M., Development, identification, and characterization of a glyphosate-tolerant soybean line. *Crop Sci.*, **35**, 1,451-1,461 (1995).
- 12) Sanger, F., Nicklen, S., Coulson, A. R., DNA sequencing with chain-terminating inhibitors. *Proc. Natl. Acad. Sci. USA*, **74**, 5,463-5,467 (1977).
- 13) Teshima, R., Akiyama, H., Okunuki, H., Sakushima, J., Goda, Y., Onodera, H., Sawada, J., Toyoda, M., Effect of GM and non-GM soybeans on the immune system of BX rats and B10A mice. *J. Food Hyg. Soc. Japan*, **41**, 188-193 (2000).
- 14) Bradford, M., A rapid and sensitive method for the quantitation of microgram quantities of protein utilizing the principle of protein-dye binding. *Anal. Biochem.*, **72**, 248-254 (1976).
- 15) Shah, D. M., Horsch, R. B., Klee, H. J., Kishore, G. M., Winter, J. A., Tumer, N. E., Hironaka, C. M., Sanders, P. R., Gasser, C. S., Aykent, S., Siegel, N. R., Rogers, S. G., Fraley, R. T., Engineering herbicide tolerance in transgenic plants. *Science*, **233**, 478-481 (1986).
- 16) Laemmli, U. K., Cleavage of structural proteins during the assembly of the head of bacteriophage T4. *Nature*, **227**, 680-685 (1970).
- 17) Millar, G., Lewendon, A., Hunter, M. G., Coggins, J. R., The cloning and expression of the *aroL* gene from *Escherichia coli* K12. Purification and complete amino acid sequence of shikimate kinase II, the *aroL*-gene product. *Biochem. J.*, **237**, 427-437 (1986).
- 18) Kishore, G. M., Brundage, L., Kolk, K., Padgette, S. R., Rochester, D., Huynh, Q. K., Della-Cioppa, G., Isolation, purification and characterization of a glyphosate tolerant mutant *E. coli* EPSP synthase. *Fed. Proc. Amer. Soc. Exp. Biol.*, **45**, 1,506 (1986).
- 19) Steinrucken, H. C., Schulz, A., Amrhein, N., Porter, C. A., Fraley, R. T., Overproduction of 5-enolpyruvylshikimate-3-phosphate synthase in a glyphosate-tolerant *Petunia hybrida* cell line. *Arch. Biochem. Biophys.*, **244**, 169-178 (1986).

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
National Food Research Institute

Mission

National Food Research Institute (NFRI) is one of independent administrative institutions (IAI) under the supervision of the Ministry of Agriculture, Forestry and Fisheries (MAFF). The mission of NFRI is to provide the society with healthy and enriched life and secure supply of safe food through conducting research. Its food research projects cover a wide range of scientific and technical fields and are aimed at developing technological systems that strengthen the agriculture, forestry, fisheries, and food processing industries.

These research projects include: the elucidation of the relation of various foods to health, namely, functional evaluation of food components; improvement of technology for analysis; development of technology to ensure food safety; improvement of food distribution and processing methods; and exploration and utilization of functions of living organisms.

Medium Term Plan for Research (2001-2005)

Elucidation of functional properties of food and development of technologies to utilize these properties

Development of methods of analysis in support of the food labeling

Development of technologies to ensure the safety and high quality of food

Development of technologies for utilizing various food materials, and advancement of food manufacturing techniques which meet the environmental needs of the society

Development of technologies for effectively using microorganisms and enzymes

Elucidation of biological functions of molecules through the latest biochemical methods and development of technologies for utilizing these molecules

Contact our Research Planning and Coordination Division for Collaboration and/or Information:

- () Joint research with private enterprises, various public corporations, and local governmental agencies;
- () Research entrusted to the Institute by private enterprises, various public corporations, and local governmental agencies, and research entrusted to other organizations or bodies by the Institute on specific topics of the Institute's projects;
- () Acceptance of trainees from private enterprises, various public organizations, and local governmental agencies, and of researchers with special status and post-doctorate fellows;
- () Cooperation with international organizations, such as the United Nations University, and exchange of researchers through fellowship programs and bilateral cooperation arrangements;
- () International/national seminars, workshops and symposia on a variety of topics including food-related research; and
- () Tour in the Institute and explanation of its research projects for visitors upon request.

Research Activities

Research Teams

Five research teams in the division of planning and coordination for needed research relating to foods

- ① Development of measures for food protection
- ② Food function and quality
- ③ Promotion of genome information analysis of useful microorganisms for fermentation industry
- ④ Creation of functionally-strengthening macromolecules based upon protein engineering, 3D analyses and bioinformatics
- ⑤ Development of micro-channel arrays and their applications including blood rheology measurement


Staphylococcus aureus

Enterohaemorrhagic Escherichia coli O157:H7

Food Function Division

The nutritional, gustatory, physical and physiological functionalities of foods and food components

- ① The physiological activity of food components to reduce serum lipid and prevent obesity
- ② Mechanism of taste response and its physiological regulation
- ③ Physicochemical properties of food and the physiological effects on humans
- ④ Research into functional food factors in vegetables and fruits and their effects
- ⑤ Mechanism of digestion and absorption, and improving effects on brain function


Lipoprotein production and hepatic fatty acid metabolism

Food Safety and Quality Division

Food safety and preservation of quality during distribution of food resources and foods

- ① Evaluation and method development of safety of food resources and foods
- ② Control of microorganisms and microbial toxins contaminating food resources and foods
- ③ Quality control of vegetables and fruits based on physiological properties
- ④ Food packaging and quality maintenance of packaged foods
- ⑤ Ecology, physiology and control of food insect pests

Aspergillus flavus: aflatoxin-producing fungus


Trichogramma chilonis Ishii, the egg parasitoid of *Coryca cephalonica* Stainton

Food Material Division

Biochemical and structural properties of food ingredients and development of new technologies for processing of agricultural products

- ① Development of new food carbohydrates and related enzymes
- ② Evaluation of quality and function of food proteins
- ③ Elucidation of structure and function of lipophilic micronutrients
- ④ Clarification of factors affecting rice quality and development of technology for quality evaluation and utilization
- ⑤ Development of new food manufacturing technology of cereals


Electron microscope of a rice grain
S: Starch granule P: Protein body

The differences between *Japanica* and *Indica* rice grains

Food Engineering Division

Food processing, distribution, nano-measurement, information technology and other technologies

- Development of advanced food processing and waste treatment system
- Development of novel dispersion, separation and reactor system
- Development of novel analyzing methodology for nano-structure and function of foods and biomaterials
- Development of technology such as soft electron, quality measurement and network application system
- Development of environmentally-friendly system for high quality food distribution


Microchannel(MC) Emulsification Mechanism

Applied Microbiology Division

Molecular analysis of food microorganisms and application of microbial functions in food processing

- Molecular biological analysis of the genes expressed in Koji mold used in Miso, Soy sauce and Sake production
- Molecular and cellular biology of baker's yeast and lactic acid bacteria
- Molecular analysis of the bacterial γ -polyglutamic acid (the sticky polymer of Natto) production
- Development of enzymatic system for the production of useful amino sugars

Metabolism of γ -polyglutamic acid, the sticky polymer of Natto

Biological Function Division

Engineering and biotechnology of biomolecules and cells

- Macromolecular structure and catalytic mechanism
- Improvement of enzyme function by gene manipulation
- Secondary metabolism and ribosome function of actinomycetes and *Bacillus subtilis*
- Cell metabolism and regulation

Tertiary structure of a protease exhibiting improved heat stability, constructed by gene replacement

Analytical Science Division


Studies on food analysis and quality evaluation

- Development of sensitive and precise analytical methods and procedures to determine food components
- Analysis of structure and dynamics of food components at molecular level in relation to their function
- Development of nondestructive methods for bio-measurement and food analysis
- Development of spectral analysis and image analysis methods for quality evaluation of foods
- Collection and preparation of data for "New Standard Tables of Food Composition in Japan"

High performance instrument for the analysis of molecular structure and dynamics

Organization

History

- 1934 Established as Rice Utilization Laboratory
- 1947 Renamed as Foods Research Institute
- 1972 Renamed as National Food Research Institute
- 1979 Moved from Tokyo to Tsukuba, the present address
- 2001 National Food Research Institute became an independent administrative institution

Access

■ By JR Joban Line



60 minutes from Ueno station in Tokyo to Ushiku station (20 minutes from Ushiku sta. to NFRI by bus)

■ By Express Bus

70 minutes from Tokyo bus terminal (Yaesu South Exit of JR Tokyo station) to NFRI

■ By Car

Approximately 60 minutes from Tokyo down town to Yatabe IC by Expressway

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II. General precautions

PCR is so powerful that even a trace amount of a template DNA can be amplified to a detectable level. It is thus critical to prevent the contamination of samples by undesired DNA, especially carryover PCR products, and by the DNases excreted from human skin that enzymatically degrade PCR samples. While ~~Part 5—Prevention of Contamination~~ of the Manual describes in detail the proper procedures to prevent contamination, the following general precautions may be useful:

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3.1.4 Equipment

- Centrifuge with a swing-bucket rotor and capable of spinning 50 ml polypropylene tube at 3,000 x g;
- Refrigerated centrifuge capable of spinning 2 ml microtubes;
- Micropipettes for volume ranges of 0.5-10, 10-100, 100-1,000 and 1,000-5,000 μl ^(Note 2);
- Water bath; and
- Vortex test tube mixer.

3.1.5 Reagents

Qiagen's DNeasy Plant Maxi kit.

using the method described in Part 6—Quantitative PCR of the Manual, that are present in the DNA extracts from the quoted and alternative extraction methods.

^(Note 2) These volume ranges are only examples. (??)

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- (15) After leaving at room temperature for 5 minutes, spin at 3,000 x g for 10 minutes at room temperature in a swing-out rotor.
- (16) Using a 100-1,000 μ l micropipette, measure the volume of the eluate and transfer it to a 2 ml sample tube. Add an equal volume of isopropanol to the tube using a 100-1,000 μ l micropipette.
- (17) After mixing by slowly turning the tube upside down 10 times, leave the solution at room temperature for five minutes.
- (18) After spinning the tube at 12,000 x g at 4°C for 15 minutes in a swing-out rotor, discard the supernatant using a 100-1,000 μ l micropipette.
- (19) Add 500 μ l 70% ethanol using a 100-1,000 μ l micropipette and tap the end of the tube with a finger until the precipitate on the bottom of the tube is completely resuspended in the solution.
- (20) After spinning at 12,000 x g for 3 minutes at 4°C in a swing-out rotor, discard the entire supernatant using a 10-100 μ l micropipette and dry the precipitate.
- (21) Add 100 μ l TE buffer, pH 8.0, using a 10-100 μ l micropipette to dissolve the precipitate. See 3.1.10 *Additional notes* for further information.
- (22) Repeat tapping the tube with a finger and spinning it until the liquid drops on the tube wall completely disappear. Incubate the solution overnight (12-24 hours) in a refrigerator.
- (23) A DNA lysate is finally produced when the solution is free from any visible insoluble matter. If any visible insoluble matter is still present even after 24 hours of incubation, use the supernatant obtained by spinning the solution at 12,000 x g for 3 minutes at 4°C and transferred to a new tube as a DNA lysate. Store the DNA lysate as well as the precipitate at -20°C or below.

3.1.6.2 Procedure B: DNA extraction using DNeasy Plant Maxi kit

- (1) Measure an appropriate amount of the ground sample into a 50 ml tube and add 10 μ l RNase (supplied in the kit) using a 0.5-10 μ l micropipette and 5 ml Buffer AP1 preheated to 65°C using a 1,000-5,000 μ l micropipette. After mixing the contents of the tube by repeatedly turning the tube upside down until no visible sample remains on the bottom of the tube, mix the solution on a vortex tube mixer.
- (2) Incubate the solution for one hour in a 65°C water bath. Every 15 minutes during incubation or 15, 30 and 45 minutes after incubation started, mix the solution by vigorously inverting the tube several times and then mix it at maximum speed on a vortex tube mixer.
- (3) After adding 1.8 ml Buffer AP2 to the tube using a 1,000-5,000 μ l micropipette and mixing the solution at maximum speed for 10 seconds on a vortex tube mixer, incubate the tube in an ice-water bath for 15 minutes.
- (4) Spin the solution at 3,000 x g for 15 minutes at room temperature in a swing-out rotor.
- (5) Aspire 4.2 ml of the supernatant using a 1,000-5,000 μ l micropipette, carefully avoiding the precipitate on the bottom of the tube or film on the liquid surface, to

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Aspirate 5 μ l of the DNA lysate solution using a 0.5-10 μ l micropipette and add TE buffer to give a total volume of 50 μ l. Scan the UV absorbance in the wavelength range from 200 to 300 nm to measure the absorbance at 230, 260 and 280 nm^(Note 3), ^(Note 4). Calculate the concentration of DNA by converting an absorbance (O.D.) of 1.0 at 260 nm to 50 ng of DNA per μ l^(Note 5). Prepare a PCR solution with a DNA concentration of 10ng/ μ l using sterilized water. If the concentration of the extracted DNA was too low for use as a PCR solution, perform the DNA extraction procedure again. If the concentration of the re-extracted DNA was still too low for use as a PCR solution, use the undiluted DNA solution as a PCR solution.

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Buffers G2, QBT, QC and QF^(Note 7).

3.2.6 Extraction procedure: DNA extraction using QIAGEN Genomic-tip 20/G

- (1) Measure an appropriate amount of the ground sample into a 50 ml tube, add 7.5 ml Buffer G2 using a 1,000-5,000 µl micropipette, and vigorously mix the solution on a vortex tube mixer.
- (2) Add another 7.5 ml of Buffer G2 using a 1,000-5,000 µl micropipette, 200 µl QIAGEN Proteinase K using a 100-1,000 µl micropipette, and 20 µl RNase A using a 10-100 µl micropipette to the tube. After mixing the contents of the tube by repeatedly turning the tube upside down until no visible sample remains on the bottom of the tube, mix the solution on a vortex tube mixer.
- (3) Incubate the solution for one hour in a 50°C water bath. Every 15 minutes during incubation or 15, 30 and 45 minutes after incubation started, mix the solution by vigorously inverting the tube several times and then mix it at maximum speed on a vortex tube mixer.
- (4) Spin the solution at 3,000 x g for 15 minutes at 4°C in a swing-out rotor.
- (5) Transfer the entire supernatant into either a new 15 ml or 50 ml tube using a 1,000-5,000 µl micropipette, carefully avoiding the precipitate on the bottom of the tube or film on the liquid surface.
- (6) Flash centrifuge the tube.
- (7) Equilibrate QIAGEN Genomic-tip 20/G with 1 ml Buffer QBT.
- (8) Load the entire supernatant onto the QIAGEN Genomic-tip 20/G in 2 ml aliquots using either a 100-1,000 or 1,000-5,000 µl micropipette, carefully avoiding the precipitate on the bottom of the tube. Allow the supernatant to move through the QIAGEN Genomic-tip 20/G by gravity flow.
- (9) Wash the column by loading 2 ml Buffer QC onto the tip using either a 100-1,000 or 1,000-5,000 µl micropipette and allowing the buffer to move through the QIAGEN Genomic-tip 20/G by gravity flow.

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3.3.5 Reagents

See Part 4—Preparation of Reagents of the Manual.

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After incubating at 60°C for 30 minutes, centrifuge at 14,000 rpm for three minutes^(Note 12).

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Remove the upper layer into a new tube. Pipette carefully to avoid drawing the intermediate layer.

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Add 100 µl TE and 2 µl RNase A (10 mg/mL) to dissolve DNA.

After incubating the solution for 30 minutes at room temperature or 37°C, add 400 µl CTAB eluent.

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Primer pairs to screen for GM maize:

- Endogenous *SSI1b* gene (Nippon Gene Cat # 315-05441, FASMAC Cat # M1-1M, or their bulk products);
- CaMV 35S promoter (Nippon Gene Cat # 317-05521, FASMAC Cat # C1-1M, or their bulk products); and
- GA21-specific sequence (Nippon Gene Cat # 312-05451, FASMAC Cat # M2-1M, or their bulk products).

Primer pairs to detect specific GM maize lines:

- Endogenous *SSI1b* gene (Nippon Gene Cat # 315-05441, FASMAC Cat # M1-1M, or their bulk products);
- GA21-specific sequence (Nippon Gene Cat # 312-05451, FASMAC Cat # M2-1M, or their bulk products);
- Bt11-specific sequence (Nippon Gene Cat # 319-05461, FASMAC Cat # M3-1M, or their bulk products);
- Event 176-specific sequence (Nippon Gene Cat # 316-05471, FASMAC Cat # M4-1M, or their bulk products); and
- MON810-specific sequence (Nippon Gene Cat # 310-05491, FASMAC Cat # M6-1M, or their bulk products).

Primer pairs to detect GM potatoes:

- Endogenous *Pss* gene (Nippon Gene Cat # 315-05701, FASMAC Cat # G5-1, their bulk products, or self-prepared);
- NewLeaf-specific sequence (Note that no commercial product of primer pairs seems to be available at present that can be used to detect this GM line. For the time being, the GM line may be detected by combination of the primer pairs for the CaMV 35S promoter and the positive control plasmid for either Roundup Resistant soybeans or maize. However, this method should be used with utmost care because the target sequence of PCR is also present in plants infected with cytomegalovirus (CMV) and many GM crops.);
- NewLeaf Plus-specific sequence (Nippon Gene Cat # 312-05711, FASMAC Cat # G1-I, their bulk products, or self-prepared);

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Table 4.1(1)
Primer pairs for soybeans*

Target gene	ID	Size	Amplified region	Remarks
Endogenous <i>Le1</i> gene	Le1-n02	118 bp	<i>Le1</i>	
Introduced <i>P-35S</i> gene	P35S-1	101 bp	<i>P-35S</i>	The results should be evaluated carefully because this primer pair also binds to the genes of non-GM plants infected with CMV.
Introduced <i>NOS-ter</i> gene	NOS ter-2	151 bp	<i>NOS-ter</i>	The results should be evaluated carefully because this primer pair also primes when the soil bacteria <i>Agrobacterium tumefaciens</i> is present.
Introduced RRS-specific gene	RRS-01	121 bp	<i>cp4 epsps</i> from <i>P. hybrida</i>	

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Table 4.1(2)
Primer pairs for maize**

Target gene	ID	Size	Amplified region	Remarks
Endogenous <i>SSIIb</i> gene	SSIIb	151 bp	<i>zSSIIb</i>	
Introduced <i>P-35S</i> gene	P35S-1	101 bp	<i>P-35S</i>	The results should be evaluated carefully because this primer pair also binds to the genes of non-GM plants infected with CMV.
Introduced <i>NOS-ter</i> gene	NOS-ter-2	151 bp	<i>NOS-ter</i>	The results should be evaluated carefully because this primer pair also primes when the soil bacteria <i>Agrobacterium tumefaciens</i> is present.
Introduced Event 176-specific gene	E176-2	100 bp	<i>cryIA(b)</i> - PEPC intron #9	
Introduced Bt11 gene	Bt11-3	127 bp	<i>adh1-1S</i> - <i>cryIA(b)</i>	
Introduced GA21 gene	GA21-3	133 bp	OTP - <i>m-epsps</i>	
Introduced T25 gene	T25-1	149 bp	<i>Pat</i> - 35S-ter	
Introduced MON810 gene	M810-2	113 bp	<i>hsp70</i> - <i>cryIA(b)</i>	

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Results file: A:\WORK_RES

Method name: A:\DEFAULT

Assay type: 260/230 & 260/280 Ratios & Conc. Background corr: [Yes] 320.0nm
Sampling device: None Pathlength: 1.0000 cm
Read average time: 0.50 sec Conc. factor: 50.000 at 260.0nm ↓ ↑

Sample ID	Net Abs 230.0nm	Net Abs 260.0nm	Net Abs 280.0nm	260.0/230.0	260.0/280.0	Dil Fact.	Conc.
1	0.0114	0.0011	0.0002	0.09683	6.49987	1.0000	0.05528
1	0.3621	0.4042	0.2297	1.11609	1.75984	50.000	1010.404
2	0.2320	0.2203	0.1258	0.94994	1.75111	50.000	550.8543
X	0.0996	0.1264	0.0707	1.26964	1.78810	50.000	316.1241
4	0.0900	0.0634	0.0360	0.70443	1.75990	50.000	158.4460
5	0.0177	0.0128	0.0069	0.72319	1.84950	50.000	✓31.94674
6	0.0233	0.0594	0.0331	2.54947	1.79373	50.000	148.3956
7	0.1781	0.2359	0.1291	1.32448	1.82650	50.000	589.6622
8	0.2203	0.2683	0.1488	1.21794	1.80242	50.000	670.6342
9) X	0.0029	0.0014	0.0017	0.46297	0.79400	50.000	✓3.38566
10) 177 X	0.0107	0.0065	0.0042	0.60688	1.52288	50.000	✓16.16356
11	0.0037	0.0042	0.0023	1.13267	1.82375	50.000	✓10.56899
12	0.0104	0.0075	0.0046	0.72348	1.64785	50.000	✓18.86953
13	0.0187	0.0197	0.0112	1.05179	1.75654	50.000	49.29562
14	0.0039	0.0063	0.0036	1.62071	1.74584	50.000	✓15.86726
15	0.0026	0.0250	0.0145	9.64179	1.72192	50.000	62.47851
16	0.0461	0.0937	0.0513	2.03400	1.82588	50.000	234.2766
17	0.0378	0.0680	0.0374	1.80035	1.81962	50.000	170.0397
19						50.000	

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BECKMAN DU-7000

Date: 03/08/11
Time: 09:54

DNA/Oligo Quant

ReadSamples RawData Method SaveClear Print Quit

Results file: A:\WORK_RES Method name: A:\DEFAULT

Assay type: 260/230 & 260/280 Ratios & Conc. Background corr: [Yes] 320.0nm
Sampling device: None Pathlength: 1.0000 cm
Read average time: 0.50 sec Conc. factor: 50.000 at 260.0nm ↓ ↑

Sample ID	Net Abs 230.0nm	Net Abs 260.0nm	Net Abs 280.0nm	260.0/230.0	260.0/280.0	Dil Fact.	Conc.
1	-0.0022	0.0004	0.0006	-0.19669	0.70187	1.0000	0.0220
② M9	0.0993	0.2581	0.1417	2.59907	1.82174	50.000	645.306
③ M10	0.0762	0.1973	0.1092	2.58765	1.80594	50.000	493.190
4 未知 DNA	0.1063	0.1492	0.0828	1.40360	1.80140	50.000	372.895
5 M125	0.1681	0.3790	0.2113	2.25469	1.79349	50.000	947.425
6 RRS	0.0593	0.1288	0.0705	2.17387	1.82817	50.000	322.036
7						50.000	

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<2. Delivering of Master Mix>

- # prepare 0.2 mL tubes for each sample, positive control, and negative control#1(no DNA).
- # Deliver 22.5 μL of Master Mix B to each 0.2 mL tube.

<3. Adding template DNA into 0.2mL tubes>

- # Add 2.5 μL of each sample into the 0.2 mL tube.
- # Mix by tapping the bottom of tube, then spin down.

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[preparation of samples]

Dilute extracted DNA with sterilized ultra pure water to a concentration of 10ng/μL

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Scheme for Developing Protocol “GM Quantification Analysis”

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


M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M

M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

**Study on the Qualitative PCR to
Detect the Presence of GMO in
Animal-derived and
Soybean-derived Products**

Hsieh, Pei-Chun

2003. 8. 4-25

Food Microbiology Division

Bureau of Food and Drug Analysis

Department of Health, Executive Yuan,

R.O.C

A. Materials

ID	Sample name	Information
1	Chicken nugget	Market
2	Baby ham	Market
3	Sausage (H)	Market
4	White sausage	Market
5	Sausage (A)	Prima Ham Company
6	Sliced ham	Prima Ham Company
7	FrankFruite sausage	Prima Ham Company
8	Ham (RH)	Prima Ham Company
9	Ham (CT)	Prima Ham Company
10	Ham (10%)	Prima Ham Company
11	Ham (15%)	Prima Ham Company
12	Ham (20%)	Prima Ham Company
13	德昌炭燒角豆乾	Market
14	德昌五香豆乾	Market
15	德昌大溪五香豆乾	Market
16	RRS powder	Positive control
17	Non-GM soybean powder	Negative control

All samples in this study content soy protein

B. Methods

1. DNA Extraction

- a. QIAgen Dneasy Maxi Kit
- b. QIAgen Genomic-tip 20/G
- c. CTAB
- d. QIAgen Dneasy Maxi Kit plus CTAB (final concentration 1%)

2. PCR Reaction

- 95°C for 10 min
- 95°C for 30 sec
- 60°C for 30 sec
- 72°C for 30 sec
- 72°C for 7 min
- 4°C forever

3. Electrophoresis

C. Results


		DNA extraction	Maxi		Genomic-tip		CTAB		Maxi+CTAB	
		PCR date	2003/8/18		2003/8/19		2003/8/20		2003/8/22	
ID	Well.	Sample name	Le1	RRS	Le1	RRS	Le1	RRS	Le1	RRS
1	1	Chicken nugget 1	+	+	+	+	-	+	+	+
	2	Chicken nugget 2	+	+	+	+	+	-	+	C&+
2	3	Baby ham 1	+	-	+	+	+	-	+	-
	4	Baby ham 2	+	-	+	-	+	-	+	-
3	5	Sausage (H) 1	+	+	+	+	+	-	+	-
	6	Sausage (H) 2	+	-	+	+	+	+	+	-
4	7	White sausage 1	-	-	-	-	-	-	-	D
	8	White sausage 2	-	-	-	-	-	-	-	-
5	9	Sausage (A) 1	+	-	+	+	+	-	+	-
	10	Sausage (A) 2	+	-	+	-	+	-	+	-
6	11	Sliced ham 1	+	-	+	+	+	-	+	-
	12	Sliced ham 2	+	-	+	+	+	-	+	-
7	13	FrankFruite sausage 1	+	-	+	+	+	+	+	+
	14	FrankFruite sausage 2	+	-	+	-	+	-	+	+
8	15	Ham (RH) 1	+	-	+	-	+	+	+	+
	16	Ham (RH) 2	+	-	+	+	+	+	+	+
9	17	Ham (CT) 1	A	-	A&B	-	A	-	A	-
	18	Ham (CT) 2	A	-	A	-	A	-	A&B	-
10	19	Ham (10%) 1	A	-	A&+	-	-	-	A	-
	20	Ham (10%) 2	A	-	-	-	-	-	A	-
11	21	Ham (15%) 1	-	-	A&B	+	-	-	A	-
	22	Ham (15%) 2	-	-	A	-	-	-	A&B	-
12	23	Ham (20%) 1	-	+	A	-	-	-	A	-
	24	Ham (20%) 2	-	-	A	-	-	-	A	-
13	25	德昌炭燒角豆乾 1	+	+	+	+	+	+	+	+
	26	德昌炭燒角豆乾 2	+	+	+	+	+	+	+	+
14	27	德昌五香豆乾 1	+	+	+	+	+	+	+	+
	28	德昌五香豆乾 2	+	+	+	+	+	+	+	+
15	29	德昌大溪五香豆乾 1	+	+	+	+	+	+	+	+
	30	德昌大溪五香豆乾 2	+	+	+	+	+	+	+	+
16	31	GM (RRS) 1	+	+	+	+	+	+	+	+
	32	GM (RRS) 2	+	+	+	+	+	+	+	+
17	33	nonGM 1	+	-	+	-	+	-	+	-
	34	nonGM 2	+	-	+	-	+	-	+	-
	35	negative(no DNA)	-	-	-	-	-	-	-	-
	36	negative(no primer)	-	-	-	-	-	-	-	-
	37	positive(Plasmid)	+	+	+	+	+	+	+	+

non-specific band length: A 160 bp, B 1 kb, C 170 bp and D 90 bp

1. Results of Qualitative PCR (DNA Extracted by QIAgen Dneasy Maxi Kit)


1.1 LE1 primer

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

1.2 RRS primer


M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

2. Results of Qualitative PCR (DNA Extracted by QIAgen Genomic-tip 20/G)


2.1 LE1 primer

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

2.2 RRS Primer

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

3. Results of Qualitative PCR (DNA Extracted by CTAB)


3.1 LE1 primer

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

3.2 RRS primer


M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

4. Results of Qualitative PCR (DNA Extracted by QIAgen Dneasy Maxi Kit plus CTAB (final concentration 1%))


4.1 LE1 primer

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

4.2 RRS primer

M 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 M M 21 22 23 24 25 26 27 28 29 30 31 32 33 34 35 36 37 M

[Reagents]

- # Sterilized ultra pure water
- # 10 x PCR Buffer II (Applied Biosystems)
- # dNTPs (Applied Biosystems)
- # MgCl₂ (Applied Biosystems)
- # Ampli Taq Gold (5U/μL) (Applied Biosystems)
- # Forward Primer (50 μM)
- # Reverse Primer (50 μM)

[Preparation of pre mix]

Prepare pre mix solution for AmpliTaq Gold according to the following table.

	Initial conc.	Pre mix conc.	Final conc.	Volume
PCR Buffer II	10 fold	3.85 fold	1 fold	900.9 μL
dNTP	2 mM	0.77 mM	0.2 mM	900.9 μL
MgCl ₂	25 mM	5.775 mM	1.5 mM	64.1 μL
total				2346 μL

[preparation of samples]

Dilute extracted DNA with sterilized ultra pure water to a concentration of 10ng/μL.

[Preparation of reaction solution]

<1. Preparation of Master Mix>

Prepare two 1.5ml tube. Mix the reagents as following to prepare the Master Mix A and B:

Master Mix A	Final conc.	for 1 tube	for 77 tubes
Sterilized ultra pure water		15.375 μL	1183.875 μL
Pre Mix (3.85 fold)	1 fold	6.5 μL	500.5 μL
Ampli Taq Gold (5U/μL)	0.025 U/μL	0.125 μL	9.625 μL
total		22 μL	1694 μL

Master Mix B	Final conc.	for 1 tube	for 37 tubes
Master Mix A		22 μL	814 μL
Le1-5' (50 μM)	0.5 μM	0.25 μL	9.25 μL
Le1-3' (50 μM)	0.5 μM	0.25 μL	9.25 μL
total		22.5 μL	832.5 μL

Master Mix B	Final conc.	for 1 tube	for 37 tubes
Master Mix A		22 μL	814 μL
RRS-5' (50 μM)	0.5 μM	0.25 μL	9.25 μL
RRS-3' (50 μM)	0.5 μM	0.25 μL	9.25 μL
total		22.5 μL	832.5 μL

After Master Mix is prepared, vortex at maximum speed for about three seconds and then spin down.

<2. Delivering of Master Mix>

prepare 0.2 mL tubes for each sample, positive control, and negative control#1(no DNA).

Deliver 22.5 μL of Master Mix B to each 0.2 mL tube.

<3. Adding template DNA into 0.2mL tubes>

Add 2.5 μL of each sample into the 0.2 mL tube.

Mix by tapping the bottom of tube, then spin down.

[Preparation of reaction solution for negative control#2(no primer)]

Prepare one 0.2 ml tube. Mix the reagents as following:

	for 1 tube
Master Mix A	22 μL
Sterilized ultra pure water	0.5 μL
Template DNA (positive control)	2.5 μL
total	25 μL

Mix by tapping the bottom of tube, then spin down.

[PCR Cycles]

Perform amplification according to the following PCR step-cycle program.

<Cycle>	<time>	
95°C	10 min	} 40cycles
95°C	30 sec	
60°C	30 sec	
72°C	30 sec	
72°C	7 min	
4°C	for ever	


No.	Sample name
CTAB 9	1 Ham (CT) 1
CTAB 9	2 Ham (CT) 2
CTAB 10	3 Ham (10%) 1
CTAB 10	4 Ham (10%) 2
CTAB 11	5 Ham (15%) 1
CTAB 11	6 Ham (15%) 2
CTAB 12	7 Ham (20%) 1
CTAB 12	8 Ham (20%) 2
Maxi 9	9 Ham (CT) 1
Maxi 9	10 Ham (CT) 2
Maxi 10	11 Ham (10%) 1
Maxi 10	12 Ham (10%) 2
Maxi 11	13 Ham (15%) 1
Maxi 11	14 Ham (15%) 2
Maxi 12	15 Ham (20%) 1
Maxi 12	16 Ham (20%) 2
M+C 9	17 Ham (CT) 1
M+C 9	18 Ham (CT) 2
M+C 10	19 Ham (10%) 1
M+C 10	20 Ham (10%) 2
M+C 11	21 Ham (15%) 1
M+C 11	22 Ham (15%) 2
M+C 12	23 Ham (20%) 1
M+C 12	24 Ham (20%) 2
G-tip 9	25 Ham (CT) 1
G-tip 9	26 Ham (CT) 2
G-tip 10	27 Ham (10%) 1
G-tip 10	28 Ham (10%) 2
G-tip 11	29 Ham (15%) 1
G-tip 11	30 Ham (15%) 2
G-tip 12	31 Ham (20%) 1
G-tip 12	32 Ham (20%) 2
	33 negative (no DNA)
	34 negative (no DNA)
	35 positive (hamid)

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

[Preparation of pre mix]

Prepare pre mix solution for AmpliTaq Gold according to the following table.

	Initial conc.	Pre mix conc.	Final conc.	Volume
PCR Buffer II	10 fold	3.85 fold	1 fold	900.9 μL
dNTP	2 mM	0.77 mM	0.2 mM	900.9
MgCl ₂	25 mM	5.775 mM	1.5 mM	541
total				2340 μL

[preparation of samples]

Dilute extracted DNA with sterilized ultra pure water to a concentration of 10ng/μL.

[Preparation of reaction solution]

<1. Preparation of Master Mix>

Prepare two 1.5ml tube. Mix the reagents as following to prepare the Master Mix A and B:

Master Mix A	Final conc.	for 1 tube	for 70 tubes
Sterilized ultra pure water		15.375 μL	1076.25 μL
Pre Mix (3.85 fold)	1 fold	6.5 μL	455 μL
Ampli Taq Gold (5U/μL)	0.025 U/μL	0.125 μL	8.75 μL
total		22 μL	1540 μL

Master Mix B	Final conc.	for 1 tube	for 33 tubes
Master Mix A		22 μL	726 μL
Le1-5' (50 μM)	0.5 μM	0.25 μL	8.25 μL
Le1-3' (50 μM)	0.5 μM	0.25 μL	8.25 μL
total		22.5 μL	742.5 μL

Master Mix B	Final conc.	for 1 tube	for 33 tubes
Master Mix A		22 μL	726 μL
RRS-5' (50 μM)	0.5 μM	0.25 μL	8.25 μL
RRS-3' (50 μM)	0.5 μM	0.25 μL	8.25 μL
total		22.5 μL	742.5 μL

After Master Mix is prepared, vortex at maximum speed for about three seconds and then spin down.

<2. Delivering of Master Mix>

prepare 0.2 mL tubes for each sample, positive control, and negative control#1(no DNA).

Deliver 22.5 μL of Master Mix B to each 0.2 mL tube.

<3. Adding template DNA into 0.2mL tubes>

Add 2.5 μL of each sample into the 0.2 mL tube.

Mix by tapping the bottom of tube, then spin down.

[Preparation of reaction solution for negative control#2(no primer)]

Prepare one 0.2 ml tube. Mix the reagents as following:

for negative control#2(no primer)		for 1 tube
Master Mix A		22 μL
Sterilized ultra pure water		0.5 μL
Template DNA (positive control)		2.5 μL
total		25 μL

Mix by tapping the bottom of tube, then spin down.

[PCR Cycles]

Perform amplification according to the following PCR step-cycle program.

<Cycle>	<time>	
95°C	10 min	} 40cycles
95°C	30 sec	
60°C	30 sec	
72°C	30 sec	
72°C	7 min	
4°C	for ever	

No.	Sample name	
CTAB 9	1 Ham (CT) 1	Korea
CTAB 9	2 Ham (CT) 2	Korea
CTAB 10	3 Ham (10%) 1	Korea
CTAB 10	4 Ham (10%) 2	Korea
CTAB 11	5 Ham (15%) 1	Korea
CTAB 11	6 Ham (15%) 2	Korea
CTAB 12	7 Ham (20%) 1	Korea
CTAB 12	8 Ham (20%) 2	Korea
x Maxi 9	9 Ham (CT) 1	Korea
x Maxi 9	10 Ham (CT) 2	Korea
Maxi 10	11 Ham (10%) 1	Korea
Maxi 10	12 Ham (10%) 2	Korea
Maxi 11	13 Ham (15%) 1	Korea
Maxi 11	14 Ham (15%) 2	Korea
Maxi 12	15 Ham (20%) 1	Korea
x Maxi 12	16 Ham (20%) 2	Korea
G-tip 9	17 Ham (CT) 1	Korea
G-tip 9	18 Ham (CT) 2	Korea
G-tip 10	19 Ham (10%) 1	Korea
G-tip 10	20 Ham (10%) 2	Korea
G-tip 11	21 Ham (15%) 1	Korea
G-tip 11	22 Ham (15%) 2	Korea
G-tip 12	23 Ham (20%) 1	Korea
G-tip 12	24 Ham (20%) 2	Korea
CTAB 1	25 Chichen nugget 1	mine
CTAB 1	26 Chichen nugget 2	mine
	27 negative(no DNA)	
	28 negative(no primer)	
	29 positive(Plasmid)	

NOTICE!

Plant Quarantine

- Prevents introduction of injurious agricultural pests and diseases from abroad
- Protects agriculture and green resources of Japan. ■

Your understanding and cooperation is solicited. ■

When you take plants or plant products to Japan, you must undergo plant quarantine,

Here are some reminders for you.

1. Prohibited items

In order to prevent the invasion of quarantine pests that may cause significant damage to agricultural products in Japan, host plants of quarantine pests, quarantine pests, and soil are prohibited to entry.

2. Plants subject to quarantine inspection

This covers practically all plants other than the import prohibited items. They must be inspected and confirmed for the freedom from quarantine pests and the conformity with quarantine requirements of Japan.

When you take plants or plant products out of Japan,

Many countries undertake plant quarantine and regulate import of plants by imposing prohibitions, restrictions, disinfection treatments and/or inspections, etc. You may consult with the Plant Quarantine Office for detailed information before your departure.

Plant Quarantine Office (TEL No.)

Yokohama Plant Protection Station 045-211-7152

Sapporo Sub-Station 011-852-1808

Shin-Chitose Airport Annex 0123-24-6154

Shiogama Sub-Station 022-362-6916

Sendai Airport Annex 022-383-4585

Niigata Sub-Station 025-244-4401

Narita Sub-Station 0476-32-6694

Haneda Airport Branch 0476-34-2352

Tokyo Sub-Station 03-5756-0229

Nagoya Plant Protection Station 03-3599-1139

Nagoya Plant Protection Station 052-651-0112

Nagoya Airport Branch 0568-28-0510

Fushiki Sub-Station 0766-44-0990

Shimizu Sub-Station 0543-52-3775

Kobe Plant Protection Station 078-331-1350

Kansai Airport Sub-Station 0724-55-1936

Osaka Sub-Station 06-6571-0801

Hiroshima Sub-Station 082-251-9593

Hiroshima Airport Annex 0848-86-8261

Sakaide Sub-Station 0877-46-4108

Moji Plant Protection Station 093-321-2601

Fukuoka Sub-Station 092-291-2504

Fukuoka Airport Branch 092-477-7575

Kagoshima Sub-Station 099-222-1046


Kagoshima Airport Branch 0995-58-2428

Naze Sub-Station 0997-52-0459

Naha Plant Protection Station 098-868-2850

Naha Airport Branch 098-857-0054

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FROM PLANT QUARANTINE — JAPAN —

Plant Protection Station
Ministry of Agriculture, Forestry and Fisheries JAPAN
<http://www.pps.go.jp>

If you carry any plants, please come to the Plant Quarantine Office before customs inspection.

- Seeds
- Bulbs and tubers
- Living plants
- Seedlings & scions
- Cutflowers
- Fresh fruits
- Fresh vegetable fruits
- Vegetables
- Grains
- Beans
- Spices
- Herbs
- Dried flowers
- Others (Soil, live insect or mushroom injurious to plants)

Arrival Passengers

Immigration


Plant Quarantine
Take plant inspection before the customs inspection.

Customs

Entry

IMPORT PROHIBITED ITEMS
No one is allowed to bring these items into Japan.
This restriction also applies to purchases from duty free shops and items imported from a third country.

ITEMS SUBJECT TO IMPORT INSPECTION
All plants other than the import prohibited must be inspected for the freedom from quarantine pests upon entry.

- ① Depending on the country of origin, underground portions of ginger, taro, greater yam, coffee, cactus, etc. and seeds for planting of watermelon, corn, pea, etc. are not allowed entry depending on the country unless accompanied with phytosanitary certificate stating that they have been inspected in the field of cultivation in exporting country.
- ② Flower bulb, fruit nursery stocks (including scions), flower bulb, potato, sweet potato, sugarcane, etc. must undergo a certain period of post-entry quarantine (about, one year) for the inspection of virus diseases.
- ③ Narcissus bulbs must be treated by hot-water.

These are typical examples.
Please call the **Plant Quarantine Office** for further information. (See reverse side for contact)
Some plants are prohibited from import or export by the Convention on International Trade in Endangered Species of Wild Fauna and Flora (Washington Convention). For particular information, call the Agriculture and Marine Product Office, Trade Control Department, Ministry of Economy, Trade and Industry (METI 03-3501-1511) or the Customs Office.

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R.O.C

A. Materials

ID	Sample name	Soy protein	Information
1	Chicken nugget	Yes	Market (Nippon Ham Company)
2	Chibini ham	Yes	Market (Marudai Foods Company)
3	Sausage (H)	Yes	Market (Maruha Company)
4	White sausage	No	Market (Maruha Company)
5	Sausage (A)	Yes	Market (Prima Ham Company)
6	Sliced ham	Yes	Market (Prima Ham Company)
7	FrankFruite sausage	Yes	Market (Prima Ham Company)
8	Ham (RH)	Yes	Market (Prima Ham Company)
9	Ham (CT)	No	Prima Ham Company
10	Ham (10%)	No	Prima Ham Company
11	Ham (15%)	No	Prima Ham Company
12	Ham (20%)	No	Prima Ham Company
13	徳昌炭角焼豆乾	Yes	Market (徳昌)
14	徳昌五香豆乾	Yes	Market (徳昌)
15	徳昌大溪五香豆乾	Yes	Market (徳昌)
16	RRS powder	Yes	Positive control
17	Non-GM soybean powder	Yes	Negative control

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A. Materials

ID	Sample name	Soy protein	Information
1	Chicken nugget	Yes	Market (Nissin Food Industry)
2	Chibini ham	Yes	Market (Nissin Food Industry)
3	Sausage (H)	Yes	Market
4	White sausage	No	Market
5	Sausage (A)	Yes	Market
6	Sliced ham	Yes	Market
7	Frankfurter sausage	Yes	Market
8	Ham (RH)	Yes	Market

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