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<p><u>1. Imaged Slice Position Accuracy and Table Indexing</u></p> <p>Laser alignment (error < 2mm)</p> <p>Localization from Scout Images (error < 2mm)</p> <p>Table Index (error <= 1mm)</p>
<p><u>2. Radiation (Dose) Profiles</u></p> <p>Radiation profiles</p>
<p><u>3. Gantry Angle Accuracy</u></p> <p>Gantry Angle</p>
<p><u>4. Film and Monitor Image Quality</u></p> <p>SMPTE printed to film</p> <p>SMPTE displayed at Operator console</p>
<p><u>5. kVp Accuracy, Waveform, and Half Value Layer</u></p> <p>A. kVp accuracy</p> <p>B. Waveform shape acceptable?</p> <p>C. Half Value Layer acceptable?</p>
<p><u>6. Radiation Safety (Texas Regulations for Control of Radiation)</u></p>
<p><u>7. CT Number Accuracy, Noise, Image Uniformity, and Artifacts</u></p> <p>7a. GE's QC Acquisitions (CT#, noise, uniformity, resolution)</p> <p>7b. CT Number of water, Noise, Uniformity, and Artifacts</p> <p>7c. CT number (water) with Image Thickness</p>
<p><u>8. Image Thickness (Slice Sensitivity Profile) :</u></p> <p>8a. Axial Image Thickness</p> <p>8b. Helical Slice Sensitivity Profiles</p>
<p><u>9. CT Number Scale Accuracy & Linearity</u></p> <p>CatPhan、GE Phantom</p>
<p><u>10. High Contrast Resolution</u></p> <p>GE QC Phantom</p> <p>Catphan with QC protocol</p>
<p><u>11. Low Contrast Detectability</u></p> <p>GE QC Phantom</p> <p>Low Contrast evaluation with Catphan rods</p> <p>Visibility of spheres with clinical protocols</p>
<p><u>10. Dosimetry of Axial Scans</u></p> <p>A. Head Mode (dose phantom in head holder) :</p> <p>B. Reproducibility, mA linearity, scan speed linearity</p> <p>C. Body Mode (dose phantom on table) :</p>

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<p><u>1. exam was sent to Stentor</u></p> <p>Slice Position Accuracy</p> <p>Table Indexing</p>
<p><u>2. CT Number Accuracy, Noise, Image Uniformity and Artifacts</u></p> <p>2a. CT number and noise</p> <p>2b. CT number with kVp</p> <p>2c. CT# Uniformity and Artifact</p> <p>2d. CT number with Img Thickness</p>
<p><u>3. Image Thickness (Slice Sensitivity Profile) :</u></p> <p>3a. Axial Image Thickness</p> <p>3b. Helical Slice Sensitivity Profile (optional)</p>
<p><u>4. CT Number Scale Accuracy</u></p> <p>CatPhan</p> <p>GE Phantom</p>
<p><u>5. High Contrast Resolution</u></p> <p>GE QC Phantom</p> <p>Catphan with clinical protocols</p> <p>Catphan with QC protocol</p>
<p><u>6. Low Contrast Detectability</u></p> <p>GE QC Phantom</p> <p>Catphan low contrast and statistical method</p> <p>Visibility of spheres with clinical protocols</p>
<p><u>7. Dose Profiles - Axial Mode</u></p> <p>Radiation profiles</p>
<p><u>8. Gantry Angle Accuracy</u></p> <p>Gantry Angle</p>
<p><u>9. Film and Monitor Image Quality</u></p> <p>SMPTE printed to film</p> <p>SMPTE displayed at Operator console</p> <p>SMPTE displayed at Operator console</p>
<p><u>10. Dosimetry of Axial Scans</u></p> <p>Head Mode :</p> <p>Body Mode :</p>

11. kVp Accuracy, Waveform, and Half Value Layer

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- Office of Communication, Education and Radiation Programs
- They apply to
 - Facilities AND
 - Accrediting bodies

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Attachment 1: ACR Mammography Accreditation Process

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MEDICAL PHYSICIST'S MAMMOGRAPHY QC TEST SUMMARY

Full-Field Digital – General Electric

Site Name: _____ Report Date: _____
 Address: _____ Survey Date: _____
 Medical Physicist's Name: _____ Signature: _____
 X-Ray Unit Manufacturer: General Electric Model: _____
 Date of Installation: _____ Room ID: _____

QC Manual Version: (check one; must use version applicable to unit tested; contact mfr if questions) 2000D 2371472-100 Rev 0, 2003
 DS 5133453-2-100 Rev 1, 2006 ESSENTIAL 5141465-2-100 Rev 1, 2006 OTHER (write in): _____

Accessory Equipment:	Manufacturer	Model	Location	QC Manual Version
Review Workstation*			<input type="checkbox"/> On-site <input type="checkbox"/> Off-site	
Laser Film Printer*			<input type="checkbox"/> On-site <input type="checkbox"/> Off-site	

*FDA recommends that only monitors and printers specifically cleared for FFDM use by FDA's Office of Device Evaluation (ODE) be used, but the use of others is also legal. See FDA's Policy Guidance Help System Modification Document #9 (page 27).

Survey Type: Mammo Eqpt Evaluation of new unit (include MQSA Rqmts for Mammo Eqpt checklist) Annual Survey

Medical Physicist's QC Tests

	Fibers	Specks	Masses	PASS/FAIL
1. Flat Field				
2. Phantom Image Quality				
Phantom IQ Test on AWS				
Phantom IQ Test on Printer				
15. Review Workstation (RWS) Tests* (for all RWS, even if located offsite)				
Overall Results ("Pass" means all tests pass; indicate "Fail" if any test fails)				

*FDA requires that all RWS comply with a QC program that is substantially the same as that recommended by the image receptor manufacturer. If the RWS has been approved by the FDA's ODE for FFDM, the FDA considers the RWS's QC manual to be "substantially the same" and you may follow it for QC. (Check with the RWS manufacturer for their system's FDA clearance status and their QC manual.) If the RWS has not been approved by the FDA's ODE for FFDM, you must follow the QC manual provided by the image receptor manufacturer. (In this case, check with the image receptor manufacturer for their required tests.)

*** YOUR MEDICAL PHYSICIST MUST SUMMARIZE HIS/HER RESULTS ON THIS FORM ***

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 FDA Policy Guidance Help System #9 – Required Medical Physicist Testing	
Equipment Evaluations on FFDM Components	
Item/Repair	MP Involvement
Bucky & detector replacement	MP must evaluate in person
Bucky (but not detector) replacement	MP must oversee
Any detector replacement or repair	MP must evaluate in person
Software modifications	MP must evaluate in person (some alternative standards otherwise)
Monitor (display) or printer replacement	Must follow FFDM mfr's QC manual

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FFDM Mfr	Model	FFDM Mfr's Printer QC Instructions
GE	2000D, DS, Essential	Follow the laser printer mfr's QC
Fischer	SenoScan	Follow the laser printer mfr's QC
Lorad	Selenia	Follow the Lorad Selenia QC Manual
Siemens	Mammomat Novation DR	Follow the laser printer mfr's QC (but conduct QC every day you print)

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Download Summary Forms from ACR Website (www.acr.org)

MEDICAL PHYSICIST'S CHECKLIST
MQSA REQUIREMENTS FOR MAMMOGRAPHY EQUIPMENT

Facility Name: _____

Unit Manufacturer: _____ Model: _____

Medical Physicist: _____ Room #: _____

Signature: _____ Date: _____

Feature	FDA Rule Section	Requirement	Applies to	Meets FDA Requirements? (if NA, please explain)
Motion of tube-image receptor assembly	3(i)	The assembly shall be capable of being fixed in any position where it is designed to operate. Once fixed in any such position, it shall not undergo unintended motion.	S-F & FFDM	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> NA
	3(ii)	This mechanism shall not fail in the event of power interruption.	S-F & FFDM	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> NA
Image receptor sizes	4(i)	Systems using screen-film image receptors shall provide, at a minimum, for operation with image receptors of 18 x 24 cm and 24 x 30 cm.	S-F	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> NA
	4(ii)	Systems using screen-film image receptors shall be equipped with moving grids matched to all image receptor sizes provided.	S-F	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> NA
	4(iii)	Systems used for magnification procedures shall be capable of operation with the grid removed from between the source and image receptor.	S-F & FFDM	<input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> NA

- Required for Equipment Evaluation report
- Addresses 900.12(b) of the FDA regulations
- Same for S-F and FFDM
- In Excel format

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	# Units*	Overall	
		Pass	Deficient
Screen-Film	14,574	88.7%	11.3%
FFDM	1711	92.7%	7.3%

**1st attempt for both initials and renewals; 2/15/03 – 7/12/06*

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Nucletron
DO NOT REMOVE THIS DOCUMENT FROM THE TRANSPORT DRUM
Source Return Document (USA) for Container type ADIC (Part No. 081.040-xx)

Receiver name: <i>Stamps</i>		Address: <i>200 2nd St NW</i>	
City: <i>Stamps</i>	State/Country: <i>PA/US</i>		
Container serial no. (on top of lid): <i>2104 05</i>	Source serial no. (on source lid): <i>1320 2102</i>		
Altitude serial no.: <i>1057</i>			
Source type: <input checked="" type="checkbox"/> mPOR V2 (GENIE 50 Beam) <input checked="" type="checkbox"/> mPOR V2 (S Beam)	<input checked="" type="checkbox"/> mPOR Classic (1.1mtr) <input type="checkbox"/> mPOR Classic (1.1mtr)		

1. Are you Nucletron certified for source exchanges?	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
2. Dose rate above the source channel (without lid) should be lower than on the side of the container.	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
3. Source is locked (lever down) and source cable is undamaged.	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
4. Length of source cable extending from unsealed tube	<input checked="" type="checkbox"/> ≤ 480mm (0.9 mm source)	<input type="checkbox"/> ≤ 480mm (1.1 mm source)
5. Source information label on container matches actual source	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
6. Container lid is closed, both nuts fastened and lid locked	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
7. Surface dose rate of the hot spot on top of the transport drum	<i>1.8 mSv/h</i>	
8. Surface dose rate of the hot spot on the side of the transport drum	<i>9.0 mSv/h</i>	
9. Dose rate for the hot spot at 1 metre from the transport drum	<i>0.3 mSv/h</i>	
10. Radioactive contamination on transport drum	<input checked="" type="checkbox"/> ≤ 0.45µCi/cm²	<input type="checkbox"/> Pass <input type="checkbox"/> Remark
11. Measurements within specification?	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
12. Transport Index (TI)	<input checked="" type="checkbox"/> 0.15TI ≤ 0.8	<i>0.3 mSv/h</i>
13. Source activity		<i>1184</i> GBq
14. Labels are correctly placed and filed in (if Yellow and address labels)	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
15. Shipping documents are completed and signed	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
16. No other reportable quality problems with source, container(s), packaging, labelling, error code 227/228, etc.?	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark
17. RMA label attached? (if yes, submit an FPR report to Nucletron)	<input checked="" type="checkbox"/> No	<input type="checkbox"/> Yes
18. Lid screwed tightly on to transport drum and red tamper inserted, source return documents enclosed (tick 'Pass' if done immediately after completing this form)	<input checked="" type="checkbox"/> Pass	<input type="checkbox"/> Remark

Nucletron Certified Source Handler: <i>M. Smith</i>	Signature: <i>M. Smith</i>	Date (YYYY-MM-DD): <i>2006-07-07</i>	Release for local and long distance shipment? <input type="checkbox"/> Yes (no unresolved remarks) <input checked="" type="checkbox"/> NO: DO NOT SHIP and call Nucletron (1)-800-4433292
-----------------------------------------------------	----------------------------	--------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

For a copy to Nucletron, fax no. (1) 410-3124196, and include the white original in the transport drum.
 Remarks place a red RMA Label on the container and on the transport drum.

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Attn: Donald L. Monroe
Radiation Control Officer

RE: Authorization to use Varian Trilogy accelerator

Dear Mr. Monroe:

We have received your letter indicating the installation of the new accelerator is complete. We have also received the initial shielding survey report performed on June 28, 2006, by Chihray Liu, Ph.D., for the Varian Trilogy accelerator, serial number 1139.

Based on the survey that you have performed, and your attestation that the procedures and operator training are adequate for the new machine, the facility is authorized to use this accelerator effective at the completion of training scheduled for July 12, 2006.

If you have questions please call us at (904) 359-6363, or FAX (904) 359-6362.

Sincerely,

A handwritten signature in black ink, appearing to read "Don K. Steiner".

Don K. Steiner
Environmental Administrator,
Radiation Machine Program

Fax: 352-846-1626

Bureau of Radiation Control • Radiation Machine Program
Post Office Box 210 • Jacksonville, FL 32231

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- 一、 American College of Radiology (ACR) , Mammography Accreditation Program , 1-800-227-6440 。
- 二、 Arkansas Department of Health and Human Services (SAR) , Division of Health , Radiation Section , Mammography Program , 1-501-661-2301 。
- 三、 Mammography Accreditation Program , Iowa Department of Public Health (SIR) , Lucas State Office Bldg. , 5th Floor , 321 East 12th Street , Des Moines, Iowa 50319 , 1-515-281-3478 。
- 四、 Texas Department of State Health Services (STX) , Mammography Accreditation Program , 1-512-834-6688 extension 2246 。

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- 一、 State of Illinois, Office of Radiation Safety, Department of Nuclear Safety, 1035 Outer Park Drive, Springfield, IL 62704, 217-785-9974, Mammography Accreditation Program。
- 二、 Iowa Department of Public Health, Iowa Department of Public Health, Lucas State Office Bldg., 5th Floor, 321 East 12th Street, Des Moines, Iowa 50319, 515-281-3478。
- 三、 State of South Carolina, Mammography Certification Program, Department of Health and Environmental Control, 2600 Bull Street, Columbia, SC 29201, 803-545-4400。

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FDA/CDRH/OCER/DMQRP

Attention: Closed Facility Notification of Records Retention

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- 一、 Mammography Accreditation Program, Iowa Department of Public Health, Lucas State Office Bldg., 5th Floor, 321 East 12th Street, Des Moines, Iowa 50319, 515-281-3478。
- 二、 State of Illinois, Office of Radiation Safety, Department of Nuclear Safety, 1035 Outer Park Drive, Springfield, IL 62704, 217-785-9974。
- 三、 State of South Carolina, Department of Health and Environmental Control, Division of Electronic Products, 2600 Bull Street, Columbia, SC 29201, 803-545-4400。

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
(Mammography Program Reporting and Information System ,MPRIS)

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- 一、 Fuji Computed Radiography Mammography Suite (FCRMS) on 07/10/06。
- 二、 GE Senographe Essential Full Field Digital Mammography (FFDM) System on 04/11/06。
- 三、 Siemens Mammomat Novation DR Full Field Digital Mammography (FFDM) System on 08/20/04
- 四、 GE Senographe DS Full Field Digital Mammography (FFDM) System on 02/19/04
- 五、 Lorad/Hologic Selenia Full Field Digital Mammography (FFDM) System on 10/2/02
- 六、 Lorad Digital Breast Imager Full Field Digital Mammography (FFDM) System on 03/15/02
- 七、 Fischer Imaging SenoScan Full Field Digital Mammography (FFDM) System on 09/25/01
- 八、 GE Senographe 2000D Full Field Digital Mammography System on 01/28/00

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DIGITAL Mammogram

GE FFDM Mammogram, Evaluation Annual Survey - [] : _____
 Evaluation required for new or moved units, major repairs (tube change, etc.)- requires MQSA checklist and check of everything.
 Annual survey requires check of all parts clinically-used.

Site: _____
 Address: _____
 Room: _____
 Manufacturer: _____
 Model : _____
 Serial Number: _____
 Install Date: _____

Physicist:
 Carl R. Keener, Ph.D., DABMP, DABR, MP0431
 D. Lloyd Goff, Ph.D., DABMP, MP0424
 José Antonio Domínguez, M.S., DABMP MP0027
 Janet Feng, Ph.D., TMP6101
 Bahadir Ozus, Ph.D., TMP6150

Texas Registration Number R-06742
 BRC# _____ expires: _____

Accreditation: ACR TDH
 Accreditation# _____ expires: _____

Lead Interpreting Physician:

Technologist:

Comments:
 Any major changes since last survey?

Notes

Pre-Survey data:
 Max kVp: _____ Max mA: _____ HF Other _____
 Target: Mo Rh W
 Focal spots (mm): 0.3 0.35 0.4 | 0.10 0.15
 Filters Mo (____µm) Rh (____µm & ____µm) Al
 Grid Type: _____ lp/cm: _____ Ratio: _____
 SID: _____cm Mag stand height: _____cm
 Receptor Size: _____ 19.2 cm x 23.0 cm
 Matrix: _____ 1920 x 2304
 Serial number: _____

Soft-copy display: _____
 Model number: _____
 Serial number: _____
 Manufacture date: _____

Hard-copy output: _____
 Model number: _____
 Serial number: _____
 Manufacture date: _____

Clinical kVp selection:
 AOP (CNT Dose Standard) _____

GE QC manual version:
 2000D 2277390-100 Rev 1, 2000 2000D 2371472-100 Rev 0, 2003
 Mobile 2371698-100 Rev 0, 2003 Senographe DS 2404641-100 Rev 2, 2004
 Seno Advantage 2391082-100 Rev 1, 2003

9. Viewing Condition Check & Setting:

Sketch RWS room if it has changed. Note all sources of light and their status..

Date:


Unit:

GE FFDM Mammography Survey - 1


Medical & Radiation Physics, Inc. (rev. 3/18/05)

5. Collimation

Check the alignment of the light field, radiation field, receptor and compression paddle for both targets. Compression paddle at 4.2 cm.

Mo

Rh

Width of penny on paddle = _____ mm

Width of penny on platform = _____ mm

$$\Delta light / radiation = edge_{radiation} - edge_{light}$$

$$\Delta radiation / receptor = (length_{whole_penny} - length_{penny_at_edge}) + \Delta light / radiation$$

SID = 66 cm

	Mo Target		Rhodium target		Pass/Fail
	mm	% of SID	mm	% of SID	
Misalignment of radiation field with light field:					
Left side					
Right side					
Total of left & right sides					
Misalignment of radiation field with receptor					
Left side					
Right side					
Chest side					
Anterior side					
Total of chest & anterior sides					
Extension of compression paddle beyond chest wall of receptor					
Paddle edge visible?					

4 & 7. Phantom and Dose

Mode	kVp	target/filter	mAs	mR/mAs	HVL
AOP-CNT					
AOP-STD					
AOP-DOSE					

Date:

Unit:

6. Evaluation of Focal Spot Performance

The system and focal spot resolutions were measured using the large and small focal spots with a high-contrast line pair phantom placed 4.5 cm above the breast support. The minimum acceptable resolution in line pairs / millimeter (lp/mm) is 13 lp/mm parallel to and 11 lp/mm perpendicular to the anode-cathode axis. The digital resolution will be poorer, but should not change between Evaluations.

Digital & film may be shot together if 2 mm Al is under film and beam is collimated. PROTECT THE RECEPTOR.

Nominal focal spot size (mm)	Resolution - bars parallel to anode-cathode axis (lp/mm)	Resolution - bars anode-cathode axis (lp/mm)	Receptor	Target	technique
0.3			digital	Mo	28 kVp, 100 mAs, 2 mm Al
0.1 (1.5x mag)			digital	Mo	28 kVp, 125 mAs, 2 mm Al
0.3			digital	Rh	28 kVp, 100 mAs, 2 mm Al
0.1 (1.5x mag)			digital	Rh	28 kVp, 125 mAs, 2 mm Al
0.3	5	5	digital-last year	Mo	NA
0.1 (1.5x mag)	8	8	digital-last year	Mo	NA
0.3	5	5	digital-last year	Rh	NA
0.1 (1.5x mag)	8	8	digital-last year	Rh	NA
0.3			film only	Mo	28 kVp, 100 mAs, 2 mm Al
0.1 (1.5x mag)			film only	Mo	28 kVp, 125 mAs, 2 mm Al
0.1 (1.8x mag)			film only	Mo	28 kVp, 150 mAs, 2 mm Al
0.3			film only	Rh	28 kVp, 100 mAs, 2 mm Al
0.1 (1.5x mag)			film only	Rh	28 kVp, 125 mAs, 2 mm Al
0.1 (1.8x mag)			film only	Rh	28 kVp, 150 mAs, 2 mm Al
minimum acceptable	13	11			

8. Artifact Evaluation & Flat Field Uniformity

Visible?	Acceptable?	Description	Pass/Fail
		Detector	
		Grid	
		Mirror/tube	
		Filter (Mo)	
		Filter (Rh)	
		Display	

Secondary Exposure Levels / Exposure to Members of the Public *Evaluation only. Use MAP phantom for scatter. PROTECT THE DETECTOR.*

- Invision 451P survey meter, SN: 0788
- Invision 451P survey meter, SN: 6020
- Cardinal Health 451P survey meter, SN: 1104
- Invision 451P survey meter, SN: 6332
-
-

kVp (target/filter) = _____ mAs per image = _____
 mA = _____ exams per week = _____

Location	Beam Direction	Exposure Level	
		mR/hr	mR/year
Control Panel	down	_____	_____
Outside barrier edge	down	_____	_____
Outside closed door	down	_____	_____
Outside wall	down	_____	_____
Background reading		_____	_____

Date: _____ Unit: _____ GE FFDM Mammography Survey - 3 Medical & Radiation Physics, Inc. (rev. 3/18/05)

Evaluation of Site's Technologist GE FFDM Quality Program

Technologist(s): _____ (list certification)

<u>QC Requirement</u>	<u>Pass / Needs Improvement / Fail</u>
1. Monitor cleanliness & viewing conditions (daily).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
2. Viewing conditions for Review Work Station (RWS) (daily).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
3. Darkroom cleanliness (daily, if applicable).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
4. Processor QC (daily, if applicable).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
5. Mobile unit quality control (after every move, if applicable).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
6. Flat field (weekly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
7. Image quality (phantom) (weekly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
8. Viewbox and viewing conditions (weekly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
9. MTF measurement (monthly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
10. AOP mode & signal-to-noise (SNR) (monthly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
11. Visual checklist (monthly - ACR only).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
12. Display device (monitor) calibration check (monthly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
13. Repeat analysis - performed, records maintained, reviewed by radiologist (quarterly).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
14. Analysis of fixer retention (quarterly, if applicable).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
15. Compression force (pressure) (semi-annually).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
16. Darkroom fog (semi-annually, if applicable).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>
17. Laser film printer QC (according to printer manufacturer recommendations).....	<input type="checkbox"/> / <input type="checkbox"/> / <input type="checkbox"/>

Comments:

Last MQSA survey?

Date: _____ Unit: _____ GE FFDM Mammography Survey - 4 Medical & Radiation Physics, Inc. (rev. 3/18/05)