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Easy Workflow & Reliability

Universal Solution for the Existing X-ray Room

The AeroDR detector is the same as an ISO 4090 compliant film cassette in size so that it will fit any existing wall-stand or table bucky tray.

Shared FPD Solution

AeroDR can be used anywhere with "the Shared FPD Solution". As soon as AeroDR is registered to any X-ray room, AeroDR will be ready to use in the X-ray room immediately.

Integrated Control Station CS-7

CS-7 can control not only the AeroDR detectors but also X-ray generators and Konica Minolta existing CR family. No need to operate the X-ray console to adjust X-ray exposure conditions.**

**5 Please contact your Konica Minolta sales representative for more details regarding connections between a X-ray generator and CS-7.

Quick Preview and Smart GUI


After exposure, a preview image immediately appears on the display of the new CS-7 console in less than two seconds. The CS-7 has a user-friendly graphic interface adding new and powerful proprietary functions. GUI design can be modified to customer preferences flexibly, succeeding the conventional console design.

● Screen Sequence

Power-saving Technology

Patient safety is of primary importance, therefore the lithium ion capacitor, the world newest technology, was adopted as a battery technology which has many advantages despite of demanding a lower power consuming panel design, which has been overcome by employing low power ICs and a power-saving control.

New Battery Technology Achieves Light-weight yet Rigid Body

The lithium ion capacitor has a charge and discharge cycle life that is tremendously longer than that of the lithium ion battery and does not markedly decrease in capacity even after it has continuously been used for many years. Therefore, it is possible to be built in to AeroDR and also friendly to the environment. In this case, the structure of the cassette case has become so simple that it is possible to significantly reduce the weight of the cassette and increase the mechanical strength of the cassette.

Reliable, Rapidly Rechargeable and Long-Life Battery

The lithium ion capacitor, which charges quickly in a battery charger or through a tethered connection, has a long charge and discharge cycle life that does not need to be replaced during the expected life cycle of the detector. If the capacitor gets exhausted in emergency, AeroDR gets over 10 images by the capacitor being recharged for only three minutes.

■ Characteristics

Battery expected life time	Above the AeroDR detector product life time
Battery charging time empty to full	Within 30 minutes (when using the AeroDR Battery Charger) Within 60 minutes (when using the dedicated wired cable)
Number of exposable images	14"x17" : 200 images / 5.5 hours 17"x17" : 173 images / 4.8 hours <small>*Under conditions that the interval between studies is five minutes and three images are captured in each study, assuming 30 seconds for each exposure to position a patient.</small>

AeroDR 的規格

WIRELESS DIGITAL RADIOGRAPHY SYSTEM AeroDR	
Model name	AeroDR P-11 (AeroDR 1417HQ) / AeroDR P-12 (AeroDR 1417S) / AeroDR P-21 (AeroDR 1717HQ)
Detection method	Indirect conversion method
Scintillator	CsI (Cesium Iodide)
Dimensions (W x D x H)	14"x17" : 363,7 x 460,2 x 15,9 mm / 17"x17" : 459,8 x 460,2 x 15,9 mm
Weight	14"x17" : 2,9 kg / 17"x17" : 3,6 kg
Pixel size	175µm
Image area size	14"x17" : 348,95x425,25mm (1,994 x2,430 pixels) / 17"x17" : 424,9x424,9mm (2,428 x2,428 pixels)
AD conversion	16 bit (65,536 gradients)
Communication	Dedicated wired ethernet connection / Wireless LAN (IEEE802.11a compliant)
WLAN encryption	Wireless encryption method : AES / Authentication method : WPA2-PSK
Cycle time	14"x17" : Approx. 8,2 seconds when connected with the dedicated wired connection / Approx. 13,3 seconds when connected with wireless LAN connection 17"x17" : Approx. 9,9 seconds when connected with the dedicated wired connection / Approx. 15,5 seconds when connected with wireless LAN connection
Dynamic range	4 digits
Battery charging time empty to full	Within 30 minutes (when using the AeroDR battery charger) Within 60 minutes (when using the dedicated wired cable)
Battery duration in standby status	Approx. 16 hours
Number of exposable images	14"x17" : 200 Images / 6,5 hours 17"x17" : 173 Images / 4,8 hours * Under conditions that the interval between studies is five minutes and three images are captured in each study, assuming 20 seconds for each exposure to position a patient.
Battery expected life time	Above the AeroDR detector product life time
Recommended storage and usage environment conditions	When operating : 10 to 30°C / 35 to 80% RH(ensure no water condensation) When not operating : -10 to 40°C / 20 to 90% RH(ensure no water condensation) In storage : -20 to 60°C* / 20 to 90% RH(ensure no water condensation) * However, performance warranty period when storing at 60°C is 6 months after packing.

- The described performance may change depending on the environment and frequency of use. (This is not a guarantee of performance.)
- The performance of battery is all performance after fully charged.

AeroDR Battery Charger	
Power	100 / 110 / 115 / 120 / 200 / 230 / 240 VAC
Weight	7,2 kg
Dimensions (W x D x H)	560 x 250 x 153mm

Control Station CS-7		
Image processing	Automatic gradation processing (G processing) , Frequency processing (F processing) , Equalization processing (E processing) , Hybrid processing (H processing) , Hybrid smooth processing (HS processing)	
Image output	Host : Up to 4 channels / Printer : Up to 2 channels	
DICOM support	Basic grayscale print management (SCU) , Storage (SCU) , Modality worklist management , Modality performed procedure step , Grayscale standard display function (Print output)	
CR / DR connections	AeroDR : Up to 4 simultaneous active detectors REGIUS cassette reader : Up to 15 units	
Main Options	Hardware options	Bar-code reader for REGIUS cassette registration , In-room sub monitor
	Software options	DICOM MWM / MPPS / DETACHED , FTP , DICOM storage output , DICOM print , X-ray generator connection , Media storage , Text annotation , Automatic body part recognition , Automatic processing parameter study

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CR CsBr(Eu)讀取機

A speedy and easy-to-operate X-ray system performing a wide range of exams with less radiation exposure.

The Fujifilm DR VELOCITY Unity fp enables radiographic exams in various positions including supine and upright. The X-ray tube and the detector both on the U-arm can be perfectly aligned, and the detector can be easily tilted for exams of angled parts such as the knee and skull. The unit has an exceptionally high image processing capability, and is fully motorized ensuring rapid set-up at exam. With the use of the Focused Phosphor Technology, excellent image quality can be obtained with less exposure of radiation.

UNCOMPROMISED IMAGE QUALITY

By a combination of **Focused Phosphor Technology** applied to the built-in detector and FUJIFILM's renowned sophisticated digital image software technologies, the **Image Intelligence™**, the FUJIFILM DR VELOCITY Unity fp offers unparalleled image quality in digital radiology.

Focused Phosphor Technology

With the use of the Focused Phosphor Technology, the X-ray exposure efficiency has been improved without changing the exposure dose.

To increase the X-ray exposure efficiency, we had to improve the Imaging Plate X-ray absorption efficiency significantly. To improve the Imaging Plate X-ray absorption efficiency, we had to increase the thickness of the phosphor layer substantially.

In practice, the Imaging Plate X-ray absorption efficiency cannot be improved by merely increasing the phosphor layer thickness for the following reasons:

- 1 sufficient stimulation light cannot penetrate deep into the phosphor layer because each phosphor particle in the phosphor layer acts as a light scattering medium, and
- 2 photo stimulated luminescence (PSL), generated deep within the phosphor layer and containing X-ray information cannot be collected on the Imaging Plate surface.

The new Focused Phosphor Technology uses a thicker phosphor plate than previous models but since the phosphor crystals of the plate are realigned in a columnar structure, it is able to permeate the stimulation light deep into the phosphor layer. The generated PSL is then successfully collected on the focused phosphor plate surface. As a result, the X-ray exposure efficiency is improved.


Image Intelligence™

"Image Intelligence™" is a set of sophisticated digital image processing software technologies incorporated in the Fujifilm DR VELOCITY Unity fp. An optimum image for examination is generated with this Image Intelligence technology.

FNC Flexible Noise Control

FNC selectively suppresses noise components while maintaining signal contrast, improving granularity in "noisy" anatomical regions.

MFP Multiple Frequency Processing

MFP is optional software that provides greater diagnostic information from a single exposure image through frequency enhancement. MFP improves visibility of both dense and peripheral tissues, simultaneously applying edge-enhancement processing to all structures in an image.

FUJIFILM DR VELOCITY Unity fp SPECIFICATIONS

Image Detector (Model: CR-IR 372 RU)

• Image size	43 x 43 cm (17 x 17 inch)
• Pixels	4280 x 4280
• Pixel size	100 microns
• Spatial frequency (Nyquist)	5.0 cy/mm
• Bit Depth	12 bit
• Preview Image	Approx. 9 seconds

Universal Arm Stand (Model: VERSO F)

• Vertical height	max. 1530 mm or more min. 500 mm
• Rotation of arm	+120/-30 degree
• Tilt of detector	+/-45 degree*
• SID	1000 - 1800 mm

*The absolute angle of detector is limited from 0 to +110 degree by software.

Generator (Model: SHF515 / SHF535 / SHF635 / SHF835)

• Output	50kW (SHF515, SHF535)/64kW (SHF635)/80kW (SHF835) (SHF835 available in some countries)
• Range of output	40 to 150kVp in 1kVp steps.
• Maximum current	640mA (SHF515, SHF535, SHF635)/800mA (SHF835)
• Console (Model: TPC 12")	12 inch Monitor and Touch screen display

X-Ray Tube (Model: E7252X / E7869X)

• Range of voltage	40 - 150kV
• Focal spot	1.2 mm (Large focus)/0.6 mm (Small focus)
• Anode heat storage capacity	300 kHU (E7252X)/600 kHU (E7869X)

Collimator (Model: 150PBL Collimator)

• Automatic collimator

Image & Information Processor (Model: CR-IR 348CL): Option

- Spatial frequency processing
- Gradation processing
- Dynamic Range compression processing
- Multi-frequency processing (option)
- Flexible noise control processing (option)
- Grid pattern removal processing (option)

Dimensions* (W x D x H)

• Universal arm stand**	2320 x 1465 x 2650 mm (91" x 58" x 104")
• Arm controller	592 x 422 x 600 mm (23" x 17" x 24")
• Generator	592 x 422 x 690 mm (23" x 17" x 27")

*This is approx. value. ** Detector is included

Weight*

• Universal arm stand**	500 kg (1102 lbs.)
• Arm controller	65 kg (143 lbs.)
• Generator	110 kg (243 lbs.)

*This is approx. value. ** Detector is included

Power Supply Requirement

• Detector	200 - 240V±10%, 2.1 2.5A, 50/60Hz (Single phase)
• Universal arm stand	230/240V±10%, 50/60Hz (Single phase)
• Generator	230/240V, 50/60Hz (Single phase) for output 50kW 400V, 50/60Hz (3-phase) for output 50kW 400V, 50/60Hz (3-phase) for output 64kW 480V, 50/60Hz (3-phase) for output 80kW

Environment Condition

• Temperature	15 - 30 degree
• Humidity	30 - 75%Rh
• Atmospheric pressure	700 - 1060 hPa

Optional Accessories


- Mobile table Carbon fiber table (Patient Capacity: 200 kg)
- Grid 10:1 or 8:1, 36 lines/cm, FID 100/140/180 cm
- Ion chamber (Model: ICX-127) 3 field
- Infrared remote controller
- Ceiling support
- Cassette holder
- Long view cassette holder
- QC phantom holder

Specifications are subject to change without notice.
All brand names or trademarks are the property of their respective owners.
All products require the regulatory approval of the importing country.
For details on their availability, contact our local representative.

FUJIFILM DR VELOCITY Unity fp (Model: CR-IR372)  0123

Dimensions

Unit mm (in.):


FUJIFILM

Dual-Side Reading Technology

The Dual-Side Imaging Plate (IP) Reading technology allows the use of a thicker phosphor layer on the transparent base, thereby increasing DQE (Detective Quantum Efficiency) by collecting the emissions from both sides of the IP.

Enhanced Image Processing

"Image Intelligence™" – a set of sophisticated digital image-processing software technologies available through the CR Console – processes image data and optimizes final output.

MFP Multi-frequency Processing*

As an optional software applicable for all types of FCR imaging, MFP is an improved version which uses frequency enhancement to provide more diagnostic data from a single exposure image, using Fujifilm's renowned Dynamic Range Control (DRC). MFP improves visibility of both dense and peripheral tissue by simultaneously applying edge enhancement processing to small and large structures within an image.

PEM Pattern Enhancement Processing for Mammography*

As an optional software specifically developed for mammographic imaging, PEM enhancement processing improves the conspicuity of micro-calcifications.

* Image processing requires the use of FCR CR Console Plus

Optimized Operation

The MAMMOASCENT AWS-c optimizes your mammography examinations with its advanced design and features.

Digital Mammography System

A digital Mammography System is a system that links the FCR PROTECT ONE and/or the FCR PROTECT CS with the SYNAPSE Mammography Viewer via the MAMMOASCENT AWS-c and/or the CR Console to optimize the viewing of any body part that may be associated with breast cancer.

FUJIFILM FCR PROTECT ONE Specifications

Standard Components:

- FCR PROTECT ONE Image Reader (Model: CR-IR 368)
- AC power cord

Other System Components (sold separately):

- CR Console Plus, MAMMOASCENT AWS-c
- Image Recorder : DRYPIX 4000/7000
- ID Card Writer
- FCR Data Management System

Supplies:

Imaging Plate:

- ST-VI (Standard): 14" x 17", 14" x 14", 10" x 12", 8" x 10", 24 x 30cm, 18 x 24cm
- HR-V (High Resolution): 24 x 30cm, 18 x 24cm
- ST-BD (Standard Dual-Side Imaging): 24 x 30cm, 18 x 24cm
- HR-BD (Dual-Side Mammography): 24 x 30cm, 18 x 24cm

IP Cassette:

- Type CC: 14" x 17", 14" x 14", 10" x 12", 8" x 10", 24 x 30cm, 18 x 24cm
- Type CH: 24 x 30cm, 18 x 24cm
- Type DS: 24 x 30cm, 18 x 24cm
- Type DM: 24 x 30cm, 18 x 24cm
- Type LC: 35.4 x 124.5cm, 35.4 x 101.7cm, 35.4 x 83.0cm, 25.2 x 58.0cm, 24.0 x 57.0cm

Time Required for IP Feed/Load:

IP Type	IP Size	Required Time
ST-VI	14" x 17" (35 x 43cm)	Approx. 60 sec.
ST-VI	14" x 14" (35 x 35cm)	Approx. 54 sec.
ST-VI	10" x 12"	Approx. 50 sec.
ST-VI	8" x 10"	Approx. 40 sec.
ST-VI	24 x 30cm	Approx. 51 sec.
ST-VI	18 x 24cm	Approx. 42 sec.
HR-V	24 x 30cm	Approx. 65 sec.
HR-V	18 x 24cm	Approx. 55 sec.
ST-BD	24 x 30cm	Approx. 85 sec.
ST-BD	18 x 24cm	Approx. 75 sec.
HR-BD	24 x 30cm	Approx. 90 sec.
HR-BD	18 x 24cm	Approx. 80 sec.

Processing Capacity

(in high-pixel density two-image output format):

IP Type	IP Size	When connected to DRYPIX 7000/CR Console Plus
ST-BD	24 x 30 cm	Approx. 42 IP/hr.
ST-BD	18 x 24 cm	Approx. 48 IP/hr.
HR-BD	24 x 30 cm	Approx. 40 IP/hr.
HR-BD	18 x 24 cm	Approx. 45 IP/hr.

Time to Print on DRYPIX 7000 (18 x 24 HR-BD & ST-BD) :
Approx. 140 sec.

Image Reading (Image output is via CR Console)

Reading Size	Standard Pixel-density		High Pixel-density	
	Pixels/mm	Number of Pixels	Pixels/mm	Number of Pixels
14" x 17" (35 x 43 cm)	5	1760 x 2140	10	3520 x 4280
14" x 14" (35 x 35 cm)	5	1760 x 1760	10	3520 x 3520
10" x 12"	6.7	1670 x 2010	10	2505 x 3015
8" x 10"	10	2000 x 2510	10	2000 x 2510
24 x 30cm (ST-VI & HR-V)	6.7	1576 x 1976	10	2364 x 2964
18 x 24cm (ST-VI & HR-V)	10	1770 x 2370	10	1770 x 2370
24 x 30cm (ST-BD & HR-BD)	10	2364 x 2964	20	4728 x 5928
18 x 24cm (ST-BD & HR-BD)	10	1770 x 2370	20	3540 x 4740

Number of Stackers: 1

Reading Gray Scale: 12 bits

Network: 10 Base T/100 Base TX

Dimensions (W x D x H): 655 x 740 x 1330mm (26" x 29" x 52")

Weight: 230kg (507 lbs.)

Power Supply Conditions:


Single phase 50-60Hz
AC120-240V ±10%
7A (max)

Environmental Conditions:

- Operating Conditions:
Temperature: 15-30°C
Humidity: 40-80%RH (No dew condensation)
- Non-operating Conditions:
Temperature: 0-45°C
Humidity: 10-90%RH (No dew condensation)

Dimensions

Unit: mm (in.)

IP Cassette with Imaging Plate

Type DM

DM Cassette with IP HR-BD for Dual-Side Mammography

Type DS

DS Cassette with IP ST-BD for Standard Dual-Side imaging

Mammography QC Program


Fujifilm's Mammography QC Program is a dedicated quality control program applicable to FUJIFILM digital mammography systems. This program enables the system to keep a stable image quality for both screening and diagnosis.

This equipment is a Class 1 laser product (0EC60825).

Specifications are subject to change without notice.
All brand names or trademarks are the property of their respective owners.
In some countries, regulatory approval may be required to import medical devices.
For the availability of these products, please contact your local sales representatives.

FUJIFILM

FUJIFILM Corporation

26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8520, JAPAN
<http://www.fujifilm.com/products/medical/>

Ref. No. XB-972ER (8K-10-03-F1120-F9711) Printed in Japan ©2009 FUJIFILM Corporation

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
Features

The Amulet is available with either an 18x24cm or 24x30cm detector. The Amulet's unique, and patented, detector technology allows for an extremely high, 50 μm , resolution image to be produced. The unit's use of "Optical Switching Technology" allows the Amulet to have an extremely low inter-exposure time of less than 15 seconds. The newly developed acquisition workstation (AWS) is based on Fujifilm's intuitive console technology and controls both image processing and X-ray control. The excellent ergonomics of the Amulet have been developed with the help of many Radiographers and patients.

Highly precise image depicting details of microcalcifications


FUJIFILM utilizes revolutionary selenium production techniques to manufacture its new flat panel detector

FUJIFILM innovative Direct Optical Switching technology

The state-of-the-art direct-conversion flat panel detector is comprised of a dual layer of a-Se, and employs the Direct Optical Switching technology instead of thin film transistors (TFTs) to directly read the image signal more efficiently with less noise. Therefore Direct Optical Switching achieves higher DQE even with a pixel size of 50 μm .


New X-ray sensor diagram

FUJIFILM exclusive amorphous selenium flat panel detector

Through vast experience in device development and vacuum deposition technology FUJIFILM developed the flat panel detector using highly purified a-Se. In addition, clinical efficiency is enhanced by employing FUJIFILM's new and proprietary Direct Optical Switching technology.

Cross-section of the sensor, and its imaging process

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The BodyTom™ was designed with patient and staff safety in mind. The core system comes standard with internal lead shielding to help safe guard against unnecessary scatter exposure. Additional protection is provided by a built-in lead shield mounted to the BodyTom™ imaging workstation.

More Standard Safety Features

System Control

- Key Lockout Feature
- Individual User Log-In
- Preset Scanning Protocols
- QA Phantom Test & Report
- Meets ACR Accreditation Guidelines
- Patient Intercom System

Dose Check

- Alert Message
- Notification Message
- Dose Displayed Prior to Scan
- Audit Trail
- Dose Reporting


Radiation Safety

- Internal Lead Shielding
- Protocols Based on Age/Weight

Dimensions

BodyTom™

Bore Diameter: 85 cm | (33.5 in)
Weight: 1587.6 kg | (3510 lbs)

BodyTom™ Overview

Portability

- Easy to Use, Built-in Drive System
- Internal Lead Shielding
- Battery Powered
- Plugs into a Standard Wall Outlet
- Wide-Angle Drive Camera
- Portable Imaging Station

Seamless Integration

- DICOM 3.1 Compliant
- Interfaces with PACS
- Wireless Communication
- Health Information System (HIS)
- Radiology Information System (RIS)
- Integrating the Healthcare Enterprise (IHE)
- Interfaces with Surgical Navigation Systems
- Dose Check

Imaging Station Technical Information

- 27" LCD | Monitor Resolution: 2560 x 1440
- Genuine Windows® 7 Professional 64 Bit
- Processor: Intel Quad Core™ Xeon Processor 2.67 GHZ
- Memory: 4 GB
- Operating System Hard Drive: 80 GB
- Secondary Storage: 1 TB
- Mounted Lead Glass Shield Option
- DVD | CD Burner
- MP3 Dock with Surround Sound System

Imaging

- 32 Slice x 1.25 mm = 4 cm Aperture
- 85 cm Gantry
- 60 cm FOV
- 1.25 mm, 2.5 mm, 5.0 mm, 10 mm Slice Thickness
- 512 x 512 Image Matrix
- Maximum Scan Length of 2 Meters
- Scout Scanning
- CT Angiography
- CT Perfusion*
- 2D, 3D, and MPR Viewing Tools

Scanning Modes

- Axial
- Helical
- Dynamic

Training

- 2 Weeks of Onsite Training Led by a NeuroLogica Clinical Applications Specialist
- All Training is CEU Accredited
- Biomedical Engineer First Responder Training Available

1 Year All Inclusive Warranty

- Covers All Parts and Labor (X-Ray Tube Included)
- Telephone Support (Within 1 Hour) 24 Hours per Day, 7 Days per Week
- On-site Personnel and Parts Within 24 Hours
- 100% Warranty on All Parts
- Required Routine Preventative Maintenance
- Software Updates (Safety or Specification)

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