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Oxford NANOPORE technologies

STEP : BASECALLING WITH GUPPY

Now we will use Guppy, our standard basecaller package, to basecall a Lambda dataset. But first we will need to download some raw data from Dropbox:

- Let's make a new directory for our Lambda data and download it:

```
ubuntu@ip-10-11-12-14:/home/ubuntu$ mkdir lambda
ubuntu@ip-10-11-12-14:/home/ubuntu$ cd lambda
ubuntu@ip-10-11-12-14:/home/ubuntu/lambda$
wget https://www.dropbox.com/s/t81ej91351d6fo/lambda.tar.gz?dl=0 -O lambda.tar.gz
```
- Let's extract this tar file:


```
ubuntu@ip-10-11-12-14:/home/ubuntu/lambda$ tar -xf lambda.tar.gz
```

This should unzip and give us 2 individual multi-read fast5 files. To double-check that the extraction completed successfully, list the contents of the directory:

```
ubuntu@ip-10-11-12-14:/home/ubuntu/lambda$ ls
FAL06504_4dd5f26403973742429ebeec317fec22caaed8da_0.fast5
FAL06504_4dd5f26403973742429ebeec317fec22caaed8da_1.fast5
```

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Overview

The workshop will be available at different Oxford Nanopore locations around the world. Oxford Nanopore experts will provide attendees (up to 20) with a nanopore sequencing overview and case study discussion.

A nanopore sequencing demonstration as well as a rapid sequencing control experiment will allow each participant to gain practical hands-on experience of running the MinION/GridION.

The workshop will cover the following:

- ✓ Comprehensive technology overview
- ✓ Library preparation in groups
- ✓ Library loading practice
- ✓ Running the MinION/GridION
- ✓ Introduction to basecalling and EPI2ME analysis workflow

Agenda

- 08:45 - Registration opens
- 09:30 - Introduction to nanopore sequencing
- 10:30 - Planning your experiments
An overview of applications, kits and protocols. Nanopore sequencing experts will also discuss hints & tips on how to maximise your sequencing runs.
- 12:00 - Networking lunch
- 13:00 - Flow cell loading practical
An Oxford Nanopore specialist will demonstrate the correct handling and loading of a flow cell. Participants will then take part in a MinION Flow Cell loading practical.
- 14:30 - Rapid control experiment
An Oxford Nanopore specialist will perform a live Rapid control experiment. Participants will then complete a Rapid barcoding experiment to gain practical hands-on experience of running the MinION/GridION.
- 16:00 - Software introduction
Our technical specialist will provide an introduction to real-time sequencing including basecalling, MinKNOW feedback and EPI2ME analysis workflows.
- 17:00 - Drinks reception
- 17:30 - Workshop closes

Data Analysis Workshop

A one day workshop with Oxford Nanopore experts dedicated to data analyses

Overview

The workshop will be available at different Oxford Nanopore locations around the world. Oxford Nanopore experts will provide attendees (up to 20) with a comprehensive data analysis overview followed by practical hands-on experience running different analysis workflows.

This workshop will give participants experience with real-time analysis tools as well as in-depth analysis workflows with command line tools.

The workshop will cover the following:

- ✓ In depth overview of EPI2ME for real-time end-to-end analysis
- ✓ Command line tutorial
- ✓ Guppy (command, parameters and folder structure)
- ✓ Overview of fast5 file format
- ✓ Quality and read length distribution plot and filtering
- ✓ Mapping to a reference
- ✓ De novo assembly and consensus correction
- ✓ Methylation detection

Agenda

- 08:45 - Registration opens
- 09:15 - Introduction to nanopore sequencing
(optional for Introduction workshop attendees)

Real-time analysis of nanopore reads

- 10:00 - Nanopore file formats
- 10:30 - Introduction to EPI2ME

Analysis workflows with command line tools

- 11:00 - Introduction to practical sessions
- 11:15 - A brief Linux tutorial
- 12:00 - Networking lunch
- 13:00 - Practical session I
Basecalling with Guppy
Generating run metrics and summaries with Nanoplot
Data filtering with Filtlong
Reference mapping using Pomoxis (Minimap2)
Visualizing Alignments with Tablet
- 14:30 - Coffee Break
- 14:40 - Practical session II
De novo assemblies with Pomoxis (Miniasm/Racon)
Error correction with Medaka
Methylation detection
- 16:30 - Workshop closes

Certificate of attendance

This certifies that

Jih-Hui Lin

Has successfully completed

**Introduction Workshop
and
Data Analysis Overview**

Signed:  . *Andy Davies*
Senior Director of Global Services Date: 29-11-2019