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
PROCEDURE FOR PROCESSING CLINICAL SAMPLES FOR LEGIONELLA CULTURE

*Concentrate by centrifugation prior to inoculating media
**Save specimen at 4°C for 5-7 days

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QUALITY CONTROL

Media Testing

Inoculating Methods

The following methods are used to test media for proper reactivity or growth characteristics. These methods are exclusive for media used for isolation and identification of the *Legionellaceae* family. The chart following this section lists media, methods of inoculation, stock isolates, and expected results when the procedures listed below are utilized.

Method 1. Suspend stock isolate from a 1 to 2 day agar culture (prepared from the working stock) into 2.0 ml sterile D.I. water at a 0.5 McFarland turbidity. Streak a .001 ml loopful for colony isolation onto the plate media for testing.

NOTE: *Candida albicans* suspension is a McFarland #3
(25-50 yeast cells per 400x field).

<u>STOCK ISOLATE LIST</u>	<u>VAMC #</u>	<u>ATCC#</u>	<u>STORAGE MEDIUM</u>
<i>Legionella pneumophila</i>	#20	43111	BCYE slant
<i>Legionella micdadei</i>	#594	VAMC	BCYE slant
<i>Candida albicans</i>	#212	14053	TSA slant
<i>Escherichia coli</i>	#210	25922	TSA slant
<i>Staphylococcus aureus</i>	#209	25923	TSA slant

Recording Results

- (1) Media information such as date tested, date prepared, and medium type are logged onto a log sheet in Quality Control log book.
- (2) Results are logged onto the log sheet, as indicated, dated, and initialed. (See Sample Log Sheet.)
- (3) All media not meeting expected results are discarded and remade. Failures are recorded in the Quality Control Log book.

File: Media Testing
G drive: Documents/Stout/Special Pathogens

>

7.Hepatic disease	1.yes 2.no 9.nil	10.Other surgery	1. cardiac 3.lung 4.brain 5.skeletomuscular 6.skin 7. BPH s/p TURB 2. No. 9. Nil 8. Sinusitis s/p op 9. Thyroid 10. Bteast 11. Hemorrhoid 12 Laryngeal Ca		
1. Etiology <i>0721=9221</i>	1. acute hepatitis A 3. Chronic hepatitis B 4. Chronic hepatitis C 5. Hepatitis D superimposed 6. Chronic nonB nonC hepatitis 7.alcoholic hepatitis 8. Hepatic failure 9.nil 10.liver cirrhosis B 11. Liver cirrhosis C 12 alcoholic liver cirrhosis 13 ascites 14 SBP 15 hepatic encephalopathy 16. Hepatorenal syndrome 17 hepatopul syndrome 18.EV bleeding 19 GV bleeding 2 No 20 polycystic liver dz 21 fatty liver 22 splenomegaly 23 hepatomegaly 25 carrier B 26 carrier C 27 Shock liver 28 liver function impairment ,cause 29 LC,cause? 30 spleen abscess 31 LC, non B non C 32 chronic liver dz	1. Years before entry	Years 2.No 9.Nil		
8.GI tract dz	1.yes 2.no 9. nil	11.CV Dz	1.Y	2.N	9.nil
1. Etiology	1. peptic ulcer 2 no 3. GU 4. DU 5 gastritis 6 diverticulitis 7 colitis 5. IBD 9 nil 10 UGI bleeding 11 EU 12 LGI bleeding 13. Angildysplasia 14 ischemia bowel dz 15 abd internal bleeding 16 malabsorption 17 irritable bowel syndrom 18.short bowel syndrome 19.gastrostomy 20 jejunostomy 21 colonostomy 22. GERD 23 Barrett's esophagitis 24 entritis	1. Etiology	1.AMI 3. Old MI 4.CAD 5.unstable angina 6. HTN 7. ISM 8.DCM 9.nil 10 CHF 11. Pul edema 12 pericarditis 2. No 13.myocarditis 14.pericardial effusion 15.vavular heart dz 16. Cardiac arrhythmia 17. AAA 18 Aortic dissection 19PAOD 20.DVT 21dyslipidemia 22. Aortic aneurysm type A 23. Aortic aneurysm type B 24.Mycotic aneurysm 25 IE 26 RHD 27. HCVD		
1. Etiology	1. peptic ulcer 2 no 3. GU 4. DU 5 gastritis 6 diverticulitis 7 colitis 5. IBD 9 nil 10 UGI bleeding 11 EU 12 LGI bleeding 13. Angildysplasia 14 ischemia bowel dz 15 abd internal bleeding 16 malabsorption 17 irritable bowel syndrom 18.short bowel syndrome 19.gastrostomy 20 jejunostomy 21 colonostomy 22. GERD 23 Barrett's esophagitis 24 entritis	12.CNS	1.Y	2.N	9.nil
1. Etiology	1. peptic ulcer 2 no 3. GU 4. DU 5 gastritis 6 diverticulitis 7 colitis 5. IBD 9 nil 10 UGI bleeding 11 EU 12 LGI bleeding 13. Angildysplasia 14 ischemia bowel dz 15 abd internal bleeding 16 malabsorption 17 irritable bowel syndrom 18.short bowel syndrome 19.gastrostomy 20 jejunostomy 21 colonostomy 22. GERD 23 Barrett's esophagitis 24 entritis	1. Etiology	1. CVA 3. Meningitis 4.encephalitis 5.brain abscess 6 ICH 7. IICP 8. EDH 9 nil 10 SDH 11 SAH 12.brain edema 13 hydrocephalus 2.no 14 senile dementia 15. SDH 16 TIA		
9.Neoplasm	1.Y 2.N 9 nil	13. renal dz	1.Y	2.N	9.nil
1.Hematology	1.AML 3.ALL 4.CML 5 CLL 6.Hodgkin's lymphoma 7. Non- Hodgkin's lymphoma 8.MM 10. MDS 99.nil 11..Preleukemia 12. Others 2.No.	1. Etiology	1.ARF 2.no 3. CRT 4. CRF 5.Uremia +H/D 6. Nephrotic syndrom 7. Nephitis 8. GM 9. Nil 10 renal cysts 11 polycystic kidney 12. Renal abscess 13 perinral abscess 14 APN 15 UTI 16 renal stone 17 uretral stone 18 U-GB stone 19. Staghorn stone 20 hydronephosis 21 BPH 22 Rhadomyolysis 23 acute epididymitis 24 prostatitis		
2.Solid tumor .	1. Brain 3.esophageal 4. Gastric 5. Intestine 6 colon 7 HCC 8. GB 10. CBD 99 nil 11. holangiocarcinoma 12. Ut Ca 13. Ovarian Ca 14 Cx Ca 15. Skin 16 Lung Ca 17 Cardiac tumor 18 parotid tumor. 2.No. 19. Du Ca 20 carcinomatosis, abd 20 lung metastasis 21 liver metastasis 22 renal tumor 23 tumor of unknown 24 pancreate Ca 25. Prostate Ca 26. NPC 27. Laryngeal Ca	14.Pul dz	1. yes	2. no	9. nil
1.Hematology	1.AML 3.ALL 4.CML 5 CLL 6.Hodgkin's lymphoma 7. Non- Hodgkin's lymphoma 8.MM 10. MDS 99.nil 11..Preleukemia 12. Others 2.No.	1. etiology	1. COPD 2.no 3. Asthma 4. Community pneumonia 5.nosocomial pnemonia 6. Aspiration pneumonia 7. Necrotizing pneumonia 8. Lung abscess 9 nil 10 viral pnemonitis 11 pep 12 fungal pneumonia 13 active pul TB 14 old pul TB s/p treatment 15. Empyema 16.pleural effusion 17pneumothroax 18 hemothroax 19acute tonsilitis 20.URI 21 ILD 22.resp failure 23.ARDS 24 acute pul edema 25 embolic lung 26 URI		


MEDIA TESTING CHART

MEDIUM	ISOLATE(S)	RESULTS (3D)
BCYE	<i>L. pneumophila</i>	1 mm colony
	<i>L. micdadei</i>	1 mm colony
Selective #2 (PAV)	<i>L. pneumophila</i>	1 mm colony/green
	<i>L. micdadei</i>	1 mm colony/blue
	<i>C. albicans</i>	No growth
	<i>E. coli</i>	0-50 colonies
	<i>S. aureus</i>	No growth
Selective #1 (PAC)	<i>L. pneumophila</i>	1 mm colony
	<i>C. albicans</i>	No growth
	<i>E. coli</i>	0-50 colonies
	<i>S. aureus</i>	No growth
DGVP	<i>L. pneumophila</i>	1 mm colony/green
	<i>L. micdadei</i>	1 mm colony/blue
	<i>E. coli</i>	0-50 colonies
	<i>S. aureus</i>	No growth

7.Hepatic disease	1.yes 2.no 9nil	10.Other surgery	1.cardiac 3.lung 4.brain 5.skeletomuscular 6.skin 7. BPH s/p TURB 2. No. 9. Nil 8. Sinusitis s/p op 9. Thyroid 10. Bteast 11. Hemorrhoid 12 Laryngeal Ca
1. Etiology 0721=9221	1.acute hepatitis A 3.Chronic hepatitis B 4. Chronic hepatitis C 5. Hepatitis D superimposed 6. Chronic nonB nonC hepatitis 7.alcoholic hepatitis 8. Hepatic failure 9.nil 10.liver cirrhosis B 11. Liver cirrhosis C 12 alcoholic liver cirrhosis 13 ascites 14 SBP 15 hepatic encephalopathy 16. Hepatorenal syndrome 17 lepatopul syndrome 18.EV bleeding 19 GV bleeding 2 No 20 polycystic liver dz 21 fatty liver 22 spleenomegaly 23 hepatomegaly 25 carrier B 26 carrier C 27 Shock liver 28 liver function impairment ,cause 29 L.C,cause? 30 spleen abscess 31 L.C, non B non C 32 chronic liver dz	1.Years before entry	Years 2.No 9.Nil
		11.CV Dz	1.Y 2.N 9.nil
		1. Etiology	1.AMI 3. Old MI 4.CAD 5.unstalbe angina 6. HTN 7. ISM 8.DCM 9.nil 10 CHF 11. Pul edema 12 pericarditis 2. No 13.mycocarditis 14.pericardial effusion 15.vavular heart dz 16. Cardiac arrhythmia 17. AAA 18 Aortic dissection 19PAOD 20.DVT 21dyslipidemia 22. Aortic aneurysm type A 23. Aortic aneurysm type B 24.Mycotic aneurysm 25 IE 26 RHD 27. HCVD
		12.CNS	1.Y 2.N 9.nil
8.GI tract dz	1.yes 2.no 9. nil	1. Etiology	1. CVA 3. Meningitis 4.encephalitis 5.brain abscess 6 ICH 7. HCP 8. EDH 9 nil 10 SDH 11 SAH 12.brain edema 13 hydrocephalus 2.no 14 senile dementia 15. SDH 16 TIA
1. Etiology	1. peptic ulcer 2 no 3. GU 4. DU 5 gastritis 6 diverticulitis 7 colitis 5. IBD 9 nil 10 UGI bleeding 11 EU 12 LGI bleeding 13. Angildysplasia 14 ischemia bowel dz 15 abd internal bleeding 16 malabsorption 17 irritable bowel syndrom 18.short bowel syndrome 19.gastrostomy 20 jejunostomy 21 colonostomy 22. GERD 23 Barrett's esophagitis 24 enrtitis	13. renal dz	1.Y 2.N 9.nil
		1. Etiology	1.AKF 2.no 3. CRF 4. CRF 5.Uremia +H/D 6. Nephrotic syndrom 7. Nephitis 8. GM 9, Nil 10 renal cysts 11 polycystic kidney 12. Renal abscess 13 perinral abscess 14 APN 15 UTI 16 renal stone 17 uretral stone 18 U-GB stone 19. Staghorne stone 20 hydronephosis 21 BPH 22 Rhadomyolysis 23 acute epididymitis 24 prostritis
9.Neoplasm	1.Y 2.N 9 nil	14.Pul dz	1. yes 2. no 9. nil
1.Hematology	1.AML 3.ALL 4.CML 5 CLL 6.Hodgkin's lymphoma 7. Non- Hodgkin's lymphoma 8.MM 10. MDS 99.nil 11..Preleukemia 12. Others 2.No.	1. etiology	1. COPD 2.no 3. Asthma 4. Community pneumona 5.nosocomial pnemonia 6. Aspiration pneumonia 7. Necrotizing pneumonia 8. Lung abscess 9 nil 10 viral pnemonitis 11 pcp 12 fungal pneumonia 13 active pul TB 14 old pul TB s/p treatment 15. Empyema 16.pleural effusion 17pneumothroax 18 hemothroax 19acure tonsilitis 20.URI 21 ILD 22.resp failure 23.ARDS 24 acute pul edema 25 embolic lung 26 URI
2.Solid tumor .	1.Brain 3.esophageal 4. Gastric 5. Intestine 6 colon 7 HCC 8. GB 10. CBD 99 nil 11. holangiocarcinoma 12. Ut Ca 13. Ovarian Ca 14 Cx Ca 15. Skin 16 Lung Ca 17 Cardiac tumor 18 parotid tumor. 2.No. 19. Du Ca 20 carcinomatosis, abd 20 lung metastasis 21 liver metastasis 22 renal tumor 23 tumor of unknown 24 pancreate Ca 25. Prostate Ca 26. NPC 27. Laryngeal Ca		

Legionella Identification Review Chart

Examine primary media plates at 3-7 days with dissecting microscope.

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LEGIONELLA IMMUNOFLUORESCENCE

Reagents and Materials (all assays)

- (1) Special Pathogens routinely purchases the following fluorescein-labeled (FA) diagnostic reagents (conjugates).
 - (a) A monoclonal DFA conjugate kit for all serogroups of *L. pneumophila* (Genetic Systems, Seattle, Washington).
 - (b) Six individual polyclonal conjugates for each *L. pneumophila* serogroup (SciMedx, Denville, New Jersey and m-Tech Monoclonal Technologies, Inc., Norcross, GA).
 - (c) A polyvalent polyclonal IFA conjugate kit for 33 *Legionella* strains, including *L. pneumophila* serogroups 1 to 10 (Meridian Diagnostics, Columbus, Ohio). *Remel*

Note: The monoclonal DFA test from Genetic Systems is the primary assay used by Special Pathogens. Ordering details for all these reagents are found in the order card file and the Inventory, Order, and Receipt file.

Definitions

- (1) Monoclonal antibody – very specific antibody that reacts with only the outer membrane of *L. pneumophila*.
- (2) Polyclonal antibody – less specific antibodies that react with multiple sites on *L. pneumophila* or other *Legionella* species.

3. List of 33 *Legionella* serogroup (22 species) Antibodies Contained in the REMEL *LEGIONELLA* Test Kit.

Serogroup	Strain	Serogroup	Strain
1	Bellingham	<i>L. gormanii</i> 1	LS-13
1	Knoxville	<i>L. micdadei</i> 1	TATLOCK
1	Philadelphia	<i>L. wadsworthii</i> 1	81-716
2	Togus 1	<i>L. feeleeii</i> 1	WO-44C-C3
3	Bloomington 2	<i>L. stainthelensi</i>	Mt. St. Helens 4
4	Los Angeles	<i>L. jordanis</i> 1	BL-540
5	Dallas 1	<i>L. anisa</i> 1	WA-316-C3
6	Chicago 2	<i>L. spiritensis</i> 1	Mt. St. Helens 9
7	Chicago 8	<i>L. hackeliae</i> 1	Lansing 2
8	Concord 3	<i>L. maceachernii</i>	PX-1-G2-E2
9	IN-23-G1-C2	<i>L. jamestowniensis</i>	JA-26-G1-E2
10	Leiden	<i>L. cherrii</i> 1	ORW
		<i>L. steigerwaltii</i> 1	SC-18-C9
		<i>L. parisiensis</i>	PF-209C-C2
		<i>L. rubrilucens</i> 1	WA-270A-C2
		<i>L. erythra</i> 1	SE-32A-C8
		<i>L. dumoffii</i> 2	NY-23
		<i>L. longbeachae</i> 1	Long Beach 4
		<i>L. longbeachae</i> 2	Tucker 1
		<i>L. bozemanii</i> 1	WIGA

LEGIONELLA IMMUNOFLUORESCENCE

2. *Legionella* Polyclonal Direct Immunofluorescence (DFA) Test Procedure

Reagents

Test reagents are purchased from SciMedx and m-Tech separately for this assay and not as a kit. The conjugates are immune rabbit sera against *L. pneumophila* or other *Legionella* species listed below.

- (1) *L. pneumophila* individual (monovalent) conjugates:

<i>L. pneumophila</i>	serogroup 1	Philadelphia 1, Knoxville 1, Bellingham strains, OLDA (m-Tech only)
<i>L. pneumophila</i>	serogroup 2	Togus 1 strain
<i>L. pneumophila</i>	serogroup 3	Bloomington 2 strain
<i>L. pneumophila</i>	serogroup 4	Los Angeles strain
<i>L. pneumophila</i>	serogroup 5	Dallas 1E strain
<i>L. pneumophila</i>	serogroup 6	Chicago 2 strain

- (2) SciMedx polyvalent I conjugate that includes all six above serogroups in one reagent.
- (3) SciMedx polyvalent II conjugate that includes serogroups 1-7 plus 6 *Legionella* species.
- (4) *Legionella* species monovalent conjugates for:
L. bozemanii
L. dumoffii
L. gormanii
L. longbeachae 1 & 2 SciMedx only)
L. jordani (SciMedx only)
L. anisa (m-Tech only)

400 Ford Road
Denville, NJ 07834 USA
(201) 625-8822 • 1 (800) 221-5598
FAX: (201) 625-8796 • TELEX: 493-0966 BIO UI

SCIMEDX
CORPORATION

MATERIALS AVAILABLE:

Legionella pneumophila FITC-labeled Rabbit Globulins, polyvalent and serogrouping.

The conjugates are prepared from immune rabbit serums produced against *L. pneumophila* serogroups 1 through 6 and the *legionella* species listed below. There are individual conjugates for each serogroup, and a polyvalent conjugate for each serogroup. The conjugates have been diluted to their optimal working dilutions with 0.25% Evans Blue diluent (which serves as a counterstain) and are distributed as Lyophilized components. The products are available as follows:

Cat. No. 1504	<i>L. pneumophila</i> serogroup 1	0.5 ml/2.0 ml	Cat. No. 150b	<i>L. bozemanii</i>	0.5 ml/2.0 ml
Cat. No. 1505	<i>L. pneumophila</i> serogroup 2	0.5 ml/2.0 ml	Cat. No. 150d	<i>L. dumoffii</i>	0.5 ml/2.0 ml
Cat. No. 1506	<i>L. pneumophila</i> serogroup 3	0.5 ml/2.0 ml	Cat. No. 150g	<i>L. gormanii</i>	0.5 ml/2.0 ml
Cat. No. 1507	<i>L. pneumophila</i> serogroup 4	0.5 ml/2.0 ml	Cat. No. 150lb	<i>L. longbeachae</i>	0.5 ml/2.0 ml
Cat. No. 1508	<i>L. pneumophila</i> serogroup 5	0.5 ml/2.0 ml	Cat. No. 150m	<i>L. micdadei</i>	0.5 ml/2.0 ml
Cat. No. 1509	<i>L. pneumophila</i> serogroup 6	0.5 ml/2.0 ml	Cat. No. 150j	<i>L. jordani</i>	0.5 ml/2.0 ml
Cat. No. 1520	<i>L. pneumophila</i> serogroups 1-6, poly	0.5 ml			
Cat. No. 1521	<i>L. pneumophila</i> serogroups 1-6, poly	2.0 ml			
Cat. No. 1530	<i>L. pneu.</i> 1-7 & b,d,g,l,m,j poly	2.0 ml			

NOTE: 0.5 ml vials provide 15 tests; 2.0 ml vials, 80 tests. Upon receipt, store the conjugates at 4-8°C. The reconstituted conjugates are stable at this temperature until expired. Check labels for specific expiration dates. Once reconstituted, the conjugates may be stored at 4-8°C for up to 90 days. Thimerosal is added as a preservative at 0.01%. The conjugates provided are at their optimum working potency for both sensitivity and specificity.

10/11/20

Fig. 1 An Approach To Disinfection


Fig. 2 Approach To Laboratory Capability For Hospitals.