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附錄(一) Laboratory directions for molecular subtyping of *Salmonella*(*E. coli* O157:H7 and *Shigella*) by Pulsed-Field Gel Electrophoresis(PFGE)

LABORATORY DIRECTIONS FOR MOLECULAR SUBTYPING Of *Salmonella*  
(*E. coli* O157:H7 and *Shigella*) By PULSED-FIELD GEL ELECTROPHORESIS (PFGE)

Monday morning, March 15, 2004

Laboratory Module IA

Preparation of PFGE Agarose Plugs from *Salmonella* Cell Suspensions

**BIOSAFETY WARNING:** These organisms may have low infectious doses and most have caused human illness. Please handle with caution.

Please read laboratory instructions on pages 1 and 2 before continuing. Treat all plasticware, glassware, pipets, spatulas, etc. that come in contact with the cell suspensions or plugs as contaminated materials and dispose of, or disinfect, according to instructions.

Please put your group initial (A, B, C or D) on all tubes, flasks, etc. that you use. You will be given plates of seven different *Salmonella* strains. Cell suspensions will be made in Cell Suspension Buffer<sup>1</sup> (CSB) from 14-16 h blood agar plates. Handle them carefully and with gloves.

1. Label seven Falcon 2054 tubes at the top with the culture numbers and your initial. Label one Falcon 2054 tube with "Blank." Add ~2 ml of Cell Suspension Buffer (CSB) to each tube.
2. Label seven 1.5-ml microcentrifuge tubes with the culture numbers and your initial (1-A<sup>3</sup>).
3. Label one well of disposable and one well of reusable PFGE plug molds with culture numbers.
4. Adjust concentration of the cell suspensions with the Dade MicroScan Turbidity Meter as described in the following steps:
  - a. Check the digital output on the instrument by inserting the CSB "blank" tube in the left position of the turbidity meter and one of the other tubes of CSB in SAMPLE position (right) of the turbidity meter. The reading should be  $0.00 \pm 0.01$ . Hold the tubes at the top.
  - b. Moisten a sterile disposable loop with sterile CSB and remove growth from the agar plates into the appropriate Falcon tube that contains 2 ml of CSB. Suspend the cells by rubbing against the wall of the tube gently so the cells are evenly dispersed and formation of aerosols is minimized.

<sup>1</sup> Formulas for reagents are on page 10.

c. Insert the tube of the cell suspension into the SAMPLE position to check the reading (desired range is 0.48 - 0.52)

d. If the reading is greater than 0.52, add additional CSB, mix well and check reading again until it is within the desired range of 0.48 - 0.52.

e. If the reading is less than 0.48, add additional growth from the plate to the cell suspension until the desired range is obtained.

5. Transfer 400  $\mu$ l (0.4 ml) of each cell suspensions to labeled microcentrifuge tubes using the 1000  $\mu$ l (1 ml) pipet and tip. Place microcentrifuge tubes in plastic holders (float) and incubate in the 37°C water bath for 5 minutes. Keep tubes containing rest of the cell suspensions on ice bath until plugs are made and they are in the shaking water bath.

6. Remove cell suspensions in 1.5 ml tubes from water bath; add 20  $\mu$ l proteinase K (20 mg/ml stock solution) to each 0.4 ml cell suspension using the 100  $\mu$ l pipet; mix by closing tubes and tapping side of tube. (The final concentration of Proteinase K [Invitrogen 25530-049] in the cell suspension is 0.5 mg/ml.)

7. Remove 1% SeaKem Gold: 1% Sodium Dodecyl Sulfate (SDS) agarose from 56°C water bath; keep agarose in beaker of water so that it will stay warm while making plugs.

Note: Do not remove melted agarose from water bath until you are ready to make the plugs; mix the agarose with the cell suspensions immediately before the agarose can cool too much.

8. Add 400  $\mu$ l (0.4 ml) melted 1% SeaKem Gold:1% SDS agarose to one of the 0.4 ml cell suspensions; mix by gently pipetting mixture up and down a few times. The 1000  $\mu$ l pipetman can be used without changing the setting.

9. Immediately, dispense part of mixture into appropriate well in both plug molds. Do not allow bubbles to form. Repeat for remaining samples. Allow plugs to solidify for 10-15 min at room temperature.

Note: When reusable plug molds (2-cm x 1-cm x 1.5-mm) are used, 2 plugs can be made from these amounts of cell suspension and agarose; when disposable plug molds (1.5 mm x 10-mm x 5-mm) are used, 4-5 plugs can be made.

10. Continue with Step 1 of **Laboratory Module IB**.

Laboratory Module IB  
Lysis of Cells in Agarose Plugs

Read instructions in this section before continuing.

1. Label seven 50 ml polypropylene screw-cap tubes with culture numbers and group initial. (One of the tubes contains ~35 ml Cell Lysis Buffer [5 ml x 7 tubes = 35 ml].)
2. Add 175  $\mu$ l Proteinase K stock solution (25  $\mu$ l x 7 tubes = 175  $\mu$ l) to the tube that contains the 35 ml Cell Lysis Buffer using 200  $\mu$ l pipet; recap tube and invert or swirl gently to mix. (The concentration of Proteinase K stock solution is 20 mg/ml. The final concentration in the Cell Lysis Buffer is 0.1 mg/ml)
3. Dispense 5 ml of Proteinase K/cell lysis buffer to each of the other six labeled 50 ml tubes. (5+ ml will be left in the original tube.)
4. If a flat edge is wanted on the plugs, trim excess agarose from top of plug with scalpel. Dispose of scalpel in red SHARPS container.
  - a. **Reusable Plug Molds:** Open mold by partially unscrewing the back and removing it. Put back of mold into plastic container that contains chlorox (5,000 ppm). Transfer plugs from mold with wide end of spatula to appropriately labeled tube. Remove tape from reusable mold and put front section of plug mold and the spatula in container with 75% Isopropyl alcohol.
  - b. **Disposable Plug Molds:** Remove white tape from bottom of molds and use tab from strip to push plugs into appropriate tube. Discard disposable mold, tab and tape in discard pan.

**Note:** Lyse both the large and small plugs of each strain in the same tube.

5. Confirm that plugs are under buffer and not on side of tube.
6. Place tubes in rack in 54°C shaking water bath incubate for 2 hours with constant agitation. Confirm that the water level is above level of lysis buffer in tubes and tubes are moving freely in the rack. Note the time the tubes went into the incubator or water bath.
7. Put bottle containing ~250 ml of sterile Ultrapure (Reagent Grade, Type 1) water into 50°C water bath.

Monday afternoon, March 15, 2004

Laboratory Module IC  
Washing of Agarose Plugs After Cell Lysis  
Read instructions in this section before continuing.

1. Remove tubes with plugs from incubator or water bath; remove caps and replace with green screened caps. Carefully pour off lysis buffer into discard pan. Gently tap bottom of tube against lab bench to force plug to bottom of tube.

**Note:** Touch top of green cap onto an absorbent paper towel so that most of the liquid during this and subsequent washing steps is removed.

2. Add 15 sterile Reagent Grade water that has been pre-heated to 50°C to each tube and screw original cap on top of green screened cap.

3. Confirm that plugs are under water and not on side of tube; return to shaking incubator or water bath (50°C). Shake tubes for 10 minutes.

4. Pour off water into sink and repeat wash steps with pre-heated water (Steps 2 and 3) one additional time.

5. Put one liter bottle of Tris-EDTA Buffer, pH 8.0 (TE) fitted with a dispenser into the 50°C water bath.

6. Pour off water from plugs, add 15 ml pre-warmed (50°C) sterile TE, mix, and shake moderately in shaking incubator or water bath (50°C) for 15 minutes.

7. Pour off TE and repeat wash step three more times (10-15 minutes each).

8. Pour off TE from last (fourth) wash step, add 10 ml TE, and place in rack (provided by instructor) for overnight refrigeration at 4°C. Confirm that plugs are in the TE and not on the side of the tube.

**Laboratory Module IIA**  
**Restriction Digestion of DNA in Agarose Plugs with Xba I**

Tuesday morning (08:45– 09:45am), March 16, 2004

Read instructions in this section before continuing.

1. Label seven 1.5-ml microcentrifuge tubes with culture numbers and group initial; label three 1.5-ml microcentrifuge tubes with H9812 (PulseNet Universal standard strain of *Salmonella* ser. Braenderup) and group initial. Label one Falcon 2057 (17-mm x 100-mm) tube with "H" and group initial.
2. Measure sterile Type 1 water and 10X H Buffer (Roche Molecular Biochemicals) into labeled Falcon 2057 tube according to the following table to make a 1:10 dilution of the buffer. Measure reagents carefully; prepare enough for 11 plug slices. Mix well.

**Note:** Wear gloves when handling 10X buffer. Keep buffer on ice.

Reagent	μl/Plug Slice	μl/11 Plug Slices
Sterile Type 1 Water	180 μl	1980 μl
H Buffer	20 μl	220 μl
<b>Total Volume</b>	<b>200 μl</b>	<b>2200 μl</b>

3. Add 200 μl diluted H Buffer to labeled 1.5-ml microcentrifuge tubes.
4. Carefully remove plug from tube containing TE with wide end of spatula and place in a sterile disposable petri dish or on the "plug cutter".
5. Cut a 2-mm-wide slice from test samples with a razor blade and transfer to the labeled 1.5-ml microcentrifuge tube containing 200 μl diluted H Buffer. Be sure plug slice is under buffer. The same spatula, razor blade or scalpel, and petri dish can be used for all samples. Replace rest of plug in original tube that contains TE Buffer.

**Notes:** The shape and size of the plug slice that is cut will depend on the size of the comb teeth that is used for casting the gel. PulseNet recommends that the combs with larger teeth (10-mm-wide teeth) be used to cast the gels because computer analysis of the gel lanes is more accurate and less tedious than analysis of gel lanes cast with combs with the smaller teeth (5.5-mm). The number of slices that can be cut from the plugs will depend on the skill of the operator, condition of the plug, and whether the slices are cut vertically or horizontally (5-mm x 10-mm plug).

6. Cut three 2-mm-wide slices of H9812 *Salmonella* ser. Braenderup standard plug and transfer to tubes of diluted H buffer. Be sure plug slices are under buffer.

7. Place tubes containing sample and control plug slices in float and incubate in a 37°C water bath for 10-15 min.

8. While plug slices are incubating in the buffer, dilute 10X H Buffer 1:10 with sterile Type 1 water according to the following table. Mix in the same Falcon 2057 tube that was used for H Buffer.

Reagent	μl/Plug Slice	μl/11 Plug Slices
Sterile Type 1 Water	175 μl	1925 μl
H Buffer	20 μl	220 μl
Enzyme (10 U/μl)	5 μl	55 μl
Total Volume	200 μl	2200 μl

9. After incubation of plug slices, remove H buffer by inserting pipet fitted with 200-250 μl tip all the way to bottom of tube and aspirating buffer. Be careful not to cut plug slice with pipet tip and that plug slice is not discarded with the tip.

10. Add<sup>2</sup> 55 μl of *Xba*I enzyme (Roche Molecular Biochemicals) to the diluted buffer; mix well.

**Note:** Keep vial of restriction enzyme on ice or in insulated storage box (-20°C) at all times. Wear gloves.

11. Add 200μl restriction enzyme mixture to each tube. Close tube and mix by tapping gently on bench top; confirm that plug slice is under enzyme mixture.

12. Place sample and control tubes in float and incubate in 37°C water bath for at least 2 hr.

---

<sup>2</sup> The instructor will dispense the *Xba*I enzyme

Laboratory Module IIB  
Casting Agarose Gel and Loading Restricted Plug Slices on the Comb

Tuesday morning (11:30am – 01:00pm), March 16, 2004

Read instructions in Steps 1-9 before continuing.

1. Place 10-well comb in 14-cm-wide gel form; use leveling platform to confirm gel is level. Confirm that front of comb holder and teeth face bottom of gel, and that the teeth of the comb touch the gel platform.
2. Remove restricted plug slices from 37°C water bath; allow to equilibrate to room temperature.
3. Remove enzyme/buffer mixture from plug slices with pipet and tip. Insert pipet fitted with 200-250 µl tip all the way to bottom of tube and aspirate buffer. Be careful not to cut plug slice with pipet tip and that plug slice is not discarded with the tip.
4. Add 200 µl 0.5X TBE to each plug slice.
5. Put comb on bench top or on top of gel form and load the plug slices on the bottom of the comb teeth; put the H9812 *Salmonella* ser. Braenderup standard plug slices on teeth 1, 5, and 10 and the sample plug slices on the remaining teeth.

**Note:** Load samples in order that is shown on the PFGE report form provided by instructors.

6. Remove excess buffer with edge of tissue; allow the plug slices to air dry for 3 minutes.
7. Position comb in gel form, confirm that plug slices are correctly aligned, and carefully pour 100 ml of melted 1% SKG agarose (55°C-60°C) into the gel form. Remove any bubbles that form with clean pipet tip. Allow gel to solidify for 30 minutes before removing the comb.
8. As soon as gel is poured, continue with **Laboratory Module IIC** below.

**Note:** Instructor will provide melted 1% SeaKem Gold (SKG) Agarose.

**Laboratory Module IIC**  
**Preparation of Electrophoresis Unit for PFGE of Restriction Digests**

Read instructions in Steps 1-4 before continuing.

1. Confirm that electrophoresis chamber is level; adjust the leveling screws on bottom of unit, if the chamber is not level. Put black gel frame in electrophoresis chamber; avoid touching the electrodes.
2. Add 2.2 L of 0.5X TBE; close cover of electrophoresis chamber.
3. Turn on power supply and pump; confirm that pump setting is -70 (the buffer flow at this setting is approximately 1 liter/min) and that buffer is circulating through the tubing.
4. Turn on cooling module and confirm that temperature setting is 14°C. (It usually takes approximately 20 minutes for the buffer to cool to 14°C).

**Laboratory Module IID**  
**Electrophoresis of Restriction Digests in PFGE Gel**

Read instructions in Steps 1-3 before continuing.

1. Unscrew and remove end gates from gel form; remove excess agarose from sides and bottom of casting platform with a tissue. Keep gel on black casting platform and carefully place gel inside casting frame in electrophoresis chamber. Close cover of chamber.
2. Select the following conditions on the CHEF Mapper electrophoresis unit:
  - Auto Algorithm
  - Low MW = 30 kb, High MW = 700 kb
  - Select default values except where noted by pressing "Enter"
  - Change run time to 16h.**
  - (Initial switch time = 2.16 s; final switch time = 63.8 s)
3. Record the milliamp (ma) value on PFGE Report form (normal values = 120-145 ma).

Laboratory Module III  
Staining and Documentation of PFGE Agarose Gel

Wednesday morning, March 17, 2004

Read instructions in this section before continuing.

1. Turn off chiller, pump, and CHEF Mapper (3 switches) when electrophoresis run is over.
- 2a. Remove gel and place it in covered plastic container that contains 40 $\mu$ l ethidium bromide in 400 ml Ultrapure water (Reagent or Type 1 water). The stock solution is 10 mg/ml; it is diluted 1:10,000 for staining.

**Notes:** Ethidium bromide is toxic and a mutagen. The diluted stain can be reused 3-5 times if it is kept in a dark bottle; discard used ethidium bromide according to the hazardous waste guidelines of your institution (See Section 10 of the PFGE Manual).

3. Place container with gel on rocker for 25-30 min.
4. Attach tubing to right port of electrophoresis chamber and drain buffer into large Erlenmeyer flask.
5. Rinse electrophoresis chamber with 1-2 liters Reagent Grade water; drain into large flask.
- 6a. Ethidium bromide - After gel has stained for 25-30 minutes, carefully pour into labeled bottle. WEAR GLOVES! Add 400 - 500 ml Reagent Grade water to gel and place on rocker to de-stain for 60-90 minutes; change water every 20-30 minutes, if possible.
7. Capture image on Gel Doc 2000 (\*.1sc file extension) according to the instruction sheet provided. Name the files according to the appropriate designation:

HKG04001-a  
HKG04001-b  
HKG04001-c  
HKG04001-d

**Note:** If background interferes with resolution, de-stain for additional 30-60 min.

8. Convert the \*.1sc file to a \*.tif file and save on diskette provided by instructor.

### Formulas for Reagents Used for PFGE<sup>3</sup>

**Cell Suspension Buffer (CSB; 100 mM Tris:100 mM EDTA, pH 8.0)**

- 10 ml of 1 M Tris, pH 8.0
- 20 ml of 0.5 M EDTA, pH 8.0
- Dilute to 100 ml with sterile Type 1 water

**Tris:EDTA Buffer (TE; 10 mM Tris:1 mM EDTA, pH 8.0)**

- 10 ml of 1 M Tris, pH 8.0
- 2 ml of 0.5 M EDTA, pH 8.0
- Dilute to 1000 ml with sterile Type 1 water

**1% SeaKem Gold:1% Sodium Dodecyl Sulfate Agarose (1% SKG:1% SDS) in TE Buffer**

- a. Weigh 0.50 (or 0.25) gram SKG into 250 ml screw-cap flask.
- b. Add 47.0 (or 23.5) ml TE Buffer; swirl gently to disperse agarose.
- c. Remove cap, cover loosely with clear film, and microwave for 30-sec; mix gently and repeat for 10-sec intervals until agarose is completely dissolved.
- d. Place flask in 56°C water bath for at least 5 minutes.
- e. Add 2.5 (or 1.25) ml of 20% SDS that has equilibrated to 56°C and mix well.
- f. Recap flask and return to 56°C water bath until ready to use.

**Cell Lysis Buffer (50 mM Tris:50 mM EDTA, pH 8.0 + 1% Sarcosine)**

- 25 ml of 1 M Tris, pH 8.0
- 50 ml of 0.5 M EDTA, pH 8.0
- 5 g of Sarcosyl (N-Lauryl-Sarcosine, Sodium salt)
- or**
- 50 ml of 10% Sarcosyl (N-Lauryl-Sarcosine, Sodium salt)
- Dilute to 500 ml with sterile Type 1 water

Add 25 µl Proteinase K stock solution (20 mg/ml) per 5 ml of cell lysis buffer just before use. The final concentration of Proteinase K is 0.1 mg/ml in the buffer.

**1% SeaKem Gold Agarose (1% SKG) in 0.5X TBE Buffer**

- a. Weigh 1.0 gram SKG into 500 ml screw-cap bottle.
- b. Add 100.0 ml 0.5 X TBE Buffer; swirl gently to disperse agarose.
- c. Remove cap, cover loosely with clear film, and microwave for 60-sec; mix gently and repeat for 15-sec intervals until agarose is completely dissolved.
- d. Place bottle in 54-58°C water bath for at least 15 minutes before pouring gel.

Use of trade names and commercial sources is for identification only and does not imply endorsement by the CDC or the U.S. Department of Health and Human Services.

---

<sup>3</sup> See Section 5a of the PulseNet PFGE Manual for detailed information.  
Laboratory Directions - PFGE