

Knowledge & Skill Development of Food Testing Laboratory Personnel from African Countries under IAFS– II



National Collateral Management Services Limited

Analysis of Metals in Food samples

- ❖ Obtain a representative sample from bulk
- ❖ Destruction of organic matter
- ❖ Separation and concentration and
- ❖ Determination

Analysis of Metals in Food samples

- ❖ Fresh foods: Macerating in a blender
- ❖ Dry products: Ground and powder is sieved
- ❖ Hard foods like chocolates:
Grating/chopping finely by hand
- ❖ Meat products: Thoroughly minced and ground in a mortar
- ❖ Fats: Melted before analysis



Analysis of Metals in Food samples

Destruction of Organic matter :

- ❖ Wet oxidation
- ❖ Dry Ashing and
- ❖ Microwave Digestion



Analysis of Metals in Food samples

❖ Wet oxidation method:

Samples are digested with Nitric acid and with one of the following on hot plate

- a) Sulfuric acid, (98%)
- b) Perchloric acid (70%)
- c) Hydrogen Peroxide (30%)



Analysis of Metals in Food samples

Dry Ash Method:

- Place the homogenized and ground sample in Muffle Furnace at 500° C
- Dissolve the ash in Hydrochloric acid (1+1)
- Dilute with water



Analysis of Metals in Food samples

Microwave digestion:

- 1g of finely ground sample treated with Nitric acid in a Teflon vessel
- Seal the vessels and digest as per the instructions
- Transfer the contents and dilute to the specified volume

7



STANDARD PREPARATION

Preparation of Arsenic standard solution

- 🧩 Take aliquot of stock solution (1000mg/l)
- 🧩 Add 5ml of HCl and 5ml of each 5% KI and L-Ascorbic acid solution
- 🧩 Allow the mixture to stand for 1 hour
- 🧩 Make up the volume with 10% HCl

8



STANDARD PREPARATION

Preparation of Selenium standard solution

- Take appropriate aliquot of stock solution (1000mg/l)
- Add equal volume of 3% Hydrochloric acid
- Heat the contents at 90°C for 20 min
- Cool and make up the volume with 3% HCl

STANDARD PREPARATION

Preparation of Mercury standard solution

- Take appropriate aliquot of stock solution (1000mg/l)
- Dilute to mark with 3% Hydrochloric acid
- Prepare different concentrations to construct calibration curve

STANDARD PREPARATION

Preparation of standard solutions of other elements (Multi element standard)


- 🧩 Acidify aliquot of stock solution with Nitric acid (HNO_3)
- 🧩 Make up the volume to obtain 2% HNO_3 in the final standard solution

Knowledge & Skill Development of Food Testing Laboratory Personnel from African Countries under IAFS– II



National Collateral Management Services Limited

Elemental Analysis in Food

 In order to cover all of the elements that should be analyzed at the appropriate ranges, several techniques are used

- ❖ As/Se/Sn–hydride generation AAS or VGA-ICP
- ❖ Hg–cold vapor technique or VGA-ICP
- ❖ All other elements at sub ppb level GFAAS or
- ❖ All other elements at ppb–ppm-% level ICP-OES

Instrumentation - AAS

- ❖ AAS is commonly used for metal analysis
- ❖ A solution of a metal compound is sprayed into a flame and vaporises to form atoms
- ❖ The metal atoms absorb light of a specific frequency, and the amount of light absorbed is a direct measure of the number of atoms of in the solution

3



Instrumentation - AAS


- ❖ Atomic Absorption Spectroscopy is a quantitative method of analysis that is applicable to many metals and few nonmetals.
- ❖ The technique was introduced in 1955 by Walsh in Australia (A. Walsh, Spectrochim. Acta, 1955, 7, 108)

4



Instrumentation - AAS

AAS



5




Instrumentation - AAS

- ❖ It is possible to measure the concentration of an absorbing species in a sample by applying the Beer-Lambert Law:

$$Abs = -\log\left(\frac{I}{I_0}\right)$$

ϵ = extinction coefficient




6




Hollow cathode lamp

Lamp



- ❖ Electron and ionic impact on cathode
- ❖ $M(s) \rightarrow M(g)$
- ❖ $M(g) \rightarrow \rightarrow \rightarrow M^*(g)$
- ❖ $M^*(g) \rightarrow M(g) + h\nu$



Thin lay of cathode material


7

Atomization – Flame mode

Two different types of flames are used

- ❖ Air-acetylene flame (2300°C)
- ❖ Usually 10-12cm long
- ❖ Nitrous oxide-acetylene flame (2700°C)
- ❖ 5cm long

8



Atomization – Flame mode


Reaction at Flame of burner:

- ❖ $M^{n+}(aq) + \text{anion}(aq) \rightarrow \text{salt}(s)$
- ❖ $\text{salt}(s) \rightarrow \text{salt}(g)$
- ❖ $\text{salt}(g) \rightarrow \text{atoms}(g)$
- ❖ $M(g) + h\nu \rightarrow M^*(g)$

9



Nebulizer-Burner



10



Nebulizer-Burner

- ❖ Nebulizer---to produce a mist or aerosol of the test solution
- ❖ Vaporizing chamber ---Fine mist is mixed with the fuel gas and the carrier gas
- ❖ Larger droplets of liquid fall out from the gas stream and discharged to waste
- ❖ Burner head ---The path length of flame is about 10 –12 cm

11



Graphite tube atomizer (GTA)




Figure 15.5 Graphite tube furnace

12



Graphite tube atomizer (GTA)

It involves mainly three steps:


- ❖ Drying or desolvation (105° C)
- ❖ Ashing or Pyrolysis (550° C)
- ❖ Atomization (depends on element)

13



Hydride generation method-For Arsenic, Selenium & Antimony


- ❖ Hydride formation of Arsenic



14




Monochromator



15



Photomultiplier tube (PMT)



16

